



Designation: E1327 – 07 (Reapproved 2018)

# Standard Test Method for Evaluation of Antimicrobial Handwash Formulations by Utilizing Fingernail Regions<sup>1</sup>

This standard is issued under the fixed designation E1327; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

## 1. Scope

1.1 This test method can be used to determine the effectiveness of antimicrobial handwashing agents (including handrubs) in the reduction of transient bacterial flora with particular emphasis on the fingernail region.

1.2 A knowledge of microbiological techniques is required for these procedures.

1.3 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety, health, and environmental practices and determine the applicability of regulatory limitations prior to use.* For more specific hazard statements, see 7.5.

1.4 *This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.*

## 2. Referenced Documents

2.1 *ASTM Standards:*<sup>2</sup>

**E1054** Test Methods for Evaluation of Inactivators of Antimicrobial Agents

**E2276** Test Method for Determining the Bacteria-Eliminating Effectiveness of Hygienic Handwash and Handrub Agents Using the Fingerpads of Adults

**E1838** Test Method for Determining the Virus-Eliminating Effectiveness of Hygienic Handwash and Handrub Agents Using the Fingerpads of Adults

<sup>1</sup> This test method is under the jurisdiction of ASTM Committee E35 on Pesticides, Antimicrobials, and Alternative Control Agents and is the direct responsibility of Subcommittee E35.15 on Antimicrobial Agents.

Current edition approved Sept. 1, 2018. Published September 2018. Originally approved in 1990. Last previous edition approved in 2012 as E1327 – 07(2012) <sup>$\epsilon$ 1</sup>. DOI: 10.1520/E1327-07R18.

<sup>2</sup> For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

## 3. Summary of Test Method

3.1 This test method, involving an improved method of recovering bacteria from hands, is used to study the effects of antimicrobial handwashes including health care personnel handwash products. The group of volunteer panelists need not refrain from using topical antimicrobials (such as deodorant soaps) before participating in the study. All subjects wash their hands with a nonantimicrobial hand soap prior to testing to remove any residual hand lotions and to lower the numbers of resident skin flora. Activity of products is measured by comparing the numbers of marker bacteria recovered from artificially contaminated fingernail regions after use of the handwashing formulations to the numbers recovered from the artificially contaminated but unwashed fingernail regions. Broth cultures of *Serratia marcescens* (a red pigmented bacterial species) and *Escherichia coli* (which produces fluorescent colonies on a special agar medium) are used as test bacteria. A spore suspension of *Bacillus subtilis* may be utilized to study (1) degree of physical removal by handwashing techniques, and (2) the recovery and precision aspects of the test method.

## 4. Significance and Use

4.1 The procedure should be used to test the degerming effectiveness of antimicrobial hand washing products used by health care personnel that are intended for frequent use, and that are intended to reduce the level of contamination acquired through contact with contaminated objects or people.

4.2 Performance of these procedures requires the knowledge of regulations pertaining to the protection of human subjects (Ref 1).<sup>3</sup>

## 5. Apparatus

5.1 *Colony Counter*—Any of several types may be used, for example, Quebec Colony Counter.

5.2 *Incubators*—One incubator capable of maintaining a temperature of  $25 \pm 2^\circ\text{C}$  (this temperature is required to ensure pigment production of *Serratia*); a second incubator

<sup>3</sup> The boldface numbers in parentheses refer to a list of references at the end of this standard.

capable of maintaining  $37 \pm 2^\circ\text{C}$  used for *E. coli* and *B. subtilis* incubation is acceptable.

5.3 *Water Bath*—Capable of maintaining temperature of  $80 \pm 2^\circ\text{C}$  for heat shocking of *B. subtilis* spores is needed.

5.4 *Sterilizer*—Any suitable steam sterilizer capable of producing the conditions of sterilization is acceptable.

5.5 *Timer*—Any stop-watch that can be read in minutes and seconds is required.

5.6 *Handwashing Sink*—A sink of sufficient size to permit panelists to wash without touching hands to sink surface or other panelists is needed.

5.6.1 *Water Faucet(s)*, to be located above the sink at a height that permits the hands to be held higher than the elbow during the washing procedure.

5.6.2 *Tap Water Temperature Regulator and Temperature Monitor*, to monitor and regulate water temperature of  $40 \pm 2^\circ\text{C}$ .

5.7 *Quad Petri plates*, 100 by 15 mm, plastic, sterile, disposable.<sup>4</sup>

5.8 *Small Petri Plates*, 60 by 15 mm, glass.

5.9 *Large Petri Plates*, 150 by 15 mm, glass.

5.10 *Tooth Brushes*:

5.10.1 *Young Size*.

5.10.2 *Battery Operated*.

5.11 *Ultraviolet Lamp*, having separate short wave and long wave bulbs.

5.12 *Germicidal Lamp Monitor Strips*.

5.13 *Inoculating Loops or Needles*, sterile.

5.14 *Plate Spreaders or Hockey Sticks*, sterile.

## 6. Reagents and Materials

6.1 *Bacteriological Pipettes*, 10.0 mL, sterile.

6.2 *Pipettors and Pipette Tips*, Eppendorf, MLA or similar types.

6.3 *Disposable Analyzer Cups*, 2 mL, plastic, not sterile.

6.4 *Sampling Solution*—Dissolve 0.4 g  $\text{KH}_2\text{PO}_4$ , 10.1 g  $\text{Na}_2\text{PO}_4$  and 1.0 g isooctylphenoxypolyethoxyethanol<sup>5</sup> in 1 L distilled water. Adjust pH to 7.8 with 0.1 N HCl or 0.1 N NaOH. Dispense in 100 mL-volumes and sterile for 20 min at  $121^\circ\text{C}$ .

6.5 *Dilution Fluid*—The sampling fluid may be used for dilutions or use Butterfields sterile phosphate buffered water (2) adjusted to pH 7.2 with suitable inactivator for the antimicrobial. Adjust pH with 0.1 N HCl or 0.1 N NaOH (see Practices E1054).

6.6 *Agar*, Tryptic soy agar or equivalent. Include the appropriate inactivator if needed.

6.7 *Agar with MUG*—Tryptic soy agar with 60 to 80  $\mu\text{g/mL}$  4-methylumbelliferyl- $\beta$ -D-glucuronide (MUG) is required.

6.8 *Test Formulations*—Directions for use of test formulation should be included if available. If these are not available, liquid antimicrobial soap formulations are tested by same routine as the nonantimicrobial control (10.5); alcoholic lotion type formulations are rubbed to dryness and then sampled for survivors (10.7).

6.9 *Nonantimicrobial Control Soap*, a liquid castile soap or other liquid soap containing no antimicrobials.

6.10 *Broth*—Tryptic soy broth or equivalent is required.

## 7. Test Organisms

7.1 *Serratia marcescens* American Type Culture Collection, ATCC No. 14756 is to be used as a marker organism. This is a strain having stable pigmentation. Grow in tryptic soy broth at  $25 \pm 2^\circ\text{C}$ .

7.2 *Escherichia coli*, ATCC No. 11229 is used as another Gram-negative marker organism. Grow in tryptic soy broth at  $35 \pm 2^\circ\text{C}$ .

7.3 *Bacillus subtilis*, ATCC No. 19659. Grow in tryptic soy broth at  $35 \pm 2^\circ\text{C}$ .

7.4 *Preparation of Spore Suspension*—Inoculate each surface of two tryptic soy agar plates (30 mL agar in 150-mm petri plates) with 1 mL of *B. subtilis* tryptic soy broth culture. Spread over the entire surface of the agar. Incubate for 5 to 10 days at  $35 \pm 2^\circ\text{C}$ . Suspend the growth in 20 mL of 0.1 % tryptone water<sup>6</sup> by rubbing the agar surface with a sterile rubber policeman. Add ethanol to the suspension to a final concentration of 80 % (wt/wt) and store in a refrigerator.

7.5 Other bacteria containing adequate markers to enable distinction from normal flora and of known safety may also be used for testing purposes. (**Warning**—The application of microorganisms to the skin may involve a health risk. Prior to applying *S. marcescens* or other bacteria to the skin, the antibiotic susceptibility profile of the strain should be determined. If the *Serratia* strain is not sensitive to Gentamicin, it should not be used. If an infection occurs, the antibiotic susceptibility profile should be made available to an attending clinician. Following the panelist's contamination and testing for the day, the panelist's hands should be decontaminated with a 70-% ethanol solution. Care should be taken to decontaminate around the fingernail regions.)

7.6 *Preparation of Marker Culture Suspension*—Inoculate a 10-mL tryptic soy broth tube with each of the test bacteria and incubate each tube at the temperature indicated to yield inocula of  $10^8$ – $10^9$  CFU/mL. When studying mixed inocula, mix equal volumes of the cultures into a sterile test tube; an equivalent

<sup>4</sup> Presterilized disposable quad plastic petri plates, the two sizes of glass petri plates and other equipment are available from most local laboratory supply houses.

<sup>5</sup> The sole source of supply of the apparatus (Triton X-100) known to the committee at this time is Rohm and Haas Co., Philadelphia, PA. If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee,<sup>1</sup> which you may attend.

<sup>6</sup> The sole source of supply of the Bacto Tryptone (Difco) water known to the committee at this time is Difco Laboratories, Detroit, MI. If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee,<sup>1</sup> which you may attend.

volume of *B. subtilis* spore suspension (that is prepared by centrifuging the alcoholic suspension and resuspending cells in water) may be added for bacterial physical removal determinations. Keep mixed suspension on ice during the day’s testing.

**8. Panelists**

8.1 Recruit a sufficient number of healthy adult human volunteers who have no clinical evidence of dermatoses, open wounds, hangnails, or other skin disorders. The number of people needed for a trial is dependent on the number of treatments within a study.

8.2 Volunteers are asked to maintain their normal use of soaps, shampoos, and so forth. They are asked to refrain from the use of acids, bases, solvents on the hands during the test period. Gloves should be provided for use where exposure to these agents is unavoidable.

**9. Experimental Design**

9.1 Each fingernail of a volunteer may be assayed separately; therefore, 10 test determinations (replicates) may be obtained from one volunteer. For the comparison of several products during a single study, a design such as a Latin Square Design may be utilized (3). For example, to compare 5 antimicrobial test products, one nonantimicrobial product and unwashed hand control (7 total variables), 7 volunteers, (or multiples of 7) should be recruited. Each person performs one testing of product or other variable on each of 7 test days, according to schedule such as the following; the numbers = day for testing that variable (see Table 1).

9.1.1 Example: Volunteer A tests Treatment 1 on Day 1, then Treatment 2 on Day 2; Volunteer B tests Treatment 2 on Day 1, Treatment 3 on Day 2, and so forth.

9.1.2 Each product or variable is tested once on each day, unless multiple numbers of volunteers are in the study.

9.1.3 The number of fingers, which are inoculated and then assayed after using the product, should be kept standard throughout. Although the number can be as high as 10, three fingers on one hand is a more convenient and cost savings approach. The ring, middle, and index fingers of the left hand have been selected for several studies; however, an operator may select the number and particular fingers to assay as long as they are held constant throughout.

**10. Procedure**

10.1 Before tests for the day, sterilize the analyzer cups by placing in suitable rack (24-well culture plates with lids are convenient) and placing the open cups under short-wave ultraviolet lamp for 15 to 30 min. To each sterile disposable

analyzer cup, add 0.9 mL of sterile diluent: set up sufficient cups only for each day’s testing.

10.2 Place 7 mL of sampling solution into each of 21 small petri plates.

10.3 Place 0.02 mL of marker culture suspension on the region surrounding the cuticle and under the fingernails of three fingers of the left hand of a volunteer. The volunteer then holds the hand in front of an electric fan for 5 min for complete drying of the suspension.

10.4 For unwashed hand determinations, proceed directly to 10.8.

10.5 When testing nonantimicrobial soap (controls), wet both hands under flowing warm tap water (40 ± 2 °C). Add 2.5 to 3.0 mL of the liquid soap to hands, rub hands together in normal washing manner for 15 s (no additional water), then rinse under the flowing water for 15 s to remove suds. Do not dry hands, proceed directly to 10.8.

10.6 For testing liquid antimicrobial soap formulations, follow the use directions on the label or follow the routine of 10.5. After washing, proceed to 10.8 without drying hands.

10.7 Alcoholic formulations are tested by placing the recommended volume on the hands and then rubbing the hands together until the alcohol has evaporated. Proceed to 10.8.

10.8 After performing the procedure for the day designated in the Latin Square Design, the technician scrubs with a toothbrush for 1 min each fingernail into a separate petri plate containing 7 mL of sampling solution.

NOTE 1—Although manual toothbrushes may be used for this purpose, greater uniformity between scrubblings may be obtained with less operator fatigue if an electric toothbrush such as the GE model TB-9 or another type is used. A brush which operates parallel with the handle is preferred because of less splashing.

10.9 After each scrubbing, the brushes are dropped into a beaker containing 70 % ethyl or isopropyl alcohol and allowed to stand for at least 10 min. The brushes are then rinsed in sterile distilled water and allowed to dry. The brushes are not sterilized.

10.10 Perform serial 10-fold dilutions. Place 0.1 mL amounts of the appropriate dilutions onto the surface of agar sections of quad plates. These drops of liquid are spread with sterile inoculation loops, needle spreaders, or hockey sticks to completely cover the quads. Allow drops to completely absorb.

10.11 Incubate inverted plates at 35 ± 2 °C for 12 to 18 h. Count the *E. coli* colony-forming units (CFU) that fluoresce under long-wave ultraviolet light. Transfer the plates to a 25 °C incubator and incubate for another day.

10.12 Count the red-pigmented *S. marcescens* CFU. Record the CFU per countable sections of the plates and convert values to the CFU obtained per finger by multiplying by the appropriate dilution factors.

10.13 Convert each CFU-per-finger determination to the log<sub>10</sub> value.

10.14 Determine the mean log<sub>10</sub> CFU per finger value. This is the mean log<sub>10</sub> value for that variable and subject for that day. These log<sub>10</sub> values are used for statistical comparisons.

**TABLE 1 Latin Square Design for Testing Seven Variables**

Day	Volunteer						
	A	B	C	D	E	F	G
1	T <sub>1</sub>	T <sub>2</sub>	T <sub>3</sub>	T <sub>4</sub>	T <sub>5</sub>	T <sub>6</sub>	T <sub>7</sub>
2	T <sub>2</sub>	T <sub>3</sub>	T <sub>4</sub>	T <sub>5</sub>	T <sub>6</sub>	T <sub>7</sub>	T <sub>1</sub>
3	T <sub>3</sub>	T <sub>4</sub>	T <sub>5</sub>	T <sub>6</sub>	T <sub>7</sub>	T <sub>1</sub>	T <sub>2</sub>
4	T <sub>4</sub>	T <sub>5</sub>	T <sub>6</sub>	T <sub>7</sub>	T <sub>1</sub>	T <sub>2</sub>	T <sub>3</sub>
5	T <sub>5</sub>	T <sub>6</sub>	T <sub>7</sub>	T <sub>1</sub>	T <sub>2</sub>	T <sub>3</sub>	T <sub>4</sub>
6	T <sub>6</sub>	T <sub>7</sub>	T <sub>1</sub>	T <sub>2</sub>	T <sub>3</sub>	T <sub>4</sub>	T <sub>5</sub>
7	T <sub>7</sub>	T <sub>1</sub>	T <sub>2</sub>	T <sub>3</sub>	T <sub>4</sub>	T <sub>5</sub>	T <sub>6</sub>