

**Designation:** E617 - 13 E617 - 18

# Standard Specification for Laboratory Weights and Precision Mass Standards<sup>1</sup>

This standard is issued under the fixed designation E617; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

#### 1. Scope

- 1.1 This specification covers weights and mass standards used in laboratories, specifically classes 000, 00, 0, 1, 2, 3, 4, 5, 6, and 7. This specification replaces National Bureau of Standards Circular 547, Section 1, which is out of print.
- 1.2 This specification and calibration method is intended for use by weight manufacturers, national metrology institutes, weight calibration laboratories, accreditation bodies, users of weights, and regulatory bodies.
  - 1.3 This specification contains the principal physical characteristics and metrological requirements for weights that are used.
  - 1.3.1 For the verification of weighing instruments;
  - 1.3.2 For the calibration of weights of a lower class of accuracy; and
  - 1.3.3 With weighing instruments.
- 1.4 Maximum Permissible Errors (formerly tolerances) and design restrictions for each class are described in order that both individual weights or sets of weights can be chosen for appropriate applications.
  - 1.4 The values stated in SI units are to be regarded as standard.
- 1.5 Weight manufacturers must be able to provide evidence that all new weights comply with specifications in this standard (e.g., (for example, material, density, magnetism, surface finish, mass values, uncertainties). Statements of compliance by calibration laboratories during subsequent calibrations must meet the requirements of ISO/IEC 17025, 5.10.4.2 and indicate on the calibration report which sections have or have not been assessed. Subsequent calibrations must meet all the requirements (including environmental parameters as shown in Table 11, of Sections 7, 8, and 9; and the requirements of ISO/IEC 17025:2005, 5.10.4.2 to make any claim of compliance to Specification E617, Maximum Permissible Errors, weight classes, or metrological traceability.
  - 1.6 The values stated in SI units are to be regarded as standard.
- 1.7 This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.

#### 2. Referenced Documents

2.1 ISO Standards:<sup>2</sup>

ISO/IEC 17025 General Requirements for the Competence of Testing and Calibration Laboratories (2005)

2.2 NIST Standards:<sup>3</sup>

NIST Handbook 143 State Weights and Measures Laboratories Program Handbook (2007)

NIST SP 811 Guide for the Use of the International System of Unit (SI) 2008 Edition

NIST SP 1038 The International System of Units (SI) – Conversion Factors for General Use (May 2006)

NISTIR 5672 Advanced Mass Calibration and Measurement Assurance Program for State Calibration Laboratories (2012)(2014)

NISTIR 6969 Selected Laboratory and Measurement Practices to Support Basic Mass Calibrations (2012)(2017)

NIST Technical Note 1297 (1994) Guidelines for Evaluating and Expressing the Uncertainty of NIST Measurement Results

<sup>&</sup>lt;sup>1</sup> This specification is under the jurisdiction of ASTM Committee E41 on Laboratory Apparatusand is the direct responsibility of Subcommittee E41.06 on Weighing Devices Laboratory Instruments and Equipment.

Current edition approved May 1, 2013 Oct. 1, 2018. Published July 2013 November 2018. Originally approved in 1978. Last previous edition approved in 2008 2013 as E617 – 97 (2008): E617 – 13. DOI: 10.1520/E0617-13. 10.1520/E0617-18.

<sup>&</sup>lt;sup>2</sup> Available from International Organization for Standardization (ISO), 1, ch. de la Voie-Creuse, CP 56, CH-1211 Geneva 20, Switzerland, http://www.iso.org.

<sup>&</sup>lt;sup>3</sup> Available from National Institute of Standards and Technology (NIST), 100 Bureau Dr., Stop 1070, Gaithersburg, MD 20899-1070, http://www.nist.gov.

2.3 OIML Standards:<sup>4</sup>

OIML D 28 Conventional Value of the Result of Weighing in Air (2004)

OIML R111-1e04 Weights of classes E1, E2, F1, F2, M1, M1-2, M2, M2-3 and M3 Part 1: Metrological and Technical Requirements (2004)

2.4 BIPM Standards:<sup>5</sup>

VIM: JCGM 200:2012 International Vocabulary of Metrology–Basic and General Concepts and Associated Terms GUM: JCGM 100:2008 Evaluation of Measurement Data–Guide to the Expression of Uncertainty in Measurement

2.5 EURAMET Standards:<sup>6</sup>

EURAMET/cg-18/V. 3.04.0 Guidelines on the Calibration of Non-Automatic Weighing Intruments (2011) Instruments (2015) 2.6 Additional Reference Documents:

CIPM-2007 Revised Formula for the Density of Moist Air, A. Picard, R. S. Davis, M. Glaser, and K. Fujii

## 3. Terminology

3.1 Definitions of Terms Specific to This Standard:

- 3.1.1 accuracy class of weights—a class of weights that meets certain metrological requirements intended to keep the errors within specified limits.
  - 3.1.2 balance—instument indicating apparent mass that is sensitive to the following forces:

$$F_g = m \cdot g$$

$$F_b = v \cdot \rho_a \cdot g = \frac{m}{\rho} \rho_a \cdot g$$

$$F_z = \mu_o \iiint_v (M + \chi H) \frac{\partial H}{\partial z} dV$$

Vertical component of the magnetic interaction between the weight and the balance or the environment, or both.

- H and M are vectors; z is the vertical cartesian coordinate. If magnetic effects are negligible, i.e. that is, the permanent magnetization (M) of the weight and the magnetic susceptibility  $(\chi)$  are sufficiently small, and the balance is calibrated with reference weights of well-known mass, the balance can be used to indicate the conventional mass,  $m_c$ , of a body under conventionally chosen conditions.
- 3.1.3 calibration (of weights)—the acts of determining the mass difference between a standard of known mass value and an "unknown" test weight or set of weights, establishing the mass value and conventional mass value of the "unknown," and of determining a quantitative estimate of the uncertainty to be assigned to the stated mass or conventional mass value of the "unknown," or both, and providing metrological traceability to the "unknown."
- 3.1.3.1 *calibration (generally)*—set of operations that establish, under specified conditions, the relationship between values of quantities indicated by a measuring instrument or measuring system, or values represented by a material measure or a reference material, and the corresponding values realized by standards.
  - 3.1.4 calibration certificate—certificate issued by calibration laboratories to document the results of a calibration.
- 3.1.5 conventional mass—conventional value of the result of weighing in air, in accordance to International Recommendation OIML D 28. For a weight taken at 20°C, the conventional mass is the mass of a reference weight of a density of 8000 kg/m³ which it balances in air of density of 1.2 kg/m³.
- 3.1.6 *correction*—mass values are traditionally expressed by two numbers, one being the nominal mass of the weight, and the second being a correction. The mass of the weight is the assigned nominal value plus the assigned correction. Positive corrections indicate that the weight embodies more mass than is indicated by the assigned nominal value. Negative corrections indicate that the weight embodies less mass than is indicated by the assigned nominal value. The correction is equivalent to the "error."
- 3.1.7 *international prototype kilogram*—the platinum-iridium cylinder maintained at the International Bureau of Weights and Measures (BIPM), at Sevres, France with an internationally accepted defined mass of 1 kg.
  - 3.1.8 *magnetism*—effect that generates an attractive or repulsive force.
- 3.1.8.1 (volume) magnetic susceptibility ( $\chi$ )—measure of the ability of a medium to modify a magnetic field. It is related to the magnetic permeability ( $\mu$ ) by the relation:  $\mu/\mu_0 = 1 + \chi$ . The quantity  $\mu/\mu_0$  is sometimes referred to as the relative permeability,  $\mu_T$ .
- 3.1.8.2 (permanent) magnetization (M)—parameter that specifies a magnetic state of material bodies such as weights, in the absenseabsence of an external magnetic field (most generally, magnetization is a vecotrvector whose magnitude and direction are

<sup>&</sup>lt;sup>4</sup> Available from Organisation Internationale de Metrologie Legale, 11 Rue Turgot, 75009 Paris, France, http://www.oiml.org.

<sup>&</sup>lt;sup>5</sup> Available from Bureau International des Poids et Mesures (BIPM), Pavilion de Breteuil, F-92312, Sèvres Cedex, France, http://www.bipm.org.

<sup>&</sup>lt;sup>6</sup> Available from Euramet, Bundesallee 100, 38116 Braunschweig, Germany, http://www.euramet.org.



not necessarily constant within the material). The magnetization of a body generates an inhomogeneous magnetic field in space and thus may produce magnetic forces on other materials.

- 3.1.9 *mass*—physical quantity, which can be ascribed to any material object and which gives a measure of its quantity of matter. The unit of mass is the kilogram.
- 3.1.10 maximum permissible errors—the maximum amount by which the sum of the conventional mass of the weight, its deviation from nominal value, and its associated uncertainty is allowed to deviate from the assigned nominal value.
- 3.1.11 metrological traceability—property of a measurement result whereby the result can be related to a reference through a documented unbroken chain of calibrations, each contributing to the measurement uncertainty. Metrological traceability requires an established calibration hierarchy. Elements for confirming metrological traceability to be an unbroken chain to an international measurement standard or a national measurement standard (IPK or NPS), shall include a documented measurement uncertainty, a documented measurement procedure, accredited technical competence, metrological traceability to the SI, and established calibration intervals (see current VIM: JCGM 200).
- 3.1.12 *reference standard*—a standard, generally of the highest metrological quality available at a given location, from which measurements made at that location are derived.
- 3.1.13 roughness parameter or R-parameter ( $R_a$  or  $R_z$ )—parameter that describes the assessed roughness profile of a sample. The letter R is indicative of the type of assessed profile, in this case R for roughness profile. The assessed profile of a sample can be in terms of different profile types: a roughness profile or R-parameter, primary profile or R-parameter, a waviness profile or R-parameter.
- 3.1.14 set of weights—a series of weights, usually presented in a case so arranged to make possible any weighing of all loads between the mass of the weight with the smallest nominal value and the sum of the masses of all weights of the series with a progression in which the mass of the smallest nominal value weight constitutes the smallest step of the series.
- 3.1.15 *temperature* (t)—in degrees Celsius, is related to the absolute thermodynamic temperature scale, called the Kelvin scale, by t = T 273.15 K.
  - 3.1.16 test weight  $(m_t)$ —weight that is to be tested according to this standard.
- 3.1.17 *tolerance test*—verification that the conventional mass of the weights and their corresponding uncertainties as tested are correct within the maximum permissible errors of the respective weight class.
- 3.1.18 *uncertainty*—non-negative parameter characterizing the dispersion of the quantity values being attributed to a measurand, based on the information used.
- 3.1.19 *units*—the units used are: (1) for mass, the milligram (mg), the gram (g) and the kilogram (kg); (2) for density, the kilogram per cubic meter (kg  $m^{-3}$ ).
- 3.1.20 *U.S. National prototype standard*—platinum-iridium kilogram identified as K20, maintained at the National Institute of Standards and Technology, with value assigned relative to the International Prototype Kilogram provides the United States access to the mass unit.
- 3.1.21 weight—material measure of mass, regulated in regard to its physical and metrological characteristics: shape, dimensions, material, surface quality, nominal value, density, magnetic properties and maximum permissible error.

Note 1—The term "weight" is also used as the physical quantity of the gravitational force of a body. From the context it is usually clear in which sense the term is used. If the sense is not clear, one may use the words "weight force" or "weight piece," depending on its meaning.

#### 3.2 Symbols:

2,		
Symbol	Unit	Definition
Α	-	represents weighing the reference
В	_	weight in a weighing cycle represents weighing the test weight in
D		a weighing cycle
C	_	correction factor for air buoyancy
D	kg	difference of balance readings
		between minimum and maximum values from eccentricity test
d	kg	scale interval
$d_1$	m	estimated distance between centers of
		weights during loading
$d_2$	m	estimated distance from the center of
$F_{b}$	N	the load receptor to one of the corners air buoyancy equal to the weight of the
ГЬ	IN	displaced air
$F_g$	N	gravitational force
$F_z^{\stackrel{\circ}{J}}$	N	magnetic force between a mass
		comparator and a weight in the vertical
	•	or z-direction
g	m s <sup>-2</sup>	gravitational acceleration

## €617 – 18

Symbol	Unit	Definition
H	A m <sup>-1</sup>	magnetizing field strength
hr	%	relative humidity
I	kg	indication of the weighing instruments
A /	les.	(scale division)
$\Delta I$	kg	indication difference of the balance, where $\Delta I = I_r - I_r$
$\Delta I_1$	ka	indication difference using an
$\Delta I_1$	kg	automatic exchange mechanism with
A /	len.	weights in first position
$\Delta l_2$	kg	indication difference using an
		automatic exchange mechanism with
A /	len.	weights in reversed position
$\Delta I_{\mathtt{S}}$	kg	change in indication of balance due to
		sensitivity weight
i	-	subscript used as an index in
•		summations
J	<del>-</del>	subscript for number of test weights or
k		number of series of measurements
M	_ A m <sup>−1</sup>	coverage factor, typically 2 or 3
IVI	AIII	permanent magnetization (see also
	les.	μ <sub>0</sub> M)
m A m	kg	mass of a rigid body (weight)
$\Delta m$	kg	mass difference, usually between test and reference weight
8 m	ka	· · · · · · · · · · · · · · · · · · ·
δm	kg	maximum permissible error on the weights
m	ka	mass, nominal value of the weight
$m_{\rm O}$	kg	(e.g. 1 kg)
m	ka	conventional mass of the weight
$m_{ m c}$ $\Delta m_{ m c}$	kg kg	conventional mass difference between
ΔIII <sub>C</sub>	ng	test weight and reference weight
A <del></del>	ka	average conventional mass difference
$\Delta m_c$	iTab Ctandon	between test weight and reference
	iTeh Standar	weight
$m_{\rm cr}$	kg	conventional mass of the reference
· · · cr		
$m_{\rm ct}$	(https://standards.	conventional mass of the test weight
m <sub>s</sub>	kg	mass of the sensitivity weight
$m_{\rm t}$		
n	Document Prev	subscript for number of measurement
		sequences
p	Pa	barometric pressure
R <sub>a</sub>	μm	mean height of roughness profile (R-
4	ASTM E617-18	parameter)
$R_{z}$	μm	maximum height of roughness profile
https://standards.iteh.a	ai/catalog/standards/sist/a1a870f8-33e8-48	8/a-81/5-1 (R-parameter) 266/astm-e61/-18
r	<del>-</del>	subscript for reference weight
S	_	subscript for sensitivity weight
s	kg	standard deviation
$s^2$	kg <sup>2</sup>	variance
T	K	thermodynamic temperature using the
		International Temperature Scale of
		1990 (ITS-90)
ΔT*	°C	initial difference between weight tem-
		perature and laboratory temperature
t	-	subscript for test weight
t	°C	temperature in degrees Celsius, where
11	lon.	t = T- 273.15 K
U	kg	uncertainty, expanded uncertainty
u	kg	uncertainty, standard uncertainty
$u_{b}$	kg	uncertainty of air buoyancy correction
u <sub>ba</sub>	kg	uncertainty of the balance
$u_{\rm c}$	kg	combined standard uncertainty
$u_{d}$	kg	uncertainty due to the display resolu-
	ka	tion of a digital balance
u <sub>E</sub>	kg kg m <sup>-3</sup>	uncertainty due to eccentricity
$u_F$	kg m	uncertainty of the formula used to cal-
11	%	culate air density
u <sub>hr</sub>		uncertainty in relative humidity
$u_{inst}$	kg	uncertainty due to instability of the ref-
"	ka	erence weight uncertainty due to magnetism
U <sub>ma</sub>	kg Pa	uncertainty due to magnetism uncertainty in barometric pressure
$u_{p}$	kg	uncertainty in barometric pressure uncertainty due to the sensitivity of the
$u_{\rm s}$	'vy	balance
$U_{t}$	°C	uncertainty in temperature
- <b>t</b>	-	



$\begin{array}{cccccccccccccccccccccccccccccccccccc$	Symbol	Unit	Definition
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$u_{w}$	kg	uncertainty due to the weighing
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	V	m <sup>3</sup>	•
$\begin{array}{c} \mu \\ \mu_0 \\ \mu_0 \\ \mu_0 \\ \end{array} \qquad \begin{array}{c} N \ A^{-2} \\ N \ A^{-2} \\ \end{array} \qquad \begin{array}{c} \text{magnetic permeability} \\ \text{magnetic constant (magnetic permeability of vacuum)}, } \mu_0 = 4\pi \times \\ 10^{-7} \ N \ A^{-2} \\ \end{array} \qquad \begin{array}{c} \mu_0 \\ M \\ \mu_r \\ V_{eff} \\ \rho \\ \end{array} \qquad \begin{array}{c} - \\ \text{relative magnetic permeability (} \mu/\mu_0 ) \\ V_{eff} \\ \rho \\ \end{array} \qquad \begin{array}{c} - \\ \text{effective degrees of freedom} \\ \text{mass of a rigid body (weight)} \\ \text{density of air as a reference value} \\ \text{equal to 1.2 kg m}^{-3} \\ \text{density of moist air} \\ \text{pal} \\ \text{pal} \\ \end{array} \qquad \begin{array}{c} kg \ m^{-3} \\ \text{density of moist air during the last} \\ \text{(previous) calibration of the reference} \\ \text{weight} \\ \end{array} \qquad \begin{array}{c} \rho_r \\ \text{density of a reference weight with} \\ \text{mass } m_r \\ \text{density of the weight being tested} \end{array}$	•		
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	μ	$NA^{-2}$	magnetic permeability
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\mu_0$	N A <sup>-2</sup>	magnetic constant (magnetic
$\begin{array}{cccccccccccccccccccccccccccccccccccc$			
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\mu_{o}M$	Т	magnetic polarization
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\mu_{r}$	_	
ρ <sub>0</sub> kg m <sup>-3</sup> density of air as a reference value equal to 1.2 kg m <sup>-3</sup> ρ <sub>a</sub> kg m <sup>-3</sup> density of moist air density of alibration of the reference weight (previous) calibration of the reference weight density of a reference weight with mass m <sub>r</sub> ρ <sub>t</sub> kg m <sup>-3</sup> density of a reference weight being tested	V <sub>eff</sub>	_	•
ρ <sub>a</sub> kg m <sup>-3</sup> density of moist air       ρ <sub>al</sub> kg m <sup>-3</sup> density of moist air during the last (previous) calibration of the reference weight       ρ <sub>r</sub> kg m <sup>-3</sup> density of a reference weight with mass m <sub>r</sub> ρ <sub>t</sub> kg m <sup>-3</sup> density of the weight being tested	ρ		
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	Ρο	kg m <sup>-3</sup>	
$\begin{array}{cccccccccccccccccccccccccccccccccccc$			
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	Pa	kg m <sup>−3</sup>	•
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\rho_{al}$	kg m <sup>-3</sup>	,
$ ho_{ m t}$ mass $m_r$ density of the weight being tested			" ,
$\rho_{t}$ kg m <sup>-3</sup> density of the weight being tested	$\rho_{r}$	kg m <sup>−3</sup>	density of a reference weight with
			mass $m_r$
χ – magnetic susceptibility	$\rho_{t}$	kg m <sup>−3</sup>	, , ,
	χ	-	magnetic susceptibility

#### 4. Maximum Permissible Errors

- 4.1 For each weight, the expanded uncertainty *U* at approximately 95 % confidence (See Section 9) of the conventional mass shall be less than or equal to one-third of the maximum permissible error given in Table 1 as defined in Section 9. Subsequent calibrations must meet all the requirements (including environmental parameters as shown in Table 11, of Sections 7, 8, and 9; and the requirements of ISO/IEC 17025:2005, 5.10.4.2 to make any claim of compliance to Specification E617, Maximum Permissible Errors, weight classes, or metrological traceability.
- 4.1.1 For each weight, the conventional mass,  $m_c$  (determined with an expanded uncertainty), shall not differ by more than the difference: maximum permissible error  $\delta m$  minus expanded uncertainty, from the nominal value of the weight,  $m_c$ :

$$m_o - (\delta m - U) \le (m_c) \le m_o + (\delta m - U) \tag{1}$$

- 4.2 Maximum permissible errors for classes 000, 00, 0, 1, 2, 3, 4, 5, 6, and 7 are given in Table 1. These maximum permissible errors apply to conventional mass values.
- Note 2—Maximum Permissible Errors for weights of denomination intermediate between those listed, the maximum permissible error shall be proportional to the values shown.
- 4.3 Maximum Permissible Errors for weights of denomination intermediate between those listed, the maximum permissible error shall be proportional to the values shown.
- 4.4 For class 000, 00, and 0 weights, which are always accompanied by certificates giving the mass values and uncertainties, the deviation from the nominal value,  $m_c m_0$ , shall be taken into account by the user.

#### 5. Physical Characteristics

- 5.1 Construction:
- 5.1.1 *Type*—Weights are divided into two types based upon the design:
- 5.1.1.1 *Type I*—These weights are of one-piece construction and contain no added adjusting material. They must be specified when weights are to be used as standards for the calibration of weights of Classes 000, 00, 0, 1, 2, and 3, and where maximum stability is required. A precise measurement of density can only be made for one-piece weights.
- 5.1.1.2 *Type II*—Weights of this type can be of any appropriate design such as screw knob, ring, or sealed plug. Adjusting material can be used as long as provided it is of a material at least as stable as the base material and is contained in such a way that it will not become separated from the weight.
- 5.1.2 Class 000, 00, and 0 shall be Type I, one piece construction. Weights with nominal values less than 1 g shall have unique shapes to differentiate the weights from one another. See Table 2. The shape of weights smaller than 1 mg shall be discussed and verified with the customer.
  - 5.1.3 Class 1, 2, 3, 4, 5, 6, and 7 may be either Type I or Type II depending on the application.
- 5.2 Design—A weight may have any shape that does not introduce features that reduce the reliability. All weights shall be free of ragged or sharp edges or ends. Both sheet metal and wire weights shall be free of cracks such as may be formed from bending.
- 5.3 Surface Area—For classes 000, 00, 0, 1, 2, 3, and 4 the surface area is not to exceed twice the area of a cylinder of equal height and diameter for weights 1 g and above. Sheet metal weights or wire weights may be used below 1 g. For Classes 5, 6, and 7 the total surface areas should be minimized to the extent possible.
  - 5.4 Material:

## **TABLE 1 Maximum Permissible Errors**

Note 1—Maximum Permissible Errors are reported in SI units, typically milligrams.

Note 2—The "grain" is the same in avoirdupois, troy, and apothecaries units of mass.

Note 3—See NIST SP 811 and NIST SP 1038 for conversion and units of measure.

Metric	Class 000	Class 00	Class 0	Class 1	Class 2	ept as noted Class 3	Class 4	Class 5	Class 6	Class
	Class 000	Ciass 00	Uiass 0	Glass I	Class 2	Ulass 3	Class 4	Ciass 5	Class b	Ciass
Metric 5000 kg					25 g	50 g	100 g	250 g	500 g	750
3000 kg					25 g 15 g	30 g	60 g	250 g 150 g	300 g	450
2000 kg					10 g	20 g	40 g	100 g	200 g	300
—1000 kg					<del>5 g</del>	10 g	<del>20 g</del>	<del>50 g</del>	<del>100 g</del>	150
1000 kg					5.0 g	10 g	20 g	50 g	100 g	150
- 500 kg					2.5 g	5 g	10 g	25 g	<del>50 g</del>	<del>75 ç</del>
500 kg					2.5 g	5.0 g	10 g	25 g	50 g	75 0
300 kg					1.5 g	3 g	6.0 g	15 g	<del>30 g</del>	75 g 45 g
300 kg					1.5 g	3.0 g	6.0 g	15 g	30 g	45 g
<del>200 kg</del>					<del>1 g</del>	2 g	4.0 g	10 g	<del>20 g</del>	<del>30 (</del>
200 kg					1.0 g	2.0 g	4.0 g	10 g	20 g	30 (
-100 kg					<del>500 mg</del>	<del>1 g</del>	2.0 g	<del>5 g</del>	<del>10 g</del>	<del>15 (</del>
100 kg					500 mg	1.0 g	2.0 g	5.0 g	<u>10 g</u>	15 (
<del>- 50 kg</del>	<del>13 mg</del>	<del>25 mg</del>	<del>63 mg</del>	<del>125 mg</del>	<del>250</del>	<del>500 mg</del>	<del>1.0 g</del>	<del>2.5 g</del>	<del>5 g</del>	<del>7.5 (</del>
50 kg	<u>13 mg</u>	25 mg	63 mg	<u>120 mg</u>	250	500 mg	<u>1.0 g</u>	<u>2.5 g</u>	<u>5.0 g</u>	7.5 <u>(</u>
<del>30 kg</del>	7.5	<del>15</del>	38	<del>75</del>	<del>150</del>	<del>300</del>	<del>600 mg</del>	<del>1.5 g</del>	<del>3 g</del>	4.5
30 kg	<u>7.5</u>	<u>15</u>	<u>38</u>	75 <del>62</del>	<u>150</u>	300	600 mg	1.5 g	<u>3.0 g</u>	4.5
<del>- 25 kg</del>	<del>6.25</del>	<del>12.5</del>	<del>31</del>	<del>62</del>	<del>125</del>	<del>250</del>	<del>500</del>	<del>1.2 g</del>	<del>2.5 g</del>	4.5
25 kg	6.25	12.5	31	<u>62</u> <del>50</del>	125	250 200	500	1.2 g	<u>2.5 g</u>	3.8
<del>- 20 kg</del>	5.0	<del>10</del>	<del>25</del>	<del>50</del>	<del>100</del>	<del>200</del>	400	<del>1.0 g</del>	<del>2 g</del>	3.8
20 kg	<u>5.0</u>	<u>10</u>	<u>25</u> <del>13</del>	<u>50</u>	100	200	400	<u>1.0 g</u>	<u>2.0 g</u>	3.8
<del>- 10 kg</del>	<del>2.5</del>	<del>5.0</del>		<del>25</del>	<del>50</del>	100	200	<del>500 mg</del>	<del>1 g</del>	2.2
10 kg	<u>2.5</u> 1.3	5.0	13	25 12	50 25	100	200	500 mg	1.0 g	2.2
5 kg		2.5	6.0	12	25	50	100	250	500 mg	1.4
3 kg	0.75	1.5	3.8	7.5	15	30	60	150	300	1.0
2 kg	0.5	1.0	2.5	5.0	10	20	40	100	200	750 r
1 kg	0.25	0.5	1.3	2.5	5.0	10	20	50	100	470
500 g	0.13	0.25	0.60	1.2	2.5	5.0	10	30	50	300
300 g	0.075	0.15	0.38	0.75	1.5	3.0	6.0	20	30	210
200 g	0.05	0.10 <del>0.05</del>	0.25	0.50 0.25	1.0 <del>0.50</del>	re 2.0	4.0	15	20 <del>10</del>	160
— <del>100 g</del>	0.025		0.13				2.0	<del>9</del> 9.0		100
<u>100 g</u> <del>50 g</del>	0.025 0.015	0.05 0.030	0.13 0.060	0.25 <del>0.12</del>	0.50 <del>0.25</del>	<u>1.0</u> <del>0.60</del>	2.0 <del>1.2</del>	9.0 5.6	10 7	100 <del>62</del>
50 g	0.015	0.030	0.060	0.12		0.60	1.2	5.6	7.0	62
30 g	0.013	0.030 0.026	0.000	0.074	$\sqrt{\frac{0.25}{0.15}}$	18 <del>0.00</del> <del>0.45</del>	<del>1.2</del> <del>0.90</del>	4.0	7.0 5	<u>62</u> 44
30 g	0.014	0.026	0.037	0.074	0.15	0.45	0.90	4.0 1204 <u>4.0</u> 66	/astr <u>5.0</u> 617	1 0 44
20 g	1000000000000000000000000000000000000	$\frac{0.025}{0.025}$	0.037	$\frac{0.074}{0.074}$ 1 a	$870\frac{0.15}{0.10}33$	0.35 a-2	$\frac{0.70}{0.70}$	1aUd. 3.0	/astm <del>3</del> 61/	$7-18 \frac{44}{33}$
20 g	0.013	0.025	0.037	0.074	0.10	0.35	0.70	3.0		33
<del>10 g</del>	0.010	0.020	0.025	0.050	0.074	0.25	0.50	2.0	3.0 2	33 21
10 g	0.010	0.020	0.025	0.050	0.074	0.25	0.50	2.0	2.0	<u>21</u> <del>13</del>
<del>-5 g</del>	0.005	0.010	0.017	0.034	0.054	0.18	0.36	<del>1.3</del>	2	13
_ 5 g	0.005	0.010	0.017	0.034	0.054	0.18	0.36	1.3	2.0	13
3 g	0.005	0.010	0.017	0.034	0.054	0.15	0.30	0.95	<u>2.0</u> 2.0	<u>13</u> 9.4
2 g	0.005	0.010	0.017	0.034	0.054	0.13	0.26	0.75	2.0	7.0
1 g	0.005	0.010	0.017	0.034	0.054	0.10	0.20	0.50	2.0	4.5
500 mg	0.002	0.003	0.005	0.010	0.025	0.080	0.16	0.38	1.0	3.0
300 mg	0.002	0.003	0.005	0.010	0.025	0.070	0.14	0.30	1.0	2.2
200 mg	0.002	0.003	0.005	0.010	0.025	0.060	0.12	0.26	1.0	1.8
100 mg	0.002	0.003	0.005	0.010	0.025	0.050	0.10	0.20	1.0	1.2
<del>- 50 mg</del>	0.002	0.003	0.005	<del>0.010</del>	0.014	0.042	0.085	<del>0.16</del>	0.50	0.80
50 mg	0.002	0.003	0.005	0.010	0.014	0.042	0.085	<u>0.16</u>		0.88
<del>- 30 mg</del>	0.002	0.003	0.005	0.010	0.014	0.038	0.075	<del>0.14</del>	0.50	0.60
30 mg	0.002	0.003	0.005	0.010	0.014	0.038	0.075	<u>0.14</u>		0.6
<del>- 20 mg</del>	0.002	0.003	0.005	0.010	0.014	0.035	0.070	<del>0.12</del>	0.50	0.50
20 mg	0.002	0.003	0.005	0.010	0.014	0.035	0.070	0.12		0.50
<del>- 10 mg</del>	0.002	0.003	0.005	0.010	0.014	0.030	0.060	0.10	0.50	0.40
10 mg	0.002	0.003	0.005	0.010	0.014	0.030	0.060	0.10	6.55	0.40
<del>- 5 mg</del>	0.002	0.003	0.005	0.010	0.014	0.028	0.055	0.080	0.20	
<u>5 mg</u>	0.002	0.003	0.005	0.010	0.014	0.028	0.055	0.080	0.00	
<del>- 3 mg</del>	0.002	0.003	0.005	0.010	0.014	0.026	0.052	0.070	0.20	
3 mg	0.002	0.003	0.005	0.010	0.014	0.026	0.052	0.070	0.00	
<del>2 mg</del>	0.002	0.003	0.005	0.010	0.014	0.025	0.050	0.060	0.20	
2 mg	0.002	0.003	0.005	0.010	0.014	0.025	0.050	0.060	0.40	
— <u>1 mg</u>	0.002	0.003	0.005	0.010	0.014	0.025	0.050	0.050	0.10	
1 mg	0.002	0.003	0.005	0.010	0.014	0.025	0.050	0.050	0.40	
<del>0.5 mg</del> 0.5 mg	<del>0.002</del> 0.002	0.003 0.003	<del>0.005</del> 0.005	<del>0.010</del> 0.010	0.014 0.014	0.025	0.050	0.050	0.10	
			0.005	0.010	0.014	0.025	0.050	0.050		



Medical					.,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,						
D. STING   D. COCC   D. COCK   D. COLID   D. D. H.   D. COCK   D. D. H.   D. COLID   D. D. H.	Denomination					±mg exce	pt as noted				
O	Metric	Class 000	Class 00	Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class 7
0.15mg 0.002 0.003 0.005 0.010  Acordapois	0.3 mg	0.002	0.003	0.005	0.010	0.014	0.025				
O.O.S.   O		0.002	0.003	0.005	0.010	0.014					
Automorphis   Class 0   Class 1   Class 2   Class 3   Class 4   Class 7   Class 7   Class 7   Pauret   mg   mg   mg   mg   mg   mg   mg   m	0.1 mg				0.010						
Power   mg		0.002	0.003								
Pound   mg   mg   mg   mg   mg   mg   mg   m	•			Class 0							
10000 ib											
10000 lb				mg	mg						
Semilar   14-9   22-9   24-9   24-9   24-9   22-9   24-9   24-9   22-9   24-9   22-9   23-9   34-9   20-9   30-9						<del>23 g</del>				•	-
SOOO   B						11_0			230 <u>y</u>		
3000   b						11 9					
3000 lb						<del>7 a</del>					
2600-lb	3000 lb					, 9	14 a	27 g			
2800 lb   46-9   9-9   48-9						<del>6 g</del>	12 g	24 g	<del>57 g</del>		170 g
2000-bb	2500 lb					· ·	-				
1900-bb	2000 lb					<del>4.5 g</del>		18 g			
1900-lb 1900-l								<u>18 g</u>			
500-lb						<del>2.3 g</del>		<del>9 g</del>			<del>68 g</del>
Solid   Soli							<u>4.5 g</u>		<u>23 g</u>		<u>68 g</u>
300 lb 200 lb 57-mg 110-mg 230-mg 1.8 g 4.5 g 9.1 g 14 g 100 lb 57-mg 110-mg 230-mg 1.8 g 4.5 g 9.1 g 14 g 100 lb 57-mg 110-mg 230-mg 1.8 g 4.5 g 9.1 g 14 g 100 lb 57-mg 110-230 450 910-mg 1.8 g 4.5 g 9.1 g 14 g 58-hb 58-hb 58-hb 59-hb 28-hb 14-hb 14-h						<del>1.7 g</del>					
200   b	500 lb						2.3 g		<u>11 g</u>	<u>23 g</u>	<u>34 g</u>
1600-lb   57-mg   140-mg   230-mg   460-mg   920-mg   2-3-g   4.5-g   6-8-g   60-lb   57   110   230   450   910-mg   2-3-g   4.5-g   6-8-g   60-lb   29   67   140   220   440   1.1-g   2.3-g   4.1-g   90-lb   1.1-g   2.1-g   2.1-								2.7 g			
100 lb				E7 ma	110 mg	220 mg					
Solib   29   57   110   220   440   11-g   23-g   41-g   30-b   30-b   11-g   23-g   41-g   30-b   17   32   68   140   260   680-mg   11-g   23-g   41-g   30-b   17   32   68   140   260   680-mg   11-g   23-g   41-g   27-g   25-b   144   28   56   140   270   680-mg   11-g   27-g   25-b   144   28   56   140   220   57-0   11-g   24-g   25-b   25-b   25-b   24   28   46   92   180   460   910-mg   2-g   20-b   22-b   23   45   91   180   460   910-mg   2-g   20-b   22-b   23   45   91   180   460   910-mg   2-g   20-b   25-b   23   45   91   180   450   910-mg   2-g   20-b   25-b					-	•					
Solb   29   57   110   230   450   1.1 g   2.3 g   4.1 g   30 b   4.1 g   30 b   4.7   34   68   446   266   680 mg   1.4 g   2.7 g   52 b   52 b   14   28   56   446   260   570   1.1 g   2.7 g   2.5 b   2.5 b   14   28   56   446   260   570   1.1 g   2.7 g   2.5 b   2.5 b   14   28   57   110   230   570   1.1 g   2.4 g   2.5 b   2.5 b   2.3 b   46   92   480   480   480   910 mg   2.9 g   20 b   20 b   12   23   45   91   180   450   910 mg   2.0 g   10 b   10 b   5.5   11   23   45   91   180   450   910 mg   2.0 g   10 b   10 b   5.5   11   23   45   91   180   450   910 mg   2.0 g   10 b   10 b   5.5   11   23   45   91   230   450   1.3 g   1.3 g   1.3 g   1.3 b   1.7   5.4   11   22   43   44   88   230   450   1.3 g   1.3 g   1.3 b   1.7   5.4   6.8   14   27   68   140   230   780 mg   5 b   2.7   5.4   11   23   45   45   110   230   780 mg   5 b   2.7   5.4   11   23   45   31   10   230   780 mg   2 b				20	<del>110</del>	230 110	<del>220</del>			2.3 g	
30 lb				29							-
30 lb				=== <del>17</del>	<del>34</del>		140				
25 lb 14 28 57 110 230 570 1.1 g 2.4 g 20-lb 12 23 45 91 180 450 910 mg 2.0 g 140-lb 182 23 45 91 180 450 910 mg 2.0 g 140-lb 185 11 23 45 91 180 450 1910 mg 2.0 g 140-lb 185 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 6 lb 2.7 5.4 11 23 45 91 10 230 780 mg 6 lb 2.7 6 lb 1.7 3.4 6.8 14 27 68 140 580 92 lb 1.2 2.3 lb 1.7 3.4 6.8 14 27 68 140 580 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.5 9.5 1 18 9.5 1 18 45 91 440 92 lb 1.5 9.5 1 18 9.5				17	32				•		2.7 g
25 lb 14 28 57 110 230 570 1.1 g 2.4 g 20-lb 12 23 45 91 180 450 910 mg 2.0 g 140-lb 182 23 45 91 180 450 910 mg 2.0 g 140-lb 185 11 23 45 91 180 450 1910 mg 2.0 g 140-lb 185 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 5.4 11 23 45 91 230 450 1.3 g 760-mg 5 lb 2.7 6 lb 2.7 5.4 11 23 45 91 10 230 780 mg 6 lb 2.7 6 lb 1.7 3.4 6.8 14 27 68 140 580 92 lb 1.2 2.3 lb 1.7 3.4 6.8 14 27 68 140 580 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.2 2.3 4.5 91 18 45 91 440 450 92 lb 1.5 9.5 1 18 9.5 1 18 45 91 440 92 lb 1.5 9.5 1 18 9.5				14	28	4 56		<del>220</del>			2.5 g
20 lb	25 lb			14	28	57	110	230	570	_	
10 lb	<del>20 lb</del>			<del>12</del>	<del>23</del>	<del>46</del>		<del>180</del>	<del>450</del>	9 <del>10 mg</del>	
10 lb	<u>20 lb</u>			12	23	45	91		<u>450</u>		
5-lb				5.5		<del>22</del>					
8-1b	<u>10 lb</u>			5.5	<u>11</u>	<u>23</u>	<u>45</u>	<u>91</u>	230	<u>450</u>	
3   b				2.7	<del>5.4</del>	11 D	22	<del>43</del>			-
3   b				2.7	$\frac{5.4}{2.4}$	11	$1 + \frac{23}{14} + \frac{23}{14}$	45	110	230	
1.0								27	69		
1.0	3 ID 2 Ib			1.7	2.4	0.8 4.6	42	27 18	45	91	
1 lb by standards itch at catalog 0.55 and 0.1.1 at 0.23 d.5 d.5 d.9 d.9 d.6 d.9 d.5 lb 0.27 0.54 d.1.1 2.2 d.3 d.5 d.5 d.23 d.6 d.6 d.5 lb 0.5 lb 0.27 0.54 d.1.1 2.3 d.5 d.5 d.6 23 d.70 0.34 b.6 d.6 d.5 lb 0.77 0.34 0.68 d.4 d.6 d.6 d.7 d.6 d.6 d.7 d.6					2/30/17			18	45	91	
1 lb by standards itch at catalog 0.55 and 0.1.1 at 0.23 d.5 d.5 d.9 d.9 d.6 d.9 d.5 lb 0.27 0.54 d.1.1 2.2 d.3 d.5 d.5 d.23 d.6 d.6 d.5 lb 0.5 lb 0.27 0.54 d.1.1 2.3 d.5 d.5 d.6 23 d.70 0.34 b.6 d.6 d.5 lb 0.77 0.34 0.68 d.4 d.6 d.6 d.7 d.6 d.6 d.7 d.6				<del>0.55</del>	111	<del>2.2</del>	<u>4.4</u>	8.8	27	4 <del>5</del>	
0.5 lb         0.27         0.54         1.1         2.3         4.5         16         23         170           0.3 lb         0.17         0.34         0.68         1.4         2.7         40         44         440           0.2 lb         0.17         0.23         0.46         0.92         4.8         8.4         9.7         91           0.2 lb         0.12         0.23         0.45         0.91         1.8         8.2         9.1         97           0.1 lb         0.055         0.11         0.23         0.45         0.91         1.8         8.2         9.1         97           0.05 lb         0.055         0.11         0.23         0.45         1.1         5.4         6.8         56           0.05 lb         0.055         0.11         0.36         0.77         3.3         4.5         33           0.05 lb         0.027         0.054         0.11         0.36         0.77         3.3         4.5         37           0.03 lb         0.021         0.034         0.068         0.32         0.59         2.4         3.2         26           0.02 lb         0.017         0.034         0.068									a0d 23 66/a	astm45.617	
0.3 lb 0.2+lb         0.17 0.2 lb         0.34 0.12         0.23 0.23         0.46 0.45 0.92         1.4 0.992         2.7 1.8 1.8 0.92         11 0.7 1.8 0.97         12 0.7 0.9 1.8 0.01         1.8 0.01         8.2 0.1 0.1 0.05         9.1 0.05         9.2 0.05         9.2 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0						1.1	2.2		<del>15</del>	<del>23</del>	
0.3 lb 0.2+lb         0.17 0.2 lb         0.34 0.12         0.23 0.23         0.46 0.45 0.92         1.4 0.992         2.7 1.8 1.8 0.92         11 0.7 1.8 0.97         12 0.7 0.9 1.8 0.01         1.8 0.01         8.2 0.1 0.1 0.05         9.1 0.05         9.2 0.05         9.2 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.05         9.2 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0				0.27	0.54	<u>1.1</u>	2.3	4.5	16	23	
0.2 lb 0.11b         0.12 0.055         0.23 0.11         0.45 0.22         0.91 0.44         1.8 1.1 1.5         8.2 6.8 6.8 59         9.1 6.6 6.8 59           0.1 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.003 lb         0.11 0.027         0.054 0.011         0.11 0.36 0.077         0.30 0.77         3.3 3.3 4.5 37         4.5 32         37 0.03 lb 0.068         0.32 0.59 0.02 lb 0.02 lb 0.02 lb 0.017         0.034 0.068         0.32 0.05 lb 0.047         0.59 0.45 0.46         2.4 3.2 2.6 0.45 0.45         3.2 4.5 3.7 3.2 2.6 0.59 0.45         2.4 3.2 2.6 0.45 0.45         3.2 4.5 3.2 2.6 0.45         2.6 0.22 0.23 0.045 0.45         1.9 1.2 1.4 1.4 1.2 1.2 1.4 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2				<del>0.17</del>	0.34	0.68	1.4	2.7		14	110
0.2 lb 0.11b         0.12 0.055         0.23 0.11         0.45 0.22         0.91 0.44         1.8 1.1 1.5         8.2 6.8 6.8 59         9.1 6.6 6.8 59           0.1 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.05 lb 0.003 lb         0.11 0.027         0.054 0.011         0.11 0.36 0.077         0.30 0.77         3.3 3.3 4.5 37         4.5 32         37 0.03 lb 0.068         0.32 0.59 0.02 lb 0.02 lb 0.02 lb 0.017         0.034 0.068         0.32 0.05 lb 0.047         0.59 0.45 0.46         2.4 3.2 2.6 0.45 0.45         3.2 4.5 3.7 3.2 2.6 0.59 0.45         2.4 3.2 2.6 0.45 0.45         3.2 4.5 3.2 2.6 0.45         2.6 0.22 0.23 0.045 0.45         1.9 1.2 1.4 1.4 1.2 1.2 1.4 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2 1.2								2.7	<u>11</u>	<u>14</u>	
0.1-lb         0.055         0.11         0.22         0.44         1.1         6.4         6.8         56           0.1 lb         0.055         0.11         0.23         0.45         1.1         5.4         6.8         59           0.05 lb         0.027         0.054         0.11         0.36         0.77         3.3         4.5         37           0.03 lb         0.017         0.034         0.068         0.32         0.59         2.0         3.2         22           0.03 lb         0.017         0.034         0.068         0.32         0.59         2.4         3.2         26           0.02 lb         0.017         0.034         0.068         0.32         0.59         2.4         3.2         26           0.02 lb         0.017         0.034         0.046         0.23         0.45         1.8         2.3         19           0.02 lb         0.011         0.034         0.045         0.23         0.45         1.9         2.3         20           0.01 lb         0.015         0.023         0.045         0.23         0.45         1.9         2.3         20           0.01 lb         0.017         0.034											
0.1 lb         0.055 b         0.11 b         0.23 b         0.45 b         1.1 b         5.4 b         6.8 b         59 b           0.05 lb         0.054 b         0.054 b         0.11 b         0.36 b         0.77 d         3.3 d         4.5 d         37 d           0.03 lb         0.017 b         0.034 b         0.068 b         0.32 b         0.59 b         2.0 d         3.2 d         22 d           0.03 lb         0.017 b         0.034 b         0.068 b         0.32 b         0.59 b         2.4 d         3.2 de         26 d           0.02 lb         0.017 b         0.034 b         0.045 b         0.23 b         0.45 d         1.9 d         2.3 d         20 d           0.01 lb         0.017 b         0.034 b         0.045 b         0.23 b         0.45 d         1.9 d         2.3 d         20 d           0.01 lb         0.012 b         0.023 b         0.034 b         0.16 b         0.34 d         1.2 d         1.4 d         12 d           0.005 lb         0.015 b         0.012 b         0.023 b         0.034 b         0.16 b         0.94 d         0.14 d         0.27 d         0.86 d         0.91 f         7.8 d           0.005 lb         0.0075 b         0.015 b									_		_
0.05 lb 0-03 lb 0-03 lb 0-02 lb 0.02 lb 0.017         0.054 0-047         0.11 0-034 0-068         0.36 0-32 0-59 0-69         0.77 2.4 2.4 2.2 2.3 2.3 2.3 2.0 2.3 2.3 2.0 2.0 2.3 2.3 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.0 2.3 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0											
0.05 lb 0-03 lb 0-03 lb 0-02 lb 0.02 lb 0.017         0.054 0-047         0.11 0-034 0-068         0.36 0-32 0-59 0-69         0.77 2.4 2.4 2.2 2.3 2.3 2.3 2.0 2.3 2.3 2.0 2.0 2.3 2.3 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.3 2.3 2.0 2.0 2.0 2.0 2.3 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0 2.0									5.4	<u>6.8</u>	<u>59</u>
0.03-lb         0.047         0.034         0.068         0.32         0.59         2.0         3.2         22           0.03 lb         0.017         0.034         0.068         0.32         0.59         2.4         3.2         26           0.02 lb         0.017         0.034         0.045         0.23         0.45         1.9         2.3         20           0.01 lb         0.012         0.023         0.034         0.16         0.34         1.2         1.4         12           0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7           0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7.8           0.003 lb         0.0075         0.015         0.024         0.11         0.22         0.64         0.91         7.8           0.003 lb         0.0075         0.015         0.024         0.11         0.22         0.64         0.91         5.4           0.002 lb         0.0075         0.015         0.024         0.01         0.22         0.64         0.91         4           0.002 lb         0.00											<del>33</del> 37
0.03 lb 0.02 lb 0.017         0.017 0.034         0.068 0.046         0.32 0.45         0.59 1.8         2.4 2.3         3.2 19         26 2.3         46 19           0.02 lb 0.01 lb 0.01 lb 0.001 lb 0.005 lb 0.005 lb 0.005 lb 0.003 lb 0.0005 l									2.5	4.5	22
0.02 lb 0.01 lb         0.017 0.012         0.034 0.023         0.045 0.034         0.23 0.16         0.45 0.34         1.9 1.2         2.3 1.4         20 1.4           0.005 lb 0.005 lb 0.005 lb         0.0075 0.005 lb         0.015 0.0075         0.015 0.015         0.024 0.014         0.14 0.27         0.86 0.91         0.91 7.8         7.8           0.003 lb 0.003 lb 0.002 lb         0.0075 0.0075         0.015 0.0015         0.024 0.011         0.11 0.022         0.64 0.91         0.91 0.50         5.4 0.91            0.002 lb 0.002 lb         0.0075 0.0075         0.015 0.0075         0.024 0.015         0.091 0.024         0.11 0.091         0.50 0.90         0.91 0.50         4.2 0.91           0.001 lb 0.0005 lb 0.0005 lb 0.0005 lb 0.0008 lb 0.0008 lb 0.0008 lb 0.0009         0.024 0.0015 0.0001         0.068 0.024 0.0024         0.15 0.024 0.068 0.068 0.15 0.068 0.15 0.068 0.15 0.068 0.15 0.36 0.90 0.90 0.90 0.90 0.90 0.90 0.90 0.9									2.4	3.2	26
0.02 lb 0.01 lb         0.017 0.012         0.034 0.023         0.045 0.034         0.23 0.16         0.45 0.34         1.9 1.2         2.3 1.4         20 1.4           0.005 lb 0.005 lb 0.005 lb         0.0075 0.005 lb         0.015 0.0075         0.015 0.015         0.024 0.014         0.14 0.27         0.86 0.91         0.91 7.8         7.8           0.003 lb 0.003 lb 0.002 lb         0.0075 0.0075         0.015 0.0015         0.024 0.011         0.11 0.022         0.64 0.91         0.91 0.50         5.4 0.91            0.002 lb 0.002 lb         0.0075 0.0075         0.015 0.0075         0.024 0.015         0.091 0.024         0.11 0.091         0.50 0.90         0.91 0.50         4.2 0.91           0.001 lb 0.0005 lb 0.0005 lb 0.0005 lb 0.0008 lb 0.0008 lb 0.0008 lb 0.0009         0.024 0.0015 0.0001         0.068 0.024 0.0024         0.15 0.024 0.068 0.068 0.15 0.068 0.15 0.068 0.15 0.068 0.15 0.36 0.90 0.90 0.90 0.90 0.90 0.90 0.90 0.9									1.8	2.3	<del>19</del>
0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7           0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7.8           0.003 lb         0.0075         0.015         0.024         0.11         0.22         0.64         0.91         5.4           0.002 lb         0.0075         0.015         0.024         0.011         0.22         0.64         0.91         5.4           0.002 lb         0.0075         0.015         0.024         0.091         0.19         0.50         0.91         4.2           0.001 lb         0.0075         0.015         0.024         0.091         0.19         0.50         0.91         4.2           0.001 lb         0.0075         0.015         0.024         0.068         0.15         0.36         0.91         2.9           0.0005 lb         0.0095         0.015         0.024         0.068         0.15         0.36         0.91         2.9           0.0004 lb         0.0095         0.023         0.07         0.15         0.3         0.9         3           0.0003 lb         0											20
0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7           0.005 lb         0.0075         0.015         0.024         0.14         0.27         0.86         0.91         7.8           0.003 lb         0.0075         0.015         0.024         0.11         0.22         0.64         0.91         5.4           0.002 lb         0.0075         0.015         0.024         0.011         0.22         0.64         0.91         5.4           0.002 lb         0.0075         0.015         0.024         0.091         0.19         0.50         0.91         4.2           0.001 lb         0.0075         0.015         0.024         0.091         0.19         0.50         0.91         4.2           0.001 lb         0.0075         0.015         0.024         0.068         0.15         0.36         0.91         2.9           0.0005 lb         0.0095         0.015         0.024         0.068         0.15         0.36         0.91         2.9           0.0004 b         0.0095         0.023         0.07         0.15         0.3         0.9         3           0.0003 lb         0.	0.01 lb			0.012	0.023	0.034	0.16	0.34	1.2	1.4	12
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	0.005 lb			0.0075	0.015	0.024	0.14	0.27	0.86	<del>0.91</del>	7
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$								0.27			<u>7.8</u>
0.002   b   0.0075   0.015   0.024   0.091   0.19   0.50   0.91   4.2											
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$											<u>5.4</u>
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$											
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$											4.2
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$											
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$								U.15	0.36	<u>0.91</u>	2.9
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$										Ð	21
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$										3	<u> </u>
0.0002 lb         0.0045         0.0091         0.023         0.045         0.045         0.045         0.095         0.20         1.2           0.0001 lb         0.0001 lb         0.009         0.023         0.07         0.15         0.3         0.9         3										J	1.5
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$										3	<u></u>
0.0001 lb 0.009 0.023 0.07 0.15 0.3 0.9 3										-	1.2
<u>0.0001 lb</u> <u>0.0045</u> <u>0.0091</u> <u>0.023</u> <u>0.041</u> <u>0.086</u> <u>0.16</u> <u>0.84</u>										3	_
	0.0001 lb			0.0045	0.0091	0.023	0.041	0.086	<u>0.16</u>		0.84



				IADLE	1 Continu	ea				
Denomination					±mg exce	pt as noted				
Metric	Class 000	Class 00	Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class
0.00005 lb			0.009	0.023	0.07	0.15	0.3	0.9	3	
0.00005 lb			0.0045	0.0091	0.014	0.036	0.073	0.14		
0.00003 lb			0.009	0.023	0.07	<del>0.15</del>	0.3	0.9	3	
0.00003 lb			0.0045	0.0091	0.014	0.032	0.064	0.11		
0.00002 lb			0.009	0.023	0.07	0.15	0.3	0.9	3	
0.00002 lb			0.0045	0.0091	0.014	0.029	0.059	0.091	_	
0.00001 lb			0.009	0.023	0.07	0.15	0.3	0.9	3	
0.00001 lb			0.0045	0.0091	0.014	0.027	0.054	0.091		
0.000005 lb			0.0045	0.0091	0.014	0.023	0.054	0.091		
0.000003 lb 0.000002 lb			0.0045	0.0091 0.0091	0.014	0.023	0.054 0.054	0.091 0.091		
0.000002 lb			0.0045	0.0091	0.014	0.023	0.054	0.091		
Avoirdupois			Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class
Ounce			mg	mg	mg	mg	mg	mg	mg	mg
<del>10 oz</del>			0.4	<del>0.7</del>	1.4	2.8	5.4	<del>19</del>	45	320
10 oz			0.35	0.70						200
8 oz			0.3	0.6	1.4 1.2	2.8 2.3	5.4 4.5	19 16	23 23	180
8 oz			027	0.54	1.1	2.3	4.5	16	23	170
<del>5 oz</del>			0.18	0.35	0.70	2.3 1.4	4.5 2.8	16 12	23 16	160
<u>5 oz</u>			$\frac{0.18}{0.14}$	0.35	0.71	1.4 1.1	2.8 2.3	<u>12</u> 9.5	11 11	130 110
4 oz				0.28	0.55	<del>1.1</del>	<del>2.3</del>	<del>9.5</del>		
<u>4 oz</u>			0.14	0.28	0.57	1.1	2.3	9.5 8.2	<u>11</u>	110
<del>3 oz</del>			0.12	0.23	0.45	0.91	<del>1.8</del>	<del>8.2</del>	<del>8.1</del>	73
3 oz			0.10	0.21	0.43	0.91	1.8 1.3	<u>8.2</u> <del>5.9</del>	5.4 5.4	
<del>2 oz</del>			0.07	0.13	0.26	0.64	<del>1.3</del>	<del>5.9</del>	<del>5.4</del>	<del>48</del>
2 oz			0.070	0.14	0.28	0.64	1.3	5.9	$\frac{5.4}{3.2}$	$\frac{70}{20}$
1 oz			0.04	0.07	0.14	0.42	0.86	3.9	<del>3.2</del>	<del>28</del>
<u>1 oz</u> <del>1/2 oz</del>			0.027 <del>0.02</del>	0.054 0.04	0.14 0.08	$\frac{0.43}{0.3}$	0.86	3.9 2.5	$\frac{3.2}{2.3}$	42 25
							0.59	<del>2.5</del>		25
1/2 oz <del>1/4 oz</del>			0.017 0.015	$\frac{0.034}{0.03}$	$\frac{0.071}{0.06}$	$\frac{0.30}{0.2}$	0.59 0.43	2.5 1.5	2.3 1.4	<u>27</u> 9.1
1/4 oz			0.013	0.023	0.034	0.20	0.43	1.6	1.4	17
1/8 oz			0.011	0.029	0.054	0.16	0.43	1.0 1.1	<del>1.4</del> <del>0.91</del>	<u>17</u> 4.3
1/8 oz			0.0075	0.015	0.024	0.16	0.31	1.1	0.91	10
1/16 oz			0.013	0.025	0.050	0.12	0.24	<del>0.73</del>	<del>0.91</del>	<u>10</u> <del>4.3</del>
1/16 oz			0.0075	0.015	0.024	0.12	0.24	0.73	0.91	6.5
<del>1/32 oz</del>			0.008	0.015	0.030	0.095	0.19	0.5	0.91	6.5 4.3
1/32 oz			0.0075	0.015	0.024	0.095	0.19	0.50		4.2
<del>1/64 oz</del>			0.006	0.012	0.024	0.077	<del>0.15</del>	0.36	<del>0.91</del>	$\frac{4.2}{4.3}$
1/64 oz			0.0075	0.015	0.023	0.077	0.15	0.36		2.8 <del>25</del>
<del>0.50 oz</del>					0.080	9.3	0.59	2.5	<del>2.3</del>	<del>25</del>
<u>0.5 oz</u>					0.071	0.30	0.59	a0d <del>1.8</del> 66/8	2.3	- 1 8 <del>27</del>
					87 <del>0.068</del> 3	e8-4 <del>0.23</del> a-8	175 <del>0.45</del> ea		astm <del>1.4</del> 617	
0.3 oz					$\frac{0.034}{0.057}$	0.23 0.19	0.45	1.8 1.4	1.4	<u>19</u> 5.9
<del>0.2 oz</del>							0.38		0.91	5.9
0.2 oz <del>0.1 oz</del>					0.034 0.050	0.19 <del>0.14</del>	0.38 0.29	<u>1.4</u> <del>0.91</del>	0.91 <del>0.91</del>	14 4.3
<del>0.1 oz</del> 0.1 oz					0.024	0.14 0.14	0.29	0.91	0.91	<del>4.3</del> 9.0
0.1 02 0.05 oz					0.024 0.050	0.14 0.11	0.29 0.23	0.91 0.64	0.91 0.91	9.0 2.0
0.05 oz					0.024	0.11	0.23	0.64	0.91	5.7
0.03 oz					0.024	0.095	0.19	0.45	0.91	2.0
0.03 oz					0.024	0.095	0.19	0.45	0.91	4.1
0.02 oz					0.023	0.077	0.18	0.4	0.91	2.0
0.02 oz					0.024	0.077	0.18	0.40	0.91	3.2
0.01 oz					0.023	0.064	0.14	0.3	0.91	2.0
0.01 oz					0.023	0.064	0.14	0.30	0.91	2.2
<del>0.005 oz</del>					0.023	0.054	0.11	0.23	**	
0.005 oz					0.023	0.054	0.11	0.23		1.5
0.003 oz					0.023	0.05	0.095	0.19	**	
0.003 oz					0.023	0.050	0.095	0.19	**	1.1
0.002 oz					0.023	0.044	0.086	0.16	**	**
0.002 oz					0.023	0.044	0.086	0.16	**	0.92
0.001 oz					0.023	0.038	0.077	0.13		
0.001 oz					0.023	0.038	0.077	0.13	**	0.66
0.0005 oz					0.023	0.031	0.064	0.11		
0.0005 oz 0.0003 oz					0.014 0.023	0.032 0.029	0.064 0.059	<u>0.11</u> 0.095	**	0.4
0.0003 oz					0.02 <del>3</del> 0.014	0.029	0.059	0.095	_	_
					0.014	0.029 0.027	0.059 0.054	0.095 0.086	**	**
0.0002.03					0.014	0.027	0.054	0.086		_
0.0002 oz					0.014	0.027 0.026	0.054 0.05	0.000 0.073	**	**
0.0002 oz							0.00	0.070		
0.0002 oz 0.0001 oz										
0.0002 oz			Class 0	Class 1	0.014 Class 2	0.026 Class 3	0.050 Class 4	0.073 Class 5	Class 6	Class



Denomination										
					±mg excep	t as noted				
Metric	Class 000	Class 00	Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class 7
	0.000 000	01400 00	01000 0	01400 1				01000 0	01400 0	0,000 7
1000					mg	mg	g & mg	4.0		-
1000 oz t					<del>160 mg</del>	<del>310 mg</del>	<del>620 mg</del>	<del>1.6 g</del>	<del>3.2 g</del>	
1000 oz t					00	310 mg	620 mg	<u>1.6 g</u>	4.5	
<del>500 oz t</del>					<del>80</del>	<del>160</del>	<del>310</del>	<del>770 mg</del>	<del>1.5 g</del>	
500 oz t						<u>160</u>	<u>310</u>	780 mg		
<del>300 oz t</del>					<del>45</del>	90	<del>190</del>	<del>470</del>	<del>940 mg</del>	
300 oz t						91 <del>62</del>	190 120	<u>470</u>		
<del>200 oz t</del>					<del>31</del>	<del>62</del>	<del>120</del>	<del>310</del>	<del>620</del>	
200 oz t						<u>62</u> <del>31</del>	120 <del>62</del>	<u>310</u>		
<del>100 oz t</del>					<del>16</del>	<del>31</del>	<del>62</del>	<del>160</del>	<del>320</del>	
100 oz t						31	62	160		
<del>50 oz t</del>					8	31 <del>16</del>	<u>62</u> <del>31</del>	160 <del>77</del>	<del>150</del>	
50 oz t						<u>16</u> 9.1	31 19	78		
30 oz t					4.6	9.1	<del>19</del>	78 47	94	
30 oz t						9.1	19 12	47		
<del>20 oz t</del>					<del>3.1</del>	<del>6.2</del>	<del>12</del>	47 35	<del>58</del>	
20 oz t							12	35		
10 oz t					<del>1.6</del>	6.2 3.1	<u>12</u> <del>6.2</del>	35 21	<del>32</del>	
10 oz t						3.1	6.2	21	02	
5 oz t					0.8	<del>3.1</del> <del>1.6</del>	6.2 3.1	12	<del>16</del>	
5 oz t					0.0	1.6	3.1		10	
3 02 t					0.46	0.91	3.1 1.9	<u>12</u> <del>8.4</del>	9.3	
					<del>0.40</del>			0. <del>4</del> 0.1	<del>3.3</del>	
3 oz t 2 oz t						0.91 <del>0.71</del>	1.9 1.4	8.4 6.4	8.0	
									0.0	
2 oz t						0.71 <del>0.45</del>	1.4 0.91	6.5 4.2	4.0	
<del>1 oz t</del>								4.2	<del>4.2</del>	
1 oz t						0.45 0.31	0.91	$\frac{4.2}{2.6}$	0.0	
0.5 oz t						0.31	0.62	2.6	<del>2.6</del>	
<u>0.5 oz t</u>						0.31 0.24 0.24 0.20	0.62	2.6		
0.3 oz t						0.24	0.49	<del>1.9</del>	<del>1.9</del>	
<u>0.3 oz t</u>						0.24	0.49 0.40	1.9 1.5		
0.2 oz t							0.40		<del>2.3</del>	
0.2 oz t						0.20 0.15	0.40	<u>1.5</u>		
0.1 oz t						0.15	0.30	0.97	<del>2.0</del>	
0.1 oz t						0.15 0.12	0.30	0.97		
0.05 oz t						0.12	0.23	0.65	<del>2.0</del>	
0.05 oz t						0.12 0.097	0.23	0.65		
0.03 oz t						0.097	$\frac{0.23}{0.19}$	0.49	<del>2.0</del>	
0.03 oz t						0.097	0.19	0.49		
0.02 oz t						0.084	0.17	0.41	<del>1.1</del>	
0.02 oz t							0.17	0.41		
0.01 oz t						<u>0.084</u> <del>0.071</del>	0.14	0.31	1.0	
						8-0.071 -8	f750.142ea		astm-e617	
414 F16655						70 <u>4.6.7.</u> 0 0		uou <u>eren</u> oon		
0.005.07.1						0.056	0.11	0.23		
0.005 oz t						<del>0.056</del> 0.056	0.11	0.23	<del>1.0</del>	
0.005 oz t						0.056	<u>0.11</u>	0.23 0.23	1.0	
0.005 oz t 0.003 oz t						0.056 0.049	0.11 0.097	0.23 0.19		
0.005 oz t <del>0.003 oz t</del> 0.003 oz t						0.056 0.049 0.049	0.11 0.097 0.097	0.23 0.23 0.19 0.19	<del>1.0</del> <del>1.0</del>	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t						0.056 0.049 0.049 0.044	0.11 0.097 0.097 0.091	0.23 0.23 0.19 0.19 0.17	1.0	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t						0.056 0.049 0.049 0.044 0.044	0.11 0.097 0.097 0.091 0.091	0.23 0.19 0.19 0.17	1.0 1.0 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t						0.056 0.049 0.049 0.044 0.044	0.11 0.097 0.097 0.091 0.091 0.078	0.23 0.23 0.19 0.19 0.17 0.17	<del>1.0</del> <del>1.0</del>	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t						0.056 0.049 0.049 0.044 0.038 0.038	0.11 0.097 0.097 0.091 0.091 0.078	0.23 0.19 0.19 0.17 0.17 0.14 0.14	1.0 1.0 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.005 oz t						0.056 0.049 0.049 0.044 0.044 0.038 0.038 0.033	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065	0.23 0.49 0.19 0.17 0.17 0.14 0.14 0.11	1.0 1.0 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t						0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11	1.0 1.0 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t						0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097	1.0 1.0 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t						0.056 0.049 0.044 0.044 0.038 0.038 0.033 0.033 0.030 0.030	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097	1.0 1.0 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0003 oz t						0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.030	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.060	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.097	1.0 1.0 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t						0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.060 0.056	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084	1.0 1.0 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t		en, at/cata.id	g/standar			0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.028	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.056 0.056 0.056	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071	1.0 1.0 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t		en, al/catalo				0.056 0.049 0.044 0.044 0.038 0.038 0.033 0.033 0.030 0.028 0.028 0.028	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.056 0.056 0.056 0.052	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.11 0.097 0.097 0.084 0.084 0.071	1.0 1.0 0.50 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t		en, al/catalo	Class 0	Class 1	Class 2	0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.028	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.060 0.056 0.056 0.052 0.052	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 Class 5	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.060 0.056 0.056 0.052 0.052	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.11 0.097 0.097 0.084 0.084 0.071	1.0 1.0 0.50 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t 0.0001 oz t		en.al/catalo		Class 1	Class 2 mg	0.056 0.049 0.044 0.044 0.038 0.038 0.033 0.033 0.030 0.028 0.028 0.028	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 9 & mg	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 0.071 Class 5 g & mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.005 oz t 0.003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.060 0.056 0.056 0.052 0.052 Class 4 g.& mg g.& mg	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.094 0.084 0.071 0.071 Class 5 g-& mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 9 & mg	0.23 0.23 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 0.071 Class 5 g & mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en, al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.030 0.028 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.060 0.056 0.056 0.052 0.052 Class 4 g.& mg g.& mg	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 0.071 Class 5 g-8 mg 0.78 g 0.39 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.030 0.028 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 g & mg 0.31 g	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.071 0.071 Class 5 g & mg 0.78 g 0.78 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en. al/catald		Class 1		0.056 0.049 0.049 0.044 0.038 0.033 0.033 0.030 0.028 0.028 0.026 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 Class 4 g & mg 0.31 g 0.16 g	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 0.071 Class 5 g-8 mg 0.78 g 0.39 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.001 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.044 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.026 Class 3 g & mg	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 g & mg g & mg 0.31 g 0.16 g 0.16 g	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 Class 5 g & mg 0.78 g 0.39 g 0.39 g 0.39 g 0.23 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight 10000 dwt 10000 dwt 5000 dwt 5000 dwt 3000 dwt		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.030 0.028 0.028 0.028 0.026 Class 3	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 9 & mg 0.31 g 0.16 g 0.16 g 91 mg	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 0.071 Class 5 g & mg 0.78 g 0.39 g 0.39 g 0.23 g 0.23 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 0.	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight 10000 dwt 5000 dwt 5000 dwt 3000 dwt 2000 dwt		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.044 0.038 0.033 0.033 0.030 0.028 0.028 0.026 Class 3  g & mg  155 78 47	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 0.052 Class 4 g. 8 mg 0.31 g 0.16 g 0.16 g 0.16 g 91 mg 62	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.097 0.084 0.071 0.071 Class 5 g-& mg 0.78 g 0.39 g 0.39 g 0.23 g 0.16 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg	
0.005 oz t 0.003 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.0005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0001 oz t 0.0001 oz t 0.0002 oz t 0.0001 oz t		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.044 0.038 0.033 0.033 0.030 0.030 0.028 0.028 0.026 Class 3 g & mg	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 Class 4 g & mg g & mg 0.31 g 0.16 g 91 mg 91 mg 62 62	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.11 0.011 0.097 0.097 0.094 0.084 0.071 Class 5 g & mg 0.78 g 0.39 g 0.39 g 0.23 g 0.23 g 0.16 g	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg 1.5 g 0.78 g 0.46 g 0.32 g	
0.005 oz t 0.003 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0003 oz t 0.0002 oz t 0.0001 oz t 0.0001 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight  10000 dwt 10000 dwt 5000 dwt 3000 dwt 2000 dwt 2000 dwt		en. al/catalo		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.028 0.028 0.026 Class 3  g & mg  155 78 47	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 Class 4 g & mg g & mg 0.31 g 0.16 g 91 mg 91 mg 62 62 31	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.097 0.084 0.084 0.071 Class 5 g-8 mg 0.78 g 0.39 g 0.23 g 0.16 g 0.16 g 78 mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 0.	
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0003 oz t 0.0000 oz t 0.0003 oz t 0.0001 oz t 0.0001 oz t 0.0001 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t		en, al/catalo		Class 1		0.056 0.049 0.049 0.044 0.044 0.038 0.033 0.033 0.030 0.028 0.028 0.026 Class 3  g & mg  155 78 47	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.052 0.052 0.052 Class 4 g & mg g & mg 0.31 g 0.16 g 91 mg 62 62 31 31	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.084 0.084 0.071 Class 5 g & mg 0.78 g 0.39 g 0.39 g 0.23 g 0.16 g 78 mg 78 mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg 1.5 g 0.78 g 0.46 g 0.32 g 0.16 g	Class 7
0.005 oz t 0.003 oz t 0.003 oz t 0.002 oz t 0.002 oz t 0.001 oz t 0.005 oz t 0.005 oz t 0.0005 oz t 0.0003 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0002 oz t 0.0001 oz t Pennyweight 10000 dwt 10000 dwt 5000 dwt 3000 dwt 2000 dwt 2000 dwt 2000 dwt 2000 dwt 2000 dwt		en.al/catak		Class 1		0.056 0.049 0.049 0.044 0.038 0.038 0.033 0.033 0.030 0.028 0.028 0.026 Class 3  g & mg  155 78 47	0.11 0.097 0.097 0.091 0.091 0.078 0.065 0.065 0.060 0.056 0.056 0.052 Class 4 g & mg g & mg 0.31 g 0.16 g 91 mg 91 mg 62 62 31	0.23 0.23 0.19 0.19 0.17 0.17 0.14 0.14 0.11 0.097 0.097 0.097 0.084 0.084 0.071 Class 5 g-8 mg 0.78 g 0.39 g 0.23 g 0.16 g 0.16 g 78 mg	1.0 1.0 0.50 0.50 0.50 0.50 0.50 0.50 Class 6 g & mg 1.5 g 0.78 g 0.46 g 0.32 g	



Denomination					±mg excep	ot as noted				
Metric	Class 000	Class 00	Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class
300 dwt							9.1	<del>28</del>	<del>56</del>	
300 dwt						<u>4.7</u>	9.1 6.2	28 21	40	
<del>200 dwt</del> 200 dwt						3.1	6.2	<del>21</del> 21	42	
100 dwt						<u>3.1</u>	<u>6.2</u> <del>3.1</del>	21 12	<del>24</del>	
100 dwt						<u>1.6</u>	3.1	12		
50 dwt							3.1 1.6	<u>12</u> <del>7.8</del>	<del>16</del>	
50 dwt						0.78	1.6 1.2	7.8 5.3		
30 dwt						0.50	<del>1.2</del>	<del>5.3</del>	11	
30 dwt 20 dwt						0.58	<u>1.2</u> <del>0.91</del>	5.3 4.2	8.4	
20 dwt						0.46	0.91	4.2	0.4	
10 dwt						<u> </u>	0.62	4.2 2.6	<del>5.2</del>	
10 dwt						0.31	0.62	2.6 1.7		
<del>5 dwt</del>							0.44	<del>1.7</del>	3.4	
5 dwt						0.22	0.44	1.7 1.3		
3 dwt						0.17	0.34	1.3	<del>2.6</del>	
3 dwt 2 dwt						<u>0.17</u>	0.34 <del>0.3</del>	<u>1.3</u> <del>0.97</del>	<del>1.9</del>	
2 dwt						0.15	0.3	0.97	1.9	
1 dwt							0.23	0.65	<del>1.3</del>	
1 dwt						0.12	0.23	0.65		
Grain			Class 0	Class 1	Class 2	Class 3	Class 4	Class 5	Class 6	Class
10000 ~=						mg	mg <del>13</del>	mg	mg	
<del>10000 gr</del> 10000 gr						6.5	13 13	<del>36</del> 36	<del>36</del> 64	
5000 gr						0.5	13 6.5	36 <del>22</del>	$\frac{64}{22}$	
5000 gr						3.2	6.5	22	32	
3000 gr							6.5 3.9 3.9 2.6	<u>22</u> <del>15</del>	32 <del>15</del>	
3000 gr						21 <u>1.9</u> S	3.9	<u>15</u> <del>11</del>	20 <del>11</del>	
<del>2000 gr</del>							<del>2.6</del>	<del>11</del>	<del>11</del>	
2000 gr						1.3	2.6 1.4	<u>11</u> <del>6.3</del>	14 6.3	
<del>1000 gr</del> 1000 gr						0.71	1.4	<del>0.3</del> 6.5	<del>6.3</del>	
500 gr						0.71	0.91	6.5 3.2	6.4 3.2	
500 gr						0.48	0.91	4.3	3.2	
300 gr							0.65	4.3 2.3	3.2 2.3	
300 gr						0.35	0.65	$\frac{3.0}{1.4}$	$\frac{2.3}{1.4}$	
<del>200 gr</del>							0.57	<del>1.4</del>		
200 gr						0.28	0.57 <del>0.4</del>	<u>2.3</u> <del>0.91</del>	<u>1.4</u> <del>0.91</del>	
<del>100 gr</del> 100 gr						<u>0.20</u>	0.40	1.5	0.91	
50 gr						8-487a-8	f75 <del>0.3</del> 2ea	a0d <del>0.91</del> 66/		
50 gr						0.15	0.30	0.97	0.91	
<del>30 gr</del>							0.25	0.71	<del>0.71</del>	
30 gr						0.12	0.25	0.71	0.91	
<del>20 gr</del>							0.21	0.58	0.58	
20 gr						<u>0.11</u>	0.21 0.17	0.58 0.42	0.91	
<del>10 gr</del> 10 gr						0.084	<del>0.17</del> 0.17		<del>0.42</del> 0.91	
<u>10 gr</u> <del>5 gr</del>						0.004	0.17	0.42 <del>0.31</del>	0.91	
5 gr						0.071	0.14	0.31	0.45	
<del>3 gr</del>							<del>0.12</del>	0.25	0.25	
<u>3 gr</u>						0.060	0.12	0.25 <del>0.22</del>	0.45	
<del>2 gr</del>						0.07:	0.11	0.22	0.22	
2 gr						0.054	0.11	$\frac{0.22}{0.17}$	$\frac{0.45}{0.17}$	
<del>1 gr</del> 1 gr						0.045	0.091	<del>0.17</del> 0.17	<del>0.17</del> 0.45	
<u>1 gr</u> <del>0.5 gr</del>						0.043	0.091 0.078	0.17 0.14	0.45 0.14	
0.5 gr						0.039	0.078	0.14	0.45	
<del>0.3 gr</del>							0.071	$\frac{0.14}{0.12}$	0.12	
0.3 gr						0.035	0.071	0.12		
0.2 gr							0.064	0.11	0.11	
0.2 gr						0.032	0.064	0.11	0.001	
<del>0.1 gr</del>						0.000	<del>0.056</del>	0.091	0.091	
<u>0.1 gr</u> <del>0.05 gr</del>						0.028	0.056 0.052	0.091 0.071	0.071	
0.05 gr						0.026	0.052	0.071	5.57	
0.03 gr							0.051	0.071	0.071	
0.03 gr						0.025	0.051	0.071		
<del>0.02 gr</del>						<del>_</del>	0.05	0.071	0.071	
0.02 gr						0.025	0.050	0.071	0.07:	
0.01 gr						0.005	0.05	<del>0.071</del>	<del>0.071</del>	
<u>0.01 gr</u>						0.025	<u>0.050</u>	<u>0.071</u>		