

Designation: D7483 - 21

# Standard Test Method for Determination of Dynamic Viscosity and Derived Kinematic Viscosity of Liquids by Oscillating Piston Viscometer<sup>1</sup>

This standard is issued under the fixed designation D7483; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon  $(\varepsilon)$  indicates an editorial change since the last revision or reapproval.

# 1. Scope\*

- 1.1 This test method covers the measurement of dynamic viscosity and derivation of kinematic viscosity of liquids, such as new and in-service lubricating oils, by means of an oscillating piston viscometer.
- 1.2 This test method is applicable to Newtonian and non-Newtonian liquids; however the precision statement was developed using Newtonian liquids.
- 1.3 The range of dynamic viscosity covered by this test method is from 0.2 mPa·s to 20 000 mPa·s (which is approximately the kinematic viscosity range of 0.2 mm²/s to 22 000 mm²/s for new oils) in the temperature range between –40 °C to 190 °C; however the precision has been determined only for new and used oils in the range of 34 mPa·s to 1150 mPa·s at 40 °C, 5.7 mPa·s to 131 mPa·s at 100 °C, and 46.5 mm²/s to 436 mm²/s at 40 °C.
- 1.4 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.
- 1.5 This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety, health, and environmental practices and determine the applicability of regulatory limitations prior to use.
- 1.6 This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.

#### 2. Referenced Documents

2.1 ASTM Standards:<sup>2</sup>

D445 Test Method for Kinematic Viscosity of Transparent and Opaque Liquids (and Calculation of Dynamic Viscosity)

D2162 Practice for Basic Calibration of Master Viscometers and Viscosity Oil Standards

D4057 Practice for Manual Sampling of Petroleum and Petroleum Products

D4177 Practice for Automatic Sampling of Petroleum and Petroleum Products

D5967 Test Method for Evaluation of Diesel Engine Oils in T-8 Diesel Engine

D6300 Practice for Determination of Precision and Bias
Data for Use in Test Methods for Petroleum Products,
Liquid Fuels, and Lubricants

D6708 Practice for Statistical Assessment and Improvement of Expected Agreement Between Two Test Methods that Purport to Measure the Same Property of a Material

D6792 Practice for Quality Management Systems in Petroleum Products, Liquid Fuels, and Lubricants Testing Laboratories - 4548 (4781510/astm-d7483-2)

2.2 ISO Standards:<sup>3</sup>

ISO/EC 17025 General Requirements for the Competence of Testing and Calibration Laboratories

2.3 NIST Standard:<sup>4</sup>

NIST Technical Note 1297 Guideline for Evaluating and Expressing the Uncertainty of NIST Measurement Results

#### 3. Terminology

- 3.1 Definitions:
- 3.1.1 *density*, *n*—mass per unit volume.

<sup>&</sup>lt;sup>1</sup> This test method is under the jurisdiction of ASTM Committee D02 on Petroleum Products, Liquid Fuels, and Lubricants and is the direct responsibility of Subcommittee D02.07 on Flow Properties.

Current edition approved Jan. 1, 2021. Published January 2021. Originally approved in 2008. Last previous edition approved in 2020 as D7483 – 20a. DOI: 10.1520/D7483-21.

<sup>&</sup>lt;sup>2</sup> For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

<sup>&</sup>lt;sup>3</sup> Available from International Organization for Standardization (ISO), 1, ch. de la Voie-Creuse, Case postale 56, CH-1211, Geneva 20, Switzerland, http://

<sup>&</sup>lt;sup>4</sup> Available from http://physics.nist.gov/ccu/Uncertainty/index.html.



FIG. 1 Viscometer with Electronics

- 3.1.2 *dynamic viscosity* ( $\eta$ ), n—the ratio between the applied shear stress and rate of shear of a liquid at a given temperature.
- 3.1.2.1 *Discussion*—It is sometimes called the coefficient of dynamic viscosity or, simply, viscosity. Thus, dynamic viscosity is a measure of the resistance to flow or to deformation of a liquid under external shear forces.
- 3.1.2.2 *Discussion*—The term dynamic viscosity can also be used in a different context to denote a frequency-dependent quantity in which shear stress and shear rate have a sinusoidal time dependence.
- 3.1.3 *kinematic viscosity*, n—the ratio of the dynamic viscosity ( $\eta$ ) to the density ( $\rho$ ) of a liquid at a given temperature.
- 3.1.3.1 *Discussion*—For gravity flow under a given hydrostatic head, the pressure head of a liquid is proportional to its density, ( $\rho$ ). Therefore, the kinematic viscosity, ( $\nu$ ), is a measure of the resistance to flow of a liquid under gravity.
- 3.1.4 rate of shear (shear rate), n—in liquid flow, the velocity gradient across the liquid.
- 3.1.5 *shear stress*, *n*—the force per unit area in the direction of the flow.
- 3.1.5.1 *Discussion*—The SI unit for shear stress is the pascal (Pa).
  - 3.2 Definitions of Terms Specific to This Standard:
- 3.2.1 oscillating piston viscometer, n—a device that measures the travel time of a piston driven electromagnetically into stationary oscillating motion through a liquid at a controlled force in order to determine the dynamic viscosity of the liquid.

# 4. Summary of Test Method

4.1 A specimen of sample is placed in the thermally controlled measurement chamber where the piston resides. The piston is driven into oscillatory motion within the measurement chamber by a controlled magnetic field. Once the sample is at the test temperature, as determined by the temperature detector, the piston is propelled repeatedly through the liquid (by the magnetic field). A shear stress (ranging from 5 Pa to 750 Pa) is imposed on the liquid under test due to the piston travel. The dynamic viscosity is determined by measuring the average travel time of the piston. The kinematic viscosity is derived by additionally measuring the ratio between the up and down

travel times. This information is then applied to a calibration curve using liquids of known viscosity to calculate the dynamic viscosity. The kinematic viscosity is derived by an externally measured density by additionally measuring the ratio between the up and down travel times. The precision and bias data for kinematic viscosity (as published in RR:D02-1755<sup>5</sup>) were derived by externally measured density and do not apply to the internal density measurement.

#### 5. Significance and Use

- 5.1 Many petroleum products, as well as non-petroleum materials, are used as lubricants for bearings, gears, compressor cylinders, hydraulic equipment, etc. Proper operation of this equipment depends upon the viscosity of these liquids.
- 5.2 Oscillating piston viscometers allow viscosity measurement of a broad range of materials including transparent, translucent and opaque liquids. The measurement principle and stainless steel construction makes the Oscillating Piston Viscometer resistant to damage and suitable for portable operations. The measurement itself is automatic and does not require an operator to time the oscillation of the piston. The electromagnetically driven piston mixes the sample while under test. The instrument requires a sample volume of less than 5 mL and typical solvent volume of less than 10 mL which minimizes cleanup effort and waste.

## 6. Apparatus

6.1 Oscillating Piston Viscometer:<sup>6,7</sup>

<sup>&</sup>lt;sup>5</sup> Supporting data have been filed at ASTM International Headquarters and may be obtained by requesting Research Report RR:D02-1755. Contact ASTM Customer Service at service@astm.org.

<sup>&</sup>lt;sup>6</sup> The Oscillating Piston Viscometer is covered by a patent. Interested parties are invited to submit information regarding the identification of an alternative to this patented item to the ASTM International headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend.

<sup>&</sup>lt;sup>7</sup> The sole sources of supply for the apparatus known to the committee at this time is Cambridge Viscosity Inc., 101 Station Landing, Medford, MA 02155 (www.cambridgeviscosity.com). If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend.

**TABLE 1 Viscosity Ranges of Oscillating Viscometer Pistons** 

| Minimum Viscosity (mPa·s ) | Maximum Viscosity (mPa·s ) | Piston Designation | Nominal Piston Diameter (mm) | Recommended Sample Volume (mL) |
|----------------------------|----------------------------|--------------------|------------------------------|--------------------------------|
| 0.02                       | 2                          | SP20               | 7.87                         | 3.2–5                          |
| 0.25                       | 5                          | SP50               | 7.83                         | 3.2-5                          |
| 0.5                        | 10                         | SP11               | 7.81                         | 3.2-5                          |
| 1                          | 20                         | SP21               | 7.76                         | 3.5–5                          |
| 2.5                        | 50                         | SP51               | 7.68                         | 3.5-5                          |
| 5                          | 100                        | SP12               | 7.62                         | 3.5-5                          |
| 10                         | 200                        | SP22               | 7.54                         | 3.5-5                          |
| 25                         | 500                        | SP52               | 7.34                         | 3.5-5                          |
| 50                         | 1000                       | SP13               | 7.21                         | 4.0-5                          |
| 100                        | 2000                       | SP23               | 6.96                         | 4.0-5                          |
| 250                        | 5000                       | SP53               | 6.27                         | 4.0-5                          |
| 500                        | 10 000                     | SP14               | 6.05                         | 4.0-5                          |
| 1000                       | 20 000                     | SP24               | 5.72                         | 4.0–5                          |

- 6.1.1 The oscillating piston viscometer (see Fig. 1) comprises a measurement chamber and calibrated piston capable of measuring the dynamic viscosity within the limits of precision given in Section 16.
- 6.1.2 *Piston*—Free moving, magnetically driven body within a Oscillating Piston Viscometer which is used for measuring the viscosity of liquids. Individual pistons are sized to measure specific viscosity ranges by varying the sensor annulus. See Table 1 for the selection of the piston according to the viscosity range.
- 6.1.3 *Measurement Chamber*—Location within Oscillating Piston Viscometer where piston motion (through the liquid under test) occurs due to an imposed electromagnetic field. See Fig. 2.
- 6.1.4 *Electronics*—Capable of controlling the electromagnetic field to propel and detect the travel time of the piston with a discrimination of  $0.01 \, \mathrm{s}$  or better and uncertainty within  $\pm 0.07 \, \%$ . The travel time is calibrated to be between  $0.4 \, \mathrm{s}$  and  $60 \, \mathrm{s}$ , at a distance of 5 mm.
- 6.1.5 *Temperature Controlled Jacket*—Sufficient for maintaining measurement chamber temperature within ±0.06 °C.
- 6.1.6 Temperature Measuring Device—Industrial platinum resistance thermometer (IPRT) or equivalent sensor with a maximum permissible error of  $\pm 0.02$  °C. It is recommended, that the temperature measuring device be verified with an independent, calibrated temperature probe at the test temperature.
  - 6.2 Temperature Regulation System:
- 6.2.1 Any liquid bath or thermoelectric means for regulating the jacket temperature.
- 6.2.2 The temperature control must be such that the temperature of the measurement chamber is held within  $\pm 0.06$  °C of the desired measurement temperature.
- 6.3 Sample Introduction Mechanism—A syringe, micropipette, or flow-through adapter for introducing between 3.2 mL and 5 mL, inclusive by pressure, into the measurement chamber.

#### 7. Reagents and Materials

7.1 Certified viscosity reference standards shall be certified by a laboratory that has been shown to meet the requirements of ISO/EC 17025 by independent assessment. Viscosity stan-

dards shall be traceable to master viscometer procedures described in Practice D2162.

- 7.2 The uncertainty of the certified viscosity reference standard shall be stated for each certified value (k = 2, 95 % confidence). See ISO/EC 17025 or NIST TN 1297.
- 7.2.1 The certified viscosity reference should have a published viscosity in accordance with Test Method D445 or equivalent means that is close to that of the liquids being tested at the test temperature. For example, if intended measurements are to be made from 5 mPa·s to 25 mPa·s at 100 °C, then a reference oil viscosity of 15 mPa·s at 100 °C would be appropriate.
- 7.3 Cleaning solvents miscible with the sample and chemically compatible with the wetted viscometer components (such as alcohols, toluene, etc.). These wetted components are typically 316L and 430 Stainless Steel.
- 2.7.4 Quality control (QC) liquid similar to 7.1, but with viscosity values internally certified as noted in 12.2.

#### 8. Sampling, Samples, and Test Units

8.1 Ensure that the sample is homogenous. Engine sampling is generally specified in the test method, for example Test Method D5967. When applicable, refer to Practice D4057 (manual) or Practice D4177 (automatic) for proper sampling techniques.

### 9. Preparation of Apparatus

- 9.1 Place the viscometer on a stable surface.
- 9.2 Select the viscosity output units. If kinematic viscosity is selected, some apparatus will internally determine density to derive kinematic viscosity. Otherwise, enter the known density and operate the unit according to the procedure in Section 13.
- 9.3 Verify calibration accuracy by testing a reference standard or QC liquid at the test temperature. Follow the procedure in Section 13.

#### 10. Calibration and Standardization

10.1 Calibrate according to manufacturer's instructions to obtain a calibration curve (using two test liquids with referenced viscosity values near, but within, the extremes of the piston range being used).




FIG. 2 Cross Sectional View of Measurement Chamber

10.2 Certified Viscosity Standards may be used as confirmatory checks on the procedure in the laboratory. This procedure is outlined in Section 13. If the dynamic viscosity result, at the calibration test point, does not agree with the certified value within the limits of precision in Section 16, each step in the procedure should be rechecked, as well as the temperature measuring device and viscometer calibration, to locate the source of error. If the source is not detected, consult the manufacturer.

#### 11. Sample Conditioning

- 11.1 Shake all new and used oil samples using the following procedure.
  - 11.1.1 Ensure cap is tight on the container.
- 11.1.2 Shake vigorously by hand for 30 s. Wait 10 s, or longer if needed, for air bubbles to dissipate.
- 11.1.3 A specimen of the sample shall be taken by pipette, pouring or pumping. Suspected nonhomogeneous samples must be conveyed for analysis promptly following the shaking and dissipation procedure of step 11.1.2.

#### 12. Quality Control/Quality Assurance (QC/QA)

- 12.1 Confirm proper performance of the instrument and the test procedure by analyzing reference oil as QC sample.
- 12.2 If suitable reference oil is not available, prepare a QC sample by replicate analyses of a batch of oil sample. Then statistically analyze the data to assign a mean value and uncertainty limit to the sample.
- 12.3 When QC/QA protocols are already established in the testing facility, these may be used to confirm the reliability of the test result.
- 12.4 When there is no QC/QA protocol established in the testing facility, guidance may be obtained from Practice D6792.

#### 13. Procedure

- 13.1 Verify or set the temperature control settings, as tested with the control standard, so the viscometer temperature reads the desired set point  $\pm 0.06$  °C while the piston is in motion.
- 13.2 Remove the piston and clean the specimen from the measurement chamber as described in the viscometer manual.
- 13.3 Load the measurement chamber with sample using the volume listed in Table 1 related to the piston size being used for the viscosity range anticipated. To minimize contamination, and if sample volume allows, pre-wet the chamber and piston with the sample material and dry wipe with a lint free cloth.
- 13.4 Load the measurement chamber with a clean piston. The piston size should be selected such that the measured viscosity is between the minimum and maximum viscosity values listed in Table 1. If the reported result is outside of this range, the measurement shall be repeated using the appropriate piston size.
- 13.5 Start the reporting software in accordance with manufacturer's instructions, which in turn will:
- 13.5.1 Allow the sample to equilibrate in the measurement chamber for at least 2 min while the piston is oscillating.
- 13.5.2 Ensure that temperature stability is within the criterion set in 6.2.
- 13.5.3 Measure the upward and downward piston travel times for each cycle and compute viscosity until the standard deviation as percent of mean, over the previous 20 dynamic viscosity computations, is less than 0.6 %.
- 13.6 Record the average result from the last 20 computations.

#### 14. Calculation and Interpretation

14.1 The calculation of dynamic viscosity and kinematic viscosity are computed and displayed automatically by the apparatus.