Designation: D5724 - 16 (Reapproved 2021)

# Standard Specification for Gouache Paints<sup>1</sup>

This standard is issued under the fixed designation D5724; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon  $(\varepsilon)$  indicates an editorial change since the last revision or reapproval.

#### 1. Scope

- 1.1 This specification establishes requirements for composition, physical properties, performance, and labeling of gouache paints.
- 1.2 This specification covers pigments, vehicles, and additives. Requirements are included for pigment identification, lightfastness, and consistency.
- 1.3 Table 1 lists some pigments meeting the lightfastness requirements in this specification. In order to identify other pigments that meet these requirements, instructions are given for test specimen preparation. Test methods for determining relative lightfastness are referenced.
- 1.4 The values stated in SI units are to be regarded as the standard. The values given in parentheses are for information only.
- 1.5 This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.

### 2. Referenced Documents alog/standards/sist/3

2.1 ASTM Standards:<sup>2</sup>

D185 Test Methods for Coarse Particles in Pigments

D279 Test Methods for Bleeding of Pigments

D476 Classification for Dry Pigmentary Titanium Dioxide Products

D823 Practices for Producing Films of Uniform Thickness of Paint, Coatings and Related Products on Test Panels

D1210 Test Method for Fineness of Dispersion of Pigment-Vehicle Systems by Hegman-Type Gage

D1535 Practice for Specifying Color by the Munsell System

D4236 Practice for Labeling Art Materials for Chronic Health Hazards

D4303 Test Methods for Lightfastness of Colorants Used in Artists' Materials

E284 Terminology of Appearance

2.2 Other Documents:

Colour Index<sup>3</sup>

#### 3. Terminology

- 3.1 Definitions:
- 3.1.1 *colour index name, n*—consists of the category (type of dye or pigment), general hue, and an assigned number given to a colorant in the Colour Index<sup>3</sup> as an international identification system.
- 3.1.1.1 *Discussion*—For example, the Colour Index Name of one phthalocyanine blue pigment is Pigment Blue 15 (PB 15).
- 3.1.2 *colour index number, n*—a five-digit number given in the Colour Index that describes the chemical constitution of a colorant.
- 3.1.2.1 *Discussion*—For example, the Colour Index Number of one phthalocyanine blue pigment is 74160.
  - 3.2 Definitions of Terms Specific to This Standard:
- 3.2.1 gouache paint, n—a pigment dispersion in a water soluble gum/resin vehicle that dries water soluble and is formulated primarily for relatively opaque and matte applications.
- 3.3 Appearance terms used in this standard are defined in Terminology E284.

#### 4. Significance and Use

- 4.1 This specification establishes quality requirements and provides a basis for common understanding among producers, distributors, and users.
- 4.2 It is not intended that all paints meeting the requirements be identical nor of uniform excellence in all respects. Variations in manufacture, not covered by this specification, may cause some artists to prefer one brand over another, either of which may be acceptable under this specification.

<sup>&</sup>lt;sup>1</sup> This specification is under the jurisdiction of ASTM Committee D01 on Paint and Related Coatings, Materials, and Applications and is the direct responsibility of Subcommittee D01.57 on Artist Paints and Related Materials.

Current edition approved Nov. 1, 2021. Published November 2021. Originally approved in 1995. Last previous edition approved in 2016 as D5724-16. DOI: 10.1520/D5724-16R21.

<sup>&</sup>lt;sup>2</sup> For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

<sup>&</sup>lt;sup>3</sup> The Society of Dyers and Colourists, *Colour Index*, 3rd ed., 5 volumes and revisions, Available from the American Association of Textile Chemists and Colorists, P.O. Box 12215, Research Triangle Park, NC 27709.



# TABLE 1 Suitable Pigments List

1,012 1 011	Note 1—Underlined information in the table and lightfastness category shall be included on every label.					
KEY						
Lightfastness Ca	0 ,					
Lightfastness I Excellent Lightfastness						
•	II Very Good Lightfastness sed in Colour Index Names:					
PB	Pigment Blue					
PBk	Pigment Black					
PBr	Pigment Brown					
PG	Pigment Green					
PO	Pigment Orange					
PR	Pigment Red					
PV	Pigment Violet					
PW PY		Pigment White				
AR	Pigment Yellow Acid Red					
BR	Basic Red					
	ns in Parenthesis:					
(CC)		s may contain up to 15 % barium sulfate for color control.				
		ain a much higher content amount of barium sulfate.				
(DL)	May darken in strong light					
(LF)	Lightfast type					
(NA)	Colour index name or number no	ot assigned				
(RS) (BS)	Red shade Blue shade					
(SM)	Sensitive to moisture					
(SS)	Sensitive to hydrogen sulfide					
(OP)	Opaque type					
Color Index	Ligthfastness Category	Common Name and Chemical Class	Color Index Number			
Name		YELLOWS				
PY 3	I Arvlic	ide Yellow IOG, with option of adding the name Hansa Yellow Light, arylide yellow	11710			
PY 6		ide Yellow, arylide yellow	11670			
PY 35		mium (hue designation), concentrated cadmium zinc sulfide (CC), (SM)	77205			
PY 37	I Cadn	mium (hue designation), concentrated cadmium sulfide (CC), (SM)	77199			
PY 42		s Yellow or Iron Oxide Yellow, synthetic hydrated iron oxide	77492			
PY 43		ow Ochre, natural hydrated iron oxide	77492			
PY 53		el Titanate Yellow, oxides of nickel, antimony and titanium	77788			
PY 65						
		ide Yellow RN, with option of adding Hanas Yellow RN, arylide yellow	11740			
PY 74 2GX70	II <u>Aryli</u> c	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)	11741			
PY 74 2GX70 PY 109	II Arylic I Isoine	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)	11741 NA			
PY 74 2GX70 PY 109 PY 110	II Arylic I Isoine I Isoine	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)  adolinone Yellow G, tetrachroloisoindolinone  adolinone Yellow R, tetrachroloisoindolinone	11741 NA 56280			
PY 74 2GX70 PY 109 PY 110 PY 139	II Arylic I Isoine I Isoine I Isoine I Isoine	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)  adolinone Yellow G, tetrachroloisoindolinone  adolinone Yellow R, tetrachroloisoindolinone  adoline Yellow, isoindoline	11741 NA 56280 NA			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)  adolinone Yellow G, tetrachroloisoindolinone  adolinone Yellow R, tetrachroloisoindolinone  adoline Yellow, isoindoline  ylide Yellow, diarylide yellow	11741 NA 56280 NA 21104			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170	II Arylic I Isoine I Isoine I Isoine I Isoine I Diary urds. iteh.ai/catalog/stand	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP)  adolinone Yellow G, tetrachroloisoindolinone  adolinone Yellow R, tetrachroloisoindolinone  adoline Yellow, isoindoline	11741 NA 56280 NA 21104			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary urds.iteh.ai/catalog/stand I Cadn	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  Traniline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide	11741 NA 56280 NA 21104 d5724-162021			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary II Diary II Cadn I Gan I Benz	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow B, tetrachroloisoindolinone indoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  Itaniline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane	11741 NA 56280 NA 21104 05724-162021 12075 77202 11780			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 PO 5 PO 20 PO 36 PO 43	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary II Cadn I Cadn I Benz I Perin	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow R, tetrachroloisoindolinone indoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  Diaryline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL)	11741 NA 56280 NA 21104 15724-162021 12075 77202 11780 71105			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary II Cadn I Cadn I Benz I Perin	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow R, tetrachroloisoindolinone indoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  Diarylide Yellow Sulfo-selenide  Traniline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide  Zimidazolone (hue designation) HL, benzimidazolane  none Orange, perinone (DL)  cle Orange, Pyrrolopyrrol	11741 NA 56280 NA 21104 05724-162021 12075 77202 11780			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73	II Arylic I Isoine I Isoine I Isoine II Diary II Dinitr I Cadn I Benz I Perin II Pyrro	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS	11741 NA 56280 NA 21104 15724-162021 12075 77202 11780 71105 NA			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalog/stand I Cadn I Benz I Perin II Pyrro	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  Traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73 PR 5 PR 9	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds.iteh.av/catalog/stand I Cadn I Benz I Perin II Pyrrc II Naph II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline idoline Yellow, diarylide yellow  ORANGES  ORANGES  draniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL	11741 NA 56280 NA 21104 15774-162021 12075 77202 11780 71105 NA			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14'	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds.iteh.ai/catalog/stand I Cadn I Benz I Perin II Pyrro II Naph II Naph II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  Itaniline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-D, naphthol AS-D	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds.iteh.ai/catalog/stand I Benz I Benz I Perin II Pyrrc II Naph II Naph II Naph II Naph II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow, isoindoline ylide Yellow, diarylide yellow ORANGES  ORANGES  Traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-D, naphthol AS-D s Red or Iron Oxide Red, synthetic iron oxide	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14'	II Arylic I Isoine I Isoine I Isoine I Isoine I Isoine I Isoine I I Diary I Cadn I Benz I Perin II Pyrro II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  Itaniline Orange, dinitraniline (SM)  mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-D, naphthol AS-D	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalog/stand I Cadn I Benz I Perin II Pyrro II Naph II Cadn II Gadn	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow R, tetrachroloisoindolinone indoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES	11741 NA 56280 NA 21104 15724-162021 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 108 PR 113 PR 122 PR 170 F3RK-70	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalo g/stand I Cadn I Benz I Perin II Pyrro II Naph II Naph II Naph II Naph II Naph II Cadn I Gadn I Gadn I Gadn I Gadn I Gadn I Gadn	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline idoline Yellow, isoindoline idoline Yellow, diarylide yellow  ORANGES  COPANGES	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalo g/stand I Cadn I Benz I Perin II Pyrrc II Naph II Naph II Naph II Naph II Cadn I Gadn I Gadn I Gadn I Saph I Hars I Cadn I Naph I Naph I Naph I Cadn I Cadn I Cadn I Cadn I Quine I Naph I Naph I Naph I Naph I Naph I Naph I Raph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow R, tetrachroloisoindolinone indolinone Yellow, isoindoline indoline Yellow, isoindoline indoline Yellow, diarylide yellow  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  ORANGES  Itini Dorange, dinitraniline (SM) inium (hue designation), concentrated cadmium sulfo-selenide iniumium (hue designation) HL, benzimidazolane inone Orange, perinone (DL) indel Orange, Pyrrolopyrrol  REDS  Inthol ITR, naphthol ITR Inthol AS-OL, naphthol AS-OL Inthol AS-D, naphthol AS-OL Inthol AS-D, naphthol AS-D Is Red or Iron Oxide Red, synthetic iron oxide inium (hue designation), concentrated cadmium-seleno sulfide (CC) inium Vermillion Red Light, Medium or Deep, cadmium mercury sulfide (CC) inacridone (hue designation), γ quinacridone Inthol AEd, naphthol Carbamide (DL) Ithol AS, naphthol AS	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 108 PR 113 PR 122 PR 170 F3RK-70	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalo g/stand I Cadn I Benz I Perin II Pyrrc II Naph II Naph II Naph II Naph II Cadn I Gadn I Gadn I Gadn I Saph I Hars I Cadn I Naph I Naph I Naph I Cadn I Cadn I Cadn I Cadn I Quine I Naph I Naph I Naph I Naph I Naph I Naph I Raph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) indolinone Yellow G, tetrachroloisoindolinone indolinone Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 PY 139 PY 170 PO 5 PO 20 PO 36 PO 43 PO 73 PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary  ards. iteh. ai/catalog/stand I Cadn I Benz I Perin II Naph II Naph II Naph II Naph II Quine II Quine II Naph II Quine II Naph II Quine	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, R, tetrachroloisoindolinone idolinone Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  Acceptable  ORANGES  ORANG	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa/ PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19 PV 14	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. ai/catalog/stand I Cadn I Benz I Perin II Pyrro II Naph II Quina II Quina	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES 43 e-92 7-002 16 83 bec/astm- traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-OL, naphthol AS-OL s Red or Iron Oxide Red, synthetic iron oxide mium (hue designation), concentrated cadmium-seleno sulfide (CC) mium Vermilion Red Light, Medium or Deep, cadmium mercury sulfide (CC) nacridone (hue designation), γ quinacridone hthol Red, naphthol AS nacridone (hue designation), γ quinacridone red  VIOLETS alt Violet, cobalt phosphate	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12467 73900			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 102 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalog/stand I Cadn I Benz I Perin II Naph II Cadn II Quina II Quina II Quina II Quina II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  O	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa/ PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19 PV 14	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalog/stand I Cadn I Benz I Perin II Naph II Cadn II Quina II Quina II Quina II Quina II Naph	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES 43 e-92 7-002 16 83 bec/astm- traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-OL, naphthol AS-OL s Red or Iron Oxide Red, synthetic iron oxide mium (hue designation), concentrated cadmium-seleno sulfide (CC) mium Vermilion Red Light, Medium or Deep, cadmium mercury sulfide (CC) nacridone (hue designation), γ quinacridone hthol Red, naphthol AS nacridone (hue designation), γ quinacridone red  VIOLETS alt Violet, cobalt phosphate	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12467 73900			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 102 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary  ards. iteh. ai/catalog/stand I Cadn I Benz I Perin II Naph II Naph II Naph II Naph II Naph II Naph II Quine II Quine II Naph II Quine II Naph II Naph II Quine II Naph II Naph II Ouine II Naph II Naph II Ouine II Naph II Naph II Naph II Ouine II Naph II Naph II Ouine II Naph II Naph II Naph II Naph II Naph II Naph II Ouine II Naph II Ouine	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  O	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa/ PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23 PB 15 PB 17:1	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary  II Dinitr I Cadn I Benz I Perin II Naph II Naph II Naph II Naph II Naph II Naph II Quina II Quina II Quina II Quina II Quina II Quina II Coba II Quina II Quina II Quina II Quina II Phtha II Dioxa	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES 43e-92a7-b02166831bcc/astm- traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-D, naphthol AS-D s Red or Iron Oxide Red, synthetic iron oxide mium (hue designation), concentrated cadmium-seleno sulfide (CC) mium Vermilion Red Light, Medium or Deep, cadmium mercury sulfide (CC) nacridone (hue designation), γ quinacridone hthol Red, naphthol AS nacridone (hue designation), γ quinacridone red  VIOLETS alt Violet, cobalt phosphate nacridone (hue designation), carbazole dioxazine  BLUES nalocyanine Blue, or Pthalo Blue, copper phthalocyanine nalocyanine Blue Lake, or Pthalo Blue Lake, trisulfonated copper phthalocyanine	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23  PB 15 PB 17:1 PB 27	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalog/stand I Cadn I Benz I Perin II Naph II Naph II Naph II Naph II Naph II Naph I Cadn I Mars I Cadn I Mars I Cadn I Mars I Cadn I Mars I Cadn I Quine I Phthe I Phthe I Phthe I Pruse	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, R, tetrachroloisoindolinone idoline Yellow, isoindoline yiide Yellow, diarylide yellow  ORANGES 43e-92a7-b02166831bcc/astm- traniline Orange, dinitraniline (SM) mium (hue designation), concentrated cadmium sulfo-selenide zimidazolone (hue designation) HL, benzimidazolane none Orange, perinone (DL) cle Orange, Pyrrolopyrrol  REDS  hthol ITR, naphthol ITR hthol AS-OL, naphthol AS-OL hthol AS-D, naphthol AS-OL s Red or Iron Oxide Red, synthetic iron oxide mium (hue designation), concentrated cadmium-seleno sulfide (CC) mium Vermilion Red Light, Medium or Deep, cadmium mercury sulfide (CC) nacridone (hue designation), γ quinacridone hthol Red, naphthol AS nacridone (hue designation), γ quinacridone red  VIOLETS alt Violet, cobalt phosphate nacridone (hue designation), carbazole dioxazine  BLUES nalocyanine Blue, or Pthalo Blue, copper phthalocyanine sain Blue, Milori Blue, alkali ferric ferrocyanide	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1 77510			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23  PB 15 PB 17:1 PB 27 PB 28	II Arylic I Isoine I Isoine I Isoine I Isoine I Isoine I Diary Irds. iteh. av catalog/stand I Benz I Perin II Pyrro II Naph II Dioxa	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, isoindoline yilide Yellow, diarylide yellow  ORANGES	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1 77510 77346			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23  PB 15 PB 17:1 PB 27	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary Irds. iteh. av catalo g/stand I Cadn I Benz I Perin II Pyrro II Naph II Ouine II Quine II Quine II Quine II Quine II Phtha II Pruss I Coba	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  OR	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1 77510			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 102 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23  PB 15 PB 17:1 PB 27 PB 28 PB 29	II Arylic I Isoine I Isoine I Isoine I Isoine I Isoine I Isoine I I Isoine I Isoine I	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow R, tetrachroloisoindolinone idoline Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow idoline i	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1 77510 77346 77007			
PY 74 2GX70 PY 109 PY 110 PY 139 PY 170 https://standa PO 5 PO 20 PO 36 PO 43 PO 73  PR 5 PR 9 PR 14' PR 101 PR 108 PR 113 PR 122 PR 170 F3RK-70 PR 188 PV 19  PV 14 PV 19 PV 23  PB 15 PB 17:1 PB 27 PB 28	II Arylic I Isoine I Isoine I Isoine I Isoine II Diary  ards. iteh. ai/catalog/stand I Cadn I Benz I Perin II Naph II Naph II Naph II Naph II Naph II Quine II Quine II Quine II Quine II Quine II Perin II Papt II Perin II Raph II Raph II Perin II Perin II Phtha II Pruss II Poba II Perins II Phtha II Pruss II Coba II Phtha II Phtha II Pruss II Coba II Phtha II Phtha II Phtha II Phtha II Pruss II Coba II Quine II Raph II Dioxa II Phtha II Phtha II Pruss II Coba II Raph II Phtha II Phtha II Phtha II Phtha II Pruss II Coba I Ultrau II Ultrau II Quina II Quina II Phtha II Phtha II Pruss II Raph II	ide Yellow 2GX70, Hansa Yellow 2GX70, arylide yellow (OP) idolinone Yellow G, tetrachroloisoindolinone idolinone Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, isoindoline ylide Yellow, diarylide yellow  ORANGES  OR	11741 NA 56280 NA 21104 12075 77202 11780 71105 NA 12490 12460 12380 77491 77202 77201 73915 12475 12467 73900 77360 73900 51319 74160 74200:1 77510 77346			

		GREENS	
<u>PG 7</u>	I	Phthalocyanine Green, or Phthalo Green, chlorinated copper phthalocyanine	74260
PG 17	I	Chromium Oxide Green, anhydrous chromium sesquioxide	77288
PG 18	I	<u>Viridian</u> , hydrous chromium sesquioxide	77289
PG 19	I	Cobalt Green, oxides of cobalt and zinc, or cobalt zincate	77335
PG 23	I	Green Earth, or Terre Verte, natural ferrous silicate containing magnesium and aluminum potassium silicates	77009
PG 36	1	Phthalocyanine Green or Phthalo Green, chlorinated and brominated copper phthalocyanine BROWNS <sup>A</sup>	e 74265
PBr 7	1	Burnt Sienna, calcined natural iron oxide	77492
PBr 7	i	Burnt Umber, calcined natural iron oxide containing manganese	77492
PBr 7	i	Raw Sienna, natural iron oxide	77492
PBr 7	i	Raw Umber, natural iron oxide containing manganese	77492
PBr 11	1	Magnesium Ferrite, synthetic iron oxide containing magnesium oxide	77495
PBr 24	1	Chrome Titanate Yellow, oxides of chrome, antimony and titanium	77310
PBr 25	1	Benzimidazolone Brown, monoazo benzimidazolone	12510
		BLACKS	
<u>PBk 1</u>	I	Jet Black, aniline black	50440
<u>PBk 6</u>	I	Lamp Black, nearly pure amorphous carbon	77266
PBk 7	I	Carbon Black, nearly pure amorphous carbon	77266
PBk 9	I	Ivory Black or Bone Black, amorphous carbon produced by charring animal bones WHITES <sup>B</sup>	77267
PW 4	1	Zinc White, zinc oxide with option of adding the name Chinese White	77947
PW 5	1	Lithopone, zinc sulfide coprecipitated with barium sulfate	77115
PW 6	I	Titanium White, titanium dioxide (rutile or anatase) with option of including some barium sulfate or zinc oxide	77891
<u>PW 7</u>	I	Zinc Sulfide, zinc sulfide	77975

<sup>&</sup>lt;sup>A</sup> Color Index Number 77491 can be used as an alternate to 77492 for PBr 7.

#### 5. Labeling Requirements

- 5.1 *Pigment(s) Identification:*
- 5.1.1 Every label shall include for each pigment contained in the paint the information underlined in Table 1, which includes the Common Name, Colour Index Name, and any additional terms necessary to identify the form of the pigment.
- 5.1.2 The complete pigment identification given in Table 1, which also includes the Colour Index Number and a simple chemical description, shall be given by the producer in an appropriate electronic version or printed publication. Manufacturers are encouraged to put this complete identification on the container label when label size permits.
- 5.1.3 The common name shall be placed on the front of the label and shall be the name of the paint except as described in 5.1.5 and 5.1.6. Other identification may be placed elsewhere on the container.
- 5.1.4 The colour index name may be spelled out in full or abbreviated depending on the size of the label. Example: Pigment Blue 15, or Pig. Blue 15 or PB 15.
- 5.1.5 Substituted Pigments—In the case of substituted pigments, except for those pigments listed in Table X5.1, the word "Hue" in equal size letters shall follow in the title, on the front of the tube, and immediately after the name of the pigment that has been simulated.
- 5.1.6 Proprietary names or optional names may be used provided the common name(s) given in Table 1 is listed with their Colour Index Names and the Lightfastness Category of the mixture somewhere on the label.
- 5.1.7 *Mixed Pigments*—Artists' paints containing more than one pigment comply with this specification if all colored pigments included in the mixture are on the suitable pigment list (Table 1) and provided the mixture itself has passed all other test requirements in this specification. The lightfastness category shall be that of the least lightfast pigment. This

lightfastness category may be changed if the mixture is tested in accordance with Test Methods D4303 and the results indicating a different category are submitted to ASTM Subcommittee D01.57 for evaluation.

- 5.1.8 Historical and Discontinued Pigments—Pigments that are either (1) primarily of a historical nature, or (2) have not been commercially manufactured for a minimum of 10 years or more, may be submitted to Subcommittee D01.57 for inclusion in Table X5.1.
- 5.1.8.1 The Common Name(s) of pigments in Table X5.1 may be used by substituted pigments without the designation of "Hue" in the title.
- 5.1.8.2 Paints using pigments listed in Table X5.1 may use the word "Genuine" in front of the title to differentiate them from substituted pigments.
  - 5.2 Provide on the Label:
  - 5.2.1 Identification of gum/resin used.
- 5.3 *Lightfastness*—The label shall contain the word "Lightfastness" followed by the appropriate rating, I, or II, as given for each pigment in Table 1, or else one of these corresponding icons (Fig. 1).
- 5.3.1 Lightfastness I pigments, when made into paint specimens as described in Section 7 and exposed, tested, and rated in accordance with Test Methods D4303, shall have a color difference ( $\Delta E^*$ ab) of 4 or less CIELAB units between the specimens measured before and after exposure.
- 5.3.2 Lightfastness II pigments, when made into paint specimens as described in Section 8 and exposed, tested, and rated in accordance with Test Methods D4303, shall have a color difference ( $\Delta E*ab$ ) of more than 4.0 but not more than 8.0 CIELAB units between the specimens measured before and after exposure.

<sup>&</sup>lt;sup>B</sup> Information on white pigments is given in Appendix X3.



FIG. 1 Lightfastness Rating

- 5.3.3 Pigments were placed in a lightfastness category on the basis of either known historical performance in art works or the ratings from four lightfastness tests conducted as described in Test Methods D4303. Results from further tests on these, or other pigments, are solicited by Subcommittee D01.57.
- 5.3.3.1 The lightfastness category of a pigment shall be changed if results from several further tests conducted in accordance with Test Methods D4303 and approved by ASTM Subcommittee D01.57, establish a different lightfastness category than the one given in Table 1.
- 5.3.3.2 Additional pigments shall be placed in Table 1 after they have been tested for lightfastness in accordance with Test Methods D4303 and the test results submitted to ASTM Subcommittee D01.57 for evaluation, provided the results demonstrate that the pigments have the lightfastness ratings required for Lightfastness I or Lightfastness II, as just described.
- 5.3.4 For information and to establish nomenclature, pigments in Lightfastness III, IV and V categories are given in Table X1.1 in Appendix X1, but are not to be used in paint conforming to this specification.
- 5.3.4.1 Lightfastness III pigments have a color difference of more than 8.0 but not more than 16.0 CIELAB units.
- 5.3.4.2 Lightfastness IV pigments have a color difference of more than 16.0 but not more than 24.0 CIELAB units.
- 5.3.4.3 Lightfastness V pigments have a color difference of more than 24.0 CIELAB units.
- 5.4 *Toxicity*—All products and labeling must conform to the Federal Hazardous Substances Act and to Practice D4236.
- 5.5 Statement of Conformance—"Conforms to ASTM Specification D5724." or "Conforms to ASTM D5724" or "Conforms to the quality requirements of ASTM D5724." This statement may be combined with other conformance statements, such as, "Conforms to the quality and health requirements of ASTM Specification D5724 and Practice D4236."
- 5.6 Address—Include on the label (1) the name and address of the manufacturer or importer, and (2) the country of manufacture.

#### 6. Quality Assurance for Gouache Paints

6.1 Conditions Not Covered in This Specification that Affect Gouache Paints:

- 6.1.1 *Substrate*—The effective pH of the paper used will affect the long-term color of the applied gouache paints.
- 6.1.2 *Environmental Conditions*—Factors such as temperature, humidity, air flow, and light conditions affect application properties, drying time, and adhesion.
- 6.1.3 *Storage*—With aging and elevated temperatures there may be a change in consistency and a discernible separation of vehicle.
  - 6.2 Vehicles—Only water soluble gums/resins shall be used.
- 6.3 *Pigments*—Pigments used in gouache paints shall be limited to those in Table 1. Their lightfastness rating shall be the numeral given in the same row.
- 6.4 Additives—Thickeners, preservatives, surfactants, and humectants may be used to achieve consistency, prevent microbe deterioration, and control application results.
- 6.5 *Inerts*—Inerts shall only be used to produce desirable working qualities.
- 6.6 Preparation of Sample—For paste and fluid paints, empty the contents of the previously unopened container onto a glass slab and mix thoroughly with a spatula to a homogeneous sample. For cake paints, take a piece of the cake on a glass slab and add water and mix until a homogeneous paint is formed.
- 6.7 Coarse Particles—Paints shall be free of oversized particles and shall form a uniform film. The maximum content of coarse particles shall be 1 % by weight, as determined by Test Methods D185.
- 6.8 Fineness of Dispersion—Determine the fineness of dispersion by Test Method D1210. For paste paint: on a glass plate, using a spatula, mix the paint with an equal volume of water until homogeneous. The maximum allowable grind reading is 1.5 mils (40 µm).
- 6.9 *Consistency*—Paints shall be smooth and easily solubilized with water to a homogeneous color.
- 6.10 Freeze-Thaw Stability—Using a freezer that has a temperature of -7°C (20°F) or lower, subject the paint to five freeze-thaw cycles. A freeze-thaw cycle shall consist of freezing the paint to a solid state (minimum of 18 h) and then thawing the paint to room temperature (minimum of 5 h). The paint shall then meet the requirements of 6.7, 6.8, and 6.9.
- 6.11 *Opacity (Hiding Power)*—Opacity rating requirements will be included in this specification as appropriate standards for opacity of paints are established. Test method described in Appendix X2 may be used to rate the opacity of a paint.
- 6.12 *Bleeding*—Bleeding rating requirements will be included in this specification as appropriate standards for bleeding of paints are established. Test Methods B of D279 may be used to determine the rate of bleeding of paints.
- 6.13 *Color Specification*—Color specification of each color by Munsell notation by Test Method D1535 may be given in an appropriate producer publication. Manufacturers are encouraged to put this specification on the container label when label size permits.

#### 7. Lightfastness Determination

7.1 If a pigment is not listed in Table 1, test specimens of a gouache paint containing the pigment shall be prepared. These test specimens shall be tested in conformance with the requirements for exposure and evaluation given in Test Methods D4303.

Note 1—A report of the results of these tests may be submitted to Subcommittee D01.57 for inclusion of the pigments in Table 1. The report shall include information on test conditions instruments used, and be accompanied by the test specimens (which will be returned).

- 7.2 Materials:
- 7.2.1 Substrates—White unlacquered paper.
- 7.2.2 Chalk Resistant Titanium Dioxide Gouache Paint.
- 7.2.3 Drawdown Bar, with 100-µm (4-mil) aperture.
- 7.3 Mixing Whites for Dilution of Colors:
- 7.3.1 Use a white containing the same gum/resin and additives as the formulation in which the pigment is incorporated if possible. The following basic composition has been found satisfactory:

	Weight,%
Gum arabic nonvolatile	6.4
Glycerine	4.6
Chalk resistant rutile titanium dioxide	66.7
(conforming to Type IV in Classification	
D476)	

- 7.4 Preparation of Test Paints:
- 7.4.1 The pigment to be tested may be milled in a soft paste consistency. If a prepared paint of known composition is available, it may be used for this test instead of preparing a standard.
- 7.4.2 Dilute the pigment paste or paint with the mixing white and draw-down on unlacquered paper until the spectro-photometric measurement of the dry paint shows from 35 to 45 % reflectance at the wavelength of maximum absorption for that pigment. The wavelength of maximum absorption is located at the point of lowest reflectance on the spectral curve between 420 and 620 nm. If using a tristimulus filter colorimeter, the lowest of the three filter readings is the region of maximum-absorption and the dilution should be adjusted so that a reading of 35 to 45 % reflectance is obtained with this filter. The diffuse white reference standard for all measurements should have an absolute reflectance between 97 and 100 %.
- 7.4.2.1 To obtain this reflectance, use the Kubelka-Munk Single Constant Method described in Appendix X4 or use a trial and error method.
- 7.4.2.2 For gouache paints, 80 g of white paint is needed since paint used in the drawdown is lost.
- 7.4.2.3 For the initial weight of white stated in 7.4.2.1 and 7.4.2.2, the weighing must be accurate to 0.05 g.
  - 7.5 Preparation of Test Panels:

- 7.5.1 Use an applicator with a 100-µm (4-mil)<sup>4</sup> aperture to make a drawdown at least 40 by 40 mm (1.5 by 1.5 in.), or other minimum size appropriate for the viewing area of the instrument used. Follow Method C in Practices D823, if a motor-driven applicator is available. Pressure of about 500 g on the applicator is recommended. For manual drawdowns follow Method E in Practices D823 when an auxiliary flatting bar is available. Return all recoverable paint to the batch to allow for repeat mixes and measurements.
- 7.5.2 When a motor-driven applicator or an auxiliary flatting bar is not available, tape the unlacquered paper on a smooth surface such as a piece of glass. Position the applicator just above the upper edge of the paper and place a pool of the diluted paint in front of it. Grasp the sides of the applicator with the fingers and pull it down the paper with a smooth uniform motion.
- 7.5.3 Prepare four specimen panels for each pigment under test. Two are used in the first lightfastness tests and two are retained in a dark place, one for visual comparisons with the exposed panels, and one in case a third test is needed to supplement results from the first two tests, as described in Test Method D4303.
- 7.5.4 Apply the test paints to the unlacquered paper as described in 7.5.1 or 7.5.2. The panels should be air dried for overnight.
- 7.5.5 Cut the test panel from a section of the drawdown that is uniformly colored. If a backing is needed to support the test panels, adhere the panels to an appropriate substrate, such as light weight posterboard.

#### 8. Exposure

8.1 Conduct exposure tests, calculate mean color difference, and assign pigments to lightfastness categories as described in Test Method D4303.<sup>5</sup>

## 9. Keywords - b02166831bee/astm-d5724-162021

9.1 gouache paints; lightfastness; quality requirements; test specimens

 $<sup>^4</sup>$  A 100- $\mu$ m (4-mil) applicator was found satisfactory to obtain complete hiding for the mixing white in 7.3.1. Deeper aperture shall be used if complete hiding is not obtained, however, too thick drawdowns may cause cracking of films during exposure.

SLightfastness categories in Table 1 and Table X1.1 were determined by Test Methods B and C in Test Methods D4303. In Test Method B using fluorescent lamps, irradiance was monitored continuously during tests by a pyranometer to a total irradiation of 1260 MJ/m². In Test Method C using an xenon arc exposure system, irradiance was monitored by a narrow bandpass radiometer at 340 nm. Calibration to a total radiation of 1260 MJ/m² was based on the solar radiation measurement in southern Florida by South Florida Test Service. The Florida data show that total spectral irradiation at 340 nm was 550.87 kJ/m² when the total irradiation reached 1260 MJ/m² in the outdoor underglass exposure rack. Accordingly, a xenon arc exposure test was carried out until the spectral irradiation at 340 nm reached 550.9 kJ/m².