INTERNATIONAL STANDARD

ISO 10048

> First edition 1991-11-01

Water quality — Determination of nitrogen — Catalytic digestion after reduction with Devarda's alloy

iTeh STANDARD PREVIEW

Qualité de l'eau — Dosage de l'azote — Minéralisation catalytique après réduction avec l'alliage de Devarda

<u>ISO 10048:1991</u> https://standards.iteh.ai/catalog/standards/sist/b835d9d9-820d-49c0-8ff2b7e62d980a9b/iso-10048-1991



Reference number ISO 10048:1991(E)

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75% of the member VIEW bodies casting a vote.

International Standard ISO 10048 was prepared by Technical Committee ISO/TC 147, Water quality, Sub-Committee SC 2, Physical, chemical, biochemical methods. https://standards.iteh.ai/catalog/standards/sist/b835d9d9-820d-49c0-8ff2-

bs://standards.iten.ai/catalog/standards/sist/b835d9d9-820d-49c0-8 b7e62d980a9b/iso-10048-1991

© ISO 1991

All rights reserved. No part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from the publisher.

International Organization for Standardization

Case Postale 56 • CH-1211 Genève 20 • Switzerland Printed in Switzerland

Water quality — Determination of nitrogen — Catalytic digestion after reduction with Devarda's alloy

1 Scope

1.1 Substance determined

This International Standard specifies a method for the determination of nitrogen present in water samples in the form of ammonium, nitrite, nitrate, and organic nitrogen compounds capable of conversion to ammonium in the conditions of the method.

1.2 Type of sample

Certain organic nitrogen compounds are known to be difficult to mineralize. Table 2 gives information on recovery of nitrogen from a variety of compounds.

2 Normative reference

The following standard contains provisions which, e of conversion method. DARP international Standard. At the time of publiof this International Standard. At the time of publication, the edition indicated was valid. All standards are subject to revision, and parties to agreements based on this International Standard are encour-

This method is applicable to the analysis of raw and 048:199 aged to investigate the possibility of applying the polluted waters. https://standards.itch.ai/catalog/standards/sism635t/recent(edition of the standard indicated below. b7e62d980a9b/iso-10(Members of IEC and ISO maintain registers of currently valid International Standards.

1.3 Range

A nitrogen content, ρ_N , of up to 200 mg/l may be determined. This range may be extended by using a smaller test portion.

1.4 Limit of detection¹⁾

Using a 50 ml test portion, the limit of detection is $\rho_N = 3$ mg/l when using the titrimetric method. Using the spectrometric method (see 7.2, note 2), the limit of detection is $\rho_N = 1$ mg/l.

1.5 Sensitivity

Using a 100 ml test portion, 1,0 ml of 0,02 mol/l hydrochloric acid is equivalent to $\rho_{\rm N} = 2,8$ mg/l.

1.6 Interferences

Falsely low results may be obtained if the mineralization procedure is over-prolonged. To minimize this effect, the procedure in 7.2 shall be carefully observed.

ISO 7150-1:1984, Water quality — Determination of ammonium — Part 1: Manual spectrometric method.

3 Principle

Reduction of oxidized nitrogen compounds to ammonium ions using Devarda's alloy. After evaporation almost to dryness, conversion of all nitrogen compounds to ammonium sulfate in the presence of concentrated sulfuric acid containing a high concentration of potassium sulfate, to raise the boiling point of the mixture, and in the presence of copper which acts as a catalyst.

Liberation of ammonia from the mixture by the addition of alkali and distillation into boric acid/indicator solution. Determination of ammonium in the distillate either by titration with standard acid or by spectrometry at a wavelength of 655 nm.

4 Reagents

During the analysis, use only reagents of recognized analytical grade, and ammonia-free water.

¹⁾ Information from Finland.

4.1 Hydrochloric acid, $\rho = 1,18$ g/ml.

4.2 Sulfuric acid, $\rho = 1.84$ g/ml.

WARNING — This reagent causes severe burns.

The highest purity sulfuric acid shall be used. Pay particular attention to the manufacturer's specification with respect to its nitrogen content and ensure that it is as low as possible.

4.3 Sodium hydroxide, approximately 300 g/l solution.

WARNING - This reagent causes severe burns.

Dissolve 320 g \pm 20 g of sodium hydroxide in about 800 ml of water. Cool to room temperature and dilute to 1 litre with water in a measuring cylinder.

Store in a polyethylene bottle.

4.4 Devarda's alloy [about 45 % (m/m) aluminium, 50 % (m/m) copper, and 5 % (m/m) zinc] powdered form.

the mark with water. Prepare this solution freshly for each batch of analyses.

4.10 Potassium standard stock solution, $\rho_N = 1.000 \text{ mg/l}.$

Dissolve 7,215 g \pm 0,001 g of potassium nitrate (KHNO₃) (previously dried at 105 °C for at least 2 h) in about 750 ml of water. Quantitatively transfer to a 1 000 ml one-mark volumetric flask and make up to the mark with water.

Store the solution in a glass bottle for not more than 2 months.

4.11 Potassium nitrate standard solution, $\rho_{\rm N} = 10 \text{ mg/l}.$

Pipette 5 ml of the standard stock solution (4.10) into a 500 ml one-mark volumetric flask and make up to the mark with water.

Store the solution in a glass bottle for not more than 1 month.

4.12 Anti-bumping granules.

iTeh STANDARD PREVIEW

Select, from commercially-available products, ma- 5 Apparatus terials as low as possible in nitrogen content Standards.itch.ai

Ordinary laboratory apparatus and the following

4.5 Potassium sulfate, (K₂SO₄).

). <u>ISO 10048:1991</u> https://standards.iteh.ai/catalog/standards/style/astrong_mantles_or_heating_block, capable of maintaining a temperature of at least 400°C. b7e62d980a9b/so-10048-1991

4.6 Boric acid/indicator solution.

4.6.1 Dissolve 0,10 g \pm 0,01 g of bromocresol green and 0,020 g \pm 0,005 g of methyl red in about 80 ml of ethanol and dilute to 100 ml with ethanol in a measuring cylinder.

4.6.2 Dissolve 20 g \pm 1 g of boric acid (H₃BO₃) in warm water. Cool to room temperature. Add 10 ml \pm 0,5 ml of indicator solution (see 4.6.1) and dilute to 1 litre with water in a measuring cylinder.

4.7 Hydrochloric acid standard volumetric solution, c(HCI) = 0.02 mol/I.

This solution shall be prepared by dilution of hydrochloric acid (4.1) followed by standardization by normal analytical procedures. Alternatively, use a commercial solution of guaranteed concentration.

4.8 Glycine solution, $\rho_{\rm N} = 1\ 000\ {\rm mg/l}$.

Dissolve 5,362 g \pm 0,002 g of glycine (H₂NCH₂COOH) in about 800 ml of water and dilute to 1 litre with water in a calibrated flask.

4.9 Glycine solution, $\rho_N = 10 \text{ mg/l}$.

Pipette 10 ml of glycine solution (4.8) into a 1 000 ml one-mark volumetric flask and make up to

5.2 Kjeldahl flasks, of capacity 100 ml to 250 ml, compatible with the heating mantles (5.1), or **digestion tubes**, of capacity 100 ml to 250 ml, fitting snugly into the heating block (5.1).

5.3 Distillation apparatus, incorporating an antisplash distillation head and a vertical condenser whose outlet may be submerged in the absorbent solution.

If the digestion vessels (5.2) are not suitable for direct attachment to the distillation apparatus, separate distillation flasks are necessary.

Alternatively, use an apparatus for steam distillation. Such apparatus is available commercially, but it shall be checked to ensure that ammonia is neither lost during distillation nor formed in the steam generated for distillation.

6 Sampling and samples

Collect laboratory samples in polyethylene or glass bottles. Analyse them as quickly as possible, or store them at between 2 °C and 5 °C until analysed.

NOTE 1 Acidification to pH 2 with sulfuric acid (4.2) may also be used as an aid to preservation, provided that

possible contamination of the acidified sample by absorption of any atmospheric ammonia is avoided.

7 Procedure

7.1 Blank test

Proceed as described in 7.2, but use about 50 ml of water instead of a test portion. Record the volume, V_2 , of hydrochloric acid (4.7) added.

7.2 Determination

WARNING — The mineralization procedure may evolve toxic sulfur dioxide gas. Hydrogen sulfide and/or hydrogen cyanide may also be liberated from polluted samples. The mineralization should therefore be carried out under an effective fume extraction system.

Pipette 50 ml of the laboratory sample (see note 2) into Kjeldahl flask (5.2). Add 4,0 ml \pm 0,1 ml of sulfuric acid (4.2), 0,20 g \pm 0,01 g of Devarda's alloy (4.4) and 2,00 g \pm 0,05 g of potassium sulfate (4.5).

Heat the distillation flask so that distillate collects at a rate of about 5 ml/min. Stop the distillation when about 30 ml have been collected. Titrate the distillate to a purple end-point, using the indicator already present in the receiving flask (4.6.2), with 0,02 mol/l hydrochloric acid (4.7) and record the volume added.

NOTES

2 The laboratory sample should have a nitrogen concentration of not more than 200 mg/l. Samples with higher concentrations than this shall be diluted with water prior to taking the 50 ml test portion.

3 Ammonium ions may be determined in the distillate by the spectrometric method described in ISO 7150-1. However, it is essential that the distillate should be collected in 1 % (V/V) hydrochloric acid instead of the boric acid/indicator solution. Clause 10 of ISO 7150-1:1984 gives details of this modification and of the modified calculation of result.

4 Where the flask is not compatible with the distillation apparatus, the contents should be transferred quantitatively to a suitable distillation flask, using a total volume of 10 ml \pm 2 ml of water and rinsing three times.

After at least 60 min, add a few anti-bumping **8** Expression of results granules (4.12) and boil the contents of the flask upder an effective fume extraction system. The volume of the contents will decrease as water is boiled away. ISO 10048:192 calculate the total nitrogen concentration, ρ_{N} , ex-

https://standards.iteh.ai/catalog/standards/s When white fumes begin to appear, insert a small/iso-1(filter funnel into the neck of the Kjeldahl flask in order to reduce evaporation. Do not allow the sample to dry. The temperature of the liquid shall not exceed 370 °C at this stage.

After white fume evolution has ended, periodically observe the mineralizate and, after it has become either colourless or light green in colour, continue heating for a further 60 min \pm 5 min.

After mineralization, allow the flask to cool to room temperature. Meanwhile, measure 20 ml \pm 2 ml of boric acid/indicator solution (4.6.2) into the receiving flask of the distillation apparatus (see note 3). Ensure that the delivery tip of the condenser is below the surface of the indicator solution.

Carefully add 10 ml \pm 1 ml of water to the mineralization flask. Then add, from a measuring cylinder, 25 ml of sodium hydroxide solution (4.3) and immediately attach the flask to the distillation apparatus (see note 3).

 $\rho_{\rm N} = \frac{V_1 - V_2}{V_0} \times c({\rm HCI}) \times 14,01 \times 1000$

pressed in milligrams per litre, using the formula

where

- V_0 is the volume, in millilitres, of the test portion (which is usually 50 ml, see 7.2);
- V_1 is the volume, in millilitres, of the standard volumetric hydrochloric acid used for titration (see 7.2);
- V_2 is the volume, in millilitres, of the standard volumetric hydrochloric acid used for the titration of the blank test (see 7.1);
- c(HCI) is the exact concentration, expressed in moles per litre, of the hydrochloric acid used for titration;
- 14,01 is the relative atomic mass of nitrogen.

Express the result either as the mass concentration of nitrogen $\rho_{\rm N}$, in milligrams per litre, or as the amount-of-substance concentration of nitrogen, $c_{\rm N}$, in micromoles per litre.

To convert $\rho_{\rm N}$ to $c_{\rm N}$, use the formula

 $c_{\rm N} = 71,4 \ \rho_{\rm N}$

8.2 Repeatability and reproducibility

Repeatability and reproducibility standard deviations have been determined as shown in table 1.

					Purine	20	82,2 - 108,01)
Table 1 — Repeatability and reproducibilitystandard deviations					Benzonitrile	20	1,5 - 8,31)
					Sodium glutamate	20	79,2 - 97,41)
C (1)	Nitrogen concen- trations	Standard deviations			4-Nitrophenol + nitrate	302)	85,2 — 97,4 ³⁾
Sample					4-Nitroaniline + nitrate	30 ²⁾	88,4 - 94,4 ³
			Repro-		Pyridine + nitrate	302)	77,6 - 87,93
		Repeatability	ducibility		Purine + nitrate	302)	82,2 - 94,33)
	PΝ	r	R			L h _ h	1
	mg/l	mg/l	mg/l		1) Range of recoveries fi laboratories.	rom between	three and eight
4-Nitrophenol	202)	0,12 — 0,48	0,36 A		2) 10 mg/l from nitrate	, 20 mg/l fra	m the organic
4-Nitroaniline	202)	0,12 — 1,98	0,78	DA	3) Results from two labo	ratories	
Pyridine	20 ²⁾	0,30 — 5,1	éstan	dar	siten al		
Purine	20 ²⁾	0,06 — 1,06	0,29)		
Sodium glutamate	202)	0,20 — 1,58	0,74	ISO 10	048:1991 Jards/sist/b835d0d0_820d_40s	0.887	
A	5,30 ³⁾	0,127	0,6577e62	1980a91	v/iso-10048-1991	0-0112-	
А	5,67 ³⁾	0,125	0,743				
В	6,994)	0,091	0,573				
В	7,984)	0,199	0,583				
В	21,8 ⁴⁾	0,37	0,93				
В	14,84)	0,45	0,92				
в	23,35)	0,93	3,18				
с	9,845)	0,57	1,93				
D	206 ⁵⁾	3,3	7,6				
1) A Pulp mill effluent, B Municipal effluent, C Surface water (river), D Industrial effluent.							
2) Data from trials carried out by 8 laboratories in Finland and Germany.							
3) Data from trials carried out by 14 laboratories in Finland.							
4) Data from trials carried out by 17 laboratories in Finland.							
5) Data from trials carried out by 12 laboratories in Germany.							
L				J			

Table 2 — Recovery of organic nitrogen

Compound

4-Nitrophenol

4-Nitroaniline

Pyridine

Nitrogen

concen-

tration

 $\rho_{\rm N}$

mg/l

20

20

20

Percentage

recovery

%

58,0 - 103,51)

90,8 - 133,81)

 $24,3 - 101,5^{1}$

9 Notes on procedure

As a check on the recovered nitrogen, carry out the procedure using the glycine solution (4.9) in the place of the laboratory sample. Similarly, check the recovery of nitrate using the nitrate solution (4.11) in the place of the laboratory sample. The nitrogen concentrations obtained for these solutions should each lie between 9,5 mg/l and 10,5 mg/l. If this is not so, scrutinize the procedure arrangement carefully in order to establish the reason for poor performance.

10 Test report

The test report shall include the following information:

- a) a reference to this International Standard;
- b) all information necessary for complete identification of the sample;
- c) details of the storage and preservation of the laboratory sample before analysis;
- d) a statement of the repeatability achieved;
- e) the results and the method of expression used;
- f) details of any other information relevant to this International Standard.

iTeh STANDARD PREVIEW (standards.iteh.ai)

<u>ISO 10048:1991</u> https://standards.iteh.ai/catalog/standards/sist/b835d9d9-820d-49c0-8ff2b7e62d980a9b/iso-10048-1991

iTeh STANDARD PREVIEW (standards.iteh.ai)

<u>ISO 10048:1991</u> https://standards.iteh.ai/catalog/standards/sist/b835d9d9-820d-49c0-8ff2b7e62d980a9b/iso-10048-1991

_

UDC 614.777:556.114:546.17:543.242

Descriptors: water, quality, chemical analysis, determination of content, nitrogen, Kjeldahl method.

Price based on 5 pages