



Designation: E 1881 – 97 (Reapproved 2002)

## Standard Guide for Cell Culture Analysis with SIMS<sup>1</sup>

This standard is issued under the fixed designation E 1881; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last approval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

### 1. Scope

1.1 This guide provides the Secondary Ion Mass Spectrometry (SIMS) analyst with a cryogenic method for analyzing individual tissue culture cells growing in vitro. This guide is suitable for frozen-hydrated and frozen-freeze-dried sample types. Included are procedures for correlating optical, laser scanning confocal and secondary electron microscopies to compliment SIMS analysis.

1.2 This guide is not suitable for cell cultures that do not attach to the substrate.

1.3 This guide is not suitable for any plastic embedded cell culture specimens.

1.4 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

### 2. Referenced Documents

2.1 *ASTM Standards:*

E 673 Terminology Related to Surface Analysis<sup>2</sup>

### 3. Terminology

3.1 *Definitions:*

3.1.1 See Terminology E 673 for definitions of terms used in SIMS.

### 4. Summary of Guide

4.1 This guide describes a cryogenic method of sample preparation for cell culture specimens for SIMS analysis. In brief, cell cultures are grown on a conducting substrate, such as silicon. When cells reach about 80 % confluency, they are fast frozen and fractured by using a sandwich method (1).<sup>3</sup> This allows freeze-fixation of cellular contents and removal of the EF-leaflet of the apical plasma membrane. Since this kind of fracture occurs in groups of cells growing together, fractured cells are easily recognized for optical, SEM and SIMS imaging.

4.2 By correlative laser scanning confocal microscopy and SIMS, the same frozen freeze-dried cell can be analyzed for organelle localization in relation to elemental content (2).

### 5. Significance and Use

5.1 The presence of cell growth medium complicates a direct analysis of cells with SIMS. Attempts to wash out the nutrient medium results in the exposure of cells to unphysiological reagents that may also alter their chemical composition. This obstacle is overcome by using a sandwich freeze-fracture method (1). This cryogenic method has provided a unique way of sampling individual cells in their native state for SIMS analysis.

5.2 The procedure described here has been successfully used for imaging Na<sup>+</sup> and K<sup>+</sup> ion transport (3), calcium alterations in stimulated cells (4,5), and localization of therapeutic drugs and isotopically labeled molecules in single cells (6). The frozen freeze-dried cells prepared according to this method have been checked for SIMS matrix effects (7). Ion image quantification has also been achieved in this sample type (8).

5.3 The procedure described here is amenable to a wide variety of cell cultures and provides a way for studying the response of individual cells for chemical alterations in the state of health and disease.

### 6. Apparatus

6.1 This guide can be used for the analysis of cell cultures with virtually any SIMS instrument.

6.2 A cold stage in the SIMS instrument is needed to analyze frozen-hydrated specimens (9).

### 7. Procedure

7.1 Cells are grown on silicon wafer pieces (approximately 1 cm<sup>2</sup> area) of any shape. Alternatively, high purity germanium wafer pieces are used for cell growth for studies involving the use of <sup>44</sup>Ca stable isotope. These substrates are nontoxic to cells and have been used for growing various cell lines (1,2,8). Sterilize the silicon or germanium pieces prior to cell seeding. After the cells reach about 80 % confluency, replace the nutrient growth medium with new medium containing 11  $\mu$ m polystyrene beads (approximately 50 000 beads per 100 mm plastic dish, see Ref (1) for details on size of the beads). These beads act as spacers during the sandwich-fracture technique. It

<sup>1</sup> This guide is under the jurisdiction of ASTM Committee E42 on Surface Analysis and is the direct responsibility of Subcommittee E 42.06 on SIMS.

Current edition approved May 10, 1997. Published July 1997.

<sup>2</sup> *Annual Book of ASTM Standards*, Vol 03.06.

<sup>3</sup> The boldface numbers in parentheses refer to a list of references at the end of this guide.