

SLOVENSKI STANDARD SIST EN 14880:2006

01-januar-2006

Dcjfý]bg_c`U_h]jbY`gbcj]`!`8c`c YjUb^Y`jgYVbcgh]`Ubcf[Ubg_Y[U`gi`ZUhU`j Ub]cbg_]\`dcjfý]bg_c`U_h]jb]\`gbcjY\`!`AYhcXU`dchYbW]caYhf]bY`h]hfUW]^Y`n idcfUVc`gj]b YjY`fDV½]cbcgY`Y_h]jbY`Y`Y_hfcXY

Surface active agents - Determination of inorganic sulphate content in anionic surface active agents - Potentiometric lead selective electrode titration method iTeh STANDARD PREVIEW

Grenzflächenaktive Stoffe - Bestimmung des Gehaltes an anorganischen Sulfaten in anionischen grenzflächenaktiven Stoffen - Potentiometrische Titration mit einer bleiselektiven Elektrode

https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93c-

dbf6cd335b53/sist-en-14880-2006

Agents de surface - Détermination de la teneur en sulfate inorganique dans les agents de surface anioniques - Méthode potentiométrique de titrage avec électrode a membrane sélective au plomb

Ta slovenski standard je istoveten z: EN 14880:2005

<u>ICS:</u>

71.100.40 Površinsko aktivna sredstva Surface active agents

SIST EN 14880:2006

en



iTeh STANDARD PREVIEW (standards.iteh.ai)

<u>SIST EN 14880:2006</u> https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93cdbfccd335b53/sist-en-14880-2006

SIST EN 14880:2006

EUROPEAN STANDARD NORME EUROPÉENNE EUROPÄISCHE NORM

EN 14880

August 2005

ICS 71.100.40

English Version

Surface active agents - Determination of inorganic sulphate content in anionic surface active agents - Potentiometric lead selective electrode titration method

Agents de surface - Détermination de la teneur en sulfate inorganique dans les agents de surface anioniques -Méthode potentiométrique de titrage avec électrode à membrane sélective au plomb Grenzflächenaktive Stoffe - Bestimmung des Gehaltes an anorganischen Sulfaten in anionischen grenzflächenaktiven Stoffen - Potentiometrische Titration mit einer bleiselektiven Elektrode

This European Standard was approved by CEN on 8 July 2005.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the Central Secretariat or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the Central Secretariat has the same status as the official versions.

CEN members are the national standards bodies of <u>Austria, Belgium, Cypr</u>us, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom, <u>abfocd335b53/sist-en-14880-2006</u>



EUROPEAN COMMITTEE FOR STANDARDIZATION COMITÉ EUROPÉEN DE NORMALISATION EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: rue de Stassart, 36 B-1050 Brussels

© 2005 CEN All rights of exploitation in any form and by any means reserved worldwide for CEN national Members. Ref. No. EN 14880:2005: E

SIST EN 14880:2006

EN 14880:2005 (E)

Contents

		Page
Foreword		
1	Scope	4
2	Normative references	4
3	Principle	4
4	Reagents	4
5	Apparatus	5
6	Sampling and preparation of the test solution	5
7	Procedure	5
8	Calculation and expression of results	7
9	Precision	7
10	Test report	7
Annex	A (informative) Titration apparatus settings ARD PREVIEW	8
Annex	B (informative) Results of inter laboratory test i.e. ai	10
Bibliog	jraphy	11
	<u>SIST EN 14880:2006</u> https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93c- dbf6cd335b53/sist-en-14880-2006	

Foreword

This European Standard (EN 14880:2005) has been prepared by Technical Committee CEN/TC 276 "Surface active agents", the secretariat of which is held by AFNOR.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by February 2006, and conflicting national standards shall be withdrawn at the latest by February 2006.

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom.

iTeh STANDARD PREVIEW (standards.iteh.ai)

<u>SIST EN 14880:2006</u> https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93cdbf6cd335b53/sist-en-14880-2006

1 Scope

This European Standard specifies a method for the determination of the content of inorganic sulphate in anionic surface active agents. It can be applied to α -olefin sulphonates and n-olefin sulphonates, alcohol sulphates, alcohol ether sulphates, alkyl benzene sulphonates and other alkyl sulphonates. This method is used also for the determination of inorganic sulphate in alkyl benzene sulphonates or other alkyl sulphonates in their acid form. This method also applies to deep-coloured samples.

NOTE Sulphate can be present as sulphuric acid, ionic salts of this acid or mixture of these.

2 Normative references

The following referenced documents are indispensable for the application of this European Standard. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN ISO 3696, Water for analytical laboratory use – Specification and test methods (ISO 3696:1987)

ISO 607, Surface active agents and detergents – Methods of sample division

3 Principle

iTeh STANDARD PREVIEW

Potentiometric titration of the test sample containing inorganic sulphate is performed in non-aqueous medium with a lead nitrate standard volumetric solution. The reaction is the following: **a**1

$$Pb(NO_3)_2 + SO_4^{2-} \rightarrow PbSO_4 \downarrow +2NO_3^{-}$$
 SIST EN 14880:2006
https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93c-

Lead sulphate precipitate is formed during the titration and the hon-aqueous medium decreases its solubility. Nitric acid is added to remove possible interference from carbonates. The endpoint is related to an increase in lead ion activity, as measured by a lead-selective electrode.

4 Reagents

During the analysis, unless otherwise specified, use only reagents of recognised analytical grade that have been checked in advance as to not interfere with the analytical results.

4.1 Water, complying with grade 3 as defined in EN ISO 3696.

NOTE If the water is purified via ion-exchange resins, ensure that no cationic or anionic species from the resins cause interference.

4.2 Nitric acid, solution $c(HNO_3) = 1 \text{ mol/l}$.

4.3 Lead nitrate standard volumetric solution, $c(Pb(NO_3)_2) = 0.01 \text{ mol/l}$.

Weigh 3,312 g of lead nitrate, previously dried at 150 °C in a vacuum oven, dissolve in water, by using a small beaker. Quantitatively transfer the solution into a 1 000 ml volumetric flask and add water up to the mark.

4.4 Lead nitrate standard volumetric solution, c(Pb(NO₃)₂) = 0,005 mol/l.

Weigh 1,656 g of lead nitrate, previously dried at 150 °C in a vacuum oven, dissolve in water, by using a small beaker. Quantitatively transfer the solution into a 1 000 ml volumetric flask and add water up to the mark.

4.5 Potassium sulphate solution, $c(K_2SO_4) = 0.01 \text{ mol/l}$.

Weigh approximately 1,74 g (*m*) of K_2SO_4 , previously dried in an oven at 180 °C, dissolve with water in a 1 000 ml volumetric flask, make up to the mark and homogenize.

Calculate the factor, *f*, of the solution, according to the Equation (1):

$$f = \frac{m}{1,7425}$$
 (1)

where

j

- *m* is the mass of potassium sulphate weighted, in grams;
- 1,742 5 is the mass, in grams, of 0,01 mole of potassium sulphate K_2SO_4 .

4.6 Propan-2-ol (C_3H_8O)

5 Apparatus

Ordinary laboratory apparatus and the following:

5.1 Automatic potentiometric titration apparatus, with drift-controlled data acquisition and dynamic titrimetric dosing equipped with a piston burette delivery system of 20 ml capacity.

- 5.2 Lead-selective electrode (Metrohm 6.0502.170¹⁾, or similar).
- 5.3 Reference electrode Ag/AgCI filled with a saturated sodium nitrate ethanolic solution.

SIST EN 14880:2006

https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93c-

6 Sampling and preparation of the test solution₈₈₀₋₂₀₀₆

6.1 Sampling

The laboratory sample shall be taken and stored in accordance with ISO 607.

6.2 Preparation of the test solution

For non-homogeneous products, it will be better to dissolve a large sample (for example 5 g) in a suitable solvent (for example propan-2-ol/water) and take an aliquot for the titration.

Exactly weigh, in a 150 ml beaker, a quantity of the sample (m_0) so as to consume a maximum of 10 ml of Pb(NO₃)₂ solution (4.3 or 4.4) during the titration (typically, the mass is in the range 0,3 g to 1 g).

7 Procedure

7.1 Standardisation of the lead nitrate standard volumetric solution

Pipette 5,00 ml of the potassium sulphate solution (4.5) into a 150 ml beaker. Add approximately 1 ml of nitric acid solution (4.2) and 100 ml of propan-2-ol (4.6).

¹⁾ Metrohm 6.0502.170 is the trade name of instruments supplied by Metrohm Itd. (CH-9101 Herisau, Switzerland). This information is given for the convenience of users of this European Standard and does not constitute an endorsement by CEN of these instruments named.

SIST EN 14880:2006

EN 14880:2005 (E)

Prepare the potentiometer (5.1) for operation, immerse the electrodes in the solution and stir by using the electromagnetic stirrer.

Titrate with the lead nitrate volumetric solution (4.3 or 4.4). Record the reagent consumption, V_2 , at the inflection point of the titration curve.

NOTE An example for instrument settings and a titration curve are given in Annex A.

Run the standardization in triplicates. The concentration of the lead nitrate standard volumetric solution, c_a , expressed in moles per litre is calculated in accordance with the following Equation (2):

$$c_a = \frac{V_1 \times c \times f}{V_2} \tag{2}$$

where

- V_1 is the volume of potassium sulphate solution (4.5), in millilitres;
- *c* is the concentration of the lead nitrate solution, in moles per litre;
- f is the factor of the potassium sulphate solution determined according to 4.5;
- V_2 is the average volume of lead nitrate standard volumetric solution used in the titration, in millilitres.

NOTE In alternative, a solution of sulphuric acid, $c(H_2SO_4) = 0.01 \text{ mol/l}$, can be used to standardize the lead nitrate solution. In this case the sulphuric acid solution can be prepared diluting in the suitable ratio a standard reagent.

7.2 Determination

SIST EN 14880:2006

https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93c-

If the expected sulphate content is lower than $(1\% m/m_3)$ use the $Pb(NO_3)_2 c(Pb(NO_3)_2) = 0,005 \text{ mol/l for the titration, otherwise use the <math>Pb(NO_3)_2$, $c(Pb(NO_3)_2) = 0,01 \text{ mol/l.}$

Exactly weigh, in a 150 ml beaker, a quantity of the sample (m_0) so as to consume a maximum of 10 ml of Pb(NO₃)₂ standard volumetric solution during the titration.

NOTE Weighing more than 1g of surfactant should be avoided because the electrical behaviour of the ion-selective electrode is strongly dependent on the surfactant concentration and could give bad potentiometric curves. Approximately 1 ml of nitric acid solution (4.2) and 100 ml of propan-2-ol (4.6) should be added.

The solution should be acidic (lower than pH 3) before titration. If 1 ml of nitric acid solution (4.2) is not sufficient to achieve this, more acid should be added.

Prepare the potentiometer (5.1) for operation, immerse the electrodes in the solution and stir by using the electromagnetic stirrer.

Titrate with the lead nitrate volumetric solution (4.3 or 4.4). Record the reagent consumption, *V*, at the inflection point of the titration curve

Clean the ion-selective electrode tip before each titration by using a soft cloth or paper together with propan-2-ol.

Calculation and expression of results 8

The content of inorganic sulphate in the sample, c_b , expressed in percent by mass of sodium sulphate (Na₂SO₄) is calculated according to the Equation (3):

$$c_b = \frac{c_a \times 142,04 \times V}{m_a \times 10} \tag{3}$$

where

- is the mass of the test sample, in grams; m_0
- is the concentration of the lead nitrate standard volumetric solution, in moles per litre; c_a
- Vis the average volume consumption of lead nitrate standard volumetric solution.

NOTE The value 142,04 is the molecular mass of Na₂SO₄ in grams.

Precision 9

9.1 Repeatability limit

iTeh STANDARD PREVIEW

The absolute difference between two independent single test results, obtained using the same method on identical test material in the same laboratory by the same operator using the same equipment within a short interval of time, will not exceed the repeatability limit, r, in more than 5 % of cases.

SIST EN 14880:2006

Precision data are given in Annex B. https://standards.iteh.ai/catalog/standards/sist/69585255-ffcc-446c-b93cdbf6cd335b53/sist-en-14880-2006

9.2 Reproducibility limit

The absolute difference between two independent single test results, obtained using the same method on identical test material in different laboratories by different operators using different equipment, will not exceed the reproducibility limit, R, in more than 5 % of cases.

Precision data are given in Annex B.

10 Test report

The test report shall include the following information:

- a) all necessary information for the complete identification of the sample;
- b) method used (a reference to this European Standard, i.e. EN 14880);
- C) test results:
- details of any operations not specified in this European Standard or in the standards to which reference is d) made, and any operations regarded as optional, as well as any incidents likely to have affected the results.