

# INTERNATIONAL STANDARD

# ISO 8988

Second edition  
1995-05-01

---

---

## Plastics — Phenolic resins — Determination of hexamethylenetetramine content — Kjeldahl method and perchloric acid method

iTeh STANDARD PREVIEW  
(standards.iteh.ai)

*Plastiques — Résines phénoliques — Détermination de la teneur en  
hexaméthylène-tétramine — Méthode Kjeldahl et méthode à l'acide*

<https://standards.iteh.ai/standards/sist/16312031-f73b-4644-a52e-43ba29cab6/iso-8988-1995>



Reference number  
ISO 8988:1995(E)

## Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

International Standard ISO 8988 was prepared by Technical Committee ISO/TC 61, *Plastics*, Subcommittee SC 12, *Thermosetting materials*.

This second edition cancels and replaces the first edition (ISO 8988:1989), of which it constitutes a minor revision.

© ISO 1995

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from the publisher.

International Organization for Standardization  
Case Postale 56 • CH-1211 Genève 20 • Switzerland

Printed in Switzerland

# Plastics — Phenolic resins — Determination of hexamethylenetetramine content — Kjeldahl method and perchloric acid method

## 1 Scope

This International Standard specifies two methods for the determination of the hexamethylenetetramine ("hexa") content of phenolic resins. The two methods are equivalent. The Kjeldahl method described in clause 3 is not applicable if there are other components containing nitrogen in the phenolic resin. The perchloric acid method described in clause 4 is only applicable if there are no other basic or acidic additives in the resin.

## 2 Normative reference

The following standard contains provisions which, through reference in this text, constitute provisions of this International Standard. At the time of publication, the edition indicated was valid. All standards are subject to revision, and parties to agreements based on this International Standard are encouraged to investigate the possibility of applying the most recent edition of the standard indicated below. Members of IEC and ISO maintain registers of currently valid International Standards.

ISO 385-1:1984, *Laboratory glassware — Burettes — Part 1: General requirements*.

## 3 Kjeldahl method

**WARNING — For safety reasons, the Kjeldahl determination must be carried out in a well ventilated hood.**

### 3.1 General

This clause specifies a method for the determination of total nitrogen, expressed as hexamethylenetetramine, in phenolic resins.

The method is applicable to hexamethylenetetramine contents of  $\geq 0,5\%$  (*m/m*).

### 3.2 Principle

The hexamethylenetetramine in a test portion is converted to ammonium bisulfate by decomposition in hot concentrated sulfuric acid and in the presence of a catalytic mixture.

The ammonium bisulfate is converted to sodium sulfate and ammonia by reaction with sodium hydroxide.

The ammonia is distilled off and collected in hydrochloric acid.

The excess hydrochloric acid is titrated with a standard volumetric solution of sodium hydroxide in the presence of an indicator.

### 3.3 Reagents

During the analysis, unless otherwise stated, use only reagents of recognized analytical grade, free of nitrogen, and only distilled water or water of equivalent purity.

**3.3.1 Sulfuric acid**, concentrated.

**3.3.2 Kjeldahl catalytic mixture**, consisting of 97 g of sodium sulfate decahydrate ( $\text{Na}_2\text{SO}_4 \cdot 10\text{H}_2\text{O}$ ), 1,5 g of copper sulfate pentahydrate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ) and 1,5 g of selenium (Se).

**3.3.3 Sodium hydroxide**, 30 % (*m/m*) solution.

**3.3.4 Hydrochloric acid**,  $c(\text{HCl}) = 0,10$  mol/l.

**3.3.5 Sodium hydroxide**, standard volumetric solution,  $c(\text{NaOH}) = 0,10 \text{ mol/l}$ .

**3.3.6 Mixed indicator**, solution.

Dissolve 60 mg of methyl red and 40 mg of methylene blue in 100 ml of ethanol.

### 3.4 Apparatus

Ordinary laboratory apparatus, plus the following:

**3.4.1 Kjeldahl flask**, capacity 250 ml or 300 ml, for the digestion.

**3.4.2 Distillation apparatus** (various models are available commercially).

**3.4.3 Burette**, capacity 50 ml, graduated in 0,1 ml steps, conforming with the requirements of ISO 385-1.

**3.4.4 Analytical balance**, accurate to 1 mg.

**3.4.5 Silicon carbide grits**, for use as anti-bumping granules.

### 3.5 Procedure

#### 3.5.1 Digestion

Into a Kjeldahl flask (3.4.1), weigh 1 g to 2 g of phenolic resin to the nearest 1 mg. Add 5 g of the catalytic mixture (3.3.2) and 25 ml of concentrated sulfuric acid (3.3.1). Heat carefully until the colour of the mixture being digested changes from black or amber to clear. When clear, increase the rate of heating until 5 min beyond the time of the colour change and, possibly, boiling. Allow the digested liquid to cool almost to room temperature, just short of solidification. Add carefully 100 ml of water and transfer the solution quantitatively into the flask of the distillation apparatus, rinsing with water. Add a few silicon carbide grits (3.4.5) to prevent bumping. Add 30 % ( $m/m$ ) NaOH solution (3.3.3) to this solution until an alkaline reaction is obtained. Then distill over the ammonia given off, together with water vapour, into a receiver containing 50 ml of hydrochloric acid (3.3.4). Continue the distillation until about 300 ml of water has been collected.

#### 3.5.2 Titration

When the distillation is completed, add a few drops of mixed indicator solution (3.3.6) to the contents of

the receiver and titrate the excess hydrochloric acid with sodium hydroxide solution (3.3.5), using the burette (3.4.3).

### 3.6 Expression of results

The hexamethylenetetramine content, expressed as a percentage by mass, is given by the formula

$$\frac{0,35(V_0 - V_1)}{m_0}$$

where

$V_0$  is the volume, in millilitres, of hydrochloric acid (3.3.4) in the receiver of the distillation apparatus;

$V_1$  is the volume, in millilitres, of sodium hydroxide solution (3.3.5) used in the back-titration;

$m_0$  is the mass, in grams, of the test portion.

### 3.7 Reproducibility

The results are reproducible to within 0,30 % ( $m/m$ ) hexamethylenetetramine.

### 3.8 Number of determinations

Carry out the determination in duplicate. If the results differ by more than 5 %, repeat the determination, again in duplicate. If not, calculate the arithmetic mean of the two individual results.

## 4 Perchloric acid method

### 4.1 General

This clause specifies a method for the determination of hexamethylenetetramine in phenolic resins by direct titration. The results of the determination may be affected by the presence of acidic or basic additives. In such cases, the use of the Kjeldahl method is recommended.

The method is applicable to hexamethylenetetramine contents of  $\geq 0,3 \%$ .

### 4.2 Principle

One of the tertiary amine groups of the hexamethylenetetramine in a test portion is determined by titration with perchloric acid.

### 4.3 Reagents

During the analysis, unless otherwise stated, use only reagents of recognized analytical grade, free of nitrogen, and only distilled water or water of equivalent purity.

**4.3.1 Hexamethylenetetramine**, dry.

**4.3.2 Acetone**.

**4.3.3 Perchloric acid**, 70 % (V/V) solution.

**WARNING — Perchloric acid is dangerous in the presence of organic matter since an explosion can occur if the perchloric acid is in excess.**

### 4.4 Apparatus

Ordinary laboratory apparatus, plus the following:

**4.4.1 Magnetic stirrer**.

**4.4.2 Automatic burette**, nominal volume at least 15 ml, graduated in 0,1 ml steps, with a stopcock made of polytetrafluoroethylene.

**4.4.3 Beakers**, capacity 100 ml.

**4.4.4 Graduated cylinder**, capacity 1 000 ml.

**4.4.5 Analytical balances**, accurate to 0,1 mg and 1 mg, respectively.

**4.4.6 pH-meter**.

### 4.5 Procedure

#### 4.5.1 Preparation and titration of solution of perchloric acid in acetone

Introduce 8 ml of perchloric acid solution (4.3.3) into the 1 000 ml graduated cylinder (4.4.4), and dilute to 1 000 ml with acetone (4.3.2).

Standardize the resulting solution with hexamethylenetetramine (4.3.1) as described below.

Weigh, to the nearest 0,1 mg, about 150 mg to 170 mg of hexamethylenetetramine (4.3.1) into a 100 ml beaker (4.4.3).

Add 30 ml to 40 ml of acetone (4.3.2), and titrate as described (4.5.2).

NOTE 1 Darkening of the solution will not affect the titration results.

The titre  $T$ , expressed in milligrams of hexamethylenetetramine per millilitre of solution, is given by the formula

$$\frac{m_1}{V_2}$$

where

$m_1$  is the mass, in milligrams, of hexamethylenetetramine;

$V_2$  is the volume, in millilitres, of perchloric acid solution needed to reduce the pH to just below zero.

#### 4.5.2 Titration

Into a 100 ml beaker weigh, to the nearest 1 mg, a quantity of phenolic resin equal to 100 times the titre determined as in 4.5.1, add 30 ml to 40 ml of acetone (4.3.2) and insert a magnetic stirring rod. Place the beaker on the magnetic stirrer (4.4.1). Insert the glass electrode of the pH-meter (4.4.6) and switch on the stirrer and the pH-meter. When the resin has dissolved, add the perchloric acid solution, prepared as in 4.5.1, slowly dropwise until the pH-value drops suddenly below zero. As the resin dissolves in acetone more rapidly than hexamethylenetetramine, the pH-reading may increase to above zero because residual hexamethylenetetramine may still be in the process of dissolving. Continue the titration until the pH-reading remains constant, slightly below zero.

### 4.6 Expression of results

The hexamethylenetetramine content, expressed as a percentage by mass, is given by the formula

$$\frac{V_2 T \times 100}{m_0}$$

where

$V_2$  is the volume, in millilitres, of perchloric acid solution used for the titration;

$T$  is the titre, expressed in milligrams of hexamethylenetetramine per millilitre, of the perchloric acid solution, as determined in 4.5.1;

$m_0$  is the mass, in milligrams, of the test portion.

### 4.7 Number of determinations

Carry out the determination in duplicate. If the results differ by more than 5 %, repeat the determination,

again in duplicate. If not, calculate the arithmetic mean of the two individual results.

## 5 Test report

The test report shall include the following information:

- a) a reference to this International Standard and the method used (Kjeldahl method or perchloric acid method);
- b) all details necessary for complete identification of the resin tested;
- c) the hexamethylenetetramine content, expressed as a percentage by mass, together with
  - the individual results,
  - the arithmetic mean;
- d) the date of the test.

# iTeh STANDARD PREVIEW (standards.iteh.ai)

ISO 8988:1995

<https://standards.iteh.ai/catalog/standards/sist/16312031-f73b-4644-a52e-43ba29cabc6/iso-8988-1995>

**iTeh STANDARD PREVIEW**  
This page intentionally left blank  
**(standards.iteh.ai)**

ISO 8988:1995

<https://standards.iteh.ai/catalog/standards/sist/16312031-f73b-4644-a52e-43ba29cab6/iso-8988-1995>

**iTeh STANDARD PREVIEW**  
**(standards.iteh.ai)**

ISO 8988:1995

<https://standards.iteh.ai/catalog/standards/sist/16312031-f73b-4644-a52e-43ba29cab6/iso-8988-1995>

---

---

**ICS 83.080.10**

**Descriptors:** plastics, thermosetting resins, phenoplasts, chemical analysis, determination of content, hexamethylenetetramine.

Price based on 4 pages

---

---