
Kakovost vode - Ugotavljanje prisotnosti in števila legionel – 2. del: Metoda neposredne membranske filtracije za vode z majhnim številom bakterij

Water quality - Detection and enumeration of Legionella - Part 2: Direct membrane filtration method for waters with low bacterial counts

**Water quality — Detection and
enumeration of *Legionella* —**

**Part 2:
Direct membrane filtration method for
waters with low bacterial counts**

Qualité de l'eau — Recherche et dénombrement des Legionella —

*Partie 2: Méthode par filtration directe sur membrane pour les eaux à
faible teneur en bactéries*



PDF disclaimer

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.

© ISO 2004

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

Contents

Page

Foreword.....	iv
1 Scope.....	1
2 Normative references	1
3 Terms and definitions.....	1
4 Safety.....	2
5 Principle	2
6 Culture media and reagents.....	2
7 Apparatus	6
8 Sampling	7
9 Procedure	7
10 Expression of results.....	9
11 Test report	9

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 11731-2 was prepared by Technical Committee ISO/TC 147, *Water quality*, Subcommittee SC 4, *Microbiological methods*.

ISO 11731 consists of the following parts, under the general title *Water quality — Detection and enumeration of Legionella*:

— *Part 2: Direct membrane filtration method for waters with low bacterial counts*

The general method will be the subject of a future Part 1 of ISO 11731.

Water quality — Detection and enumeration of *Legionella* —

Part 2:

Direct membrane filtration method for waters with low bacterial counts

WARNING — Persons using this part of ISO 11731 should be familiar with normal laboratory practice. This part of ISO 11731 does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to ensure compliance with any national regulatory conditions.

1 Scope

This part of ISO 11731 describes a monitoring method for the isolation and enumeration of *Legionella* organisms in water intended for human use (e.g. hot and cold water, water used for washing), for human consumption and for treated bathing waters (e.g. swimming pools). It is especially suitable for waters expected to contain low numbers of *Legionella*. As the growth of *Legionella* may be inhibited by overgrowth of other bacterial colonies on the membrane, the method is only suitable for waters containing low bacterial counts.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 3696:1987, *Water for analytical laboratory use — Specification and test methods*

ISO 8199:—¹⁾, *Water quality — General guidance on the enumeration of micro-organisms by culture*

ISO 11731:1998, *Water quality — Detection and enumeration of Legionella*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

3.1

Legionella

genus of Gram-negative bacteria normally capable of growth in no less than 2 days on buffered charcoal yeast extract agar containing L-cysteine and iron(III), and forming colonies, often white, purple to blue or lime green in colour

NOTE Some species fluoresce under long wavelength UV light. The colonies have a ground-glass appearance when viewed with a low power stereomicroscope. Growth does not occur in the absence of L-cysteine with the exception of *L. oakridgensis* and *L. spiritensis*. *L. oakridgensis* and *L. spiritensis* require L-cysteine and iron for primary isolation but can grow weakly in the absence of added L-cysteine thereafter.

1) To be published. (Revision of ISO 8199:1988)

4 Safety

The reagents used in this part of ISO 11731 should be subject to assessment in accordance with control substances hazardous to health.

Legionella species can be handled safely by experienced microbiologists on the open bench in a conventional microbiology laboratory conforming to containment level 2. Infection is caused by inhalation of the organism and it is advisable therefore to assess all techniques for their ability to produce aerosols. If in any doubt, carry out the work in a safety cabinet.

5 Principle

5.1 General

Bacteria, including *Legionella* organisms, in the water sample are concentrated by membrane filtration. After filtration, the filter is treated with acid buffer added directly into the funnel to reduce the growth of non-*Legionella* organisms. The filter is subsequently transferred onto a plate of agar medium selective for *Legionella* and incubated. Samples expected to contain sufficient numbers of *Legionella* need not be subjected to concentration prior to culture (9.1).

5.2 Enumeration

After incubation, morphologically characteristic colonies which form on the selective medium are regarded as presumptive *Legionella*.

5.3 Confirmation

Presumptive colonies are confirmed as *Legionella* organisms by subculture to demonstrate their growth requirement for L-cysteine and iron. Further biochemical and serological tests are needed for species identification. Species identification may not be considered necessary for routine monitoring but is indispensable in outbreak situations.

NOTE *L. pneumophila* serogroup 1 is the causative agent of most legionellosis cases and is therefore considered the most "critical" type of *Legionella* to be found in the water system. Since increasing numbers of cases of legionellosis caused by other serogroups of *L. pneumophila* and other *Legionella* species are being described, even the presence of other *Legionella* species in water is considered a potential risk.

6 Culture media and reagents

6.1 General

Use chemicals of analytical grade in the preparation of media and reagents unless otherwise stated. Alternatively, use commercially available dehydrated media and reagents. Prepare the media according to the manufacturer's instruction and add freshly prepared (or thaw the stored material at room temperature prior to use) selective agents or growth supplements at the concentrations recommended. Prepare media using glass distilled water or water of equivalent quality complying with ISO 3696:1987, Grade 3. Other grades of chemicals may be used providing they can be shown to produce the same results.