

Designation: B 472 - 04a

Standard Specification for Nickel Alloy Billets and Bars for Reforging¹

This standard is issued under the fixed designation B 472; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

This standard has been approved for use by agencies of the Department of Defense.

1. Scope*

- 1.1 This specification covers UNS N06002, UNS N06030, UNS N06035, UNS N06022, UNS N06200, UNS N06230, UNS N06600, UNS N06617, UNS N06625, UNS N08020, UNS N08026, UNS N08024, UNS N08120, UNS N08926, UNS N08367, UNS N10242, UNS N10276, UNS N10665, UNS N10675, UNS N12160, UNS R20033, UNS N06059, UNS N06686, UNS N10629, UNS N08031, UNS N06045, UNS N06025, and UNS R30556* billets and bars for reforging.
- 1.2 The values stated in inch-pound units are to be regarded as the standard. The values given in parentheses are for information only.
- 1.3 This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to become familiar with all hazards including those identified in the appropriate Material Safety Data Sheet for this product/material as provided by the manufacturer, to establish appropriate safety and health practices, and determine the applicability of regulatory limitations prior to use.

2. Referenced Documents

- 2.1 ASTM Standards: 2
- A 262 Practices for Detecting Susceptibility to Intergranular Attack in Austenitic Stainless Steels
- B 880 Specification for General Requirements for Chemical Check Analysis Limits for Nickel, Nickel Alloys and Cobalt Alloys
- E 1473 Test Methods for Chemical Analysis of Nickel, Cobalt, and High-Temperature Alloys

3. Terminology

- 3.1 Definitions of Terms Specific to This Standard:
- 3.1.1 *billets and bars*, *n*—terms billets and bars as used in this specification shall be understood as billets and bars for reforging.

4. Ordering Information

- 4.1 It is the responsibility of the purchaser to specify all requirements that are necessary for material ordered under this specification. Examples of such requirements include, but are not limited to, the following:
 - 4.1.1 Quantity (weight or number of pieces),
 - 4.1.2 Name of material or UNS number,
 - 4.1.3 Form (bar or billet),
 - 4.1.4 Dimensions,
 - 4.1.5 ASTM designation and year of issue,
 - 4.1.6 Inspection (12.1),
- 4.1.7 Certification—State if certification or a report of test results is required (Section 14),
 - 4.1.8 Supplementary requirements, if any, and
 - 4.1.9 If possible, the intended end use.

Note 1—A typical ordering description is as follows: 10 000 lb (4536 kg), UNS N08020, forging bar, 41/4 in. (107.95 mm) round, Specification B 472.

5. Materials and Manufacture

- 5.1 The products shall be hot worked from ingots by rolling, forging, extruding, hammering, or pressing.
- 5.2 The products may be conditioned by chipping, grinding, or machining to remove injurious surface defects provided the depth of conditioning does not exceed that which will affect the surface condition or dimensions of the article to be forged from the bar or billet.

6. Chemical Composition

- 6.1 The material shall conform to the requirements as to chemical composition prescribed in Table 1.
- 6.2 If a product (check) analysis is performed by the purchaser, the material shall conform to the requirements specified in Table 1 subject to the permissible tolerances in Specification B 880.

¹ This specification is under the jurisdiction of ASTM Committee B02 on Nonferrous Metals and Alloys and is the direct responsibility of Subcommittee B02.07 on Refined Nickel and Cobalt and Their Alloys.

Current edition approved Oct. 1, 2004. Published October 2004. Originally approved in 1968. Last previous edition approved in 2004 as B 472 - 04.

^{*} New designation established in accordance with ASTM E 527 and SAE J1086, Practice for Numbering Metals and Alloys (UNS).

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.



TABLE 1 Chemical Requirements

Element —		Composition, %							
		UNS N08026	0 UNS N	108020	UNS N08024	UNS N08367	UNS N	108926	UNS R20033
Carbon, max		0.03	0.07		0.03	0.030	0.020		0.015
Manganese, m		1.00	2.00		1.00	2.00	2.00		2.0
Phosphorus, m	nax	0.03	0.045		0.035	0.040	0.03		0.02
Sulfur, max		0.03	0.035		0.035	0.030	0.01		0.01
Silicon, max		0.50 1.00			0.50	1.00	0.50		0.50
Nickel		33.00–37.20 32.00–38			35.00-40.00	23.50-25.50	24.00-		30.0–33.0
Chromium				-21.00	22.50–25.00	20.00-22.00	19.00-		31.0–35.0
Molybdenum				3.00	3.50-5.00	6.00-7.00	6.0–7.		0.50-2.0
Copper		2.00-4.00	3.00-4		0.50-1.50	0.75 max	0.5–1.	5	0.30-1.20
Columbium (NI	b) +		8 × ca	ırbon–1.00	0.15-0.35				
tantalum		0.40.040				0.40.005	0.45.0		0.05.000
Nitrogen Iron ^A		0.10-0.16				0.18–0.25	0.15-0		0.35-0.60
Iron''		remainder	remair	naer	remainder	remainder	balanc	e	balance
Element						nposition, %			
Element		UNS N06030	UNS N06022	UNS N06200	UNS N10276	UNS N10665	UNS N10675	UNS N06002	UNS N06230
0									
Carbon, max		0.03	0.015	0.010	0.010	0.02	0.01	0.05-0.15	
Manganese, max		1.5	0.50	0.50	1.0	1.0	3.0	1.00	0.30-1.00
Phosphorous, ma	X	0.04	0.02	0.025	0.04	0.04	0.030	0.04	0.030
Sulfur, max		0.02	0.02	0.010	0.03	0.03	0.010	0.03	0.015
Silicon, max		0.8	0.08	0.08	0.08	0.10	0.10	1.00	0.25-0.75
Nickel Chromium		remain				remainder	65.0 min	remainde	
Molybdenum		28.0-31 4.0-6.0				1.0 max 26.0-30.0	1.0-3.0 27.0-32.0	20.5-23.0 8.0-10.0	1.0-3.00
•		4.0-6.0 1.0-2.4		5 15.0-17. 1.3-1.9					
Copper Columbium		0.30-1.					0.20		•••
Nb) + tantalum		0.30-1.	50						
ron		13.0-17	7.0 2.0-6.0	3.0 max	4.0-7.0	_ 2.0 max	1.0-3.0	17.0-20.0	3.0 max
Cobalt, max		5.0	2.5	2.0	2.5	1.0	3.0	0.5-2.5	5.0
Tungsten		1.5-4.0			3.0-4.5		3.0 max	0.3-2.3	13.0-15.0
Vanadium, max			0.35		0.35		0.20		10.0 10.0
Titanium, max		(ht	tn c //			ifah oi	0.2		
Zirconium, max			UD5.//	Stan(iai us.	Ittii.a	0.10		
Columbium (Nb)			I				0.20 max		
Tantalum			Danasi	1000	4 Duck		0.20 max		
Nickel +					ı Erev	1ew	94.0-98.0		
Molybdenum									
Aluminum, max				0.50			0.50		0.20-0.50
Lanthanum					470 "0 4				0.005-0.05
Boron				<u>ASIM E</u>	<u>84 / 2 - 04a</u>				0.015 max
	ards.iteh.a	ai/catalog/s	standards/si	st/bfca74b	Composition,	3%9d1b-447b	2a359b41	/astm-b4	472-04a
lement	UNS	UNS	UNS	UNS	UNS	UNS	UNS		UNS
	N12160	R30556	N06625	N06600	N10242	N08120	N06617		N06035
arbon	0.15 max	0.05-0.15	0.10 max	0.15 max	0.03	0.02-0.10	0.05-0.15		0.050 max
anganese	1.5 max	0.50-2.00	0.50 max	1.0 max	0.80 max	1.5 max	1.0 max		0.50 max
	0.030 max	0.04 max	0.015 max		0.030 max	0.040 max			0.030 max
							0.015 ma	ιX	0.015 max
ulfur	0.015 max	0.015 max	0.015 max	0.015 max		0.03 max			0.60 max
ulfur licon	0.015 max 2.4-3.0	0.20-0.80	0.015 max 0.50 max	0.50 max	0.80 max	1.0 max	1.0 max		
ulfur Ilicon ickel	0.015 max 2.4-3.0 remainder ^A	0.20-0.80 19.0-22.5	0.015 max 0.50 max 58.0 min ^A	0.50 max 72.0 min ^A	0.80 max remainder ^A	1.0 max 35.0-39.0	1.0 max 44.5 min ²	4	remainder ^A
ulfur ilicon ickel hromium	0.015 max 2.4-3.0 remainder ^A 26.0-30.0	0.20-0.80 19.0-22.5 21.0-23.0	0.015 max 0.50 max 58.0 min ^A 20.0-23.0	0.50 max 72.0 min ^A 14.0-17.0	0.80 max remainder ^A 7.0-9.0	1.0 max 35.0-39.0 23.0-27.0	1.0 max 44.5 min ² 20.0-24.0	4	remainder ^A 32.25-34.25
ulfur ilicon ickel hromium olybdenum	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0	0.015 max 0.50 max 58.0 min ^A 20.0-23.0 8.0-10.0	0.50 max 72.0 min ^A 14.0-17.0 	0.80 max remainder ^A 7.0-9.0 24.0-26.0	1.0 max 35.0-39.0 23.0-27.0 2.50 max	1.0 max 44.5 min ² 20.0-24.0 8.0-10.0	4	remainder ^A 32.25-34.25 7.60-9.00
ulfur ilicon ickel hromium olybdenum opper	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0	0.015 max 0.50 max 58.0 min ^A 20.0-23.0 8.0-10.0	0.50 max 72.0 min ^A 14.0-17.0 0.5 max	0.80 max remainder ^A 7.0-9.0 24.0-26.0 	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max	1.0 max 44.5 min/ 20.0-24.0 8.0-10.0 0.5 max	4	remainder ^A 32.25-34.25 7.60-9.00 0.30 max
ulfur ilicon ickel hromium olybdenum opper olumbium	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0	0.015 max 0.50 max 58.0 min ^A 20.0-23.0 8.0-10.0	0.50 max 72.0 min ^A 14.0-17.0 	0.80 max remainder ^A 7.0-9.0 24.0-26.0	1.0 max 35.0-39.0 23.0-27.0 2.50 max	1.0 max 44.5 min ² 20.0-24.0 8.0-10.0	4	remainder ^A 32.25-34.25 7.60-9.00
ulfur licon ckel nromium olybdenum opper olumbium lb) + tantalum	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 	0.015 max 0.50 max 58.0 min ^A 20.0-23.0 8.0-10.0 3.15-4.15	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 	0.80 max remainder ^A 7.0-9.0 24.0-26.0	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9	1.0 max 44.5 min ² 20.0-24.0 8.0-10.0 0.5 max 	4	remainder ^A 32.25-34.25 7.60-9.00 0.30 max
ulfur licon ckel nromium olybdenum opper olumbium lb) + tantalum trogen	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 	0.80 max remainder ^A 7.0-9.0 24.0-26.0	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 	4	remainder ^A 32.25-34.25 7.60-9.00 0.30 max
ulfur licon ckel hromium olybdenum opper olumbium lb) + tantalum trogen	0.015 max 2.4-3.0 remainder ⁴ 26.0-30.0 1.0 max 3.5 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ⁴	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max	4	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max
ulfur licon ickel hromium olybdenum opper olumbium lb) + tantalum itrogen on	0.015 max 2.4-3.0 remainder ⁴ 26.0-30.0 1.0 max 3.5 max 27.0-33.0	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ⁴ 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min-	.15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00
ulfur licon ickel hromium olybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten	0.015 max 2.4-3.0 remainder ⁴ 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0 	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 	4 15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max
ulfur licon ickel hromium olybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 	.15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
ulfur licon ickel hromium olybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max tanium	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0 	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 	4 15.0 max	remainder ^A 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
ulfur ilicon ickel hromium olybdenum opper olumbium Ib) + tantalum itrogen on obalt, max ungsten anadium, max tanium irconium	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 0.001-0.10	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 0.6 max	4 15.0 max	remainder ^A 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
ulfur licon ickel hromium olybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max tanium rconium olumbium	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80 	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 0.001-0.10 0.30 max	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 0.4 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 0.6 max 	4 15.0 max	remainder ^A 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
ulfur ilicon ickel hromium lolybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max ttanium irconium olumbium antalum	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80 	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 0.001-0.10 0.30 max 0.30-1.25	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 0.4 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 0.20 max	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 0.6 max 	4 15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
hosphorous ulfur iliticon ickel hromium lolybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max itanium irconium olumbium antalum luminum, max	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80 1.0 max	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 0.001-0.10 0.30 max 0.30-1.25 0.10-0.50	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 0.4 max 0.4	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 0.20 max 	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 0.6 max 	4 15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20
ulfur ilicon ickel hromium lolybdenum opper olumbium lb) + tantalum itrogen on obalt, max ungsten anadium, max ttanium irconium olumbium antalum	0.015 max 2.4-3.0 remainder ^A 26.0-30.0 1.0 max 3.5 max 27.0-33.0 1.0 max 0.20-0.80 	0.20-0.80 19.0-22.5 21.0-23.0 2.5-4.0 0.10-0.30 remainder ^A 16.0-21.0 2.0-3.5 0.001-0.10 0.30 max 0.30-1.25	0.015 max 0.50 max 58.0 min ⁴ 20.0-23.0 8.0-10.0 3.15-4.15 5.0 max 0.4 max 	0.50 max 72.0 min ^A 14.0-17.0 0.5 max 6.0-10.0	0.80 max remainder ^A 7.0-9.0 24.0-26.0 2.0 max 1.00 max	1.0 max 35.0-39.0 23.0-27.0 2.50 max 0.50 max 0.4-0.9 0.15-0.30 remainder ^A 3.0 2.50 max 0.20 max	1.0 max 44.5 min' 20.0-24.0 8.0-10.0 0.5 max 3.0 max 10.0 min- 0.6 max 	4 .15.0 max	remainder ⁴ 32.25-34.25 7.60-9.00 0.30 max 2.00 max 1.00 0.60 max 0.20