



SLOVENSKI STANDARD
SIST ISO 11264:2006

01-september-2006

Učinkovitost določanja herbicidov v tleh s pomočjo HPLC metode z UV-detekcijo

Soil quality -- Determination of herbicides -- Method using HPLC with UV-detection

iTeh STANDARD PREVIEW
(standards.iteh.ai)

Qualité du sol -- Dosage des herbicides -- Méthode par CLHP avec détection par UV

[SIST ISO 11264:2006](https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-40c4de15cc17/sist-iso-11264-2006)

Ta slovenski standard je istoveten z: **ISO 11264:2005**

ICS:

13.080.10 Sestava in lastnosti tal Chemical characteristics of soils

SIST ISO 11264:2006

en

iTeh STANDARD PREVIEW
(standards.iteh.ai)

SIST ISO 11264:2006

<https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006>

INTERNATIONAL
STANDARD

ISO
11264

First edition
2005-05-01

**Soil quality — Determination of
herbicides — Method using HPLC with
UV-detection**

*Qualité du sol — Dosage des herbicides — Méthode par CLHP avec
détection par UV*

iTeh STANDARD PREVIEW
(standards.iteh.ai)

[SIST ISO 11264:2006](https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006)

<https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006>



Reference number
ISO 11264:2005(E)

© ISO 2005

ISO 11264:2005(E)**PDF disclaimer**

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.

iTeh STANDARD PREVIEW (standards.iteh.ai)

[SIST ISO 11264:2006](https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006)

<https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006>

© ISO 2005

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

Contents

Page

Foreword.....	iv
1 Scope.....	1
2 Normative reference	1
3 Principle	1
4 Reagents	1
5 Apparatus.....	4
6 Procedure.....	5
6.1 Sample preparation.....	5
6.2 Sample preservation and pretreatment	5
6.3 Addition of water to the sample	5
6.4 Preparation of soil extract.....	6
6.5 HPLC-determination	6
7 Calculation	9
7.1 Calculation of the mass concentration ρ_{EX} of the substance in extract solution E2 according to Equation (2).....	9
7.2 Calculation of the content w of the substance in the soil extract solution E2 according to Equation (3)	10
7.3 Expression of results.....	10
8 Test report.....	10
9 Accuracy	10
Annex A (informative) Chromatograms.....	11
Annex B (normative) Validation	13
Annex C (informative) Additional compounds which were tested with this method (see 6.4 and 6.5)	16
Bibliography	17

ISO 11264:2005(E)**Foreword**

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 11264 was prepared by Technical Committee ISO/TC 190, *Soil quality*, Subcommittee SC 3, *Chemical methods and soil characteristics*.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

SIST ISO 11264:2006

<https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-bf41-46c4de15cc17/sist-iso-11264-2006>

Soil quality — Determination of herbicides — Method using HPLC with UV-detection

1 Scope

This International Standard specifies a high-performance liquid chromatography (HPLC) method for qualitative and quantitative determination of herbicides of various substance classes in soils. This method covers triazines including their related metabolites, phenyl urea compounds and other herbicides. Compounds are identified and quantified with UV-detection.

The limit of detection for triazines and phenyl urea compounds is ~ 0,01 mg/kg dry matter. It is dependent upon both the compound and the soil matrix.

2 Normative reference

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 10381-1, *Soil quality — Sampling — Part 1: Guidance on the design of sampling programmes*

ISO 10381-2, *Soil quality — Sampling — Part 2: Guidance on sampling techniques*

ISO 11465, *Soil quality — Determination of dry matter and water content on a mass basis — Gravimetric method*

ISO 14507, *Soil quality — Pretreatment of samples for determination of organic contaminants*

3 Principle

The field-moist soil sample is extracted with a mixture of acetone and water (2:1). After addition of NaCl and dichloromethane or petroleum, the isolated organic phase is concentrated and transferred to a acetonitrile/water mixture and without further clean-up is analysed using RP-HPLC, gradient of acetonitrile and water, with UV-detection. Results can be confirmed using diode array UV spectra, GC-MS, GC-NPD or GC-AED (some may need derivatisation).

4 Reagents

All reagents shall be of known analytical grade. The purity of the reagents used shall be checked by running a blank determination as described in 6.5. If the blank value is unreasonably high, i.e. more than 10 % of the lowest value of interest, find the cause through a step-by-step examination of the whole procedure. For measurements at the limit of determination, even reagents suitable for residue analyses may not fulfil this criterion. In this case, sufficient blank determinations shall be incorporated in each series of samples.

4.1 Water, for residue analysis, normally tap water (drinking water) is suitable.

4.2 Acetone, for residue analysis.

ISO 11264:2005(E)

- 4.3 Sodium chloride.**
- 4.4 Dichloromethane**, for residue analysis.
- 4.5 Petroleum ether**, boiling range 40 °C to 60 °C, for residue analysis.
- 4.6 Sodium sulfate**, anhydrous, for residue analysis.
- 4.7 Acetonitrile**, HPLC grade.
- 4.8 Water**, HPLC grade.
- 4.9 2-propanol**, for chromatography.
- 4.10 Acetonitrile (4.7)/water (4.8) mixture**, 1 + 1 (*V + V*) for dissolving the sample.
- 4.11 Methanol**, HPLC grade.
- 4.12 Standard substances** (including Chemical Abstracts Number).

Alachlor [015972-60-8], Atrazine [001912-24-9], Atrazine-desethyl [06190-65-4], Atrazine-desisopropyl [01007-28-9], Bromacil [000314-40-9], Chlortoluron [015545-48-9], Chloroxuron [001982-47-4], Cyanazine [021725-46-2], Dichlobenil [001194-65-6], Diuron [000330-54-1], Ethofumesate [026225-79-6], Hexazinon [051235-04-2], Isoproturon [034123-59-6], Metazachlor [061729-08-2], Metamitron [041394-05-2], Metabenzthiazuron [018691-97-9], Metobromuron [003060-89-7], Metolachlor [051218-45-2], Metoxuron [019937-59-8], Monuron [000150-68-5], Pendimethalin [040487-42-1], Propazine [000139-40-2], Propyzamide [023950-58-5], Sebuthylazine [00728-69-3], Simazine [000122-34-9], Terbutryn [000886-50-0], Terbutylazine [005915-41-3].

4.13 Standard solutions for HPLC determination

4.13.1 Stock solutions, $\rho_{11} = 1 \text{ mg/ml}$

Weigh 10,0 mg of each standard substance (4.12) and place in individual 10 ml measuring flasks (5.13). Dissolve with 2-propanol (4.9) or, in the case of sebuthylazine and simazine, with small portions of acetone (4.2) or, in the case of metamitron, with methanol (4.11). Then make up to the mark with 2-propanol (4.9).

4.13.2 Intermediate dilutions, $\rho_{12} = 100 \text{ }\mu\text{g/ml}$

Pipette 1,0 ml of each of the stock solutions (4.13.1) into an individual 10 ml measuring flask (5.13) and make up to the mark with 2-propanol (4.9).

4.13.3 Working solutions I, $\rho_{13} = 10 \text{ }\mu\text{g/ml}$

NOTE Working solutions are made to establish the retention time. Different concentrations are made because the analytes have different detector responses.

Pipette 1,0 ml of each of the intermediate dilutions (4.13.2) into an individual 10 ml measuring flask (5.13) and make up to the mark with acetonitrile/water mixture (4.10).

4.13.4 Working solutions II, $\rho_{14} = 1 \text{ }\mu\text{g/ml}$

Transfer, with a microlitre syringe (5.14), 100 μl of each of the intermediate dilutions (4.13.2) into a 10 ml measuring flask (5.13) and make up to the mark with acetonitrile/water mixture (4.10).

4.13.5 Working solutions III, $\rho_{15} = 0,1 \mu\text{g/ml}$

Transfer, with a microlitre syringe (5.14), 100 μl of working solution I (4.13.3) into a 10 ml measuring flask (5.13) and make up to the mark with acetonitrile/water mixture (4.10).

4.13.6 Mixed solutions I, $\rho_{16} = 5 \mu\text{g/ml}$ to 100 $\mu\text{g/ml}$

According to Table 1, transfer with a microlitre syringe (5.14), between 50 μl and 1 000 μl of the stock solutions (4.13.1) into a 10 ml measuring flask (5.13) and make up to the mark with 2-propanol (4.9).

4.13.7 Mixed solutions II, $\rho_{17} = 0,5 \mu\text{g/ml}$ to 10 $\mu\text{g/ml}$

Pipette 1,0 ml of the mixed solution I (4.13.6) into a 10 ml measuring flask (5.13) and make up to the mark with acetonitrile/water mixture (4.10) (see Table 1).

4.13.8 Mixed solutions III, $\rho_{18} = 0,05 \mu\text{g/ml}$ to 1 $\mu\text{g/ml}$

Pipette 1,0 ml of the mixed solution II (4.13.7) into a 10 ml measuring flask (5.13) and make up to the mark with acetonitrile/water mixture (4.10).

Table 1 — Concentration of compounds of mixed solutions I (4.13.6), II (4.13.7) and III (4.13.8)

No.	Compound ^a	Volume ^b μl	Mixed solution I $\mu\text{g/ml}$	Mixed solution II $\mu\text{g/ml}$	Mixed solution III $\mu\text{g/ml}$
1	Atrazine-desisopropyl	200	20	2	0,2
2	Metamitron	500	50	5	0,5
3	Atrazine-desethyl	200	20	2	0,2
4	Hexazinon	200	20	2	0,2
5	Metoxuron	200	20	2	0,2
6	Bromacil	400	40	4	0,4
7	Simazine	100	10	1	0,1
8	Monuron	100	10	1	0,1
9	Cyanazine	100	10	1	0,1
10	Methabenzthiazuron	100	10	1	0,1
11	Chlortoluron	100	10	1	0,1
12	Atrazine	100	10	1	0,1
13	Isoproturon	100	10	1	0,1
14	Diuron	200	20	2	0,2
15	Metobromuron	100	10	1	0,1
16	Metazachlor	500	50	5	0,5
17	Sebutylazine	100	10	1	0,1
18	Propazine	100	10	1	0,1
19	Dichlobenil	200	20	2	0,2
20	Terbutylazine	100	10	1	0,1
21	Chloroxuron	100	10	1	0,1
22	Propyzamide	100	10	1	0,1
23	Terbutryn	50	5	0,5	0,05
24	Ethofumesate	500	50	5	0,5
25	Metolachlor	500	50	5	0,5
26	Alachlor	1 000	100	10	1
27	Pendimethalin	100	10	1	0,1

^a Compounds are in elution order using the conditions stated in 6.5. For additional compounds, see Annex C.

^b Volume to be taken from the stock solution (4.13.1) to prepare 10 ml of the mixed solution I (4.13.6).

ISO 11264:2005(E)

5 Apparatus

5.1 Glass bottles, with screw cap and sealing ring of polytetrafluor ethylene, nominal capacity 250 ml and 1 000 ml.

5.2 Extraction apparatus (see Figure 1), consisting of the following.

5.2.1 Measuring cylinder, of 1 000 ml nominal volume with ground neck and ground-glass stopper.

5.2.2 Ground-glass insert (wash-bottle principle) with adjustable U-tube, the bent end of which can be placed at the interface between the two phases.




Figure 1 — Extraction apparatus

5.3 Shaking apparatus, head-over-head or horizontal. Complete mixing of the soil and extractant should be obtained.

5.4 Vacuum rotating evaporator, with water bath and vacuum control.

5.5 Pear-shaped flask, of 10 ml nominal volume, with graduations at 0,5 ml and 1,0 ml and a ground-glass stopper.

5.6 Membrane filter, porosity at least 0,45 µm, solvent-resistant.

5.7 Test tube, of 2 ml nominal volume, sealable, with polytetrafluor ethylene disc.

5.8 HPLC apparatus, consisting of high-pressure pumps with gradient programmer, injection valve of 20 µl to 50 µl, column thermostat, solvent de-gassing device.

- 5.9 UV-detector**, with variable wavelengths or diode array detector (DAD).
- 5.10 Electronic recording device** (integrator or laboratory data system).
- 5.11 Separation column**, length 250 mm, inner diameter 4 mm, and pre-column holders (see 6.5.1).
- 5.12 Pre-column**, length 4 mm, inner diameter 4 mm (see 6.5.1)
- 5.13 Measuring flasks**, 10 ml nominal volume, with glass stopper.
- 5.14 Microlitre syringe**, 50 µl to 1 000 µl.
- 5.15 Ultrasonic bath or vortex mixer**.

6 Procedure

6.1 Sample preparation

Obtain representative field-moist soil samples in accordance with ISO 10381-1, using sampling apparatus in accordance with ISO 10381-2.

6.2 Sample preservation and pretreatment

Field-moist samples shall be pretreated in accordance with ISO 14507 as soon as possible (no chemical drying).

Store the samples in a dark place at a temperature below – 20 °C. Herbicides may be subject to microbial conversion under certain conditions. It is recommended that samples be frozen, if they are stored for more than 2 days.

<https://standards.iteh.ai/catalog/standards/sist/cf07b03c-288f-40a8-b411-46c4de15cc17/sist-iso-11264-2006>

Determine the water and dry matter contents of a subpart of the sample in accordance with ISO 11465.

6.3 Addition of water to the sample

Add to the soil mass a specific mass of water calculated according to Equation (1) (each prepared sample will have a water mass of 100 g).

$$m_w = 100 - \frac{m_s \cdot w_w}{100} \quad (1)$$

where

m_w is the mass of added water, in g;

m_s is the mass of sample, in g;

w_w is the mass fraction of water, in %.

NOTE Addition of water is important to have a fixed ratio (water acetone dichloromethane of 1:2:1,5). Under these conditions, the organic phase will be the upper layer.