
**Implants for surgery — Retrieval and
analysis of surgical implants —**

**Part 1:
Retrieval and handling**

Implants chirurgicaux — Retrait et analyse des implants chirurgicaux —

Partie 1: Retrait et manipulation

**iTeh STANDARD PREVIEW
(standards.iteh.ai)**

ISO 12891-1:2011

<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>



iTeh STANDARD PREVIEW
(standards.iteh.ai)

ISO 12891-1:2011

<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>



COPYRIGHT PROTECTED DOCUMENT

© ISO 2011

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

Contents

Page

| | |
|--|-----------|
| Foreword | iv |
| Introduction..... | v |
| 1 Scope | 1 |
| 2 Terms and definitions | 1 |
| 3 Method..... | 2 |
| 3.1 Obtaining the clinical history of the implant and patient..... | 2 |
| 3.2 Pre-explantation checks and examinations..... | 2 |
| 3.3 Collecting the surgical implant..... | 2 |
| 3.4 Collecting the tissue and fluid samples..... | 3 |
| 3.5 Photographic record of the explantation..... | 3 |
| 3.6 Containing and labelling the retrieved surgical implant, tissues and fluids for future identification | 4 |
| 3.7 Cleaning the retrieved surgical implant..... | 4 |
| 3.8 Decontaminating the retrieved surgical implant..... | 7 |
| 3.9 Packaging the retrieved surgical implant, tissues and fluids for shipment..... | 7 |
| 3.10 Use of coolant materials..... | 8 |
| 3.11 Labelling of the packing materials..... | 8 |
| 3.12 Documentation to be supplied with retrieved surgical implants..... | 8 |
| 3.13 Unpacking following shipment..... | 9 |
| 3.14 Cleaning and decontamination following shipment..... | 9 |
| 3.15 Documentation to be maintained during examination, analysis and storage | 9 |
| 4 Analysis of retrieved surrounding tissues and fluids..... | 10 |
| 5 Infection control | 10 |
| 5.1 General | 10 |
| 5.2 Work practices..... | 10 |
| 5.3 Personal protective equipment..... | 10 |
| 5.4 Maintenance of the worksite | 11 |
| 5.5 Human waste disposal..... | 13 |
| 5.6 Special practices | 13 |
| Annex A (informative) Suggested minimum information to be obtained for retrieved surgical implants..... | 14 |
| Annex B (informative) Generic procedures for the decontamination of surgical implants | 16 |
| Annex C (informative) Analyses to be performed on retrieved tissue samples and fluids..... | 22 |
| Bibliography..... | 23 |

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 12891-1 was prepared by Technical Committee ISO/TC 150, *Implants for surgery*.

This second edition cancels and replaces the first edition (ISO 12891-1:1998), which has been technically revised.

ISO 12891 consists of the following parts, under the general title *Implants for surgery — Retrieval and analysis of surgical implants*:

- ISO 12891-1:2011
<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>
- *Part 1: Retrieval and handling*
 - *Part 2: Analysis of retrieved metallic surgical implants*
 - *Part 3: Analysis of retrieved polymeric surgical implants*
 - *Part 4: Analysis of retrieved ceramic surgical implants*

Introduction

The investigation of retrieved surgical implants, adjacent tissues and associated fluids can be undertaken to:

- determine the cause of a clinical complication or surgical implant failure;
- improve knowledge of surgical implant performance and safety;
- improve knowledge of the interactions of surgical implants and the human body;
- develop materials with improved biocompatibility and implants with improved functional longevity.

This International Standard specifies methods for the retrieval, handling and analysis of surgical implants and associated tissue samples and fluids which are removed from patients during retrieval surgery or post-mortem. ISO 12891-2 to ISO 12891-4 specify methods for the detailed analysis of specific types of surgical implants, in which protocols are provided for the collection of data and examinations for metallic, polymeric and ceramic surgical implants in relation to their typical applications. For particular investigation programmes, additional, more specific, protocols can be required. If special analytical techniques are employed, the appropriate handling procedures need to be specified.

The purpose of this International Standard is to:

- specify a method for the retrieval of surgical implants, which is intended to prevent damage to the implants, the associated tissues, and fluids;
- ensure that retrieved materials are handled safely and decontaminated correctly, and that the risk of transmission of infectious diseases is minimized;
- ensure that the retrieval process is properly documented;
- allow comparisons between investigation results from different sources.

Many variables are involved when undertaking the retrieval of surgical implants. The retrieval can be for the routine replacement of a pacemaker battery or it can be for the revision of a defective surgical implant. The retrieval can be from a living patient or it can be a post-mortem study. The retrieval can involve the removal of a single surgical implant or multiple components as, for example, in the case of hip replacements or certain fracture fixation or spinal devices. In addition to the retrieval of the surgical implant, associated tissues and fluids might also need to be removed. The retrieval can involve a wide variety of personnel such as surgeons, nurses, other hospital staff, the surgical implant manufacturer, the investigator, and the shipping service. Finally, the type of analysis to be performed can vary and can include visual, chemical, histological and microbiological studies and the eventual analysis can have an impact on the retrieval process. These variables make it impossible to specify a single method which has to be followed in all retrieval cases. For this reason, certain requirements listed in this part of ISO 12891 might only be applicable in certain circumstances and for this reason some of the requirements are prefaced with statements such as “If applicable” or “Whenever possible”.

This International Standard presents a methodology for the systematic retrieval of surgical implants. In particular, it focuses on the practical requirements. In addition to these requirements, there are legal and ethical considerations which might need to be taken into account. These considerations include matters relating to the ownership of the implant, the obtaining of the patient's consent before the implant is retrieved, the patient's right to confidentiality and the need to protect the patient's safety, health and litigation rights throughout. For a detailed consideration of these issues, appropriate advice can be sought.

NOTE The methods specified in this International Standard can also be applicable to the retrieval and analysis of surgical implants in animal studies.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

ISO 12891-1:2011

<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>

Implants for surgery — Retrieval and analysis of surgical implants —

Part 1: Retrieval and handling

1 Scope

This part of ISO 12891 specifies the method to be followed for the retrieval and handling of surgical implants and associated tissues and fluids. In particular, it specifies the essential steps to be followed for the safe and proper obtaining of the clinical history, pre-explantation checks and examinations, collection, labelling, cleaning, decontamination, documentation, packing and shipping. This part of ISO 12891 also provides guidance on infection control.

NOTE National or other regulations, which can be more stringent, can apply.

This part of ISO 12891 does not apply in cases of explantation where there is no intention to collect retrieval data. However, many clauses give useful information which can apply in these cases also.

2 Terms and definitions

[ISO 12891-1:2011](https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011)

<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>

For the purposes of this document, the following terms and definitions apply.

2.1

absorbent

material capable of absorbing liquids

NOTE Absorbent material can be either particulate or non-particulate.

2.2

contamination

unintentional addition or modification, including exposure to a potentially infectious agent

2.3

infectious waste

waste containing or suspected to contain human pathogenic microbiological agents

2.4

outer shipping container

outermost container in which the package is finally shipped

2.5

primary container

tube, envelope, or other impermeable container which holds the retrieved material to be shipped

2.6

secondary container

container into which the primary container is placed

3 Method

3.1 Obtaining the clinical history of the implant and patient

Whenever possible, the clinical history of the patient and the surgical implant shall be obtained and recorded. This clinical history shall include at least the following, if available:

- the name or identification number of the patient as permitted by the applicable national regulations;
- the original diagnosis which resulted in the use of the surgical implant;
- an X-ray of the surgical implant *in situ* taken after the insertion operation;
- the patient's activity level including the ability to perform work, sports and recreational activities;
- the patient's medical history relevant to the surgical implant, including the hospital or clinic at which the surgical implant was implanted;
- information on the patient's experience with the surgical implant just before surgical implant retrieval;
- the date of retrieval;
- the hospital or clinic at which the surgical implant was retrieved.

The information obtained should be treated as confidential.

NOTE 1 Annex A gives an example of the information to be obtained. The annex can be modified as necessary, e.g. for special investigations.

NOTE 2 In any analysis of an explanted surgical implant, it is advantageous to have as much relevant knowledge of the clinical history as possible.

3.2 Pre-explantation checks and examinations

Whenever possible, data which might be lost post-explantation should be collected prior to explantation.

Whenever possible, a functional check of the implant involving an objective measurement shall be performed before explantation surgery.

NOTE A functional check assists in the understanding of the post-explantation behaviour.

If applicable, electronic or other data associated with the surgical implant shall be collected before the implant is explanted. These data should be provided to the evaluator of the surgical implant.

Where appropriate and justifiable, taking into account the need for patient safety, non-invasive examinations of the implantation site with the implant *in situ* shall be performed before the implant is explanted. Such examinations may include X-ray, computed axial tomography scan or magnetic resonance imaging.

3.3 Collecting the surgical implant

Taking into account the need for patient safety, the surgical implant shall be retrieved in a manner which causes as little damage as possible to both the surgical implant and the surrounding tissues. As far as possible, functional surfaces, e.g. bearing surfaces of joint prostheses, mechanical connections, e.g. hinges, joints, screws, and fracture surfaces of broken surgical implants, shall be protected during and after explantation.

Fragments and debris, which can provide valuable information, shall also be retrieved.

Retrieved surgical implants should be handled with care either by hand or using appropriate instruments.

The following shall be documented:

- a) the position or orientation or state of the retrieved components, if there is more than one retrieved component and if the position or orientation or state is abnormal;
- b) the location and type of damage, if damage occurs during explantation.

NOTE For proper scientific examination, it is advantageous for the surgical implant to be maintained in a state as close as possible to that in which it existed at the time of retrieval.

3.4 Collecting the tissue and fluid samples

Taking into account the need for patient safety, if tissue and/or fluid samples are to be collected for analysis, then these shall be retrieved in a manner which causes as little damage as possible to both the surgical implant and the tissues.

For microbiological investigation, swabs, tissue and/or fluid samples shall be taken from a location adjacent to the implant as soon as possible after the surgical implant has been exposed. Where and how the specimens are taken shall be recorded.

NOTE Special culturing techniques can be required to reveal unusual organisms. Sampling for immunological investigations requires expert advice and can call for special procedures.

For histological examination, tissue samples shall be taken from a location adjacent to the implant and/or from other relevant sites (e.g. lymph nodes or any tissue with abnormal appearance).

When possible, tissue samples for histological examination shall include portions extending into the normal tissue.

The site of the tissue excision and the orientation of the tissue relative to the surgical implant shall be indicated and recorded. Where possible, the proximal end of the tissue shall be marked (e.g. with a suture). Where necessary, the original length of the tissue shall be maintained (e.g. with plastic muscle biopsy clamps or by pinning the tissue to a corkboard or by other means, which avoid contact with metal which could corrode).

The tissue samples for histological examination shall be transferred as early as possible to an appropriate fixative or other media. The type of fixative used and the time between excision and placement in the fixative or media shall be documented. The tissue sample shall be treated in a routine manner as required for histological examination, unless a special method is needed for special investigations.

If appropriate, the media used to preserve tissue attached to a retrieved surgical implant shall be selected so as not to affect the surgical implant. When it is not possible to preserve the tissues without affecting the retrieved surgical implant, the portions of the retrieved surgical implant to be analysed should be determined and the tissue preserved accordingly.

Fluids obtained by aspiration shall be appropriately preserved for examination unless a special method is needed for special investigations. The preservation method should be chosen taking into account the intended analysis.

In post-mortem studies, histological examination of remote tissues, e.g. liver and kidney, should also be performed, if there is a need to assess toxicity in these locations.

3.5 Photographic record of the explantation

Where appropriate, photographic records shall be made of the surgical implant *in situ*, of the surgical site and of the explanted surgical implant and any associated tissue specimens.

Where appropriate, the orientation of all removed surgical implant components in relation to each other and their placement in relation to the body and associated excised material shall be recorded. If not self-explanatory, the proximal end and the orientation in the transverse plane of the implant shall be marked and documented. Any observed abnormalities in the appearance or condition of the device shall be recorded.

NOTE As appropriate, either the surgical implant and the tissues themselves or the support upon which they are placed (e.g. a corkboard) can be marked.

3.6 Containing and labelling the retrieved surgical implant, tissues and fluids for future identification

Immediately after collection or photography, all surgical implants, associated tissues and fluids, which are retrieved for analysis, shall be placed in appropriate containers that can be sealed in such a way that any subsequent opening of the containers can be detected. Suitable containers include envelopes, bags, jars, pots and boxes.

NOTE Adhesive tape is normally used to seal the containers.

For example, when the container is an envelope, the flap of the envelope shall be taped so that the tape covers both the flap and the envelope itself. When the container is a bag, the opening of the bag shall be sealed with tape. Similarly, when the container is a jar, pot or box, the juncture of the lid and the container shall be taped. The retriever shall place his initials across the tape. In this way, opening the container tears the tape and disturbs the initials.

Immediately after containing the surgical implants, tissues and fluids, all containers shall be labelled to ensure their precise identification at some later date. The label shall contain at least the following information:

- a) an accurate description of the contents of the container (e.g. vascular graft, type XYZ);
- b) name or initials of retriever;
- c) date, time and place of retrieval;
- d) name or identification number of patient, if available, as permitted by the applicable national regulations;
- e) the container number or identifier, if there is more than one container;
- f) the orientation of each component relative to the others, if there is more than one retrieved component and if the orientation is abnormal.

The labels used shall be of a non-removable type (labels that tear when someone tries to remove them).

3.7 Cleaning the retrieved surgical implant

All surgical implants which are retrieved for analysis shall be cleaned before decontamination, unless otherwise specifically instructed. Cleaning can be performed off-site.

Retrieved surgical explants shall be cleaned as follows.

The retrieved surgical implant shall be thoroughly rinsed under running water, but not scrubbed, to remove all biological contaminants, unless such contaminants are important to the analysis. Adherent tissues considered important to the analysis shall be treated as a tissue sample (see 3.4). Loosely adherent material of possible interest should be preserved before the surgical implant is rinsed.

The retrieved surgical implant shall in addition be cleaned as recommended by the manufacturer. If the manufacturer cannot be contacted, or is unable to supply a means for cleaning the surgical implant, the method chosen shall be that given in Table 1 or any other method which has been shown to be effective while preserving the integrity of the implant, e.g. peracetic acid.

All solutions to be used in the cleaning of retrieved surgical implants shall be prepared at the time of cleaning, and shall not be stored in the laboratory for future use. Proteolytic enzyme solutions and ultrasonic bath solutions shall be disposed of according to the manufacturer's instruction, or may be decontaminated using a chemical disinfectant, and discarded in the sanitary sewer. Depending on local sewage company requirements and on the recommendation of the manufacturer, chemical cleaning agents may need to be neutralized before discarding into the sanitary sewer.

Any biological debris removed from the surgical implants shall be decontaminated by autoclaving, or disinfected via a chemical disinfectant, before disposal (see 5.5).

NOTE Table 1 presents general recommendations and is intended for use only when a manufacturer's recommendations cannot be obtained.

When chemical cleaning agents and/or an ultrasonic bath are used, cleaning should be performed inside a class II, type B (see Reference [2]), biological safety cabinet, which should be exhausted to the outside. In cases where there is tissue in-growth present, a proteolytic enzyme solution may be used in conjunction with ultrasonic cleaning, but only when no histological investigation is planned. Retrieved surgical implants which are too large to be placed in an ultrasonic bath shall be sprayed or surface-wiped with an appropriate chemical cleaning agent, or disinfected according to the ultrasonic bath solutions manufacturer. Such surgical implants should be cleaned in a biological safety cabinet of the class and type described above or in an isolated and well-ventilated area in the laboratory. Proper protective precautions, as specified in Clause 5, should be followed. Disposable swabs, brushes and wipes may be used to remove visible debris from such implants, in conjunction with an appropriate chemical agent.

iTeh STANDARD PREVIEW (standards.iteh.ai)

[ISO 12891-1:2011](https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011)

<https://standards.iteh.ai/catalog/standards/sist/7980eb97-5ffa-4efa-921f-0685ff46dde7/iso-12891-1-2011>

Table 1 — Generic recommendation for cleaning and decontaminating explanted surgical implants

| Device or implant ^a disinfected | Cleaning method ^b | Decontamination method ^{bc} |
|---|---|--|
| Cardiac pacemaker housing | Proteolytic enzyme solution or 70 % to 80 % isopropanol | Ethylene oxide gas or 70 % to 80 % ethanol or 3 % stabilized hydrogen peroxide |
| Leads | 70 % to 80 % ethanol or 70 % to 80 % isopropanol | 70 % to 80 % ethanol or 70 % to 80 % isopropanol |
| Cardiac valve: Mechanical valves | Proteolytic enzyme solution at or below room temperature with subsequent ultrasonic treatment | Ethylene oxide gas |
| Xenografts | Proteolytic enzyme solution | Ethylene oxide gas or buffered, alkaline 2 % solution of glutaraldehyde |
| Allografts | Broad-spectrum antibiotic solution | Ethylene oxide gas or buffered, alkaline 2 % solution of glutaraldehyde |
| Vascular grafts, biologic | 2 % buffered alkaline glutaraldehyde | Buffered, alkaline 2 % glutaraldehyde |
| Vascular grafts, synthetic | Proteolytic enzyme or 3 % stabilized hydrogen peroxide solution with subsequent ultrasonic treatment | Ethylene oxide gas or buffered, alkaline 2 % solution of glutaraldehyde or 4 % formaldehyde solution |
| Intra-aortic balloons and other temporary cardiac-assist implants | Peracetic acid ^d with subsequent ultrasonic treatment or solution sodium hypochlorite (500 mg/l to 600 mg/l) | Ethylene oxide gas or 70 % aqueous solutions of ethanol or isopropanol |
| Breast implants | Intense water rinse, proteolytic enzyme solution with subsequent ultrasonic treatment | 2 % glutaraldehyde, 4 % formaldehyde or ethylene oxide gas |
| Hydrocephalus shunts | Proteolytic enzyme solution at or below room temperature with subsequent ultrasonic treatment | Buffered, alkaline 2 % glutaraldehyde, ethylene oxide gas or 4 % formaldehyde ^e |
| Vascular port and peritoneal access implants | Sodium hypochlorite solution (50 mg/l to 60 mg/l) or 3 % hydrogen peroxide | Buffered, alkaline 2 % solution of glutaraldehyde or 70 % ethanol or isopropanol with 0,2 % glutaraldehyde |
| Intraocular lenses (HEMA) | | Ethylene oxide gas |
| Silicone elastomeric and polymeric implant components | 70 % to 80 % aqueous ethanol or isopropanol with subsequent ultrasonic treatment or sodium hypochlorite solution (50 mg/l to 60 mg/l) or 3 % hydrogen peroxide | Ethylene oxide gas or buffered, alkaline 2% solution of glutaraldehyde |
| Polymeric implant components (PMMA, PE-UHMW) | Proteolytic enzyme solution, with ultrasonic treatment, or sodium hypochlorite solution (50 mg/l to 60 mg/l) or 3 % hydrogen peroxide | Buffered, alkaline 2 % solution of glutaraldehyde or ethylene oxide gas |
| Metallic implant components | Intense water rinse, 70 % to 80 % aqueous ethanol or isopropanol with subsequent ultrasonic treatment or proteolytic enzyme or sodium hypochlorite solution (50 mg/l to 60 mg/l) or 3 % hydrogen peroxide | Steam autoclave or ethylene oxide |
| Ceramic implant components | Proteolytic enzyme solution, with ultrasonic treatment or sodium hypochlorite solution (50 mg/l to 60 mg/l) or 3 % hydrogen peroxide | Buffered, alkaline 2 % solution of glutaraldehyde or ethylene oxide gas |

^a When tissues are to be preserved, methods such as glutaraldehyde fixation may be used.
^b Percentages are volume fractions.
^c For disinfecting, a soaking time of 2 h to 3 h is sufficient. However, a 24 h contact time may be used to provide an extra margin of safety.
^d **WARNING — Peracetic acid is an explosive; it should be used with caution and stored in an explosion-proof refrigerator.**
^e KOH ($c = 4 \text{ mol/l}$) shall be used for final disposition of central nervous system explants.