



Designation: D4641 – 94 (Reapproved 2006)

Standard Practice for Calculation of Pore Size Distributions of Catalysts from Nitrogen Desorption Isotherms¹

This standard is issued under the fixed designation D4641; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This practice covers the calculation of pore size distributions for catalysts and catalyst carriers from nitrogen desorption isotherms. The computational procedure is particularly useful for determining how the pore volume is distributed in catalyst samples containing pores whose sizes range from approximately 1.5 to 100 nm (15 to 1000 Å) in radius. It should be used with caution when applied to isotherms for samples containing pores both within this size range and pores larger than 100 nm (1000 Å) in radius. In such instances the isotherms rise steeply near $P/P_o = 1$ and the total pore volume cannot be well defined. The calculations should be begun at a point on the isotherm near saturation preferably in a region near $P/P_o = 0.99$, establishing an upper limit on the pore size distribution range to be studied. Simplifications are necessary regarding pore shape. A cylindrical pore model is assumed, and the method treats the pores as non-intersecting, open-ended capillaries which are assumed to function independently of each other during the adsorption or desorption of nitrogen.

NOTE 1—This practice is designed primarily for manual computation and a few simplifications have been made for this purpose. For computer computation, the simplified expressions may be replaced by exact expressions.

1.2 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

2. Referenced Documents

2.1 *ASTM Standards:*²

D3766 Terminology Relating to Catalysts and Catalysis

¹ This practice is under the jurisdiction of ASTM Committee D32 on Catalysts and is the direct responsibility of Subcommittee D32.01 on Physical-Chemical Properties.

Current edition approved Oct. 1, 2006. Published October 2006. Originally approved in 1987. Last previous edition approved in 1999 as D4641-94(1999)¹. DOI: 10.1520/D4641-94R06.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

D4222 Test Method for Determination of Nitrogen Adsorption and Desorption Isotherms of Catalysts and Catalyst Carriers by Static Volumetric Measurements

3. Terminology

3.1 *Definitions*—Consult Terminology D3766.

3.2 *Symbols*:

$P_4(i)$	=	pressure after equilibration during desorption, torr.
$P_0(i)$	=	liquid nitrogen vapor pressure, torr.
V_{de}	=	see 12.4.10 and 12.5 in Test Method D4222.
$r_k(i)$	=	radius of inner core calculated from Kelvin equation, Å.
T	=	boiling point of nitrogen, K.
V_L	=	liquid nitrogen molar volume at T , cm ³ /mole.
γ	=	liquid nitrogen surface tension at T , mN/m.
$T(i)$	=	average thickness of the nitrogen film adsorbed on the pore walls, Å.
$r_p(i)$	=	radius of cylindrical pore given by $r_k(i) + t(i)$, Å.
Q	=	volume correction factor defined as $(\bar{r}_p/\bar{r}_k)^2$.
$\Delta V_T(i)$	=	decrease in the amount of nitrogen adsorbed caused by a lowering in relative pressure, mm ³ /g.
$\Delta V_f(i)$	=	volume of liquid nitrogen desorbed from pore walls during thinning of the film, mm ³ /g.
$\Delta V_k(i)$	=	liquid volume of the inner core in which capillary condensation of the nitrogen occurs, mm ³ /g.
$\Delta V_p(i)$	=	liquid volume contained in a group of pores having mean radius \bar{r}_p , mm ³ /g.
$\Sigma \Delta V_p$	=	cumulative pore volume, mm ³ /g.
$\Delta S_p(i)$	=	area of the pore walls of a cylinder having volume ΔV_p , m ² /g.

4. Summary of Practice

4.1 The pore size distribution is determined by analyzing the desorption data of the nitrogen isotherm. The nitrogen uptake is caused by the multilayer adsorption of a film of nitrogen on the pore walls and by capillary condensation of the nitrogen in the “inner core” regions of the pores. The relative pressure at which filling of the core occurs for a given pore size

by capillary condensation is predicted from the Kelvin equation (1).³ During desorption, thinning of the multilayer film adsorbed on the pore walls occurs in pores which have previously lost their capillary condensate. Corrections for film thinning are determined by a procedure involving the surface area and radius of the film which becomes exposed as desorption proceeds. In principle, the computational procedure can be applied to either the adsorption branch or desorption branch of the nitrogen isotherm. Unless the presence of ink-bottle shaped pores is suggested by an abrupt closure of the desorption branch on the adsorption branch, the distribution curve derived from the desorption data is preferred, and is described in this procedure. The computational method is essentially the procedure developed by Barrett, Joyner, and Halenda (2) except for the incorporation of a few simplifications.

NOTE 2—In cases where it has been established that the adsorption branch of the nitrogen isotherm is to be analyzed, the procedure proposed by Cranston and Inkley (3) is recommended.

5. Significance and Use

5.1 Pore volume distribution curves obtained from nitrogen sorption isotherms provide one of the best means of characterizing the pore structure in porous catalysts, provided that the limitations of the method are kept in mind. Used in conjunction with the BET treatment for surface area determination (4), these methods provide an indispensable means for studying the structure associated with pores usually important in catalysts. This practice is particularly useful in studying changes in a series of closely related samples caused by treatments, such as heat, compression, or extrusion often used in catalyst manufacturing. Pore volume distribution curves can often provide valuable information during mechanistic studies dealing with catalyst deactivation.

6. Computational Procedure

6.1 This procedure requires the use of a series of experimentally measured relative pressures [$P_4(i)/P_0(i)$] and the corresponding quantities of nitrogen gas adsorbed [V_{de}] expressed in units of cm^3 STP/g. The experimental data required in the use of this procedure can be measured by following the steps outlined in Test Method D4222. Inspect the nitrogen sorption isotherm in the region above $P/P_0 = 0.95$. If the solid contains no pores larger than 100 nm (1000 Å) radius, the isotherm remains nearly horizontal over a range of P/P_0 approaching unity and it is a simple matter to select a starting relative pressure within this region, establishing an upper limit on the pore size range to be studied. If pores larger than 100 nm (1000 Å) are present however, the isotherm rises rapidly near $P/P_0 = 1$ and the total pore volume cannot be well defined. This limiting adsorption can then be identified reliably only if the temperature is very carefully controlled and there are no “cold spots” in the apparatus (which lead to bulk condensation of the gas and a false measure of the adsorption in the volumetric method). Selecting the starting relative pressure for

the computational procedure is then made more difficult. In most cases a starting relative pressure of 0.99 will be suitable, which corresponds to an upper limit on pore size of 100 nm (1000 Å) in radius. If necessary, interpolate the values of V_{de} to determine the quantity of nitrogen gas adsorbed at the chosen starting relative pressure.

6.2 The procedure requires numerous arithmetical steps which can best be carried out with the aid of a work sheet. An example (4) of a form found useful in the calculations is provided in Table 1. List in descending order the experimentally determined relative pressures [$P_4(i)/P_0(i)$] in Column 1, beginning with the value chosen as the starting relative pressure. Generally, values below a relative pressure of 0.25 will not be required in the calculations. Convert the uptake values into a liquid volume (mm^3/g) by multiplying the value of V_{de} in cm^3 STP/g with the conversion factor 1.5468 derived from $V_L = 34.67 \text{ cm}^3/\text{mole}$. List in Column 9 the corresponding quantities of nitrogen adsorbed.

6.3 For each relative pressure, calculate a value for the radius of the core, r_k , by means of the Kelvin equation,

$$RT \ln (P_4/P_0) = - \frac{2\gamma V_L}{r_k} \quad (1)$$

given in the form

$$r_k (\text{Å}) = - \frac{9.574}{\ln (P_4/P_0)} \quad (2)$$

with $T = 77.35 \text{ K}$; $\gamma = 8.88 \text{ mN/m}$; and $V_L = 34.67 \text{ cm}^3/\text{mole}$. List the values in Column 2. For each successive decrement in relative pressure, calculate \bar{r}_k , the mean of the values of r_k for the present and previous pressures, and list these mean values in Column 3.

6.4 The average thickness, t , of the multilayer film of nitrogen adsorbed on the walls of the pores at each relative pressure is used to calculate the amount of nitrogen desorbed from the film in pores which have lost their capillary condensate. For each relative pressure, calculate a value for the film thickness from the expression (5)

$$t (\text{Å}) = \left[\frac{13.99}{0.034 - \log (P_4/P_0)} \right]^{\frac{1}{2}} \quad (3)$$

and list the values in Column 4. For each successive decrement in relative pressure, calculate the differences in the values of t , and list these differences as Δt in Column 5.

6.5 Since a cylindrical pore model is assumed, the radius of the pore, r_p , is given by addition of the core radius, r_k , and the film thickness value, t . Add the values in Column 2 to the corresponding values in Column 4 and record the results in Column 6 as r_p . For each successive decrement in relative pressure, calculate \bar{r}_p , the mean of the pore radii, r_p , for the present and previous pressures, and record these values in Column 7.

6.6 Compute the quantity $(\bar{r}_p/\bar{r}_k)^2$ from the values listed in Columns 7 and 3. This quantity will be used later to correct the core volume to the volume for each group of pores. The core volume is the region within the pore that fills by capillary condensation of the nitrogen. List the computed values in Column 8 as the volume correction factor, Q .

³ The boldface numbers in parentheses refer to a list of references at the end of this standard.

TABLE 1 Pore Distribution Computational Work Sheet

Sample Identification _____

Date _____

1 P_4/P_0	2 r_k Å	3 \bar{r}_k Å	4 t Å	5 Δt Å	6 r_p Å	7 \bar{r}_p Å	8 Q	9 V_{de} mm ³ /g	10 ΔV_T mm ³ /g	11 ΔV_f mm ³ /g	12 ΔV_k mm ³ /g	13 ΔV_p mm ³ /g	14 ΔS_p m ² /g	15 $\Sigma \Delta S_p$ m ² /g	16 $\Sigma \Delta V_p$ mm ³ /g
		—		—		—	—								
										[xcirc]					

Column 2: $r_k = \frac{-9.574}{\ln(P_4/P_0)}$

Column 11: $\Delta V_f = 0.085 \cdot \Delta t \cdot \Sigma \Delta S_p$ (preceding line)

Column 13: $\Delta V_p = \Delta V_k \times Q$

Column 4: $t = \left[\frac{13.99}{0.034 - \log(P_4/P_0)} \right]^{1/2}$

Column 12: $\Delta V_k = \Delta V_T - \Delta V_f$

Column 14: $\Delta S_p = 20 \times (\Delta V_p / \bar{r}_p)$

iTech Standards
<https://standards.iteh.ai>
 Document Preview

NOTE 3—For a cylindrical pore \bar{r}_p and \bar{r}_k are related to Q by the exact expression:

$$Q = \left[\bar{r}_p / (\bar{r}_k + \Delta t) \right]^2 \quad (4)$$

For $\bar{r}_k > 30$ Å, $\Delta t < 1$ % \bar{r}_k . Simplifying Q by eliminating Δt gives $(\bar{r}_p / \bar{r}_k)^2$.

6.7 The amounts of nitrogen desorbed for each successive decrement in relative pressure are calculated by progressive subtraction of the values of the amounts of nitrogen adsorbed, V_{de} , listed in Column 9 from the succeeding one. Compute these differences and list the values in Column 10 as, ΔV_T . Each value of ΔV_T except for the initial one in line 2 contains contributions from the amounts of nitrogen given up by loss of capillary condensate and by thinning of the nitrogen film adsorbed on the walls of pores which have previously released their capillary condensate. The initial value of ΔV_T is due solely to the amount of nitrogen contributed from loss of capillary condensate, since it is assumed that at the highest relative pressure all of the pores are completely filled with nitrogen, and that no thinning of the nitrogen film occurs for the first decrement in relative pressure.

6.8 In completing the calculation to obtain a value for the pore volume, ΔV_p , corresponding to each mean pore size, \bar{r}_p , of a group of pores, it will be necessary to work through the calculation for each line before proceeding to the next line. ΔV_f in Column 11 is the amount of nitrogen given up during desorption from thinning of the nitrogen film adsorbed on pore

walls. For line 2, ΔV_f equals zero and the value of ΔV_T is assigned to ΔV_k , the volume of the inner core which fills by capillary condensation of the nitrogen. Multiply the value of ΔV_k by the corresponding volume correction factor, Q , listed in Column 8 to obtain ΔV_p . List the value in Column 13 of line 2.

6.9 Calculate the surface area of the pore walls contained in volume ΔV_p as follows:

$$\Delta S_p \text{ (m}^2\text{/g)} = 20 \times \frac{\Delta V_p}{\bar{r}_p} \quad (5)$$

Using the value of ΔV_p from Column 13 and the corresponding value of \bar{r}_p from Column 7, compute a value for ΔS_p and list it in Column 14. A value for the total surface area of the pores that have become exposed is obtained by summation of the value for ΔS_p with other ΔS_p values in all preceding lines of Column 14. List the value of total area in Column 15 as $\Sigma \Delta S_p$. A value for the cumulative pore volume is obtained by summation of the value ΔV_p with other ΔV_p values in all preceding lines of Column 13. List the value of the cumulative pore volume in Column 16 as $\Sigma \Delta V_p$.

NOTE 4—The expression relating the surface area of the pore walls contained in volume ΔV_p is as follows:

$$\Delta S_p \text{ (m}^2\text{/g)} = 2 \times \frac{\Delta V_p}{\bar{r}_p} \times 10^4 \quad (6)$$

with ΔV_p in cubic centimetres and \bar{r}_p in angstroms. Converting cubic