



Designation: F 2212 – 02 (Reapproved 2007)<sup>ε1</sup>

## Standard Guide for Characterization of Type I Collagen as Starting Material for Surgical Implants and Substrates for Tissue-Engineered Medical Products (TEMPs)<sup>1</sup>

This standard is issued under the fixed designation F 2212; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

<sup>ε1</sup> NOTE—Formatting and grammar were corrected editorially throughout in April 2007.

### INTRODUCTION

Collagen-based medical devices are becoming more prevalent, especially in the area of soft tissue augmentation. The use of collagen in surgery dates back to the late 1800s, with the use of catgut sutures, human cadaveric skin, and fascia. More recently, collagen has been used in hemostatic sponges, dermal equivalents, injectables for soft tissue augmentation, as a matrix for cell-based products and as a vehicle for drug delivery. It is because of the versatility of collagen in medical applications that specific characterizations should be performed as a way to compare materials.

### 1. Scope

1.1 This guide for characterizing collagen-containing biomaterials is intended to provide characteristics, properties, and test methods for use by producers, manufacturers, and researchers to more clearly identify the specific collagen materials used. With greater than 20 types of collagen and the different properties of each, a single document would be cumbersome. This guide will focus on the characterization of Type I collagen, which is the most abundant collagen in mammals, especially in skin and bone. Collagen isolated from these sources may contain other types of collagen, for example, Type III and Type V. This guide does not provide specific parameters for any collagen product or mix of products or the acceptability of those products for the intended use. The collagen may be from any source, including, but not limited to, animal or cadaveric sources, human cell culture, or recombinant sources. The biological, immunological, or toxicological properties of the collagen may vary, depending on the source material. The properties of the collagen prepared from each of the above sources must be thoroughly investigated, as the changes in the collagen properties as a function of source materials is not thoroughly understood. This guide is intended to focus on purified Type I collagen as a starting material for surgical implants and substrates for Tissue-Engineered Medical

Products (TEMPs); some methods may not be applicable for gelatin nor for tissue implants. This guide may serve as a template for characterization of other types of collagen.

1.2 The biological response to collagen in soft tissue has been well documented by a history of clinical use (1, 2)<sup>2</sup> and laboratory studies (3, 4, 5, 21). Biocompatibility and appropriateness of use for a specific application(s) is the responsibility of the device manufacturer.

1.3 The following precautionary caveat pertains only to the test method portion, Section 5, of this guide. *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory requirements prior to use.*

### 2. Referenced Documents

2.1 *ASTM Standards:*<sup>3</sup>

**E 1298** Guide for Determination of Purity, Impurities, and Contaminants in Biological Drug Products

**F 619** Practice for Extraction of Medical Plastics

**F 720** Practice for Testing Guinea Pigs for Contact Allergens: Guinea Pig Maximization Test

<sup>1</sup> This guide is under the jurisdiction of ASTM Committee F04 on Medical and Surgical Materials and Devices and is the direct responsibility of Subcommittee F04.42 on Biomaterials and Biomolecules for TEMPs.

Current edition approved Feb. 1, 2007. Published February 2007. Originally approved in 2002. Last previous edition approved in 2002 as F 2212 – 02.

<sup>2</sup> The boldface numbers in parentheses refer to the list of references at the end of this standard.

<sup>3</sup> For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

- F 748** Practice for Selecting Generic Biological Test Methods for Materials and Devices
- F 749** Practice for Evaluating Material Extracts by Intracutaneous Injection in the Rabbit
- F 756** Practice for Assessment of Hemolytic Properties of Materials
- F 763** Practice for Short-Term Screening of Implant Materials
- F 813** Practice for Direct Contact Cell Culture Evaluation of Materials for Medical Devices
- F 895** Test Method for Agar Diffusion Cell Culture Screening for Cytotoxicity
- F 981** Practice for Assessment of Compatibility of Biomaterials for Surgical Implants with Respect to Effect of Materials on Muscle and Bone
- F 1251** Terminology Relating to Polymeric Biomaterials in Medical and Surgical Devices
- F 1439** Guide for Performance of Lifetime Bioassay for the Tumorigenic Potential of Implant Materials
- F 1903** Practice for Testing For Biological Responses to Particles *in vitro*
- F 1904** Practice for Testing the Biological Responses to Particles *in vivo*
- F 1905** Practice For Selecting Tests for Determining the Propensity of Materials to Cause Immunotoxicity
- F 1906** Practice for Evaluation of Immune Responses In Biocompatibility Testing Using ELISA Tests, Lymphocyte Proliferation, and Cell Migration
- F 1983** Practice for Assessment of Compatibility of Absorbable/Resorbable Biomaterials for Implant Applications
- F 2148** Practice for Evaluation of Delayed Contact Hypersensitivity Using the Murine Local Lymph Node Assay (LLNA)
- 2.2 ISO Standards:**<sup>4</sup>
- ISO 10993-1** Biological Evaluation of Medical Devices—Part 1: Evaluation and Testing
- ISO 10993-3—Part 3** Tests for Genotoxicity, Carcinogenicity and Reproductive Toxicity
- ISO 10993-9—Part 9** Framework for Identification and Quantification of Potential Degradation Products
- ISO 10993-10** Biological Evaluation of Medical Devices—Part 10: Tests for Irritation and Delayed-Type Hypersensitivity
- ISO 10993-17—Part 17** Methods for Establishment of Allowable Limits for Leachable Substances Using Health-Based Risk Assessment
- ISO 13408-1: 1998** Aseptic Processing of Health Care Products—Part 1: General Requirements
- 2.3 EN (European Norm) Documents:**<sup>5</sup>
- EN 12442-1** Animal Tissues and their Derivatives Utilized in the Manufacture of Medical Devices—Part 1: Analysis and Management of Risk
- EN 12442-2—Part 2** Controls on Sourcing, Collection and Handling
- EN 12442-3—Part 3** Validation of the Elimination and/or Inactivation of Virus and Transmissible Agents
- 2.4 U. S. Pharmacopeia Documents:**<sup>6</sup>
- United States Pharmacopeia (USP), Edition XXIV (24)**
- USP 24/NF 19** Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin
- 2.5 Code of Federal Regulations:**<sup>7</sup>
- Code of Federal Regulations**, Title 21, Part 820
- Federal Register Vol. 43**, No. 141, Friday, July 21, 1978
- Human Cells, Tissues and Cellular and Tissue-Based Products**, Establishment Registration and Listing. 21 CFR Parts 207, 807, and 1271
- Federal Register/Vol. 66**, No. 13, Jan 19, 2001/Rules and Regulations, page 5447
- Suitability Determination for Donors of Human Cell and Tissue-based Products**, 21 CFR 1271 Part C, Proposed Rule
- Current Good Tissue Practice for Manufacturers of Human Cellular and Tissue-Based Products**, Inspection and Enforcement. Proposed Rule. Federal Register/Vol. 66, No. 5/January 8, 2001/Proposed Rules, pages 1552-1559
- Guidance for Screening and Testing of Donors of Human Tissue Intended for Transplantation**, Availability. Federal Register/Vol. 62, No. 145/July 29, 1997/NoticesDraft Guidance for Preclinical and Clinical Investigations of Urethral Bulking Agents used in the Treatment of Urinary Incontinence. November 29, 1995. (ODE/DRARD/ULDB), Document No. 850
- Guidance for Industry and for FDA Reviewers**, Medical Devices Containing Materials Derived from Animal Sources (Except for *In Vitro* Diagnostic Devices), November 6, 1998, U.S. Department of Health and Human Services, Food and Drug Administration, Center for Devices and Radiological Health
- CFR 610.13(b)**, Rabbit Pyrogen Assay
- 2.6 ICH Documents:**<sup>8</sup>
- International Conference on Harmonization (1997)** Guidance for Industry M3 Nonclinical Safety Studies for the Conduct of Human Clinical Trials for Pharmaceuticals 62 FR 62922
- International Conference on Harmonization (1996)** Guideline for Industry S2A Specific Aspects of Regulatory Genotoxicity Tests for Pharmaceuticals. 61 FR 18199
- International Conference on Harmonization (1997)** Guidance for Industry S2B Genotoxicity: A Standard Battery for Genotoxicity Testing of Pharmaceuticals 62 FR 62472
- International Conference on Harmonization (1994)** Guideline for Industry S5A Detection of Toxicity to Reproduction for Medicinal Products. 59 FR 48746

<sup>4</sup> Available from International Organization for Standardization (ISO), 1 rue de Varembe, Case postale 56, CH-1211, Geneva 20, Switzerland.

<sup>5</sup> Available from European Committee for Standardization, CEN Management Centre 36, rue de Stassart B-1050 Brussels, Belgium.

<sup>6</sup> Available from U.S. Pharmacopeia (USP), 12601 Twinbrook Pkwy., Rockville, MD 20852.

<sup>7</sup> Available from U.S. Government Printing Office Superintendent of Documents, 732 N. Capitol St., NW, Mail Stop: SDE, Washington, DC 20401.

<sup>8</sup> ICH Secretariat, c/o IFPMA, 30 rue de St-Jean, P.O. Box 758, 1211 Geneva 13, Switzerland.

**International Conference on Harmonization (1996)** Guidance for Industry S5B Detection of Toxicity to Reproduction for Medicinal Products: Addendum on Toxicity to Male Fertility. 61 FR 15360

**International Conference on Harmonization (1996)** Guideline for Industry S1A The Need for Long-term Rodent Carcinogenicity Studies of Pharmaceuticals. 61 FR 8153

**International Conference on Harmonization (1998)** Guidance for Industry S1B Testing for Carcinogenicity of Pharmaceuticals. 63 FR 8983

**International Conference on Harmonization (1995)** Guideline for Industry S1C Dose Selection for Carcinogenicity Studies of Pharmaceuticals. 60 FR 11278

**International Conference on Harmonization (1997) S1C(R)** Guidance for Industry Addendum to Dose Selection for Carcinogenicity Studies of Pharmaceuticals: Addition of a Limit Dose and Related Notes. 62 FR 64259

**International Conference on Harmonization (ICH) Q1A ICH** Harmonized Tripartite Guidance for Stability Testing of New Drug Substances and Products (September 23, 1994)

**U.S. Food and Drug Administration (FDA and Committee for Proprietary Medicinal Products (CPMP), 1998** International Conference on Harmonization (ICH), Quality of Biotechnological Products: Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin, Consensus Guideline ICH Viral Safety Document: Step 5

2.7 *FDA Documents*.<sup>9</sup>

**FDA** Guideline on Validation of the Limulus Amebocyte Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products and Healthcare Products, DHHS, December 1987

**U.S. Food and Drug Administration (FDA) Center for Biologics Evaluation and Research (CBER), 1993** Points to Consider in the Characterization of Cell Lines Used to produce Biologicals

**U.S. Food and Drug Administration (FDA) Center for Biologics Evaluation and Research (CBER), 1997** Points to Consider in the Manufacture and Testing of Monoclonal Antibody Products for Human Use, 94D-0259

**FDA** Interim Guidance for Human and Veterinary Drug Products and Biologicals, Kinetic LAL techniques, DHHS, July 15, 1991

2.8 *AAMI Documents*.<sup>10</sup>

**ANSI/AAMI/ISO 11737-1: 1995** Sterilization of Medical Devices—Microbiological Methods—Part 1: Estimation of Bioburden on Product

**ANSI/AAMI/ISO 11737-2: 1998** Sterilization of Medical Devices—Microbiological Methods—Part 2: Tests of Sterility Performed in the Validation of a Sterilization Process

**AAMI TIR No. 19-1998** Guidance for ANSI/AAMI/ISO 10993-7: 1995, Biological Evaluation of Medical

Devices—Part 7: Ethylene Oxide Sterilization Residuals  
**AAMI/ISO 14160-1998** Sterilization of Single-Use Medical Devices Incorporating Materials of Animal Origin—Validation and Routine Control of Sterilization by Liquid Chemical Sterilants

**AAMI ST67/CDV-2: 1999** Sterilization of Medical Devices—Requirements for Products Labeled “Sterile”

2.9 *Other References*:

Draft Guidance for Preclinical and Clinical Investigations of Urethral Bulking Agents Used in the Treatment of Urinary Incontinence, November 29, 1995. (ODE/DRARD/ULDB), Document No. 850.<sup>11</sup>

### 3. Terminology

3.1 *Definitions*:

3.1.1 *adventitious agents, n*—an unintentionally introduced microbiological or other infectious contaminant. In the production of TEMPs, these agents may be unintentionally introduced into the process stream or the final product, or both.

3.1.2 *biocompatibility, n*—a material may be considered biocompatible if the materials perform with an appropriate host response in a specific application (22).

3.1.3 *collagen, n*—Type I collagen is a member of a family of structural proteins found in animals. Type I collagen is part of the fibrillar group of collagens. It derives from the COL1A1 and COL1A2 genes, which express the alpha chains of the collagen. All collagens have a unique triple helical structure configuration of three polypeptide units known as alpha-chains. Proper alignment of the alpha chains of the collagen molecule requires a highly complex enzymatic and chemical interaction *in vivo*. As such, preparation of the collagen by alternate methods may result in improperly aligned alpha chains and, putatively, increase the immunogenicity of the collagen. Collagen is high in glycine, L-alanine, L-proline, and 4-hydroxyproline, low in sulfur, and contains no L-tryptophan. Natural, fibrillar Type I collagen is normally soluble in dilute acids and alkalis. When heated (for example, above approximately 40°C), collagen is denatured to single alpha chains (gelatin). At each end of the chains are short non-helical domains called telopeptides, which are removed in some collagen preparations. Through non-covalent interactions with sites on adjacent helices, fibrillogenesis is achieved. Subsequently, non-reducible cross links are formed. Type I collagen can be associated with Type III and Type V collagen and also with the other non-collagenous proteins like elastin and other structural molecules like glycosaminoglycans and complex lipoproteins and glycoproteins.

3.1.4 *degradation, n*—change in chemical, physical, or molecular structure or appearance (that is, gross morphology) of material.

3.1.5 *endotoxin, n*—a high-molecular-weight lipopolysaccharide (LPS) complex associated with the cell wall of gram-negative bacteria that is pyrogenic in humans.

3.1.5.1 *Discussion*—Though endotoxins are pyrogens, not all pyrogens are endotoxins.

<sup>9</sup> Available from the Food and Drug Administration, 5600 Fishers Ln., Rockville, MD 20857.

<sup>10</sup> Association for the Advancement of Medical Instrumentation, 1110 N. Glebe Rd., Suite 220, Arlington, VA 22201-4795.

<sup>11</sup> Available from the FDA, 5600 Fishers Ln., Rockville, MD 20857. <http://www.fda.gov/cdrh/ode/oderp850.html>.

3.1.6 *microorganisms, n*—bacteria, fungi, yeast, mold, viruses, and other infectious agents.

3.1.6.1 *Discussion*—Not all microorganisms are infectious or pathogenic.

3.1.7 *solubility, n*—a measure of the extent to which the material can be dissolved.

3.1.7.1 *Discussion*—In the context of collagen, refers to the dissociation of the fibrillar aggregates of collagen molecules into a solution. Native Type I collagen which is soluble in dilute acids, but not soluble in neutral pH conditions is termed “insoluble” or “acid-soluble,” while simple aggregates of non-fibrillar collagen soluble in neutral salt solutions are termed “neutral salt-soluble.” Post translational surface charge modifications may alter the solubility of collagen in neutral pH condition.

3.1.8 *sterilization, n*—the destruction or removal of all microorganisms in or about an object.

3.1.8.1 *Discussion*—Examples are by chemical agents, electron beam, gamma irradiation, ultraviolet (UV) exposure, or filtration.

3.1.9 *suspension, n*—the dispersion of a solid through a liquid with a particle size large enough to be detected by purely optical means.

#### 4. Significance and Use

4.1 The objective of this guide is to provide guidance in the characterization of Type I collagen as a starting material for surgical implants and substrates for TEMPs. This guide contains a listing of physical and chemical parameters that are directly related to the function of collagen. This guide can be used as an aid in the selection and characterization of the appropriate collagen starting material for the specific use. Not all tests or parameters are applicable to all uses of collagen.

4.2 The collagen covered by this guide may be used in a broad range of forms or applications, for example (but not limited to) medical devices, tissue-engineered medical products (TEMPs) or cell-, drug-, or DNA-delivery devices for implantation. The use of collagen in a practical application should be based, among other factors, on biocompatibility and physical test data. Recommendations in this guide should not be interpreted as a guarantee of clinical success in any tissue-engineered medical product or drug delivery application.

4.3 The following general areas should be considered when determining if the collagen supplied satisfies requirements for use in TEMPs. These are source of collagen, chemical and physical characterization and testing, and impurities profile.

#### 5. Chemical and Physical Characterizations

5.1 These methods are suggested assays; however, other validated assay methods may be used. Selection of assay systems will vary, depending on the configuration of the collagen (that is, soluble or insoluble). The user should ensure that the method selected is reliable and commonly accepted in protein chemistry. A review of collagen materials may be found in Li, 2000 (7), while a review of the collagen family of proteins may be found in Refs (9-13 and 20). When selecting an appropriate test method, the user should note that impurities in highly purified collagen fibers can be as low as 1 to 2 %, while test method uncertainties can be as high as 10 to 15 %

between samples. For soluble collagen, the following represents a non-inclusive list of assay systems available: Sodium Dodecyl Sulfate Polyacrylamide Gel Electrophoresis (SDS-PAGE); peptide mapping; and physico-chemical analysis. A similar list for insoluble collagen may include, but not be limited to, assay methods for hexosamine (that is, detection of glycoproteins), lipid, total sugar, desmosine (that is, elastin), and amino acid composition (that is, collagen composition profile; non-collagenous amino acids). Additionally, methods such as transmission electron microscopy may be helpful in characterizing the collagen.

5.2 The concentration of collagen should be expressed in mass/volume or mass/mass. Colorimetric assays or amino acid analysis for hydroxyproline are commonly used methods to measure collagen content.

5.3 Amino acid analysis will provide information on the composition of the amino acids of collagen (that is, the amino acids must be within the range of published data for highly purified collagen preparations, generally in the form of acid soluble). The collagen must not contain amino acids that are not in collagen, such as tryptophan. Amino acid analysis is routinely performed on hydrolyzed collagens by reverse phase High Performance Liquid Chromatography (HPLC). Amino acid analysis should contain no cysteine (sulfur-containing amino acid), low in tyrosine, and no tryptophan. This method can also be used for the quantification of hydroxyproline. There are other methods available for amino acid analysis.

5.4 Purity of soluble collagen can be analyzed by SDS-PAGE either on the collagen directly or after digestion of the collagen by purified bacterial collagenase to detect any remaining proteins by SDS-PAGE.

5.5 *Elastin Assay*—Elastin can be a component of the impurities in an insoluble collagen preparation. One method to assay for elastin, although other methods are available, involves the detection of desmosine (6). These impurities can be detected by Western blots, ELISAs, and other types of assays.

5.6 Peptide mapping is one possible method to identify Type I collagen. The most commonly used peptide mapping method utilizes Cyanogen Bromide (CNBr) digest. The digest can be analyzed by SDS-PAGE or HPLC.

5.7 *Impurities Profile*—The term impurity relates to the presence of extraneous substances and materials in the collagen introduced during the purification process. These impurities can be detected by Western blots, ELISAs, and other types of assays. The user is also directed to Guide E 1298 for additional information. If there is a concern for the presence of processing aids or other impurities associated with the collagen, they should be addressed with the supplier. The major impurities of concern include, but are not limited to the following: endotoxins, glycosaminoglycans, elastin, lipids, improperly aligned collagen molecules, host cell contaminants, cell culture contaminants, heavy metals, bioburden, viruses, transmissible spongiform encephalopathy (TSE) agents, cross-linking, and enzymatic agents. While not considered an impurity, Type III collagen may also be associated with Type I collagen. At minimum, any protein impurity of greater than 1 % in the final collagen preparation should be identified and quantified.

5.8 *Endotoxin Content*—Endotoxin contamination is difficult to prevent because it is ubiquitous in nature, stable and small enough to pass through sterilizing filters (0.22 µm). Endotoxin tests for collagen include the gel clot, endpoint assay and the kinetic assay. The gel clot test is the simplest and easiest of the *Limulus* ameobocyte lysate (LAL) test methods, although much less sensitive than the kinetic assay. The quantitative kinetic assay, which measures the amount of time required to reach a predetermined optical density, is the most sensitive (Food and Drug Administration, Guideline on Validation of the *Limulus* Ameobocyte Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products and Healthcare Products.) Each new lot of reagents should be validated. The endotoxin level in collagen will ultimately be critical to its use in biomedical applications where there are regulatory limits to the amount of endotoxin that can be implanted into humans. Relevant **FDA guidance for allowable levels and information regarding validation of endotoxin assays should be consulted if human trials are contemplated (Interim Guidance for Human and Veterinary Drug Products and Biologicals)**. The user is also directed to CFR 610.13(b) for information pertaining to the rabbit pyrogen assay.

5.9 *Heavy Metal Content by the USP Method*—This test is provided to demonstrate that the content of heavy metal impurities does not exceed a limit in the individual product specification. This method is based on <231> Heavy Metals, 1st and 6th Supplement USP-NF. Substances that typically respond to this test are lead, mercury, bismuth, arsenic, antimony, tin, cadmium, silver, copper, and molybdenum. Under the specified test conditions, the limit is determined by a concomitant visual comparison of metals that are colored by sulfide ion with a control prepared from a Standard Lead Solution. Additional heavy metal contaminants may be present due to processing. If necessary, the user may detect these contaminants by various methods, including not limited to, spectrographic, chromatographic, and flame atomic absorption.

5.10 *Microbiological Safety*—The presence of bacteria, yeast, and mold are also impurities that can arise in an biological sample. User will validate sterilization and characterize effect(s) on product. The presence of bacteria may also contribute to the presence of endotoxins. The following Microbiological Tests in USP 24 are of particular relevance: Microbial Limit Tests <61>, Sterility Tests <71>, Sterilization and sterility assurance of compendial articles <12211>, and the Biological Tests and Assays: Bacterial Endotoxins Tests <85>. The user should also consider other relevant standards, such as, but not limited to, Association for the Advancement of Medical Instrumentation (AAMI) standards and international standards, of which the following are examples: **ANSI/AAMI/ISO 11737-1: 1995 ; ANSI/AAMI/ISO 11737-2: 1998; and ISO 13408-1: 1998**. The collagen is first dissolved in a sterile, aqueous solution, then filtered using sterile techniques through a 0.45-µm membrane filter. The filters are subsequently incubated on Tryptic Soya Agar to determine the presence of bacteria, and on Sabouraud Dextrose Agar to determine the presence of yeast and mold. If collagen products are intended

to serve as a barrier to microorganisms, this function will need to be validated with specific experiments.

5.11 Carbohydrate analysis of collagens can be carried out by classical gas-liquid chromatographic methods or spectrophotometric methods. If novel source, then glycosylation of proteins may need to be considered.

5.12 Trypsin Susceptibility will detect that portion of collagen that has been denatured from the purification steps such as acid and base treatment, solvent treatment, and so forth. Trypsin will digest that portion of the collagen and can be measured by assaying the hydroxyproline content of the supernatant. Triple helical collagen is resistant to digestion by most proteases. Susceptibility to trypsin or other appropriate proteases is determined by exposing the collagen to the enzyme and assaying the digest for degradation. There are several methods for this test.

5.13 Differential Scanning Calorimetry (DSC) determines dissociation temperature of collagens in fibrils, as well as detecting microfibrils and denatured collagen at lower melting temperatures.

5.14 Viscosity is more applicable to gels or suspensions but may be useful with collagen configured in forms such as, but not limited to, pastes or films (8). Viscosity of collagen-based materials depends on a number of factors, but is not limited to the following: solution or dispersion/suspension, concentration, temperature, operating condition, and so forth. It is not feasible to determine the viscosity of films. This is a routine test performed with a viscometer (not a u-tube). The user must clearly state the conditions of the test.

5.15 Transmission electron microscopy may be used to show the quality of collagen fibers. Unraveling or changes in banding will be obvious.

5.16 DNA sequence data on recombinant or transgenic source cells: Verify sequence data for expression gene, that is, COL1A1 or COL1A2.

5.17 The collagen material shall have specifications for an extensive set of chemical and physical properties such as, but not limited to, those listed in **Table 1**. The table represents methods which may or may not be appropriate for characterizing a particular collagen sample. Not all the methods listed may be required to characterize the sample and the specificity and sensitivity vary among the methods listed. The user should be familiar with the limitations of the appropriate test methods.

## 6. Product Development Considerations

6.1 *Storage Conditions/Shelf Life Stability of Collagen*—For collagen, the most relevant stability-indicating parameters are those related to the functionality of the polymer. Depending on what function the collagen will have in the final formulation, parameters such as viscosity (apparent and intrinsic) and biological activity may also be considered. Storage conditions are of importance, especially for collagen solutions. **International Conference on Harmonization (ICH) guidance documents should be consulted for information on stability testing of pharmaceuticals (that is, ICH Q1A ICH Harmonized Tripartite Guideline for Stability Testing of New Drug Substances and Products)**.

**TABLE 1 Characterization Methods for Type I Collagen**

| Characterization Method   | Applicable to        |
|---|----------------------|
| <b>Chemical</b>   |                      |
| Appearance  | Soluble or Insoluble |
| Concentration   | Soluble or Insoluble |
| Purity  | Soluble or Insoluble |
| Amino acid analysis   | Soluble or Insoluble |
| Peptide mapping   | Soluble or Insoluble |
| Impurities profile, includes Heavy Metal Analysis               | Soluble or Insoluble |
| Carbohydrate analysis   | Soluble or Insoluble |
| Trypsin resistance  | Soluble or Insoluble |
| pH of implantable   | Soluble or Insoluble |
| Additives (cross-linkers, lubricants, drugs, sterilents)        | Soluble or Insoluble |
| <b>Physical</b>   |                      |
| Shrink Temperature (DSC)  | Insoluble            |
| Viscosity   | Mainly soluble       |
| TEM   | Insoluble            |
| SDS-PAGE  | Soluble or Insoluble |
| Moisture Content (5 to 20 %), dependent on storage environment  | Insoluble            |
| Electron Micrograph (native banded 640 Å structure for fibrils) | Insoluble            |
| <b>Biochemical</b>  |                      |
| Endotoxin level   | Soluble or Insoluble |
| Bioburden   | Soluble or Insoluble |
| % Type I collagen/Total Protein                                 | Soluble or Insoluble |
| % Other Types Collagen and List of which Types present          | Soluble or Insoluble |
| Total DNA (ppm or %)  | Soluble or Insoluble |
| Total Lipid   | Soluble or Insoluble |
| % native collagen (by trypsin resistance, circular Dichroism)   | Soluble or Insoluble |

Abbreviation in Table:

DSC = Differential Scanning Calorimetry

TEM = Transmission Electron Microscopy

SDS-PAGE = Sodium Dodecyl Sulfate Polyacrylamide Gel Electrophoresis

**6.2 Sterilization Method, if Applicable, and Effects of Sterilization on Product**—The user should verify that the sterilization method does not adversely effect the collagen end product. Collagen can be sterilized by gamma irradiation (with subsequent degradation of the collagen) or by ethylene oxide (with potential residuals), or prepared using aseptic processing steps. Solutions of collagen may be (1) filter sterilized if the viscosity of the collagen solution permits; (2) gamma-irradiated with a resulting loss in viscosity (molecular weight). Selection of the method of sterilization will depend upon the viscosity or molecular weight needs of the final application. Use of ethylene oxide will also require testing for residuals for ethylene oxide. The reader should refer to the relevant standards regarding the sterilization of healthcare products by radiation, steam and ethylene oxide gas, such as **AAMI TIR No. 19-1998**; **AAMI/ISO 14160-1998**; and **AAMI ST67/CDV-2: 1999**.

**6.3 Sourcing**—The criteria to consider for safe sourcing include appropriate human or animal donor selection and the tissue collection procedures to assure that the source material is unlikely to contain TSE infectivity. Additional information can be obtained from the following documents: **EN 12442-1**, **EN 12442-2**, **EN 12442-3**; **21 CFR Parts 207, 807, 820, and 1271 Part C**; **Federal Register Vol. 43**; **Federal Register Vol. 62**; **Federal Register Vol. 66, No. 5, January 8, 2001, Pages 1552-59**; **Federal Register, Vol. 66, No. 13, January 19, 2001, Page 5447**; **ISO 13408, 1998**. Additional documents may be

available. The user should consult the most current version of the document. The collagen can be isolated from tissues or cell cultures by any method, including, but not limited to, extraction by dilute acids or dilute salt solutions or by enzymatic digestion of the tissue (**14-17**). The user should be aware that even though Type III collagen is less abundant, it is often associated with Type I, except in bones and tendons. Type V collagen is also associated with Type I.

**6.4 Viral and Transmissible Spongiform Encephalopathy (TSE) Agent Inactivation**—Viruses and TSE agents can be introduced into a product as a result of raw materials sourcing or through adventitious means. Appropriate measures should be taken so that the resultant product is free from viruses and TSE agents. For further guidance on viral or TSE clearance, or both, the user is directed to the references throughout this guide as well as **USP 24/NF 19 <1050>**, and other pertinent references, as appropriate. Additional information may be found in the following FDA Guidance Document, FDA points to consider and International Conference on Harmonization (ICH) documents: **Guidance for Industry and for FDA Reviewers: Medical Devices Containing Materials Derived from Animal Sources (Except for *In Vitro* Diagnostic Devices)**; **U.S. Food and Drug Administration (FDA and Committee for Proprietary Medicinal Products (CPMP), 1998**; **“International Conference on Harmonization (ICH), Quality of Biotechnological Products: Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin.” Consensus Guideline ICH Viral Safety Document: Step 5**; **U.S. Food and Drug Administration (FDA) Center for Biologics Evaluation and Research (CBER), 1993**; **“Points to Consider in the Characterization of Cell Lines Used to Produce Biologicals; U.S. Food and Drug Administration (FDA) Center for Biologics Evaluation and Research (CBER), 1997**; **“Points to Consider in the Manufacture and Testing of Monoclonal Antibody Products for Human Use.” 94D-0259**.

**6.4.1 Viral Clearance**—To reduce or eliminate potential infectivity, the user should exclude sources of raw materials from humans or animals likely to be infected with known viral pathogens. Viral clearance methods can include, but not be limited to, methods such as high or low pH, urea treatment, chemical treatments, and filtration. However, even these harsh treatments may not ensure complete viral inactivation. Viral clearance should be demonstrated by an appropriately validated viral clearance study protocol. The user should verify that the viral clearance procedure is compatible with the starting material or the configured end product. For human tissue sources for manufacturing into collagen, the observance of Good Tissue Practices should be considered.

**6.4.2 TSE Clearance**—Due diligence should be made to the sourcing of raw materials, process design to remove potential TSE agents and treatments to inactivate TSE agents for those products which can withstand the harsh treatments required to inactivate TSE agents. The user is referred to the **“Meeting Report, International Workshop on Clearance of TSE Agents from Blood Products and Implanted Tissues,” (23)** and the FDA Guidance Document **“Guidance for Industry and for FDA Reviewers: Medical Devices Containing Materials Derived from Animal Sources (Except for *In Vitro* Diagnostic Devices)”**