



SLOVENSKI STANDARD
oSIST prEN 14031:2020
01-april-2020

Izpostavljenost na delovnem mestu - Kvantitativno določevanje lebdečih endotoksinov

Workplace exposure - Quantitative measurement of airborne endotoxins

Exposition am Arbeitsplatz - Quantitative Messung von luftgetragenen Endotoxinen

Exposition sur les lieux de travail - Mesure quantitative des endotoxines aéroportées

ITeH STANDARD PREVIEW
(standards.iteh.ai)

Ta slovenski standard je istoveten z: prEN 14031

[ksIST FprEN 14031:2021](https://standards.iteh.ai/catalog/standards/sist/8521eef9-9e37-46d1-8c5f-d2f985c7444b/ksist-fpren-14031-2021)

<https://standards.iteh.ai/catalog/standards/sist/8521eef9-9e37-46d1-8c5f-d2f985c7444b/ksist-fpren-14031-2021>

ICS:

13.040.30 Kakovost zraka na delovnem mestu Workplace atmospheres
mestu

oSIST prEN 14031:2020

en,fr,de

iTeh STANDARD PREVIEW
(standards.iteh.ai)

[kSIST FprEN 14031:2021](#)

<https://standards.iteh.ai/catalog/standards/sist/8521eef9-9e37-46d1-8c5f-d2f985c7444b/ksist-fpren-14031-2021>

EUROPEAN STANDARD
NORME EUROPÉENNE
EUROPÄISCHE NORM

DRAFT
prEN 14031

March 2020

ICS 13.040.30

Will supersede EN 14031:2003

English Version

Workplace exposure - Quantitative measurement of airborne endotoxins

Exposition sur les lieux de travail - Mesure quantitative des endotoxines aéroportées

Exposition am Arbeitsplatz - Quantitative Messung von luftgetragenen Endotoxinen

This draft European Standard is submitted to CEN members for parallel enquiry. It has been drawn up by the Technical Committee CEN/TC 137.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

<https://standards.iteh.ai/catalog/standards/sist/8521ee9-9e37-46d1-8c5f-d2985c744414/sist-fren-14031-2021>

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

Warning : This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

Contents	Page
European foreword.....	3
Introduction	4
1 Scope	5
2 Normative references.....	5
3 Terms and definitions	5
4 Requirements	6
4.1 Operator skills	6
4.2 Sampling	7
4.2.1 General.....	7
4.2.2 Equipment.....	7
4.2.3 Filters	7
4.3 Sample transport.....	7
4.4 Sample storage at the laboratory	7
4.5 Glassware	7
4.6 Extraction liquids, water and detergents.....	8
4.7 Test for contamination.....	8
4.8 Sampling documentation.....	8
5 Extraction	9
5.1 General.....	9
5.2 Equipment.....	9
5.3 Procedure	9
5.4 Storage conditions	9
5.5 Documentation of extraction.....	10
6 Analytical method.....	10
6.1 General.....	10
6.2 Equipment.....	10
6.3 Procedure	11
6.3.1 Measurement	11
6.3.2 Calibration.....	11
6.3.3 Calculation.....	11
6.4 Validation	11
6.5 Documentation of analysis	11
7 Expression of results	12
8 Precision	12
9 Test report	12
Annex A (informative) Properties of endotoxins and sources of occupational exposure	14
Annex B (informative) Principle of LAL-assays	16
Annex C (informative) Deposition of endotoxins on inner sampler walls.....	18
Bibliography.....	19

European foreword

This document (prEN 14031:2020) has been prepared by Technical Committee CEN/TC 137 “Assessment of workplace exposure to chemical and biological agents”, the secretariat of which is held by DIN.

This document is currently submitted to the CEN Enquiry.

This document will supersede EN 14031:2003.

The major technical changes between this European Standard and the previous edition are as follows:

- a) document title adjusted to the wording used in the Scope;
- b) terms and definitions taken over from EN 13098:2019, where appropriate;
- c) Clause 4 restructured;
- d) requirements for sample storage at the laboratory rewritten;
- e) recommendations for extraction of samples revised;
- f) storage conditions described more precisely;
- g) further requirements on documentation of extraction added;
- h) new Table A.1 on examples of different occupational environments where endotoxins are present added;
- i) whole document editorially revised.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

oSIST prEN 14031:2020
<https://standards.iteh.ai/catalog/standards/sist/8521ee9-9e37-46d1-8c5f-d2f985c7444b/ksist-pr-en-14031-2021>

Introduction

The term 'endotoxin' refers to lipopolysaccharides that are present in the outer membrane of the cell walls of Gram-negative bacteria. These lipopolysaccharides are a class of pure lipid carbohydrate molecules (free of protein and other cell wall components) that are held responsible for most of the biological properties characteristic of bacterial endotoxins.

Endotoxins play an important role in the development of organic dust-related symptoms (for example, inflammatory reactions in the airways and/or systemic reactions) among workers exposed via inhalation.

To investigate ill health resulting from occupational exposure to airborne endotoxins it is important therefore to measure their presence in workplace air accurately. At present, different methods are used to measure airborne endotoxins. Standardization with respect to sampling, transportation, extraction, analytical methods and storage of samples or extracts is important in order to obtain accurate and comparable results and reduce uncertainties in exposure assessment.

iTeh STANDARD PREVIEW (standards.iteh.ai)

[kSIST FprEN 14031:2021](https://standards.iteh.ai/catalog/standards/sist/8521eef9-9e37-46d1-8c5f-d2f985c7444b/ksist-fpren-14031-2021)

<https://standards.iteh.ai/catalog/standards/sist/8521eef9-9e37-46d1-8c5f-d2f985c7444b/ksist-fpren-14031-2021>

1 Scope

This document specifies methods for the quantitative measurement of airborne endotoxins and gives general requirements for sampling on filters, transportation, storage as well as the analysis of samples.

This document provides also guidelines for the assessment of workplace exposure to airborne endotoxins.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 481, *Workplace atmospheres — Size fraction definitions for measurement of airborne particles*

EN 1540, *Workplace exposure — Terminology*

EN ISO 13137, *Workplace atmospheres — Pumps for personal sampling of chemical and biological agents - Requirements and test methods (ISO 13137)*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in EN 1540 and the following apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

— IEC Electropedia: available at <http://www.electropedia.org/>

— ISO Online browsing platform: available at <https://www.iso.org/obp>

Note 1 to entry: In particular, the following terms used in this document are defined in EN 1540: aerosol, biological agent, bioaerosol, exposure (by inhalation), inhalable fraction, measuring procedure, personal sampling, static sampling and workplace.

3.1

control standard endotoxin

CSE

standard that is traceable to the reference standard endotoxin (RSE)

3.2

endotoxin

constituent of the external membrane of Gram-negative bacteria (lipopolysaccharide), consisting of a complex lipid, lipid A, which is covalently bound to a polysaccharide

Note 1 to entry “Free endotoxin” is liberated after cell death and by budding from living cells. Lipid A is the active (toxic) part and is a potent pro-inflammatory substance and can induce febrile, bronchial and other symptoms in exposed workers. The composition and the toxicity of endotoxin differ between bacteria species

Note 2 to entry: A brief overview with respect to physical and chemical properties of endotoxins and sources of occupational exposure is given in Annex A.

[SOURCE: EN 13098:2019, 3.9 – modified, Note 2 to entry added]

prEN 14031:2020 (E)**3.3****endotoxinfree liquid**

water or solution containing no measurable endotoxin

Note 1 to entry Endotoxinfree and pyrogenfree are used interchangeable in the literature

Note 2 to entry: In practise, less than 0,2 endotoxin units per millilitre is often used as threshold value for endotoxinfree.

3.4**endotoxin unit**

EU

unit standardized against the defined reference material, reference standard endotoxin

[SOURCE: EN 13098:2019, 3.10]

3.5**Limulus Amoebocyte Lysate**

LAL

enzymes extracted from the blood cells of the horse shoe crab (*Limulus polyphemus*) that are activated by endotoxin and other molecules (glucans etc.)

[SOURCE: EN 13098:2019, 3.16]

3.6**Limulus Amoebocyte Lysate-assay (standards.iteh.ai)**

LAL-assay

functional assay to measure endotoxin concentrations

Note 1 to entry The method is highly sensitive and based on the activation of a clotting enzyme present in the lysate of hemolymph of the horseshoe crab (*Limulus polyphemus*).

3.7**lipopolysaccharide**

LPS

water-soluble and stable molecule composed of lipid and polysaccharide present in Gram-negative bacteria

Note 1 to entry The lipid part of LPS is termed 'lipid A' and is essential for the toxic properties of LPS

Note 2 to entry The terms 'endotoxin' and 'lipopolysaccharide' are often used interchangeably in scientific literature regarding occupational health.

3.8**reference standard endotoxin**

RSE

purified lipopolysaccharide from *Escherichia coli* that serves as an international reference standard

4 Requirements**4.1 Operator skills**

The operator shall be trained in aseptic techniques to avoid contamination of the sample during any measurement phase. This includes knowledge about sampling equipment and how to perform the sampling.

4.2 Sampling

4.2.1 General

Measurement of airborne endotoxins should preferably be performed by sampling of the inhalable fraction of aerosols. Aerosols shall be sampled on a filter using a pump to draw air through the filter.

Personal sampling shall be used for risk evaluation including comparison with exposure guidelines.

Static sampling can be used to identify sources of exposure and to measure background concentrations of airborne endotoxins.

NOTE Other comparable sampling methods can be used for evaluation of other specific parameters, for example, to measure endotoxins in different size distributions or particle size fractions, or for short- or long-term exposure measurement.

4.2.2 Equipment

The pumps used for personal sampling and static sampling shall fulfil the requirements specified in EN ISO 13137. Personal samplers shall collect the inhalable fraction of aerosols as specified in EN 481.

The filter holder or cassette should be rendered endotoxinfree. Heat-resistant reusable samplers shall be treated as indicated in 4.3. Non heat-resistant reusable samplers shall be cleaned thoroughly so that contamination of the samples is avoided. Cleaning agents shall be removed by subsequent washing in endotoxinfree water. Disposable samplers do not need to be sterilized when it is assumed that the manufacturing process is endotoxinfree. Possible contamination shall be controlled by using blank filters (see 4.7).

4.2.3 Filters

Binderfree glass fibre filters should be used for sampling of airborne endotoxins because these filters perform very well regarding their collection efficiency, blank values and recovery. If any other filter is used the performance of these parameters shall be comparable.

NOTE Several filter materials (for example, polyvinylchloride, polytetrafluoroethylene and polycarbonate) are commonly used but some materials adsorb endotoxins and thereby give false low values.

4.3 Sample transport

After sampling, filters shall be prepared for transport by sealing in a container or sampler cassette. The sample shall be transported in conditions that prevent large variations in temperature and the gaining of moisture. If transport time is above 24 h, or there is risk of large temperature variations, actions to prevent moisture shall be taken, either by using an isolating container or by including a dehumidifier.

Conditions during transport shall be documented, including transport time and temperature log.

If the temperature is above 30 °C the filters shall be placed in cooler conditions (for example, by using a cooling box).

4.4 Sample storage at the laboratory

If the sample is not extracted within two weeks after arrival at the laboratory the sample shall be stored in conditions that prevent large variations in temperature and gain of moisture. Alternatively, for storage of more than two weeks they should be frozen at a temperature of about -20 °C or below, because endotoxins are very stable when frozen. However, samples shall not be frozen and thawed more than once, as this can affect the detectable endotoxin content of them.

4.5 Glassware

Glassware shall be rendered endotoxinfree, for example, by heating to at least 180 °C for about 4 h.

prEN 14031:2020 (E)**4.6 Extraction liquids, water and detergents**

All liquids used shall be endotoxinfree.

4.7 Test for contamination

To test for contamination a series of at least two blank filters shall be included in each sampling session. Except for the actual sampling, blank filters shall be treated equal to the other filters that are used for sampling.

4.8 Sampling documentation

The sampling operations carried out shall be documented to obtain comparable and reliable concentration values of endotoxins.

The sampling documentation shall include at least the following information:

- a) name of the organization and person performing the sampling;
- b) date of sampling;
- c) purpose of sampling;
- d) a unique identifier code for the sample and sampling device (for example, sampling pump, used);
- e) name and address of the company where the sampling was carried out, or a unique identifier to preserve confidentiality;
- f) workplace description including bioaerosol generating activities;
- g) type and name of sampler used;
- h) type and name of filter used;
- i) type of sampling (personal or static) and in case of personal sampling: name and function of the person who was involved in the sampling, or a unique identifier to preserve confidentiality; or in case of static sampling, placement of sampling equipment;
- j) positioning of sampling equipment;
- k) location of sampling inlet and orientation relative to air movement, where applicable;
- l) start and end time of sampling and its duration;
- m) flow rate, in litres per minute, at least before and after sampling period;
- n) sampled volume;
- o) sampled health-related fraction (see EN 481);
- p) details about sample transport and sample storage;
- q) environmental conditions during sampling (temperature, relative humidity, outdoor weather conditions);
- r) other relevant observations.

5 Extraction

5.1 General

Samples should not be extracted unless they can be analysed immediately. If this is not possible, extracts may be stored refrigerated overnight. For long-term storage, extracts may be frozen as specified in 5.4.

5.2 and 5.3 specify the equipment to be used in a recommended method for the extraction of endotoxins from filters as well as the extraction procedure. If other equipment and methods are used they shall be demonstrably equivalent in performance to this recommended method.

NOTE For deposition of endotoxins in inner sampler walls see Annex C.

5.2 Equipment

5.2.1 Laboratory shaker, able to shake the sample vigorously in extraction liquid.

5.2.2 Vortex shaker

5.2.3 Centrifuge, able to centrifuge at $1\ 000 \times G$

5.2.4 Tubes or containers, made of glass or polypropylene, endotoxinfree

5.3 Procedure

iTeh STANDARD PREVIEW

The loaded filter shall be transferred aseptically into a tube or container. Extraction shall be performed in endotoxinfree water using a method validated to maximize extraction.

NOTE 1 Due to lack of data no recommendation can be given for use of detergents in the extraction procedure.

<https://standards.iteh.ai/catalog/standards/sist/8521ee9-9e37-46d1-8c5f-4f8550441403/pr-en-14031-2020>

For filters with a diameter of less than 50 mm a minimum of 5 ml of an extraction liquid should be used, and a minimum of 10 ml for filters with a diameter of 50 mm or more. In order to obtain maximum extraction of endotoxin into suspension, samples are rocked/shaken for about 1 h at a speed of at least $2\ 000\ \text{min}^{-1}$ using the laboratory shaker at room temperature.

Before extraction the pH value should be measured in a subsample of the extraction liquid to ensure compatibility with the performance of the LAL assay and to comply with the recommendations given by the manufacturer of the LAL assay kit used.

NOTE 2 The pH value of bioaerol extracts is normally in the range from 6 to 8 which usually complies with the recommendations given by manufacturers of LAL assay kits.

Samples shall be centrifuged for about 15 min at about $1\ 000 \times G$. Supernatant shall be collected, well mixed and divided into aliquots of extracts.

5.4 Storage conditions

If supernatants of extracts need to be stored they shall either be stored overnight at about $4\ ^\circ\text{C}$ to $8\ ^\circ\text{C}$ or be stored frozen at about $-20\ ^\circ\text{C}$ or below for long-term storage. Extracts can be stored frozen in tightly sealed tubes for longer time periods (several years). An extract shall not be repeatedly frozen and thawed as this can affect the measurable endotoxin concentration in the fluid. Thawed samples should be mixed using a Vortex shaker before analysed.

Results from the analysis of previously frozen samples shall not be compared with results from samples analysed immediately after sampling. It shall be documented that analysis has been done on previously frozen samples.