
Blato, obdelani biološki odpadki in tla - Določevanje polikloriranih bifenilov (PCB) s plinsko kromatografijo z masno selektivnim detektorjem (GC/MS) in s plinsko kromatografijo z detektorjem z zajetjem elektronov (GC/ECD)

Soil, treated biowaste and sludge - Determination of polychlorinated biphenyls (PCB) by gas chromatography with mass selective detection (GC-MS) and gas chromatography with electron-capture detection (GC-ECD)

Boden, behandelter Bioabfall und Schlamm - Bestimmung von polychlorierten Biphenylen (PCB) mittels Gaschromatographie mit Massenspektrometrie-Kopplung (GC-MS) und Gaschromatographie mit Elektroneneinfangdetektion (GC-ECD)

Sols, biodéchets traités et boues - Dosage des polychlorobiphényles (PCBs) par chromatographie en phase gazeuse-spectrométrie gazeuse couplée avec un détecteur de masse (CG-SM) ou un détecteur par capture d'électrons (CG-ECD)

Ta slovenski standard je istoveten z: prEN 16167

ICS:

13.030.20	Tekoči odpadki. Blato	Liquid wastes. Sludge
13.080.10	Kemijske značilnosti tal	Chemical characteristics of soils

oSIST prEN 16167:2017

en,fr,de

EUROPEAN STANDARD
NORME EUROPÉENNE
EUROPÄISCHE NORM

DRAFT
prEN 16167

October 2017

ICS 13.030.01; 13.080.10

Will supersede EN 16167:2012

English Version

**Soil, treated biowaste and sludge - Determination of
polychlorinated biphenyls (PCB) by gas chromatography
with mass selective detection (GC-MS) and gas
chromatography with electron-capture detection (GC-
ECD)**

Schlamm, behandelter Bioabfall und Boden -
Bestimmung von polychlorierten Biphenylen (PCB)
mittels Gaschromatographie mit
massenspektrometrischer Detektion (GC-MS) und
Gaschromatographie mit Elektroneneinfangdetektion
(GC-ECD)

This draft European Standard is submitted to CEN members for enquiry. It has been drawn up by the Technical Committee CEN/TC 444.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

Warning : This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents

Page

European foreword.....	4
Introduction	5
1 Scope	6
2 Normative references	6
3 Terms and definitions	7
4 Principle	8
5 Interferences	8
5.1 Interference with sampling and extraction	8
5.2 Interference with GC	8
6 Safety remarks	9
7 Reagents	9
8 Apparatus.....	15
9 Sample storage and preservation	16
9.1 Sample storage	16
9.2 Sample pretreatment.....	16
10 Procedure.....	17
10.1 Blank test	17
10.2 Extraction.....	17
10.2.1 General.....	17
10.2.2 Extraction procedure 1: Samples using acetone/petroleum ether and agitation	18
10.2.3 Extraction procedure 2: Samples using Soxhlet.....	19
10.2.4 Extraction procedure 3: Samples using acetone/petroleum ether/sodium chloride and agitation	19
10.3 Concentration.....	19
10.4 Clean-up of the extract	20
10.4.1 General.....	20
10.4.2 Clean-up A – Aluminium oxide	21
10.4.3 Clean-up B – Silica gel	21
10.4.4 Clean-up C – Gel permeation chromatography.....	21
10.4.5 Clean-up D – Florisil®	22
10.4.6 Clean-up E – silica H ₂ SO ₄ /silica NaOH	22
10.4.7 Clean-up F – Benzenesulfonic acid/sulfuric acid	22
10.4.8 Clean-up G – DMF/ <i>n</i> -hexane partitioning.....	22
10.4.9 Clean-up H – concentrated sulfuric acid	22
10.4.10 Clean-up I – TBA sulfite reagent	23
10.4.11 Clean-up J – Clean-up using pyrogenic copper	23
10.4.12 Clean up K – AgNO ₃ /silica.....	23
10.5 Addition of the injection standard.....	23
10.6 Gas chromatographic analysis (GC).....	23
10.6.1 General.....	23
10.6.2 Setting the gas chromatograph	23
10.7 Mass spectrometry (MS)	24

10.7.1	Mass spectrometric conditions.....	24
10.7.2	Calibration of the method using an internal standard.....	25
10.7.3	Measurement	26
10.7.4	Identification.....	26
10.7.5	Check on method performance	26
10.7.6	Calculation.....	27
10.8	Electron capture detection (ECD)	28
10.8.1	General	28
10.8.2	ECD conditions.....	28
10.8.3	Calibration of the method using internal standards.....	28
10.8.4	Measurement	28
10.8.5	Identification.....	28
10.8.6	Check on ECD method performance	29
10.8.7	Calculation.....	30
11	Performance characteristics	30
12	Precision	30
13	Test report	30
Annex A	(informative) Repeatability and reproducibility data.....	31
A.1	Materials used in the interlaboratory comparison study.....	31
A.2	Interlaboratory comparison results.....	31
Annex B	(informative) Examples for retention times of PCBs.....	33
Annex C	(informative) Calculation method for the estimation of total PCB content.....	34
Bibliography	39


 (https://standards.iteh.ai)
 Document Preview

SIST EN 16167:2018

<https://standards.iteh.ai/catalog/standards/sist/85068b6c-e437-4510-baa4-0ebf6fe6ba53/sist-en-16167-2018>

prEN 16167:2017 (E)**European foreword**

This document (prEN 16167:2017) has been prepared by Technical Committee CEN/TC 444 “Test methods for environmental characterization of solid matrices”, the secretariat of which is held by NEN.

This document is currently submitted to the CEN Enquiry.

This document will supersede EN 16167:2012.

This document has been prepared under a mandate given to CEN by the European Commission and the European Free Trade Association.

The preparation of this document by CEN is based on a mandate by the European Commission (Mandate M/330), which assigned the development of standards on sampling and analytical methods for hygienic and biological parameters as well as inorganic and organic determinants, aiming to make these standards applicable to sludge, treated biowaste and soil as far as this is technically feasible.

iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

[SIST EN 16167:2018](https://standards.iteh.ai/catalog/standards/sist/85068b6c-e437-4510-baa4-0ebf6fe6ba53/sist-en-16167-2018)

<https://standards.iteh.ai/catalog/standards/sist/85068b6c-e437-4510-baa4-0ebf6fe6ba53/sist-en-16167-2018>

Introduction

Polychlorinated biphenyls (PCB) have been widely used as additives in industrial applications where chemical stability has been required. This stability on the other hand creates environmental problems when PCBs are eventually released into the environment. Since some of these PCB compounds are highly toxic, their presence in the environment (air, water, soil, sediment and waste) is regularly monitored and controlled. At present determination of PCB is carried out in these matrices in most of the routine laboratories following the preceding steps for sampling, pretreatment, extraction, clean-up by measurement of specific PCB by means of gas chromatography in combination with mass spectrometric detection (GC-MS) or gas chromatography with electron capture detector (GC-ECD).

This European Standard was developed in the European project 'HORIZONTAL'. It is the result of a desk study "3-12 PCB" and aims at evaluation of the latest developments in assessing PCBs in sludge, soil, treated biowaste and neighbouring fields. Taken into account the different matrices and possible interfering compounds, this European Standard does not contain one single possible way of working. Several choices are possible, in particular relating to clean-up. Detection with both MS-detection and ECD-detection is possible. Two different extraction procedures are described and 11 clean-up procedures. The use of internal and injection standards is described in order to have an internal check on choice of the extraction and clean-up procedure. The method is as far as possible in agreement with the method described for PAHs (see CEN/TS 16181). It has been tested for ruggedness.

This European Standard is applicable and validated for several types of matrices as indicated in Table 1 (see also Annex A for the results of the validation).

Table 1 — Matrices for which this European Standard is applicable and validated

Matrix	Materials used for validation
Sludge	Municipal sewage sludge
Biowaste	Compost
Soil	Sandy soil

WARNING — Persons using this European Standard should be familiar with usual laboratory practice. This European Standard does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to ensure compliance with any national regulatory conditions.

IMPORTANT — It is absolutely essential that tests conducted according to this European Standard be carried out by suitably trained staff.

1 Scope

This draft European Standard specifies a method for quantitative determination of seven selected polychlorinated biphenyls (PCB28, PCB52, PCB101, PCB118, PCB138, PCB153 and PCB180) in sludge, treated biowaste and soil using GC-MS and GC-ECD (see Table 2).

Table 2 — Target analytes of this European Standard

Target analyte		CAS-RN ^a
PCB28	2,4,4'-trichlorobiphenyl	7012-37-5
PCB52	2,2',5,5'-tetrachlorobiphenyl	35693-99-3
PCB101	2,2',4,5,5'-pentachlorobiphenyl	37680-37-2
PCB118	2,3',4,4',5-pentachlorobiphenyl	31508-00-6
PCB138	2,2',3,4,4',5'-hexachlorobiphenyl	35056-28-2
PCB153	2,2',4,4',5,5'-hexachlorobiphenyl	35065-27-1
PCB180	2,2',3,4,4',5,5'-heptachlorobiphenyl	35065-29-3
^a CAS-RN Chemical Abstracts Service Registry Number.		

The limit of detection depends on the determinants, the equipment used, the quality of chemicals used for the extraction of the sample and the clean-up of the extract.

Under the conditions specified in this European Standard, limit of application of 1 µg/kg (expressed as dry matter) can be achieved.

Sludge and treated biowaste may differ in properties and also in the expected contamination levels of PCBs and presence of interfering substances. These differences make it impossible to describe one general procedure. This European Standard contains decision tables based on the properties of the sample and the extraction and clean-up procedure to be used.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 15934, *Sludge, treated biowaste, soil and waste — Calculation of dry matter fraction after determination of dry residue or water content*

EN 16179, *Sludge, treated biowaste and soil — Guidance for sample pretreatment*

EN ISO 5667-15, *Water quality — Sampling — Part 15: Guidance on the preservation and handling of sludge and sediment samples (ISO 5667-15)*

EN ISO 16720, *Soil quality — Pretreatment of samples by freeze-drying for subsequent analysis (ISO 16720)*

EN ISO 22892, *Soil quality — Guidelines for the identification of target compounds by gas chromatography and mass spectrometry (ISO 22892)*

ISO 8466-1, *Water quality — Calibration and evaluation of analytical methods and estimation of performance characteristics — Part 1: Statistical evaluation of the linear calibration function*

ISO 18512, *Soil quality — Guidance on long and short term storage of soil samples*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

- IEC Electropedia: available at <http://www.electropedia.org/>
- ISO Online browsing platform: available at <http://www.iso.org/obp>

3.1

polychlorinated biphenyl

PCB

biphenyl substituted by one to ten chlorine atoms

[SOURCE: EN 15308:2016, 3.1]

3.2

congener

member of the same kind, class or group of chemicals, e.g. anyone of the two hundred and nine individual PCB

Note 1 to entry: The IUPAC congener numbers are for easy identification; they do not represent the order of chromatographic elution.

[SOURCE: EN 15308:2016, 3.2]

3.3

critical pair

pair of congeners that will be separated to a predefined degree (e.g. $R = 0,5$) to ensure chromatographic separation meets minimum quality criteria

[SOURCE: EN 15308:2016, 3.6]

3.4

internal standard

$^{13}\text{C}_{12}$ -labelled PCB or other PCB that are unlikely to be present in samples added to the sample before extraction and used for quantification of PCB content

[SOURCE: EN 15308:2016, 3.4, modified – "waste samples" is replaced here with "samples".]

3.5

injection standard

$^{13}\text{C}_{12}$ -labelled PCB or other PCB that is unlikely to be present in samples added to the sample extract before injection into the gas chromatograph, to monitor variability of instrument response and the recovery of the internal standards

[SOURCE: EN 15308:2016, 3.5, modified – "waste samples" is replaced here with "samples".]

4 Principle

Due to the horizontal character of this European Standard, different procedures for different steps (modules) are allowed. Which modules should be used depends on the sample. A recommendation is given in this European Standard. Performance criteria are described and it is the responsibility of the laboratories applying this European Standard to show that these criteria are met. Using of spiking standards (internal standards) allows an overall check on the efficiency of a specific combination of modules for a specific sample. But it does not necessarily give the information upon the extensive extraction efficiency of the native PCB bonded to the matrix.

After pretreatment according to the methods referred to in 9.2, the test sample is extracted with a suitable solvent.

The extract is concentrated by evaporation: If necessary, interfering compounds are removed by a clean-up method suitable for the specific matrix. The eluate is concentrated by evaporation.

The extract is analysed by gas chromatography with either mass spectrometric (GC-MS) or electron capture detection (GC-ECD). Applying GC-MS the various compounds are separated using a capillary column with a stationary phase of low polarity. In case of GC-ECD, extracts are analysed using two columns of different polarity (see 8.2.1).

PCBs are identified and quantified by comparison of relative retention times and relative peak heights (or peak areas) with respect to internal standards added. The efficiency of the procedure depends on the composition of the matrix that is investigated.

5 Interferences

5.1 Interference with sampling and extraction

Use sampling containers of materials (preferably of steel, aluminium or glass) that do not change the sample during the contact time. Avoid plastics and other organic materials during sampling, sample storage or extraction. Keep the samples from direct sunlight and prolonged exposure to light.

During storage of the samples, losses of PCBs may occur due to adsorption on the walls of the containers. The extent of the losses depends on the storage time.

5.2 Interference with GC

Substances that co-elute with the target PCB may interfere with the determination. These interferences may lead to incompletely resolved signals and may, depending on their magnitude, affect accuracy and precision of the analytical results. Peak overlap does not allow an interpretation of the result. Asymmetric peaks and peaks being broader than the corresponding peaks of the reference substance suggest interferences.

Chromatographic separation between the following pairs can be critical. The critical pair PCB28 and PCB31 is used for selection of the capillary column (see 8.2.2). If molecular mass differences are present, quantification can be made by mass selective detection. If not or using ECD, the specific PCB is reported as the sum of all PCBs present in the peak. Typically, the concentrations of the co-eluting congeners compared to those of the target congeners are low. When incomplete resolution is encountered, peak integration shall be checked and, when necessary, corrected.

- PCB28 – PCB31
- PCB52 – PCB73
- PCB101 – PCB89 / PCB90

- PCB118 – PCB106
- PCB138 – PCB164 / PCB163

Presence of considerable amounts of mineral oil in the sample may interfere with the quantification of PCB in GC-MS. In presence of mineral oil, GC-ECD may be preferred or mineral oil can be removed using clean-up procedure G (see 10.4.8) using DMF/*n*-hexane.

Presence of tetrachlorobenzyltoluene (TCBT)-mixtures may disturb the determination of the PCB with GC-ECD.

6 Safety remarks

PCBs are highly toxic and shall be handled with extreme care. Avoid contact with solid materials, solvent extracts and solutions of standard PCB with the body shall not be allowed to occur. It is strongly advised that standard solutions are prepared centrally in suitably equipped laboratories or are purchased from suppliers specialized in their preparation.

Solvent solutions containing PCB shall be disposed of in a manner approved for disposal of toxic wastes.

For the handling of hexane precautions shall be taken because of its neurotoxic properties.

National regulations shall be followed with respect to all hazards associated with this method.

7 Reagents

All reagents shall be of recognized analytical grade. The purity of the reagents used shall be checked by running a blank test as described in 10.1. The blank shall be less than 50 % of the lowest reporting limit.

7.1 Reagents for extraction:

7.1.1 Acetone, $(\text{CH}_3)_2\text{CO}$.

7.1.2 *n*-heptane, C_7H_{16} .

7.1.3 Petroleum ether, boiling range 40 °C to 60 °C.

Hexane-like solvents with a boiling range between 40 °C and 89 °C are allowed.

7.1.4 Anhydrous sodium sulfate, Na_2SO_4 :

The anhydrous sodium sulfate shall be kept carefully sealed.

7.1.5 Distilled water or water of equivalent quality, H_2O .

7.1.6 Sodium chloride, NaCl , anhydrous.

7.1.7 Keeper substance. High boiling compound, i.e. octane, nonane.

7.2 Reagents for clean-up:

7.2.1 Clean-up A using aluminium oxide:

7.2.1.1 Aluminium oxide, Al_2O_3 .

Basic or neutral, specific surface 200 m²/g, activity Super I according to Brockmann.

prEN 16167:2017 (E)**7.2.1.2 Deactivated aluminium oxide**

Deactivated with approximately 10 % water.

Add approximately 10 g of water (7.1.5) to 90 g of aluminium oxide (7.2.1.1). Shake until all lumps have disappeared. Allow the aluminium oxide to condition before use for some 16 h, sealed from the air, use it for maximum two weeks.

NOTE The activity depends on the water content. It can be necessary to adjust the water content.

7.2.2 Clean-up B using silica gel 60 for column chromatography:

7.2.2.1 Silica gel 60, particle size 63 μm to 200 μm .

7.2.2.2 Silica gel 60, water content: mass fraction $w(\text{H}_2\text{O}) = 10\%$.

Silica gel 60 (7.2.2.1), heated for at least 3 h at 450 °C, cooled down in a desiccator and stored containing magnesium perchlorate or a suitable drying agent. Before use heat at least for 5 h at 130 °C in a drying oven. Then allow cooling in a desiccator and add 10 % water (mass fraction) in a flask. Shake for 5 min intensively by hand until all lumps have disappeared and then for 2 h in a shaking device. Store the deactivated silica gel in the absence of air, use it for maximum of two weeks.

7.2.3 Clean-up C using gel permeation chromatography (GPC):

7.2.3.1 Bio-Beads® S-X3¹⁾.

7.2.3.2 Ethyl acetate, $\text{C}_4\text{H}_8\text{O}_2$.

7.2.3.3 Cyclohexane, C_6H_{12} .

Preparation of GPC, for example: Put 50 g Bio-Beads® S-X3 (7.2.3.1) into a 500 ml Erlenmeyer flask and add 300 ml elution mixture made up of cyclohexane (7.2.3.3) and ethyl acetate (7.2.3.2) 1:1 (volume) in order to allow the beads to swell; after swirling for a short time until no lumps are left, maintain the flask closed for 24 h. Drain the slurry into the chromatography tube for GPC. After approximately three days, push in the plungers of the column so that a filling level of approximately 35 cm is obtained. To further compress the gel, pump approximately 2 l of elution mixture through the column at a flow rate of 5 ml \cdot min⁻¹ and push in the plungers to obtain a filling level of approximately 33 cm.

7.2.4 Clean-up D using Florisil®²⁾

7.2.4.1 Florisil®, baked 2 h at 600 °C. Particle size 150 μm to 750 μm .

7.2.4.2 Iso-octane, C_8H_{18} .

1) Bio-Beads® is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by CEN of this product. Equivalent products may be used if they can be shown to lead to the same results.

2) Florisil® is a trade name for a prepared diatomaceous substance, mainly consisting of anhydrous magnesium silicate. This information is given for the convenience of users of this European Standard and does not constitute an endorsement by CEN of this product. Equivalent products may be used if they can be shown to lead to the same results.