

SLOVENSKI STANDARD oSIST prEN 17411:2019

01-september-2019

Gnojila - Določevanje perklorata v mineralnih gnojilih s tekočinsko kromatografijo in z detekcijo na osnovi tandemske masne spektrometrije (LC-MS/MS)

Fertilizers - Determination of perchlorate in mineral fertilizers by liquid chromatography and tandem mass spectrometry detection (LC-MS/MS)

Düngemittel - Bestimmung von Perchlorat in mineralischen Düngemitteln - Bestimmung mittels Flüssig-Chromatographie und Tandem-Massenspektrometrischer Detektion (LC-MS/MS)

Engrais — Dosage du perchlorate dans les engrais minéraux par chromatographie liquide couplée à la spectrométrie de masse (LCMS/MS)

Ta slovenski standard je istoveten z: prEN 17411

ICS:

65.080 Gnojila Fertilizers

71.040.50 Fizikalnokemijske analitske Physicochemical methods of

metode analysis

oSIST prEN 17411:2019 en,fr,de

oSIST prEN 17411:2019

iTeh STANDARD PREVIEW (standards.iteh.ai)

SIST EN 17411:2021

https://standards.iteh.ai/catalog/standards/sist/25d94cec-57fa-4157-99c1-f39e10e8d2bd/sist-en-17411-2021

EUROPEAN STANDARD NORME EUROPÉENNE EUROPÄISCHE NORM

DRAFT prEN 17411

July 2019

ICS 65.080

English Version

Fertilizers - Determination of perchlorate in mineral fertilizers by liquid chromatography and tandem mass spectrometry detection (LC-MS/MS)

Engrais - Dosage du perchlorate dans les engrais minéraux par chromatographie liquide couplée à la spectrométrie de masse (LCMS/MS) Düngemittel - Bestimmung von Perchlorat in mineralischen Düngemitteln - Bestimmung mittels Flüssig-Chromatographie und Tandem-Massenspektrometrischer Detektion (LC-MS/MS)

This draft European Standard is submitted to CEN members for enquiry. It has been drawn up by the Technical Committee CEN/TC 260.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

Warning: This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION COMITÉ EUROPÉEN DE NORMALISATION EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

Con	Contents	
Euro	opean foreword	eword
1	Scope	4
2	Normative references	4
3	Terms and definitions	4
4	Principle	4
5	Reagents	4
6	Apparatus and equipment	5
7 7.1 7.2 7.3	Sampling and sample preparation Matrix calibration curve	5 5
8	Calculation	7
9 9.1 9.2 9.3	Inter-laboratory testRepeatabilityReproducibilityReproducibilityReproducibilityReproducibility	7 7 7
10	Additional remarks Standards itah ai)	8
11	Test report	8
Anne	ex A (informative) Results of the inter-laboratory test	9
Bibliography 230340000 3703 1008 2303 2303 2003 2003 2003 2003 2003 2		

European foreword

This document (prEN 17411:2019) has been prepared by Technical Committee CEN/TC 260 "Fertilizers and liming materials", the secretariat of which is held by DIN.

This document is currently submitted to the CEN Enquiry.

iTeh STANDARD PREVIEW (standards.iteh.ai)

<u>SIST EN 17411:2021</u> https://standards.iteh.ai/catalog/standards/sist/25d94cec-57fa-4157-99c1 f39e10e8d2bd/sist-en-17411-2021

1 Scope

This document specifies a method for the determination of traces of perchlorate by liquid chromatography and tandem mass spectrometry detection (LC-MS/MS). The method is applicable to mineral fertilizers.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 12944-1, Fertilizers and liming materials and soil improvers – Vocabulary – Part 1: General terms

EN 12944-2, Fertilizers and liming materials and soil improvers – Vocabulary – Part 2: Terms relating to fertilizers

EN ISO 3696, Water for analytical laboratory use – Specification and test methods (ISO 3696)

EN ISO 14820-2, Fertilizers and liming materials – Sampling and sample preparation – Part 2: Sample preparation

3 Terms and definitions

For the purposes of this document, the terms and definitions given in EN 12944-1 and EN 12944-2 apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

- IEC Electropedia: available at http://www.electropedia.org/
- ISO Online browsing platform: available at http://www.iso.org/obp

4 Principle

The perchlorate content of the fertilizer sample is determined by liquid chromatography with tandem mass spectrometry detection after water extraction. The sample is dissolved, filtered and directly analysed by LC-MS/MS. Quantification is performed using matrix calibration curve.

5 Reagents

Use only reagents of recognized analytical grade.

- **5.1 Water**, grade 3 according to EN ISO 3696.
- **5.2 Blank solution**, free of perchlorate similar to the samples.

Weigh 5 g of potassium nitrate in 250 ml eluent (5.4).

5.3 Perchlorate standard solutions, $\rho(\text{ClO}_{4}^{-}) = 1 \text{ mg/l}$ and 10 mg/l.

Prepare the standard perchlorate stock solutions from a certified perchlorate solution.

5.4 Eluent, used for LC-MS/MS (mobile phase), water+acetonitrile 1+1 (v+v), acetic acid 0,1 %.

6 Apparatus and equipment

Usual laboratory equipment, glassware, and the following:

6.1 LC-MS/MS instrumentation:

6.1.1 Column.

For perchlorate determination, it is recommended to use a KP-RPX 250 \times 4.6, 5 μ m column, from K' (PRIME) Technologies¹. Similarly performing columns may also be used.

6.1.2 Equipment.

Use a mass spectrometer equipped with electrospray ionization (ESI) in the negative mode. The mass transitions for the perchlorate ion are: (-) 99/83, (-) 101/85.

- **6.2 Analytical balance**, with an accuracy of \pm 0,1 mg.
- **6.3 Graduated flasks**, capacity 1 000 ml, 250 ml, 50 ml and 10 ml volumetric flasks.
- **6.4 Nylon filters**, with 0,22 μm filter membrane.
- **6.5 Autosampler vials**, suitable for LC auto-samplers.

7 Procedure Teh STANDARD PREVIEW

7.1 Sampling and sample preparation

Sampling is not part of this European Standard. A recommended sampling method is given in EN ISO 14820-1 [1].

Sample preparation shall be performed according to EN ISO 14820-2.7fa-4157-9901-

7.2 Matrix calibration curve

The calibration curve is prepared in matrix to correct for matrix interferences.

For the calibration solution, take 5 ml of acetonitrile and transfer quantitatively to a 10 ml volumetric flask (6.3). Add an adequate volume of perchlorate standard solution (5.3) as given in Table 1 and bring to volume with the blank solution.

 $^{^{1}~}$ This column – KP-RPX 250 × 4.6, 5 μm column, from K' (PRIME) Technologies is an example of a suitable instrument available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by CEN of this product. Equivalent products may be used if they can be shown to lead to the same results.

Calibration Standard solution Standard solution Mass concentration of solution calibration solution $\rho(\text{ClO}_4^-) = 1 \text{ mg/l}$ $\rho(\text{ClO}_{4}^{-}) = 10 \text{ mg/l}$ mg/l (5.3)(5.3)μl μl 1 50 0,005 2 100 0,01 3 20 0,2 4 50 0,05 5 100 0,1 200 0,2 6

Table 1 — Preparation of calibration solutions

Prepare a calibration curve by injecting the six calibration solutions (see Table 1) according to the operation procedure of the available instrument (6.1) at the beginning of every day of the analysis.

The calibration curve is described by the following formula:

$$y = mx + b ag{1}$$

Calibration levels lower than 0,005 mg/l may result in not linear response when matrix calibration is used.

7.3 Determination of perchlorate in test solutions

- **7.3.1** With the analytical balance (6.2) weigh 5 g of ground sample.
- **7.3.2** Transfer the sample to a 250 ml volumetric flask (6.3) and add 150 ml 0,1 % acetic acid (5.4) to dissolve the sample. Bring to volume with 0,1 % acetic acid (5.4). Stir manually to homogenize (solution 1).
- **7.3.3** Take an aliquot of 25 ml of solution 1 and transfer quantitatively to a 50 ml volumetric flask (6.3). Add acetonitrile to the mark and stir manually to homogenize (solution 2).
- NOTE Solution 1 (sample diluted with 1% acetic acid) is diluted with acetonitrile to obtain solution 2 (solution 2 = 25 ml solution 1 + 25 ml acetonitrile).
- **7.3.4** Depending on the concentration of perchlorates in the sample, solution 2 may be diluted again using blank solution (5.2).
- **7.3.5** Filter the solution using a nylon filter (6.4). Discard the first filtrate and then collect the filtrate in a vial (6.5).
- **7.3.6** Proceed with the analysis according to the operation procedure of the available instrument (6.1) using the same instrument and the same conditions as used for the preparation of the calibration curve.

8 Calculation

Calculate the the mass fraction of ClO_4^- in the fertilizer, in mg/kg, by using the following Formula (2):

$$w(\text{ClO}_4^-) = \frac{\rho \times D \times V}{m} \tag{2}$$

where

 $w(ClO_4^-)$ is the mass fraction of ClO_4^- in the fertilizer, in mg/kg;

 ρ is the mass concentration of the reading of the sample from the calibration curve, in mg/l;

D is the dilution used (here: 2);

V is the volume of the sample, in ml (here: 250 ml);

m is the mass of the test portion, in g (here: 5 g).

9 Precision

9.1 Inter-laboratory test

An inter-laboratory test was carried out in 2018 with seven participating laboratories analysing four different samples. Repeatability and reproducibility were calculated according to ISO 5725-2 [2]. Each laboratory carried out three determinations on each sample. The statistical results are explained in Table A.1.

The values derived from the inter-laboratory test might not be applicable to concentration ranges and matrices other than those given in Annex A.

9.2 Repeatability

The absolute difference between two independent test results obtained with the same method on identical test material, in the same laboratory, by the same operator, using the same equipment within a short interval of time, will in no more than 5 % of the cases exceed the values of *r* given in Table 2.

9.3 Reproducibility

The absolute difference between two single test results obtained with the same method on identical test material, in different laboratories, by different operators, using different equipment, will in no more than 5 % of the cases exceed the values of *R* given in Table 2.

Sample	Mean value mg/kg	r mg/kg	R mg/kg
1	21,0	3,9	6,6
2	5,8	0,7	4,1
3	34,7	4,4	18,3
4	13,7	1,1	6,6

Table 2 — Mean values, repeatability and reproducibility limits

10 Additional remarks

- When samples with high perchlorate concentrations have been measured, there is an increased risk for carry-over. Wash the column properly after detecting high concentrations. It is recommended to apply a sample dilution.
- Prolonged use of the column may result in an increase of retention time. Retention time shall be checked at regular intervals using a sample with known concentration.

11 Test report

The test report shall contain at least the following information:

- a) all information necessary for the complete identification of the sample;
- b) test method used with reference to this draft European Standard, prEN 17411;
- c) test results obtained expressed as mass fraction of perchlorate in the fertilizer;
- d) date of sampling and sampling procedure (if known);
- e) date when the analysis was finished;
- f) whether the requirement of the repeatability limit has been fulfilled;
- g) all operating details not specified in this document, or regarded as optional, together with details of any incidents occurred when performing the method, which might have influenced the test result(s).

https://standards.iteh.ai/catalog/standards/sist/25d94cec-57fa-4157-99c1f39e10e8d2bd/sist-en-17411-2021