

## SLOVENSKI STANDARD oSIST prEN 17134-2:2022

01-september-2022

Tekstilije in tekstilni izdelki - Določevanje biocidnega dodatka - 2. del: Konzervansi na osnovi klorofenola, metoda z uporabo plinske kromatografije

Textiles and textile products - Determination of biocide additives – Part 2, Chlorophenol-based preservatives, method using gas chromatography

Textilien und textile Erzeugnisse - Bestimmung von Biozid-Zusatzstoffen - Teil 2: Konservierungsmittel auf Chlorphenolbasis, Verfahren mittels Gaschromatographie

Textiles et produits textiles - Détermination des additifs biocides - Partie 2 : Conservateurs à base de chlorophénol, méthode par chromatographie en phase gazeuse

Ta slovenski standard je istoveten z: prEN 17134-2

ICS:

59.080.01 Tekstilije na splošno Textiles in general

71.040.50 Fizikalnokemijske analitske Physicochemical methods of

metode analysis

oSIST prEN 17134-2:2022 en,fr,de

oSIST prEN 17134-2:2022

## iTeh Standards (https://standards.iteh.ai) Document Preview

SIST EN 17134-2:2023

## EUROPEAN STANDARD NORME EUROPÉENNE EUROPÄISCHE NORM

# **DRAFT prEN 17134-2**

July 2022

ICS 59.080.01

Will supersede EN 17134:2019

#### **English Version**

# Textiles and textile products - Determination of biocide additives - Part 2, Chlorophenol-based preservatives, method using gas chromatography

Textilien und textile Erzeugnisse - Bestimmung von Biozid-Zusatzstoffen - Teil 2: Konservierungsmittel auf Chlorphenolbasis, Verfahren mittels Gaschromatographie

This draft European Standard is submitted to CEN members for enquiry. It has been drawn up by the Technical Committee CEN/TC 248.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

**Warning**: This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION COMITÉ EUROPÉEN DE NORMALISATION EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

### Contents

		Page
Europ	ean foreword	4
1	Scope	6
2	Normative references	6
3	Terms and definitions	6
1	Abbreviations	6
5	Principle	7
6	Apparatus	
5.1	Gas chromatograph (GC), with mass selective detector (MSD)	
5.2	Analytical balance, weighing with an accuracy of 0,1 mg	
5.3	Heating block, sand bath or drying cabinet, suitable for a temperature of 90 °C ± 1 in the KOH solution	
5.4	Gas-tight glass vials, for example headspace vials, (20 ml)	
<b>5.5</b>	GC vials, about 2 ml	
5.6	Pasteur pipettes, graduated pipettes, suitable autopipettes	
5.7	Vortex shaker	
5.8	Horizontal shaker, capable of at least 200 min <sup>-1</sup>	
7	Reagents	
7.1	General Andreas	
7.2	Chlorophenol mixture Stock solution	
7.3	Tetrachloroguaiacol (tetrachloro-o-methoxyphenol) CAS number: 2539-17-5	
7.4	Isotope labelled internal standards	
7.5	n-Hexane, for residue analysis.	
7.6	Potassium hydroxide solution (KOH), aqueous solution (1 mol/l)	9
7.7	Tetrachloroguiacol (TCG) solution, at a concentration of 1 μg/ml in acetonitrile (internal standard for injection control)	0
<b>7.8</b> tps	Stock solution isotope-labelled internal standards in KOH solution, each 0,1 mg/m	
7.9	Extraction solution (containing isotope-labelled internal standards)	
7.10	Potassium carbonate (K <sub>2</sub> CO <sub>3</sub> ), aqueous solution (0,1 mol/l)	
7.11	Acetic anhydride (C <sub>4</sub> H <sub>6</sub> O <sub>3</sub> ), for analysis.	
7.12	Water, grade 3, according to ISO 3696.	
7.13	Acetone	10
7.14	Acetonitrile	
7.15	Working solution calibration level 1, for daily calibration	10
7.16	Working solution calibration level 2, for daily calibration	10
7.17	Working solution calibration level 3, for daily calibration	
7.18	Working solution calibration level 4, for daily calibration	
7.19	Working solution calibration level 5, for daily calibration.	10
3	Sampling	10
9	Procedure	
9.1	General	
9.2	Extraction with KOH	
9.3	Extractive acetylation	
9.4	Derivatization of the chlorophenol mixtures and the TCG standards (calibration fo	r

9.5	Gas chromatography (GC)	11
10	Expression of results	12
10.1	General	12
10.2	Calculation of the individual chlorphenols in the test specimen	12
10.3	Reliability of the procedure	
11	Test report	13
Annex	A (informative) Chromatographic test conditions	14
<b>A.1</b>	Preliminary remarks	14
<b>A.2</b>	Gas chromatography with mass selective detector (GC-MSD)	14
Annex	B (informative) Procedure characteristics	16
<b>B.1</b>	General	16
<b>B.2</b>	Ring test results	16
<b>B.2.1</b>	Results of the first ring test	16
B.2.2	Results of the second ring test from one real test specimen and two spikes of th test specimen, with 4 laboratories	
B.3	Influence of standard solutions	25
<b>B.4</b>	Influence of internal standards on contents of pentachlorophenol and other chlorophenols	31
Annex	c C (informative) Verification method for free mono- and dichlorophenol	41
<b>C.1</b>	General	41
<b>C.2</b>	Procedure (https://standards.iteh.ai)	41

#### **Document Preview**

SIST EN 17134-2:2023

#### **European foreword**

This document (prEN 17134-2:2021) has been prepared by Technical Committee CEN/TC 248 "Textiles and textile products", the secretariat of which is held by BSI.

This document is currently submitted to the CEN Enquiry.

This document will supersede EN 17134:2019.

## iTeh Standards (https://standards.iteh.ai) Document Preview

SIST EN 17134-2:2023

#### Introduction

In Europe, according to REGULATION (EU) 2019/1021[1] of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants pentachlorophenol (PCP) and its salts and esters as constituents of articles are prohibited. According to COMMISSION DELEGATED REGULATION (EU) 2021/277 of 16 December 2020 amending Annex I to Regulation (EU) 2019/1021 of the European Parliament and of the Council on persistent organic pollutants as regards pentachlorophenol and its salts and esters, articles containing PCP in concentrations equal or lower than 5 mg/kg are allowed.

Further chlorinated phenols are restricted by voluntary specifications (ecolabel criteria, industry initiatives and standards).

This document specifies a method in which chlorophenols (CP) are acetylated before chromatographic detection and the amount of detected chlorophenol acetate is quantified by a correction with an internal standard.

**WARNING -** The use of this document involves hazardous materials. It does not purport to address all of the safety or environmental problems associated with its use. It is the responsibility of users of this document to take appropriate measures to ensure the safety and health of personnel and the environment prior to application of the document, and fulfil statutory and regulatory requirements for this purpose.

iTeh Standards (https://standards.iteh.ai) Document Preview

SIST EN 17134-2:2023

#### 1 Scope

This document specifies a method to determine the content of free pentachlorophenol (PCP), and tetrachlorophenol- (TeCP), trichlorophenol- (TriCP), dichlorophenol- (DiCP), monochlorophenol- (MoCP) isomers. The procedure also releases chlorophenols from salts and esters. The method is applicable to natural and synthetic textiles, including coated fabrics and clothing components (e.g. buttons, zippers, etc.).

#### 2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

<std>ISO 3696, Water for analytical laboratory use — Specification and test methods</std>

#### 3 Terms and definitions

No terms and definitions are listed in this document.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

- ISO Online browsing platform: available at https://www.iso.org/obp
- IEC Electropedia: available at https://www.electropedia.org/

#### 4 Abbreviations

The following abbreviations are used in this document for chlorophenols:

СР	chlorophenols Document Preview			
MoCP	monochlorophenol			
DiCP://standardichlorophenobg/standards/sist/bf44eb9b-194d-4412-b459-daf0fecd11c2/sist-en-17134-2-2023				
TriCP	trichlorophenol			
TeCP	tetrachlorophenol			
PCP	pentachlorophenol			
2-MoCP-d4	2-chlorophenol-d₄			
2,3-DiCP-d3	$2,3$ -dichlorophenol- $d_3$			
TeCP-13C6	$2,3,4,6$ -tetrachlorophenol- $^{13}C_6$			
PCP-C13	$pentachlorophenol \hbox{$^{-13}$C}_6$			
TCG	tetrachloroguaiacol (= tetrachloro-o-methoxyphenol)			

The following abbreviations are used in this document for the designation of textile fibres:

CO cotton

PES polyester

CV viscose

EL elastane

LI linen

polyamide

#### 5 Principle

PA

The sample is extracted with a potassium hydroxide solution (KOH solution) for 16 h at 90 ° C.

During subsequent extractive derivatization with n-hexane and acetic acid anhydride, CPs are acetylated and the chlorinated acetates are analyzed using gas chromatography with a mass-selective detector (MSD). The quantitative determination is made by correction with  $^{13}$ C and  $^{2}$ H labelled internal standards. For pentachlorophenol, the  $^{13}$ C derivative of PCP is used, for tetrachlorophenols the  $^{13}$ C derivative of 2,3,4,6-tetrachlorophenol, for tri- and dichlorophenols 2,3-dichlorophenol- $^{13}$ C and for monochlorophenols 2-chlorophenol- $^{13}$ C serves as an internal standard for injection control.

#### 6 Apparatus

- **6.1 Gas chromatograph (GC),** with mass selective detector (MSD).
- **6.2** Analytical balance, weighing with an accuracy of 0,1 mg.
- **6.3 Heating block, sand bath or drying cabinet,** suitable for a temperature of 90 °C ± 1 °C in the KOH solution.
- $\textbf{6.4 Gas-tight glass vials,} for example headspace vials, (20 ml) \\ \textbf{1.} \\ \textbf{2.} \\$
- **6.5 GC vials,** about 2 ml.
- 6.6 Pasteur pipettes, graduated pipettes, suitable autopipettes
- 6.7 Vortex shaker
- 6.8 Horizontal shaker, capable of at least 200 min<sup>-1</sup>
- 7 Reagents
- 7.1 General

Unless otherwise stated, chemicals of analytical grade shall be used.

#### 7.2 Chlorophenol mixture Stock solution

A mixture of chlorophenols which, at a concentration of 50 µg/ml in a suitable solvent (for example

#### Acetonitrile) contains the following isomers:

2-chlorophenol CAS number: 95-57-8

3-chlorophenol CAS number: 108-43-0

4-chlorophenol CAS number: 106-48-9

2,3-dichlorophenol CAS number: 576-24-9

2,4-dichlorophenol CAS number: 120-83-2

2,5-dichlorophenol CAS number: 583-78-8

2,6-dichlorophenol CAS number: 87-65-0

3,4-dichlorophenol CAS number: 95-77-2

3,5-dichlorophenol CAS number: 591-35-5

2,3,4-trichlorophenol CAS number: 15950-66-0

2,3,5-trichlorophenol CAS number: 933-78-8

2,3,6-trichlorophenol CAS number: 933-75-5

2,4,5-trichlorophenol CAS number: 95-95-4

2,4,6-trichlorophenol CAS number: 88-06-2 dands itch ai)

3,4,5-trichlorophenol CAS number 609-19-8

2,3,4,5-tetrachlorophenol CAS number 4901-51-3

2,3,4,6-tetrachlorophenol CAS number: 58-90-2

2,3,5,6-tetrachlorphenol CAS number: 935-95-5

pentachlorophenol CAS number: 87-86-5

NOTE This mixture of chlorophenols is available from laboratory chemical suppliers.

#### 7.3 Tetrachloroguaiacol (tetrachloro-o-methoxyphenol) CAS number: 2539-17-5

#### 7.4 Isotope labelled internal standards

2-chlorophenol-d<sub>4</sub> CAS number: 93951-73-6

2,3-dichlorophenol-d<sub>3</sub> CAS number: 93951-74-7

2,3,4,6-tetrachlorophenol-<sup>13</sup>C<sub>6</sub> CAS number:1246820-81-4

pentachlorophenol-13C<sub>6</sub> CAS number: 85380-74-1

Instead of the isotope-labelled derivatives of tetra-, di- and monochlorophenol mentioned, other isotope-labelled analogues of the same chlorination levels, e.g. 4-chlorophenol-d<sub>4</sub>, may also be used as internal standards.

#### **7.5** *n***-Hexane,** for residue analysis.

#### **7.6 Potassium hydroxide solution (KOH),** aqueous solution (1 mol/l).

Weigh 56,1 g KOH into a 1 l volumetric flask and dissolve with 100 ml water (7.12) (caution: heat generation!). After cooling, fill up to 1 l with water (7.12).

## 7.7 Tetrachloroguiacol (TCG) solution, at a concentration of 1 $\mu$ g/ml in acetonitrile (internal standard for injection control)

Dissolve 10 mg TCG in 100 ml acetonitrile (7.14). Take 1,0 ml from this stock solution and dilute to 100 ml with acetonitrile (7.14).

#### 7.8 Stock solution isotope-labelled internal standards in KOH solution, each 0,1 mg/ml

Dissolve 10 mg of each isotope-labelled internal standard (7.4) in 100 ml KOH solution (7.6).

#### 7.9 Extraction solution (containing isotope-labelled internal standards)

1,0 ml stock solution isotope-labelled internal standards (7.8) is filled up with KOH solution (7.6) in a 1 l volumetric flask, the concentration of each of the isotope-labelled internal standards (7.4) is 0,1  $\mu$ g/ml .

The concentration of isotope-labelled internal standards shall be adapted to the respective calibration range of the chlorophenols (see solutions 7.15 to 7.19). If a more sensitive measuring instrument is used and the calibration is lowered, then the concentration of isotope-labelled internal standards in the KOH solution should also be reduced.

#### **7.10 Potassium carbonate** ( $K_2CO_3$ ), aqueous solution (0,1 mol/l).

Weigh 13,82 g K<sub>2</sub>CO<sub>3</sub> into a 1 l volumetric flask, dissolve with water (7.12) and fill up to 1 l.

#### 7.11 Acetic anhydride (C<sub>4</sub>H<sub>6</sub>O<sub>3</sub>), for analysis.

**7.12 Water,** grade 3, according to ISO 3696.

#### 7.13 Acetone

#### 7.14 Acetonitrile

#### 7.15 Working solution calibration level 1, for daily calibration

30  $\mu$ l stock solution chlorophenol mixture (7.2) is filled up with acetonitrile (7.14) in a 10 ml volumetric flask; the concentration of chlorophenols is 0,15  $\mu$ g/ml.

#### 7.16 Working solution calibration level 2, for daily calibration

400  $\mu$ l stock solution chlorophenol mixture (7.2) is filled up with acetonitrile (7.14) in a 10 ml volumetric flask; the concentration of chlorophenols is 2,0  $\mu$ g/ml.

#### **7.17 Working solution calibration level 3,** for daily calibration

 $800 \,\mu l$  stock solution chlorophenol mixture (7.2) are filled up with acetonitrile (7.14) in a 10 ml volumetric flask; the concentration of chlorophenols is  $4 \,\mu g/ml$ .

#### 7.18 Working solution calibration level 4, for daily calibration

600  $\mu$ l stock solution chlorophenol mixture (7.2) is filled up with acetonitrile (7.14) in a 2 ml volumetric flask; the concentration of chlorophenols is 15  $\mu$ g/ml.

#### 7.19 Working solution calibration level 5, for daily calibration.

1 000  $\mu$ l stock solution chlorophenol mixture (7.2) is filled up with acetonitrile (7.14) in a 2 ml volumetric flask; the concentration of chlorophenols is 25  $\mu$ g/ml.

The concentrations of the calibration solutions (7.15 to 7.19) are examples and should be adjusted according to the specified limit values and the linearity of the device.

#### 8 Sampling

From the material to be examined, a representative test specimen shall be taken and cut into pieces of about  $0.5 \text{ cm} \times 0.5 \text{ cm}$ .

#### 9 Procedure

#### 9.1 General

If it is only required to test for free mono- and dichlorophenols, the procedure described in Annex C shall be applied.

#### 9.2 Extraction with KOH

Place approximately 1 g (minimum sample weight 0,2 g) of the cut test specimen (to the nearest 0,01 g) in the reaction vessel (6.4). Add 10 ml extraction solution (7.9). After closing, the vessel is transferred to a heating apparatus (6.3) and left for 16 h  $\pm$  15 min at (90  $\pm$  1) °C. The temperature in the reaction vessel is checked with extraction solution. If the test specimen to be examined is too voluminous to be