

TECHNICAL SPECIFICATION
SPÉCIFICATION TECHNIQUE
TECHNISCHE SPEZIFIKATION

CEN/TS 17784-2

April 2022

ICS 65.080

English Version

Organo-mineral fertilizers - Identification of complexing agents - Part 2: Method using high-performance liquid chromatography (HPLC)

Engrais organo-minéraux - Identification des agents complexants - Partie 2 : Méthode par chromatographie liquide à haute performance (HPLC)

Organisch-mineralische Düngemittel - Identifizierung von Komplexbildnern - Teil 2: Verfahren mittels Hochleistungs-Flüssigkeitschromatographie (HPLC)

This Technical Specification (CEN/TS) was approved by CEN on 21 February 2022 for provisional application.

The period of validity of this CEN/TS is limited initially to three years. After two years the members of CEN will be requested to submit their comments, particularly on the question whether the CEN/TS can be converted into a European Standard.

CEN members are required to announce the existence of this CEN/TS in the same way as for an EN and to make the CEN/TS available promptly at national level in an appropriate form. It is permissible to keep conflicting national standards in force (in parallel to the CEN/TS) until the final decision about the possible conversion of the CEN/TS into an EN is reached.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

<https://standards.iteh.ai/catalog/standards/sist/0c961d07-0d8a-4f9d-be15-3917d601fad6/sist-ts-cen-ts-17784-2-2023>



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

Contents	Page
European foreword	3
Introduction	4
1 Scope	5
2 Normative references	5
3 Terms and definitions	5
4 Principle	5
5 Interferences	5
6 Reagents	6
7 Apparatus	7
8 Sampling and sample preparation	7
9 Procedure	7
10 Calculation of the heptagluconic acid content and expression of the results	9
11 Test report	9
Annex A (informative) Chromatograms of the standard and a typical sample solution	10
Annex B (informative) Absorption spectra of the HGA	12
Annex C (informative) Calculation of the molar ratio HGA:Metal	13
Annex D (informative) Complete names of chelating agents	14
Bibliography	15

European foreword

This document (CEN/TS 17784-2:2022) has been prepared by Technical Committee CEN/TC 260 “Fertilizers and liming materials”, the secretariat of which is held by DIN.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN shall not be held responsible for identifying any or all such patent rights.

This document has been prepared under a Standardization Request given to CEN by the European Commission and the European Free Trade Association.

Any feedback and questions on this document should be directed to the users’ national standards body. A complete listing of these bodies can be found on the CEN website.

According to the CEN/CENELEC Internal Regulations, the national standards organisations of the following countries are bound to announce this Technical Specification: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

[SIST-TS CEN/TS 17784-2:2023](https://standards.iteh.ai/catalog/standards/sist/0c961d07-0d8a-4f9d-be15-3917d601fadb/sist-ts-cen-ts-17784-2-2023)

<https://standards.iteh.ai/catalog/standards/sist/0c961d07-0d8a-4f9d-be15-3917d601fadb/sist-ts-cen-ts-17784-2-2023>

Introduction

Micronutrients are considered to be, in plant nutrition, a number of elements known to be needed in small amounts for proper plant growth and development. The most common are Iron (Fe), Manganese (Mn), Molybdenum (Mo), Copper (Cu), Zinc (Zn) and Boron (B).

If an organo-mineral fertilizer contains a substance, or one of the substances in the mixture, which is intended to enhance the long term availability to plants of micronutrients in the EU fertilizing product, that substance can be either a chelating agent or a complexing agent.

The incorporation of heptagluconic acid as complexing agent in organo-mineral fertilizers is intended to enhance the long term availability to plants of micronutrients in such EU fertilizing products.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

[SIST-TS CEN/TS 17784-2:2023](https://standards.iteh.ai/catalog/standards/sist/0c961d07-0d8a-4f9d-be15-3917d601fad6/sist-ts-cen-ts-17784-2-2023)

<https://standards.iteh.ai/catalog/standards/sist/0c961d07-0d8a-4f9d-be15-3917d601fad6/sist-ts-cen-ts-17784-2-2023>

1 Scope

This document specifies a chromatographic method which allows the identification of heptagluconic acid (HGA) in organo-mineral fertilizers containing heptagluconic acid metal complexes.

NOTE For the complete names of the chelating agents mentioned in this document, see Annex D.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 12944-1, *Fertilizers and liming materials — Vocabulary — Part 1: General terms*

EN 12944-2, *Fertilizers and liming materials — Vocabulary — Part 2: Terms relating to fertilizers*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in EN 12944-1 and EN 12944-2 and the following apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

— ISO Online browsing platform: available at <https://www.iso.org/obp>

— IEC Electropedia: available at <https://www.electropedia.org/>

3.1

complexing agent

organic substance forming a flat or steric structure with one di- or tri-valent transition metal cation (zinc (Zn), copper (Cu), iron (Fe), manganese (Mn) or cobalt (Co))

4 Principle

The method is based on demetallation with phosphoric acid of the micronutrient HGA complex present in an aqueous solution of the sample.

The complexing agent is then identified and determined by high-performance liquid chromatography.

The separation is carried out on an NH₂ phase bonded to silica column and an aqueous solution of phosphoric acid and acetonitrile as eluent.

The detection is based on UV photometry at 210 nm.

WARNING — Users of this document should be familiar with normal laboratory practice. This document does not purport to address all of the safety issues, if any, associated with its use. It is the responsibility of the user to establish appropriate health and safety practices and to ensure compliance with any national regulatory conditions.

IMPORTANT — It is absolutely essential that tests conducted according to this document are carried out by suitably trained staff.

5 Interferences

Current knowledge on eventual interferences from other substances is summarized in the list below:

CEN/TS 17784-2:2022 (E)

- a) High concentrations of phosphate in the sample solution can create a large peak preventing the identification/determination of HGA.
- b) High concentrations of chloride, sulfate and nitrate do not interfere in the identification/determination of the complexing agent.
- c) The presence of the chelates of EDDHSA, [o,o] EDDHA, [o,o] EDDHMA, EDTA, DTPA, CDTA, HEEDTA, IDHA as well as the corresponding chelating agent do not interfere since they are separated from HGA.

These substances can be detected in the chromatogram by the appearance of a peak at larger retention times. Therefore, the presence of these kinds of substances shall be taken into account when successive injections are scheduled.

- d) The presence of gluconic acid does interfere in the determination of the complexing agent.
- e) The presence of aspartic acid, humic substances and lignosulfonic acid can interfere with the identification/determination of HGA.

6 Reagents

All reagents shall be of recognized analytical grade.

6.1 Water

All water used should conform to EN ISO 3696, be degassed and be free of organic contaminants.

6.2 Sample preparation solvent.

Add to 800 ml of water, 2 ml of *ortho*-phosphoric acid 85 % and 25 ml of methanol in a 1 l volumetric flask. Dilute to the mark with water and homogenize.

6.3 HGA stock solution.

The mass concentration (ρ) of this solution is $\rho(\text{HGA}) = 19\,893 \text{ mg/l}$.

This solution shall be freshly prepared daily, because of the formation of the corresponding lactone if it is let standing for a long period of time.

Weigh to the nearest 0,1 mg about 2 500 mg of the heptagluconic acid, sodium salt dihydrate (CAS # 10094-62-9) 99 % (mass concentration), add 50 ml of water in a 100 ml volumetric flask. After dissolution, dilute to the mark with water and homogenize.

6.4 Eluent A: *ortho*-phosphoric acid.

The molar concentration (c) of this solution is $c(\text{H}_3\text{PO}_4) = 30 \text{ mmol/l}$ and methanol.

Add to 800 ml of water, 2 ml of *ortho*-phosphoric acid 85 % (mass concentration) and 25 ml of methanol high-performance liquid chromatography (HPLC) grade in a 1 l volumetric flask. Dilute to the mark with water and homogenize. Before use, filter the solution through a 0,45 μm membrane filter (7.5).

6.5 Eluent B: acetonitrile (HPLC-grade).

7 Apparatus

Usual laboratory equipment, glassware, and the following:

7.1 Magnetic stirrer.

7.2 Chromatograph.

Equipped with:

- an isocratic pump delivering the eluent at a flow rate of 1 ml/min;
- an injection valve with a 20 µl injection loop or equivalent;
- a NH₂ column¹; internal diameter: 4,6 mm; column length: 250 mm; dp = 5 µm;
- a NH₂ guard-column (recommended);
- a UV-VIS detector with a 210 nm-filter;
- an integrator.

7.3 Chromatographic conditions.

According to Table 1.

Table 1 — Chromatographic conditions

Flow rate	Eluent A (6.4)	Eluent B (6.5)
	%	%
1 ml/min	75	25

7.4 Balance.

Balance, with a maximum permissible error of ±0,1 mg.

7.5 Membrane filters.

Micro membrane filters resistant to aqueous solutions, with porosity of respectively 0,45 µm and 0,2 µm.

8 Sampling and sample preparation

Sampling and sample preparation are not part of the method specified in this document.

Recommended sampling methods are given in EN 1482-1 and, for sample preparation, in EN 1482-2.

9 Procedure

9.1 Preparation of the HGA-metal complex sample solution.

The mass of the test portion to be used to prepare the sample solution is dependent on the declared metal content of the product.

¹ PhenoSphere™ NH₂ 80A 5 µm 250 mm x 4,6 mm or equivalent is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by CEN of this product.

CEN/TS 17784-2:2022 (E)

Weigh into a 150 ml beaker, approximately the amount of sample specified in Table 2, to the nearest 0,1 mg:

Table 2 — Amount of sample

Declared metal content % (mass fraction)	Mass of test portion mg
10 ≤ to 15	300
5 ≤ to < 10	500
< 5	1 000

Add 50 ml of sample preparation solvent (6.2) and dissolve it with a magnetic stirrer (7.1) during 5 min. Make up to volume in a 100 ml volumetric flask with sample preparation solvent (6.2). Let the solution stand overnight in darkness to allow the metal phosphate to form.

9.2 Preparation of the calibration solutions.

Pipette a volume (V in ml) (see Table 3) of the HGA stock solution (6.3) in six 100 ml volumetric flasks respectively. Make up to volume with the sample preparation solution (6.2) and homogenize. Let the solution stand overnight.

Table 3 — Composition of the calibration solutions

Solution	V ml	Concentration of HGA mg HGA/l
1	1	199
2	2	398
3	6	1 194
4	8	1 591
5	10	1 989
6	16	3 183

NOTE The molecular mass of heptagluconic acid, sodium salt dihydrate corresponds to 284 g/mol, whereas HGA has a molecular mass of 226 g/mol.

9.3 Chromatographic analysis.

Immediately before injection, all calibration and sample solutions shall be filtered through a 0,2 μ m membrane filter (7.5).

After stabilization of the chromatographic conditions (7.3), inject the calibration solutions (9.2) into the chromatographic system (7.2).

The major peak obtained corresponds to heptagluconic acid.

NOTE 1 Since the calibration solutions are not freshly prepared (see 9.2), two defined peaks can appear in the chromatograms, one tentatively assigned to the lactone and the other corresponding to the heptagluconic acid.

Adjust the attenuation on the integrator, in order to obtain a suitable range for the HGA peak from the standard solution. A typical chromatogram is given in Figure A.1. Measure the retention time.