



SLOVENSKI STANDARD
oSIST prEN 17430:2022

01-september-2022

Kemična razkužila in antiseptiki - Higieno razkuževanje rok z drgnjenjem z virucidnim sredstvom - Preskusna metoda in zahteve (faza 2, stopnja 2)

Chemical disinfectants and antiseptics - Hygienic handrub virucidal - Test method and requirements (phase 2, step 2)

Chemische Desinfektionsmittel und Antiseptika - Viruzide hygienische Händedesinfektion - Prüfverfahren und Anforderungen (Phase 2, Stufe 2)

Désinfectants chimiques et antiseptiques - Traitement hygiénique de mains par frictions virucide - Méthode d'essai et prescriptions (phase 2, étape 2)

Ta slovenski standard je istoveten z: prEN 17430

ICS:

11.080.20 Dezinfektanti in antiseptiki Disinfectants and antiseptics

oSIST prEN 17430:2022

en,fr,de

EUROPEAN STANDARD
NORME EUROPÉENNE
EUROPÄISCHE NORM

DRAFT
prEN 17430

August 2022

ICS 11.080.20

English Version

Chemical disinfectants and antiseptics - Hygienic handrub virucidal - Test method and requirements (phase 2, step 2)

Désinfectants chimiques et antiseptiques ; Traitement
hygiénique de mains par frictions virucide ; Méthode
d'essai et prescriptions (phase 2, étape 2)

Chemische Desinfektionsmittel und Antiseptika -
Viruzide hygienische Händedesinfektion -
Prüfverfahren und Anforderungen (Phase 2, Stufe 2)

This draft European Standard is submitted to CEN members for enquiry. It has been drawn up by the Technical Committee CEN/TC 216.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

Warning : This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

Contents	Page
European foreword.....	3
1 Scope	4
2 Normative references	4
3 Terms and definitions	4
4 Requirements	5
5 Test methods	5
5.1 Principle.....	5
5.2 Materials and reagents.....	5
5.2.1 Test organism.....	5
5.2.2 Culture media and reagents.....	6
5.3 Apparatus and glassware.....	8
5.3.1 General.....	8
5.3.2 Usual microbiological laboratory equipment and, in particular, the following:.....	8
5.4 Preparation of test organism suspensions and product test solutions.....	9
5.4.1 Test virus suspension (test organisms suspensions, contamination fluid).....	9
5.4.2 Cell lines.....	10
5.5 Cell culture preparation for virucidal testing.....	10
5.6 Infectivity assay.....	10
5.6.1 Quantal tests (end point titration).....	10
5.6.2 Plaque assay.....	10
5.6.3 Product test solutions.....	11
5.7 Procedure for assessing the virucidal activity of the product on volunteers' hands.....	11
5.7.1 General.....	11
5.7.2 Test procedure with volunteers.....	12
6 Experimental data and calculation	14
6.1 Protocol of the results.....	14
6.2 Calculation of virus titre.....	14
6.3 Calculation of PFU.....	14
6.4 Calculation of the lg reduction R (lg prevalue minus lg postvalue).....	14
6.5 Verification of the methodology — Test validation.....	14
6.6 Statistical evaluation (significance testing), expression of results and precision.....	14
6.7 Conclusion.....	15
6.8 Test report.....	15
Annex A (normative) Standard handrub procedure	17
Annex B (informative) Calculation of the viral infectivity titre	18
B.1 Example of TCID ₅₀ determination by the Spearman-Kärber method.....	18
B.2 Calculation of plaque-forming units (PFU).....	18
Annex C (informative) Quality control of soft soap	20
Annex D (informative) Examples of reporting of results and significance testing	21
Annex E (normative) Test for non-inferiority	28
Bibliography	30

European foreword

This document (prEN 17430:2022) has been prepared by Technical Committee CEN/TC 216 “Chemical disinfectants and antiseptics”, the secretariat of which is held by AFNOR.

This document is currently submitted to CEN-Enquiry.

iTeh STANDARD PREVIEW
(standards.iteh.ai)

[oSIST prEN 17430:2022](https://standards.iteh.ai/catalog/standards/sist/606577ae-e7f0-40aa-9abd-42b5ba4792eb/osist-pren-17430-2022)

<https://standards.iteh.ai/catalog/standards/sist/606577ae-e7f0-40aa-9abd-42b5ba4792eb/osist-pren-17430-2022>

prEN 17430:2022 (E)

1 Scope

This document specifies a test method simulating practical conditions for establishing whether a product for hygienic handrub reduces the release of virus contamination on hands when rubbed onto the artificially contaminated hands of volunteers.

NOTE 1 Attention is drawn to the fact that tests on human volunteers are the subject of legal provisions in certain European countries/regions.

This document applies to products for hygienic handrub for use in areas and situations where disinfection is medically indicated. Such indications occur in patient care, for example:

- in hospitals, in community medical facilities and in dental institutions;
- in clinics of schools, of kindergardens and of nursing homes;

and can occur in the workplace and in the home. It can also include services such as laundries and kitchens supplying products directly for the patient.

EN 14885 specifies in detail the relationship of the various tests to one another and to “use recommendations”.

NOTE 2 This method corresponds to a phase 2, step 2 test.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 12353, *Chemical disinfectants and antiseptics — Preservation of test organisms used for the determination of bactericidal (including Legionella), mycobactericidal, sporicidal, fungicidal and virucidal (including bacteriophages) activity*

EN 14885, *Chemical disinfectants and antiseptics — Application of European Standards for chemical disinfectants and antiseptics*

EN 14476, *Chemical disinfectants and antiseptics — Virucidal quantitative suspension test for chemical disinfectants and antiseptics used in human medicine — Test method and requirements (phase 2/step 1)*

EN 16777, *Chemical disinfectants and antiseptics - Quantitative non-porous surface test without mechanical action for the evaluation of virucidal activity of chemical disinfectants used in the medical area - Test method and requirements (phase 2/step 2)*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in EN 14885 apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

- IEC Electropedia: available at <http://www.electropedia.org/>
- ISO Online browsing platform: available at <http://www.iso.org/obp>

4 Requirements

When tested in accordance with Clause 5, the mean reduction of the test organism murine norovirus strain S99 Berlin achieved by the hygienic handrub with the product under test shall be at least not inferior to that achieved by a specified reference hygienic handrub (a volume fraction of 70 % concentration of ethanol).

Table 1 — Minimum and additional test conditions

	Hygienic handrub
Test virus	<p>Virucidal activity^a</p> <p>Murine norovirus</p> <p>Limited spectrum virucidal activity^b</p> <p>Murine norovirus</p>
Test temperature	room temperature of between 21,5 °C ± 3,5 °C
Contact time	according to the manufacturer's recommendation, but between 30 s and 120 s
Additional conditions ^d	Further contact time(s)
<p>^a To claim the virucidal activity the product shall pass standards EN 14476 with <i>Poliovirus</i>, <i>Adenovirus</i> and <i>Murine Norovirus</i> and EN 16777 with <i>Adenovirus</i> and <i>Murine Norovirus</i>.</p> <p>^b To claim limited spectrum virucidal activity the product shall pass EN 14476 with <i>Adenovirus</i> and <i>Murine Norovirus</i>.</p>	

5 Test methods

5.1 Principle

Hands of volunteers are artificially contaminated with test organisms. The number of test organisms released from their fingertips into sampling fluids is assessed before and after the hygienic handrub. The ratio of the two resulting values (virus titres) represents a measure for the antimicrobial activity of the product tested. The necessary precision is achieved by repeating the test on 18 to 22 volunteers. To compensate for extraneous influences it is compared with the reduction obtained by a reference handrub which is performed with the same volunteers, on the same day and under comparable environmental conditions.

5.2 Materials and reagents

5.2.1 Test organism

Murine norovirus strain S99 Berlin¹

¹ Murine norovirus can be obtained from Friedrich-Loeffler-Institut Bundesinstitut für Tiergesundheit, Hauptsitz Insel Riems Südufer 10, 17493 Greifswald-Insel Riems; phone: +49 38351 7-0; fax: +49 38351 7-121. <http://www.fli.de/>. This

prEN 17430:2022 (E)

This test organism has been specifically chosen to meet health and safety guidance and ethical committee considerations. According to German Ordinance on Biological Substances (BiostoffV)/TRBA 462 Murine Norovirus is classified in risk group 1.

5.2.2 Culture media and reagents**5.2.2.1 General**

All weights of chemical substances given in this document refer to the anhydrous salts. Hydrated forms may be used as an alternative, but the weights required shall be adjusted to allow for consequent molecular weight differences.

The reagents shall be of analytical grade and/or appropriate for microbiological purposes. They shall be free from substances that are toxic or inhibitory to the test organisms. To improve reproducibility, it is recommended that commercially available dehydrated material is used for the preparation of culture media. The manufacturer's instructions relating to the preparation of these products should be rigorously followed. For each culture medium and reagent, a time limitation for use should be fixed.

5.2.2.2 Water

The water shall be freshly glass-distilled water and not demineralized water. If distilled water of adequate quality is not available, water for injections (see bibliographic reference [3]) may be used.

Sterilize in the autoclave [5.3.2.1 a)]. Sterilization is not necessary if the water is used e.g. for preparation of culture media and subsequently sterilized.

If the water is sterilized during sterilization of the reagents, this is not necessary.

5.2.2.3 Phosphate buffered saline (PBS)

Sodium chloride (NaCl)	8,00 g
Potassium chloride (KCl)	0,20 g
Disodium hydrogen phosphate, 12-hydrate ($\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$)	2,89 g
Potassium phosphate, monobasic (KH_2PO_4)	0,20 g
Water (5.2.2.2)	to 1 000,0 ml

5.2.2.4 Neutral Red (1:1000 solution) or crystal violet solution

Prepare neutral red (e.g. Sigma N7005²) stock solution at 0,1 mg/ml in water (5.2.2.2).

Prepare 0,1 % crystal violet in 20 % ethanol.

Filter through a 0,40 µm pore size filter and store 4 °C in the dark.

5.2.2.5 Fetal calf serum (FCS)

FCS shall be certified free of viruses and mycoplasma. Extraneous viruses and mycoplasma may interfere with cell and virus growth resulting in false results.

For RAW 264.7 cells, special FCS shall be used due to the cells' high sensitivity to endotoxins.

information is given for the convenience of users of this document and does not constitute an endorsement by CEN of the product named.

² Sigma N7005 is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by CEN of this product.

5.2.2.6 Trichloroacetic acid (10 % solution) (TCA)

Dissolve 10 g of TCA crystals in 80 ml of water (5.2.2.2), then adjust the volume to 100 ml with water. Stir to complete solution.

5.2.2.7 Hard water

For the preparation of 1 l of hard water, the procedure is as follows:

- prepare solution A: dissolve 19,84 g magnesium chloride (MgCl_2) and 46,24 g calcium chloride (CaCl_2) in water (5.2.2.2) and dilute to 1 000 ml. Sterilize by membrane filtration (5.3.2.13) or in the autoclave [5.3.2.1 a)]. Autoclaving – if used – may cause a loss of liquid. In this case make up to 1 000 ml with water (5.2.2.2) under aseptic conditions. Store the solution in the refrigerator (5.3.2.17) for no longer than one month;
- prepare solution B: dissolve 35,02 g sodium bicarbonate (NaHCO_3) in water (5.2.2.2) and dilute to 1000 ml. Sterilize by membrane filtration (5.3.2.13). Store the solution in the refrigerator (5.3.2.17) for no longer than one week;
- place 600 ml to 700 ml of water (5.2.2.2) in a 1 000 ml volumetric flask (5.3.2.10) and add 6,0 ml (5.3.2.9) of solution A, then 8,0 ml of solution B. Mix and dilute to 1 000 ml with water (5.2.2.2). The pH (5.3.2.4) of the hard water shall be $7,0 \pm 0,2$. If necessary, adjust the pH by using a solution of approximately 40 g/l (about 1 mol/l) of sodium hydroxide (NaOH) or approximately 36,5 g/l (about 1 mol/l) of hydrochloric acid (HCl).

The hard water shall be freshly prepared under aseptic conditions and used within 12 h.

NOTE When preparing the product test solutions (5.6.3), the addition of the product to the hard water produces different final water hardness in each test tube. In any case, the final hardness in the test tube expressed as calcium carbonate (CaCO_3) is lower than 375 mg/l.

5.2.2.8 Growth and maintenance media

Dulbecco's Modified Eagle's Medium (DMEM) or equivalent, supplemented with appropriate concentration of heat inactivated (e.g. 56 °C for 30 min) and mycoplasma-free fetal calf serum FCS (5.2.2.5), antibiotics, and other growth factors as needed shall be used. Additional supplements like glucose, glutamine, and sodium pyruvate shall be indicated in the test report.

- a) A *growth medium* for cell multiplication is supplemented with 10 % FCS (5.2.2.5). Add 10 parts of FCS to 90 parts of DMEM.
- b) A *maintenance medium* to maintain the cell culture metabolism without stimulation of cell proliferation is supplemented with 2 % FCS. Add 2 parts of FCS (5.2.2.5) to 98 parts of DMEM.

Other media may be used if appropriate for certain cell lines.

See also bibliographic reference [4]. See EN 12353 for a detailed description.

NOTE Materials and reagents for cell culture can be purchased from biological supply companies.

5.2.2.9 Diluted soft soap

Linseed oil	50 parts by weight
Potassium hydroxide	9,5 parts by weight
Ethanol (min. 95 %)	7 parts by weight

prEN 17430:2022 (E)

Hot distilled water ($75\text{ °C} \pm 5\text{ °C}$) as needed

Prepare a solution of 9,5 parts potassium hydroxyde in 15 parts water (5.2.2.2) and add 50 parts linseed oil. Heat up to approximately 70 °C while constantly stirring. Add the ethanol and continue heating while stirring until the saponification process is completed and a sample dissolves clearly in water and almost clearly in alcohol. The weight of the soft soap is then brought up to 100 parts by addition of water (5.2.2.2), heated up to $75\text{ °C} \pm 5\text{ °C}$ to dilute the soft soap. Take 200 g of the soft soap, fill up to 1 000 g with water (5.2.2.2) and sterilize in the autoclave [5.3.2.1 a)]. The pH of the final diluted soft soap shall range between 10,0 and 11,0.

For quality control of the soft soap, see Annex C.

5.2.2.10 Ethanol as reference handrub [volume fraction of 70 % at $21,5\text{ °C} \pm 3,5\text{ °C}$]

Fill 547,2 g ethanol with a purity of min volume fraction of 99,5 % (determined by gas chromatography; density 0,790) in a 1 000 ml flask equipped with a glass stopper on the weighing platform of a scale (precision 0,1 g). Add 322,7 g water (5.2.2.2). Close the flask with the matching glass stopper and shake the contents of the flask thoroughly. During mixing the solution warms up about to 30 °C . After cooling down to $21,5\text{ °C} \pm 3,5\text{ °C}$, it will give a volume of approximately 980 ml.

NOTE This solution can be kept indefinitely at approximately room temperature if protected from light, provided that the flask is closed with a ground-in glass stopper to avoid any evaporation.

5.3 Apparatus and glassware**5.3.1 General**

Sterilize all glassware and parts of the apparatus that will come into contact with the culture media and reagents or the sample, except those which are supplied sterile, by one of the following methods:

a) by moist heat, in the autoclave [5.3.2.1 a)];

b) by dry heat, in the hot air oven [5.3.2.1 b)].

5.3.2 Usual microbiological laboratory equipment³ and, in particular, the following:**5.3.2.1 Apparatus for sterilization (moist and dry heat).**

a) For moist heat sterilization, an autoclave capable of being maintained at $(121_0^{+3})\text{ °C}$ for a minimum holding time of 15 min;

b) for dry heat sterilization, a hot air oven capable of being maintained at $(180_0^{+5})\text{ °C}$ for a minimum holding time of 30 min, at $(170_0^{+5})\text{ °C}$ for a minimum holding time of 1 h or at $(160_0^{+5})\text{ °C}$ for a minimum holding time of 2 h.

5.3.2.2 Water baths, capable of being controlled at $21,5\text{ °C} \pm 3,5\text{ °C}$.

5.3.2.3 Inverted microscope for reading cell cultures microscopically.

5.3.2.4 pH meter, having an inaccuracy of calibration of 0,1 pH units at $21,5\text{ °C} \pm 3,5\text{ °C}$.

5.3.2.5 Stopwatch

³ Disposable sterile equipment is an acceptable alternative to reusable glassware.

- 5.3.2.6 Deep freezer** (-20°C and -70°C or less).
- 5.3.2.7 Electromechanical agitator** e.g. Vortex[®] mixer⁴.
- 5.3.2.8 Containers:** sterile test tubes, culture bottles or flasks of suitable capacity.
- 5.3.2.9 Graduated pipettes** of nominal capacities 10 ml, 1 ml and 0,1 ml. Calibrated automatic pipettes, with disposable tips, may be used.
- 5.3.2.10 Volumetric flasks**, calibrated at $21,5^{\circ}\text{C} \pm 3,5^{\circ}\text{C}$.
- 5.3.2.11 Sterile microtitre plates**, six or eight well plates for cell cultures.
- 5.3.2.12 flasks** for cell cultures.
- 5.3.2.13 Membrane filtration apparatus** for filtration of media, $0,2\ \mu\text{m}$ pore size.
- 5.3.2.14 CO₂ incubator** (95 % air, 5 % CO₂), capable of being controlled at either $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$, or at $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for incubation of cell cultures.
- 5.3.2.15 Biological safety cabinet**, class II.
- 5.3.2.16 Centrifuge** ($400\ g_{\text{N}}$ to $1000\ g_{\text{N}}$).
- 5.3.2.17 Refrigerator**, capable of being controlled at 2°C to 8°C .

5.4 Preparation of test organism suspensions and product test solutions

5.4.1 Test virus suspension (test organisms suspensions, contamination fluid)

The test organisms and their stock cultures shall be prepared and kept in accordance with EN 12353.

The stock virus suspension is multiplied in an appropriate cell line that produces high titres of infectious viruses. The cell debris is separated by centrifugation (5.3.2.16) ($400\ g_{\text{N}}$ for 15 min). This preparation is called "test virus suspension".

It is suggested that the minimum titre of the virus suspension - determined by a quantal test (5.6.1) or by plaque test (5.6.2) - is at least $10^8\ \text{TCID}_{50}/\text{ml}$. In any case, it shall be sufficiently high to at least $10^4\ \text{TCID}_{50}/\text{ml}$ may be recovered from the hands.

In exceptional cases the test virus suspension may be concentrated by appropriate methods (e.g. ultracentrifugation).

The test virus suspension is kept in small volumes below -70°C or preferably at -196°C under nitrogen. Due to safety reasons, and – in some cases – to limit the possibility of genetic mutations, only 10 passages from the original seed (e.g. virus from culture collection) are allowed.

The test virus suspension is used undiluted for the test procedure (5.7).

The titre of the virus suspension is determined as described (see 6.2 and 6.3).

⁴ Vortex[®] is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by CEN of this product.

prEN 17430:2022 (E)**5.4.2 Cell lines**

Murine norovirus is multiplied in RAW 264.7 cells (ATCC TIB-71) or other cell lines of appropriate sensitivity.

5.5 Cell culture preparation for virucidal testing

Cells used as suspension in quantal tests shall be added to the dilutions of the test mixture in such a density to enable the formation of a monolayer in at least two days in the cell control (see 5.6.1.1).

Cell monolayers shall be > 90 % confluent before inoculation. Cell lines and culture media will be selected in accordance with the test virus (see 5.2.1).

5.6 Infectivity assay**5.6.1 Quantal tests (end point titration)****5.6.1.1 Virus titration on cells in suspension**

Dilute the test virus suspension (5.4.1) by tenfold dilution steps with a volume of 0,5 ml or 0,1 ml of test virus suspension plus an appropriate volume of maintenance medium (see 5.2.2.8). Pipettes shall be changed after each dilution step.

Transfer 0,1 ml of each dilution into six or eight wells of a microtitre plate (5.3.2.11), beginning with the highest dilution. Add 0,1 ml of cell culture in such a density to enable the formation of a monolayer in at least 2 days in the cell control. The last row of six or eight wells (5.3.2.11) shall not receive any test virus suspension and serve as the cell control.

After the appropriate incubation time, according to the virus type, the viral cytopathic effect is read using an inverted microscope (see 5.3.2.3).

5.6.1.2 Virus titration on cell monolayers

Transfer 0,1 ml of each dilution (see 5.6.1.1) into six or eight wells (5.3.2.11) of a microtitre plate containing a confluent (>90 %) cell monolayer. The last row of six or eight wells will receive 0,2 ml of maintenance medium (see 5.2.2.8) and will serve as cell control. If the cells were inoculated without any medium after 1 h of incubation at 37 °C, 0,1 ml of maintenance medium is added to each well.

The plates are incubated at 37 °C for the appropriate incubation time. The viral cytopathic effect is read using an inverted microscope (see 5.3.2.3).

5.6.2 Plaque assay

Plastic tray wells (surface diameter 30 mm to 35 mm) with confluent cell monolayers are washed once with phosphate buffered saline (PBS) (5.2.2.3) and inoculated with 0,2 ml of serial dilutions of virus in DMEM + 2 % FCS (see 5.2.2.8). Three wells are generally used per dilution. After adsorption period of 1 h at 37 °C, during which the cell monolayers are kept moist by tilting the dishes every 8 min to 10 min, the inoculum is removed and the cell monolayers are washed once with PBS. Subsequently, the wells are overlaid with 3 ml of a mixture consisting of 2 % melted agarose or another appropriate semisolid medium and 2times concentrated DMEM with 4 % FCS. The cultures are incubated for 5 days to 6 days at 37 °C in a CO₂ incubator (see 5.3.2.14). Plaques can be counted after addition of 2 ml of a second overlay with the same composition of the first and also containing 5 % of a 1:1 000 solution of neutral red and further incubation (in the dark) at 37 °C for 24 h to 48 h in a CO₂ incubator (see 5.3.2.14). Counting can be performed also after addition of crystal violet. The cell monolayers are fixed by adding 2 ml of 10 % trichloroacetic acid (TCA) (5.2.2.6) over the agar overlay for 10 min to 15 min at 21,5 °C ± 3,5 °C. The agar overlay is then removed and 2 ml of 0,1 % crystal violet in 20 % ethanol are added. After 10 min to 15 min at 21,5 °C ± 3,5 °C, the wells are extensively washed with water and the plaques (white spots) are counted.