



Standard Test Method for Mercurous Nitrate Test for Copper Alloys¹

This standard is issued under the fixed designation B154; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reappraisal. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reappraisal.

This standard has been approved for use by agencies of the Department of Defense.

1. Scope*

1.1 This test method describes the technique for conducting the mercurous nitrate test for residual stresses in wrought copper alloy mill products.

~~NOTE 1—For any particular copper alloy, reference should be made to the material specification.~~

~~NOTE 1—For any particular copper alloy, reference should be made to the material specification.~~

~~NOTE 2—Test Method B858 may be considered as a possible alternative test method which does not involve the use of mercury.~~

~~NOTE 3—This test method is considered historically reliable for determining the potential state of residual stress in copper alloys, but not promoted for use due to the hazards relating to mercury use and environmentally appropriate disposal.~~

1.2 *Units*—The values stated in SI units are to be regarded as the standard. The values given in parentheses are for information only.

1.3 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determines the applicability of regulatory limitations prior to use.* For specific precautionary and hazard statements see Sections 1, 6, and 7. (~~Warning—Mercury is a definite health hazard in use and disposal.~~)has been designated by EPA and many state agencies as a hazardous material that can cause central nervous system, kidney and liver damage. Mercury, or its vapor, may be hazardous to health and corrosive to materials. Caution should be taken when handling mercury and mercury containing products. See the applicable product Material Safety Data Sheet (MSDS) for details and EPA's website – <http://www.epa.gov/mercury/faq.htm> – for additional information. Users should be aware that selling mercury and/or mercury containing products into your state may be prohibited by state law.)

2. Referenced Documents

2.1 *ASTM Standards:*²

[B846 Terminology for Copper and Copper Alloys](#)

[B858 Test Method for Ammonia Vapor Test for Determining Susceptibility to Stress Corrosion Cracking in Copper Alloys](#)

[D1193 Specification for Reagent Water](#)

3. Terminology

3.1 For terms related to copper and copper alloys, refer to Terminology [B846](#).

4. Summary of Test Method

4.1 The prepared test specimen is completely immersed in the mercurous nitrate test solution for 30 min at ambient temperature. Upon removal from the solution, the test specimen is ~~washed~~wiped and immediately examined visually for cracks.

5. Significance and Use

5.1 This test method is an accelerated test for detecting the presence of residual (internal) stresses that might result in failure of individual parts in storage or in service due to stress corrosion cracking.

5.2 This test method is not intended for use on assemblies or parts under applied stress. If used for that purpose, the results shall be for information only and not a cause for rejection of the assembly, its component parts, or the original mill product.

6. Reagents and Materials

6.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents conform to the specifications of the Committee on Analytical Reagent of the American Chemical Society where such

*A Summary of Changes section appears at the end of this standard

specifications are available.³ Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

6.2 *Purity of Water*—Unless otherwise indicated, references to water shall be understood to mean Type IV reagent water or better, as defined of Specification **D1193**.

6.3 *Mercurous Nitrate Solution*—The solution shall be an aqueous mercurous nitrate solution containing 10 g of mercurous nitrate solution (HgNO_3) and 10 mL of nitric acid (HNO_3) (sp gr 1.42) per litre of solution.

6.4 *Preparation*—The aqueous mercurous nitrate solution shall be prepared by either of the following procedures, A or B. Used solutions may be replenished as described in 6.5.

6.4.1 *Procedure A*—Dissolve 11.4 g of $\text{HgNO}_3 \cdot 2\text{H}_2\text{O}$ or 10.7 g of $\text{HgNO}_3 \cdot \text{H}_2\text{O}$ in approximately 40 mL of distilled water acidified with 10 mL of HNO_3 (sp gr 1.42). After the crystals are completely dissolved, dilute the solution with water to 1000 mL. (**Warning**—The mercurous nitrate crystals are obtainable in both the monohydrate and dihydrate form and should be handled with caution because of their highly toxic effects.) (**Warning**—When weighing crystals, the weight of the water of crystallization should be taken into consideration. The mercurous nitrate crystals are photosensitive and when they have turned yellow are difficult to dissolve.) (**Warning**—Care should be exercised when handling and mixing chemicals. Qualified personnel using appropriate chemical-laboratory techniques should only do the handling and mixing.)

6.4.2 *Procedure B*—Dissolve 76 g of mercury in 114 mL of diluted HNO_3 (1 part water to 1 part HNO_3) (sp gr 1.42). Carefully dilute with distilled water to 1000 mL. This provides a concentration of 100 g of HgNO_3 after a slight loss due to heating. Add the water in small portions while stirring to prevent local overdilution. This gradual dilution, together with the excess acid, will prevent precipitation of basic salts of mercury. Dilute 100 mL of this solution (10 %) with 7 mL of HNO_3 (sp gr 1.42) and 893 mL of water. (**Warning**—Mercury is a definite health hazard and therefore equipment for the detection and removal of mercury vapor produced in volatilization is recommended. The use of rubber gloves in testing is advisable.)

6.5 *Replenishment of Solution*—The spent solution may be reclaimed by replenishing the mercurous nitrate solution, to a 1 volume percent concentration, as follows:

6.5.1 Measure 50 mL of the spent HgNO_3 solution in a graduated cylinder.

6.5.2 Transfer to an Erlenmeyer flask, and add 10 mL of HNO_3 (1 + 1).

6.5.3 Add slowly 1 weight per volume percent potassium permanganate (KMnO_4) solution from a buret with a constant shaking until there is an excess as indicated by the pink color, which persists for several minutes.

6.5.4 Add iron (II) sulfate (FeSO_4) crystals until the solution, when shaken, becomes clear. Then titrate the solution with 0.1 N potassium thiocyanate (KCNS) solution to the appearance of a reddish brown color. Repeat this procedure with 50 mL of a standard 1 weight per volume percent of HgNO_3 solution.

6.5.5 The ratio, *R*, of the number of millilitres of KCNS solution required to titrate the spent solution, to the number of millilitres required to titrate the standard solution, determines the number of millilitres, *X*, of 10 volume percent HgNO_3 in 3 volume percent HNO_3 solution required to replenish 1 L of spent solution. Values of *R* and *X* for a litre volume are given in **Table 1**.

<https://standards.iteh.ai/catalog/standards/sist/6a3dec89-e0c4-44a6-9a17-cd15ea548f62/astm-b154-12>

7. Hazards

7.1 **Warning**—Mercury is a definite health hazard in use and disposal.

7.2 *Suggested Mercurous Nitrate Disposal:*

7.2.1 To mercurous nitrate solutions add sodium hydroxide (NaOH) to pH 10 to 11.

7.2.2 Filter precipitated mercury and other heavy metals.

7.2.3 Though the filtrate is low in free mercurous or mercuric ions, it must be further treated before disposal.

7.2.4 To each litre of filtrate, add two drops (0.1 cm³) of 24 volume percent ammonium sulfide ($(\text{NH}_4)_2\text{S}$).

7.2.5 After the second filtering, the filtrate may be discarded.

NOTE 4—If heating is used in either of the previous procedures, the container should be covered with a watch glass to prevent loss of HNO_3 and water to the atmosphere. After solution is complete, use a small volume of retained dilution water to rinse the watch glass into the container.

7.2.5.1 Monitor the filtrate to assure it meets appropriate health safety standards, or is disposed of properly.

7.2.6 The precipitates should be collected and stored with the mercury-contaminated test samples and sold to a licensed mercury disposal-service.

¹ This test method is under the jurisdiction of ASTM Committee **B05** on Copper and Copper Alloys and is the direct responsibility of Subcommittee **B05.06** on Methods of Test.

Current edition approved Oct. 1, 2005/Oct. 1, 2012. Published October 2005/November 2012. Originally approved in 1941. Last previous edition approved in 2004/2005 as **B154-01**/**B154-05**. DOI: [10.1520/B0154-05](https://doi.org/10.1520/B0154-05); [10.1520/B0154-12](https://doi.org/10.1520/B0154-12).

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

³ *Reagent Chemicals, American Chemical Society Specifications*, American Chemical Society, Washington, DC. For Suggestions on the testing of reagents not listed by the American Chemical Society, see *Analar Standards for Laboratory Chemicals*, BDH Ltd., Poole, Dorset, U.K., and the *United States Pharmacopeia and National Formulary*, U.S. Pharmacopeial Convention, Inc. (USPC), Rockville, MD.