



Designation: C1307 – 15

Standard Test Method for Plutonium Assay by Plutonium (III) Diode Array Spectrophotometry¹

This standard is issued under the fixed designation C1307; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This test method describes the determination of total plutonium as plutonium(III) in nitrate and chloride solutions. The technique is applicable to solutions of plutonium dioxide powders and pellets (Test Methods C697), nuclear grade mixed oxides (Test Methods C698), plutonium metal (Test Methods C758), and plutonium nitrate solutions (Test Methods C759). Solid samples are dissolved using the appropriate dissolution techniques described in Practice C1168. The use of this technique for other plutonium-bearing materials has been reported (1-5), but final determination of applicability must be made by the user. The applicable concentration range for plutonium sample solutions is 10–200 g Pu/L.²

1.2 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

1.3 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

2. Referenced Documents

2.1 ASTM Standards:³

C697 Test Methods for Chemical, Mass Spectrometric, and Spectrochemical Analysis of Nuclear-Grade Plutonium Dioxide Powders and Pellets

C698 Test Methods for Chemical, Mass Spectrometric, and

Spectrochemical Analysis of Nuclear-Grade Mixed Oxides ((U, Pu)O₂)

C757 Specification for Nuclear-Grade Plutonium Dioxide Powder, Sinterable

C758 Test Methods for Chemical, Mass Spectrometric, Spectrochemical, Nuclear, and Radiochemical Analysis of Nuclear-Grade Plutonium Metal

C759 Test Methods for Chemical, Mass Spectrometric, Spectrochemical, Nuclear, and Radiochemical Analysis of Nuclear-Grade Plutonium Nitrate Solutions

C833 Specification for Sintered (Uranium-Plutonium) Dioxide Pellets

C859 Terminology Relating to Nuclear Materials

C1168 Practice for Preparation and Dissolution of Plutonium Materials for Analysis

3. Terminology

3.1 Definitions:

3.1.1 For definitions of terms related to nuclear materials, refer to Terminology C859.

4. Summary of Method

4.1 In a diode array spectrophotometric measurement, as in a conventional spectrophotometric measurement, the substance being determined absorbs light at frequencies characteristic of that substance. The amount of light absorbed at each wavelength is directly proportional to the concentration of the species of interest. The absorption is a function of the oxidation state and the complexation obtained in the solution matrix selected for measurement. Beer's Law permits quantifying the species of interest relative to a traceable standard when both solutions are measured under the same conditions. The array of photosensitive diodes permits the (virtually) simultaneous collection of spectral information over the entire range of the instrument, for example, 190–820 nm (or any selected portion of the range). An entire absorption spectrum can be obtained in 0.1 s; however, optimum precision is obtained from multiple spectra collected over a 4-s period.

4.2 Reduction to plutonium(III) is accomplished by the addition of a measured quantity of reductant solution to the sample aliquant.

4.2.1 For nitrate solutions, ferrous sulfamate is the recommended reductant. Aliquants (1 mL or less) of the sample

¹ This test method is under the jurisdiction of ASTM Committee C26 on Nuclear Fuel Cycle and is the direct responsibility of Subcommittee C26.05 on Methods of Test.

Current edition approved Jan. 1, 2015. Published January 2015. Originally approved in 1995. Last previous edition approved in 2014 as C1307 – 14. DOI: 10.1520/C1307-15.

² For solid samples, select the sample size and dissolved solution weight to yield sample solutions in the 10–30 g Pu/L range. With special preparation and spectral analysis techniques, the method has been applied to nitrate solutions in the 0.1–10 g Pu/L range.

³ For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

solution are diluted with 10 mL of a ferrous reductant/matrix solution to 1 g Pu/L, and measured.

4.2.2 For chloride solutions, ascorbic acid is the recommended reductant. Aliquants of the sample solution, each containing 50–100 mg of plutonium, are diluted with 2 mL of zirconium solution to complex fluoride ions, 2 mL ascorbic acid reductant solution, and 1.0 M HCl to a total volume of 25 mL, yielding 2–4 g Pu/L solutions for measurement.

4.3 Plutonium concentration is determined from light absorption measurements taken on the sample solution in the blue-green region from 516 to 640 nm where a broad doublet band is observed. Spectral quantifying software capable of fitting the sample spectrum with spectral information from standard solutions is used to calculate the plutonium concentration. Both commercially available (6) and custom-designed (7-12) spectral fitting software have been developed which may be used for plutonium measurements. The users of this procedure are responsible for selecting or customizing, or both, the spectral fitting (and instrument control) software that best meets their individual measurement methodology and needs. Software selection will dictate many of the procedural specifics not included in this procedure. This procedure is intended to address key measurement requirements and to allow users discretion in establishing appropriate procedural details and technique variations. The software package selected should include a feature that indicates the quality of spectral fit, thereby providing information on the measurement reliability and the presence of interferences that absorb light or otherwise alter the plutonium(III) spectrum without requiring supplemental measurements.

5. Significance and Use

5.1 This test method is designed to determine whether a given material meets the purchaser's specification for plutonium content.

6. Interferences

6.1 Materials meeting the applicable material specifications of the ASTM standard for which this procedure was developed, when dissolved and diluted without introduction of interfering contaminants as described in Practice C1168, will contain no interfering elements or species.

NOTE 1—Fluoride, if present, would interfere if the zirconium, routinely added to the sample solution aliquant for the chloride matrix, were omitted from the procedure. Zirconium may be added to the nitrate matrix. Ferrous-Reductant Solution to handle fluorides if present. Zirconium, when used, should be added to all samples, blanks, and standards to obtain a consistent matrix. Refer to Specifications C833 and C757.

6.2 Interferences are caused by: (1) materials that absorb light in the region of the plutonium absorption, (2) undissolved solids that cause light scattering, (3) strong oxidizing or complexing agents that prevent complete reduction of the plutonium to the plutonium(III) oxidation state, and (4) anions that shift the spectrum.

6.2.1 Absorption of light in the region of interest by another species is a potential interference. Identification of potentially interfering species and inclusion of their spectra in the spectral curve fitting process will significantly reduce their effect. At a

minimum, sample measurements should be flagged when the higher than normal fitting error occurs, resulting from the presence of unidentified absorbing species. Enhancement of the spectral curve fitting capabilities of the DAS can be achieved by taking double derivatives of the spectrum collected. The spectral curve fitting software of the DAS is then used to quantitate the mathematically manipulated spectrum.

NOTE 2—Care must be taken in the choice of the preprocessing methods (derivatives, mean centering, autoscaling, or channel selection, or combinations thereof) as these may affect the robustness of the final model, particularly with regard to unknown interferences. Use of residual analysis will not always detect unknown interferences and results will vary depending on the preprocessing methods and models employed.

6.2.2 This spectrophotometric assay method should not be used on turbid (cloudy) solutions or solutions containing undissolved material. In addition to visual or turbidity meter measurements, or both, the presence of undissolved solids may be identified by the resulting shifts in the spectral baseline and by elevated spectral fitting errors.

NOTE 3—Plutonium oxides, mixed oxides, and plutonium metals meeting the material specifications for which this test method is intended, will dissolve when procedures in Practice C1168 are followed. Failure to achieve dissolution is an indication that the material does not meet the specifications, and the application of this test method for plutonium assay must be verified by the user. The user and customer are cautioned: when undissolved solids that persist after exhaustive dissolution efforts are to be removed by filtration through filter paper or other inert material of appropriate porosity, the subsequent plutonium assay measurements require close scrutiny. While filtration of undissolved solids may permit the reliable measurement of the concentration of plutonium in the filtrate, the resulting analysis may not be representative of their source material. Solids may indicate incomplete dissolution of the plutonium in the sample material, not necessarily a plutonium-free refractory residue. When this technique is utilized in support of preprocessing operations, process solutions containing solids may be an indication of incomplete dissolution of the plutonium-bearing material being processed or of an error in process operations. In addition to process control considerations, the undissolved solids may represent accountability and criticality control problems.

6.2.3 Strong oxidizing agents and complexing agents in sufficient concentration to prevent complete reduction typically are not present in plutonium nitrate samples. Appreciable concentrations of fluoride and sulfate anions have been found to interfere. The concentration of hydrofluoric acid, added to catalyze dissolution of oxides, may be removed by evaporation prior to measurement to ensure that the zirconium effectively complexes the traces of fluoride ion. Changes in the plutonium spectrum from incomplete reduction due to oxidizing agents and shifts in the spectrum due to complexing agents are also indicated by increases in the spectral curve fitting error.

6.2.4 Anion identity and concentration will shift the location and alter the shape of the absorption curve. The system calibration must include the anion shift effect by encompassing the expected range of anion identities and concentrations or by using appropriate spectral fitting features that identify and correct for the effect.

6.3 A study was conducted at the Los Alamos National Laboratory to determine the immunity of the Pu(III) spectrophotometric assay method to a diverse species of potential interferences. The elements studied were element numbers 1, 9, 11–13, 17, 19, 22–31, 35, 42, 44–46, 48, 50, 53, 57, 58, 60,

62, 73, 74, 76, 77, 79, 83, 90, 92, 93, and 95. Potential interference from nitrate, phosphate, sulfate, and oxalic acid is also documented (13).

7. Apparatus

7.1 *Diode Array Spectrophotometer (DAS)*—Wavelength range 190–820 nm; wavelength accuracy ± 2 nm; wavelength reproducibility ± 0.05 nm; full dynamic range 0.0022 to 3.3; photometric accuracy at 1 AU with a NBS 931 filter at 512 nm is ± 0.005 AU; baseline flatness < 0.0013 AU; noise at 500 nm is 0.0002 AU RMS; stray light measured with a Hoya 056 filter at 220 nm < 0.05 %;⁴

7.2 *Analytical Balance*—Readability of 0.1 mg; linearity 0.1 mg over any 10 g range and 0.2 mg over 160 g full scale.

7.3 *Solution Density Meter*—Readability of 0.1 mg/mL; precision of 0.3 mg/mL; linearity and accuracy 0.5 mg/mL over the range 0 to 2.0 g/mL.

7.4 *Adjustable, Fixed-Volume Pipetters*—Calibrated to deliver the desired range of volumes for sample and matrix-reductant solutions.

8. Reagents and Materials

8.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society where such specifications are available.⁵ Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

8.2 *Purity of Water*—Unless otherwise indicated, references to water shall be understood to mean distilled or deionized water.

8.3 *Ascorbic Acid-Reductant Solution* ($C_6H_8O_6$, aminoguanidine bicarbonate ($CH_6N_4H_2CO_3$), 0.4 M in each reagent)—Prepare fresh daily by dissolving 7 g of ascorbic acid and 5.5 g aminoguanidine bicarbonate in 80 mL of 1 M HCl, then dilute to a final volume of 100-mL 1 M HCl.

NOTE 4—The ascorbic acid is stabilized by the addition of aminoguanidine (Guanylhiazine: $HN:C(NH_2)NHNH_2$). The stabilized reductant solution has been found to be effective when ascorbic acid stability problems are encountered.

⁴ The optical specifications listed in 7.1 have been found to provide satisfactory results. In addition to these specifications, an acceptable spectrometer system should also provide multicomponent spectral quantitating software, computer control, and fiber optic capabilities. Systems which meet these specifications include, but are not limited to, the Hewlett Packard 8451A and 8452A, the Agilent 8453, and the TIDAS II. Fiber Optic Coupler, Fibers, and 1, 3, and 4-cm Flow Cell/Pump System have been developed and reported on (2, 3, 5). Specifications listed are those of the HP-8451A and the HP-8452A. Although the wavelength range for plutonium requires only a fraction of the 190–820 nm range specified (plutonium absorption spectrum is measured over the 520–634 nm region) spectrophotometers with significantly smaller range would be of little general use to the purchaser.

⁵ *Reagent Chemicals, American Chemical Society Specifications*, American Chemical Society, Washington, DC. For suggestions on the testing of reagents not listed by the American Chemical Society, see *Analar Standards for Laboratory Chemicals*, BDH Ltd., Poole, Dorset, U.K., and the *United States Pharmacopeia and National Formulary*, U.S. Pharmacopeial Convention, Inc. (USPC), Rockville, MD.

8.4 *Ferrous-Reductant Solution* (ferrous sulfamate, 0.05 M; sulfamic acid, 0.25 M; nitric acid, 1.0 M)—Prepare fresh weekly by adding 12 mL of freshly prepared ferrous sulfamate (2 M) to 90 mL of sulfamic acid (1.5 M). Stir, then add 175 mL of nitric acid (3.0 M) and dilute to 500 mL with water.

8.5 *Ferrous Sulfamate* ($Fe(NH_2SO_3)_2$, 2.0 M)—Prepare fresh for the preparation of the ferrous-reductant solution. Add 220 g of solid sulfamic acid to 450 mL of water, stir, and heat at 70–80°C until dissolved. Continue stirring and heating, while adding approximately 0.5-g portions of iron metal powder (Fe^0) until 56 g of iron have been dissolved in the heated sulfamic acid. Filter the solution while hot; allow to cool; then dilute to a final volume of 50 mL.

NOTE 5—The dissolution of the sulfamic acid need not be quantitative before beginning the addition of the iron powder. Excessive heating beyond the time required to achieve the dissolution of the sulfamic acid/iron powder or at temperatures above 80°C will cause excessive decomposition of the sulfamic acid

8.6 *Hydrochloric Acid* (HCl, 12 M)—Concentrated, sp gr 1.19.

8.7 *Hydrochloric Acid* (HCl, 1.0 M)—Add 84 mL of hydrochloric acid (sp gr 1.19) to approximately 500 mL of water. Stir, then dilute to a final volume of 1 L.

8.8 *Nitric Acid* (HNO_3 , 15.8 M)—Concentrated, sp gr 1.42.

8.9 *Nitric Acid* (1.0 M)—Add 63 mL of nitric acid (sp gr 1.42) to approximately 500 mL of water. Stir, then dilute to a final volume of 1 L.

8.10 *Nitric Acid* (3.0 M)—Add 190 mL of nitric acid (sp gr 1.42) to approximately 500 mL of water. Stir, then dilute to a final volume of 1 L.

8.11 *Plutonium Standard Solutions*—Prepare standards traceable to a national measurement system, which cover the range of concentrations over which sample measurements will be performed.

8.12 *Sulfamic Acid* (NH_2SO_3H , 1.5 M)—Dissolve 145 g of solid sulfamic acid in 900 mL of water with stirring. Filter, then dilute with water to a final volume of 1 L.

8.13 *Zirconium Reagent* ($ZrOCl_2 \cdot 8H_2O$, 0.75 M)—Dissolve 120.5 g zirconium chloride octahydrate in 450 mL of 1.0 M HCl; dilute to a final volume of 500 mL of 1.0 M HCl.

9. Calibration and Standardization of Instrument

9.1 Calibrate the system prior to each use. To calibrate, prepare several aliquants of at least two different plutonium standard solutions in the same concentration range as the samples to be measured in accordance with the preparation procedure described in Section 10. At least one of the standard solutions prepared will be used independently to ensure the accuracy of the calibration and should not be used in generating the calibration curve.

9.2 Following spectral referencing, measure each of the aliquants from one or more of the standard solutions. Quantitate each of the resulting spectra using appropriate software fitting routines; establish the calibration curve; and test the