
**Surface active agents — Microbiology
— Microbiological test methods for
liquid hand dishwashing**

*Agents de surface — Microbiologie — Méthodes d'essai
microbiologique pour les détergents liquides de lavage de vaisselle
à la main*

iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 21703:2019

<https://standards.iteh.ai/catalog/standards/iso/40930c74-850c-4bc3-89df-88e21febc947/iso-21703-2019>



iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 21703:2019

<https://standards.iteh.ai/catalog/standards/iso/40930c74-850c-4bc3-89df-88e21febc947/iso-21703-2019>



COPYRIGHT PROTECTED DOCUMENT

© ISO 2019

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
CP 401 • Ch. de Blandonnet 8
CH-1214 Vernier, Geneva
Phone: +41 22 749 01 11
Fax: +41 22 749 09 47
Email: copyright@iso.org
Website: www.iso.org

Published in Switzerland

Contents

Page

Foreword	v
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Principle	2
5 Diluents and culture media	2
5.1 General	2
5.2 Neutralizing diluent and diluents	2
5.2.1 General	2
5.2.2 Neutralizing diluent: Fluid casein digest — soy lecithin — polysorbate 20 medium (SCDLP 20 broth)	3
5.2.3 Diluent for the bacterial suspension (Tryptone sodium chloride solution)	3
5.3 Culture medium for counting	3
5.3.1 General	3
5.3.2 Soybean-casein digest agar medium (SCDA) or tryptic soy agar (TSA)	3
5.3.3 Other medium for counting	4
5.4 Culture medium for detection	4
5.4.1 General	4
5.4.2 Enrichment broth: Eugon LT 100 broth	5
5.5 Culture media for isolation and identification	5
5.5.1 MacConkey agar medium	5
5.5.2 Levine eosin-methylene blue agar medium (EMB agar medium)	6
5.5.3 Cetrimide agar medium	6
5.5.4 Pseudomonas agar medium for detection of pyocyanin (Pseudomonas agar P)	7
6 Apparatus and glassware	7
7 Strains of microorganisms	7
8 Handling of the product and laboratory samples	8
9 Procedure	8
9.1 General recommendation	8
9.2 Preparation of the initial suspension	8
9.3 Enumeration methods	8
9.3.1 Dilutions for enumeration	8
9.3.2 Plate-count methods	8
10 Method for detection and identification	11
10.1 Enrichment	11
10.2 Isolation	11
10.3 Detection of aerobic mesophilic bacteria	11
10.4 Detection and Identification of <i>Escherichia coli</i>	11
10.4.1 General	11
10.4.2 Gram's stain	11
10.4.3 Culture on Levine eosin-methylene blue agar medium (EMB agar medium)	11
10.4.4 Expression of results	12
10.5 Detection and identification of <i>Pseudomonas aeruginosa</i>	12
10.5.1 General	12
10.5.2 Gram's stain	12
10.5.3 Oxidase test	12
10.5.4 Culture on Pseudomonas agar medium for detection of pyocyanin	12
10.5.5 Expression of results	12
11 Neutralization of the antimicrobial properties of the product	13
11.1 General	13

11.2	Preparation of inoculum.....	13
11.3	Suitability of counting methods.....	13
11.3.1	Principle.....	13
11.3.2	Suitability test of the pour-plate method	13
11.3.3	Suitability of the surface spread method.....	13
11.3.4	Suitability of the membrane filtration method.....	14
11.4	Suitability of the detection method.....	14
11.4.1	Procedure.....	14
11.4.2	Interpretation.....	14
12	Test report.....	15
Annex A (informative)	Neutralizers of antimicrobial activity of preservatives and rinsing liquids.....	16
Annex B (informative)	Other culture media.....	17
Bibliography		19

iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 21703:2019

<https://standards.iteh.ai/catalog/standards/iso/40930c74-850c-4bc3-89df-88e21febc947/iso-21703-2019>