
**Animal and vegetable fats and oils —
Determination of polycyclic aromatic
hydrocarbons by on-line donor-acceptor
complex chromatography and HPLC with
fluorescence detection**

*Corps gras d'origines animale et végétale — Détermination de la teneur
en hydrocarbures aromatiques polycycliques par chromatographie de
complexe donneur-accepteur et CLHP avec détection par fluorescence*

Sample Document

get full document from standards.iteh.ai



PDF disclaimer

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.

Sample Document

get full document from standards.iteh.ai



COPYRIGHT PROTECTED DOCUMENT

© ISO 2009

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

Contents

Page

Foreword.....	iv
Introduction	v
1 Scope	1
2 Normative references	1
3 Terms and definitions.....	1
4 Principle	2
5 Reagents, materials and standards	3
6 Apparatus	5
7 Sampling and preparation of the test sample.....	5
8 Sample preparation	6
9 Procedure	7
9.1 HPLC analysis	7
9.2 Identification of PAHs.....	7
10 Calculation of individual PAHs.....	7
11 Method validation data	8
11.1 Accuracy.....	8
11.2 Within-laboratory precision	8
11.3 Recovery.....	8
11.4 Dynamic range	8
11.5 Limit of quantification	8
12 Precision.....	9
12.1 International collaborative trial	9
12.2 Repeatability.....	9
12.3 Reproducibility.....	9
13 Test report	9
Annex A (informative) Example of the individual parts of an HPLC system.....	10
Annex B (informative) Example of the operating conditions of the pumps of the HPLC system.....	11
Annex C (informative) Example of the operating conditions of the column thermostat and the detector of the HPLC system	12
Annex D (informative) Example of the autosampler programme of the HPLC system	13
Annex E (informative) Tubing connections of the HPLC system	14
Annex F (informative) Chromatogram of a standard calibration sample	15
Annex G (informative) Determination precision for a sunflower oil, range 0,1 µg/kg to 3,5 µg/kg	16
Annex H (informative) Determination precision for an olive oil, range 0,1 µg/kg to 3,5 µg/kg.....	17
Annex I (informative) Determination precision for a soybean oil, range 0,1 µg/kg to 3,5 µg/kg	18
Annex J (informative) Determination precision for a coconut oil, range 0,1 µg/kg to 3,5 µg/kg.....	19
Annex K (informative) Interlaboratory test	20
Bibliography	22

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 22959 was prepared by Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 11, *Animal and vegetable fats and oils*.

Sample Document

get full document from standards.iteh.ai

Introduction

Polycyclic aromatic hydrocarbons (PAHs) are formed during pyrolytic processes such as the incomplete combustion of organic substances or have a petrogenic origin (mineral oils). Edible fats and oils may be contaminated by environmental pollution and/or processing steps prior to refining. The presence of PAHs in fats and oils is a health concern due to their carcinogenicity. Different levels of PAHs have been observed in crude edible oils. Refining of the oils (deodorization, bleaching, charcoal treatment) under the appropriate conditions reduces the content of the individual PAHs to the microgram per kilogram level. The known methods of analysis of PAHs in edible fats and oils include complex and laborious extraction and clean-up procedures to isolate the low levels of PAHs present.

With the donor-acceptor complex-chromatography (DACC) technique, PAHs can be extracted from different matrices. PAHs are electron donors (π -electrons) and the strong interaction of the PAHs with an electron acceptor stationary phase results in retention of the PAHs and elution of (the bulk of) the other components of the oil. This International Standard specifies an automated on-line method for the determination of PAHs in edible oils and fats, which can easily be applied as a routine analysis. The method consists of an LC-LC coupling of a clean-up DACC column to an analytical column for the separation. PAHs are quantified by fluorescence detection.

Compared to older techniques, this automated on-line method significantly reduces the amount of solvent used and saves considerable time. The DACC column clean-up is fast and is carried out during the HPLC run of the previous sample. The total analysis time for one sample is approximately 90 min, compared with the traditional methods which require 8 h to 10 h. Moreover, the system can run 24 h/day. Finally, losses of volatile PAHs during solvent evaporation, for example, are eliminated. The quantification limits of 0,1 $\mu\text{g/kg}$ of the individual PAHs have been retained with the DACC method, which automatically corrects for possibly incomplete recoveries because the calibration samples are subjected to the same treatment as the samples to be analysed. The system uses conventional HPLC instrumentation.

Sample Document

get full document from standards.iteh.ai

Animal and vegetable fats and oils — Determination of polycyclic aromatic hydrocarbons by on-line donor-acceptor complex chromatography and HPLC with fluorescence detection

1 Scope

This International Standard specifies a high performance liquid chromatographic (HPLC) procedure for the determination of polycyclic aromatic hydrocarbons (PAHs) in edible fats and oils.

The method has been validated for coconut (CN), olive (OV), sunflower (SF), and soybean (BO) oil, and is possibly applicable to other oils, dependent on the determination of appropriate parameters.

The lowest level of quantification for the PAHs is 0,1 µg/kg. The lowest possible amount of each PAH which can be distinguished from the baseline noise has not been determined. The validated concentration range of the method is 0,1 µg/kg to 3,5 µg/kg for each individual PAH. For samples containing (light) PAH contents > 3,5 µg/kg, dilution to bring the contents into the validated range is possible. It is also possible to adjust the range of the calibration curves. However, ranges exceeding 3,5 µg/kg have not been validated.

PAHs which can be determined by this method are: anthracene, phenanthrene, fluoranthene, pyrene, chrysene, benzo[a]anthracene, benzo[e]pyrene, benzo[a]pyrene, perylene, benzo[ghi]perylene, anthanthrene, dibenzo[a,h]anthracene, coronene, indeno[1,2,3-cd]pyrene, benzo[a]fluoranthene, benzo[b]fluoranthene, benzo[k]fluoranthene.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 661, *Animal and vegetable fats and oils — Preparation of test sample*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

3.1

polycyclic aromatic hydrocarbon

PAH

compound that contains two or more condensed (fused) aromatic hydrocarbon rings and whose content can be determined according to the method specified in this International Standard

3.2

light polycyclic aromatic hydrocarbon

compound with two to four condensed (fused) aromatic hydrocarbon rings

EXAMPLES

Compound	CAS No.	Compound	CAS No.	Compound	CAS No.
acenaphthene	83-32-9	benzo[e]pyrene	192-97-2	naphthalene	50-32-8
acenaphthylene	208-96-8	chrysene	218-01-9	phenanthrene	85-01-8
anthracene	120-12-7	fluoranthene	206-44-0	pyrene	129-00-0
benzo[a]anthracene (1,2-benzanthracene)	56-55-3	fluorene	86-73-7		

3.3

heavy polycyclic aromatic hydrocarbon

compound with five and more condensed (fused) aromatic hydrocarbon rings

EXAMPLES

Compound	CAS No.	Compound	CAS No.	Compound	CAS No.
benzo[a]pyrene (1,2-benzopyrene)	50-32-8	benzo[k]fluoranthene	207-08-9	dibenzo[a,h]anthracene (1,2,5,6-dibenzoanthracene)	53-70-3
benzo[a]fluoranthene	203-33-8	benzo[ghi]perylene (1,12-benzoperylene)	191-24-2	indeno[1,2,3-cd]pyrene	193-39-5
benzo[b]fluoranthene	205-99-2	coronene	191-07-1	perylene	198-55-0

3.4

polycyclic aromatic hydrocarbon content

mass fraction of a polycyclic aromatic hydrocarbon or polycyclic aromatic hydrocarbon mixture in a matrix

EXAMPLES Individual polycyclic aromatic hydrocarbon content; light polycyclic aromatic hydrocarbon content; heavy polycyclic aromatic hydrocarbon content.

NOTE The content is expressed as a mass fraction in micrograms per kilogram.

4 Principle

The PAHs in edible oils are determined by on-line coupling of donor-acceptor complex chromatography (DACC) and HPLC with fluorescence detection. The oil samples are eluted over a column with a modified stationary phase (DACC column) which will act as an electron acceptor. This column will retain the PAHs (electron donors) by π - π interactions. After elution of the oil, the PAHs are transferred on-line to the analytical reversed phase column. The individual PAHs are detected at different wavelengths. The retention times of the PAHs are used to identify the individual compounds. The levels of the PAHs in the oil samples are calculated by external calibration.

5 Reagents, materials and standards

WARNING — The method requires harmful reagents. Respect normal laboratory safety regulations. All PAHs are suspected carcinogenic compounds. Therefore, it is essential that the preparation of the stock solutions, the standard dilutions and the samples of the calibration curve (5.3) are performed by preference in a class-2 laboratory. Furthermore, a laboratory coat and safety gloves are essential for these steps. Contaminated tissues and gloves shall be collected in a plastic bag and removed after sealing the bag.

5.1 Reagents.

5.1.1 Acetonitrile, HPLC grade, mass fraction $w[\text{C}_2\text{H}_3\text{N}] > 99,9 \%$.

5.1.2 Ethyl acetate, HPLC grade, mass fraction $w[\text{C}_4\text{H}_8\text{O}_2] > 99,8 \%$.

5.1.3 2-Propanol, HPLC grade, mass fraction $w[\text{C}_3\text{H}_8\text{O}] > 99,9 \%$.

5.1.4 Toluene, HPLC grade, mass fraction $w[\text{C}_7\text{H}_8] > 99,9 \%$.

5.1.5 Water, HPLC grade.

5.2 Standards.¹⁾

5.2.1 Anthracene, mass fraction $w[\text{C}_{14}\text{H}_{10}] > 99 \%$.

5.2.2 Phenanthrene, mass fraction $w[\text{C}_{14}\text{H}_{10}] > 99 \%$.

5.2.3 Fluoranthene, mass fraction $w[\text{C}_{16}\text{H}_{10}] > 99 \%$.

5.2.4 Pyrene, mass fraction $w[\text{C}_{16}\text{H}_{10}] > 99 \%$.

5.2.5 Chrysene, mass fraction $w[\text{C}_{18}\text{H}_{12}] > 99 \%$.

5.2.6 Benzo[a]anthracene (1,2-Benzoanthracene), mass fraction $w[\text{C}_{18}\text{H}_{12}] > 99 \%$.

5.2.7 Benzo[e]pyrene, mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.8 Benzo[a]pyrene (1,2-Benzopyrene), mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.9 Perylene, mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.10 Benzo[ghi]perylene (1,12-Benzoperylene), mass fraction $w[\text{C}_{22}\text{H}_{12}] > 99 \%$.

5.2.11 Anthanthrene, mass fraction $w[\text{C}_{22}\text{H}_{12}] > 99 \%$.

5.2.12 Dibenzo[a,h]anthracene (1,2,5,6-Dibenzoanthracene), mass fraction $w[\text{C}_{22}\text{H}_{14}] > 99 \%$.

5.2.13 Coronene, mass fraction $w[\text{C}_{24}\text{H}_{12}] > 99 \%$.

5.2.14 Indeno[1,2,3-cd]pyrene, mass fraction $w[\text{C}_{22}\text{H}_{12}] > 99 \%$.

5.2.15 Benzo[a]fluoranthene, mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.16 Benzo[b]fluoranthene, mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.17 Benzo[k]fluoranthene, mass fraction $w[\text{C}_{20}\text{H}_{12}] > 99 \%$.

5.2.18 BCR certified reference material 458, coconut oil with 6 PAHs.

1) IRMM (<http://www.irmm.jrc.be>) and Sigma-Aldrich (<http://www.sigmaaldrich.com>) are suitable suppliers. This information is given for the convenience of users of this International Standard and does not constitute an endorsement by ISO of products so supplied. Products from other sources may be used if they can be shown to lead to the same results.