

---

---

**Milk-based infant formula powders —  
Quantification of whey protein content  
by sodium dodecyl sulfate-capillary  
gel electrophoresis (SDS-CGE)**

*Formules infantiles en poudre à base de lait — Quantification de la  
teneur en protéine de lactosérum par électrophorèse capillaire sur gel  
contenant du dodécylsulfate de sodium (SDS-CGE)*

*iteh Standards*  
(<https://standards.iteh.ai>)  
Document Preview

ISO 23293:2020

<https://standards.iteh.ai/catalog/standards/iso/112485af-b68b-49ad-99d7-f3ff48939cfb/iso-23293-2020>



Reference numbers  
ISO 23293:2020(E)  
IDF 247:2020(E)

© ISO and IDF 2020

# iTeh Standards (<https://standards.iteh.ai>) Document Preview

ISO 23293:2020

<https://standards.iteh.ai/catalog/standards/iso/112485af-b68b-49ad-99d7-f3ff48939cfb/iso-23293-2020>



## **COPYRIGHT PROTECTED DOCUMENT**

© ISO and IDF 2020

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office  
CP 401 • Ch. de Blandonnet 8  
CH-1214 Vernier, Geneva  
Phone: +41 22 749 01 11  
Fax: +41 22 749 09 47  
Email: [copyright@iso.org](mailto:copyright@iso.org)  
Website: [www.iso.org](http://www.iso.org)

Published in Switzerland

International Dairy Federation  
Silver Building • Bd Auguste Reyers 70/B  
B-1030 Brussels  
Phone: +32 2 325 67 40  
Fax: +32 2 325 67 41  
Email: [info@fil-idf.org](mailto:info@fil-idf.org)  
Website: [www.fil-idf.org](http://www.fil-idf.org)

# Contents

Page

|                                                                  |           |
|------------------------------------------------------------------|-----------|
| <b>Foreword</b>                                                  | <b>iv</b> |
| <b>1 Scope</b>                                                   | <b>1</b>  |
| <b>2 Normative references</b>                                    | <b>1</b>  |
| <b>3 Terms and definitions</b>                                   | <b>1</b>  |
| <b>4 Principle</b>                                               | <b>1</b>  |
| <b>5 Reagents and materials</b>                                  | <b>1</b>  |
| <b>6 Apparatus</b>                                               | <b>2</b>  |
| <b>7 Procedure</b>                                               | <b>3</b>  |
| 7.1 Sample preparation                                           | 3         |
| 7.2 CGE analysis                                                 | 3         |
| 7.3 Electropherogram processing                                  | 4         |
| 7.3.1 Electropherogram integration                               | 4         |
| 7.3.2 Identification of the casein region                        | 4         |
| <b>8 System suitability</b>                                      | <b>4</b>  |
| <b>9 Calculations and expression of results</b>                  | <b>5</b>  |
| 9.1 Calculation of whey protein content                          | 5         |
| 9.2 Expression of results                                        | 5         |
| <b>10 Precision</b>                                              | <b>5</b>  |
| 10.1 Interlaboratory tests                                       | 5         |
| 10.2 Infant formulas with a whey protein content of 40 % or more | 6         |
| 10.2.1 Repeatability                                             | 6         |
| 10.2.2 Reproducibility                                           | 6         |
| 10.3 Infant formulas with a whey protein content below 40 %      | 6         |
| 10.3.1 Repeatability                                             | 6         |
| 10.3.2 Reproducibility                                           | 6         |
| <b>11 Test report</b>                                            | <b>6</b>  |
| <b>Annex A (informative) Example of a CGE analysis</b>           | <b>7</b>  |
| <b>Annex B (informative) Precision data</b>                      | <b>10</b> |
| <b>Bibliography</b>                                              | <b>11</b> |