
**Cigarettes — Determination of
ammonia in cigarette mainstream
smoke with an intense smoking
regime using ion chromatography**

*Cigarettes — Dosage de l'ammoniac dans le courant principal
de la fumée de cigarette avec un régime de fumage intense par
chromatographie par échange d'ions*

iteh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 23920:2020

<https://standards.iteh.ai/catalog/standards/iso/017c1158-c5c6-440f-a55c-b18c73c1d92d/iso-23920-2020>



iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 23920:2020

<https://standards.iteh.ai/catalog/standards/iso/017c1158-c5c6-440f-a55c-b18c73c1d92d/iso-23920-2020>



COPYRIGHT PROTECTED DOCUMENT

© ISO 2020

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
CP 401 • Ch. de Blandonnet 8
CH-1214 Vernier, Geneva
Phone: +41 22 749 01 11
Email: copyright@iso.org
Website: www.iso.org

Published in Switzerland

Contents

Page

Foreword	iv
Introduction	v
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Principle	1
5 Apparatus	2
6 Reagents	2
7 Preparation	3
7.1 General	3
7.2 Preparation of solutions	3
7.2.1 Sulphuric acid, 0,01 M — Impinger solution (trapping system 1)	3
7.2.2 Sulphuric acid, 0,1 M — Solution C (Ion Chromatography Eluent)	3
7.2.3 MSA 0,003 M — Solution A (Ion chromatography eluent)	3
7.2.4 Sulphuric acid, 0,01 M — Ammonium standards preparation	3
7.2.5 Hydrochloric acid, 0,05 M — Solution for GFP impregnation (trapping system 2)	3
7.2.6 Hydrochloric acid, 0,01 M — Extraction solution (trapping system 2)	3
7.3 Preparation of standards	3
7.3.1 Primary ammonium stock solution	3
7.3.2 Calibration standards	4
8 Sampling	4
9 Tobacco product preparation	4
10 Sample generation — Smoking of cigarettes	4
10.1 General	4
10.2 Smoking machine setup	4
10.2.1 General	4
10.2.2 Trapping system 1 — GFP and impinger traps	4
10.2.3 Trapping system 2 — Impregnated GFP	5
10.3 Smoking	6
11 Sample analysis	6
11.1 Preparation of sample	6
11.1.1 Extraction of GFP — Trapping system 1	6
11.1.2 Extraction of GFP — Trapping system 2	6
11.2 Determination	6
11.2.1 Ion chromatography operating conditions	6
11.2.2 Calibration	7
11.2.3 Calculation	8
12 Repeatability and reproducibility	8
12.1 General	8
12.2 Results of the 2015 collaborative study	9
13 Report	9
Annex A (informative) Examples of chromatograms	10
Bibliography	11