
**Sterile hypodermic syringes for
single use —**

**Part 3:
Auto-disabled syringes for fixed-dose
immunization**

Seringues hypodermiques stériles, non réutilisables —

Partie 3: Seringues autobloquantes pour vaccination à dose fixe

iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 7886-3:2020

<https://standards.iteh.ai/catalog/standards/iso/1d7c598a-4ffe-409b-99a2-a45cbe7a87b1/iso-7886-3-2020>



iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 7886-3:2020

<https://standards.iteh.ai/catalog/standards/iso/1d7c598a-4ffe-409b-99a2-a45cbe7a87b1/iso-7886-3-2020>



COPYRIGHT PROTECTED DOCUMENT

© ISO 2020

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
CP 401 • Ch. de Blandonnet 8
CH-1214 Vernier, Geneva
Phone: +41 22 749 01 11
Fax: +41 22 749 09 47
Email: copyright@iso.org
Website: www.iso.org

Published in Switzerland

Contents

Page

Foreword	v
Introduction	vi
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Nomenclature	2
5 General requirements	3
6 Extraneous matter	4
6.1 General	4
6.2 Limits for acidity or alkalinity	4
6.3 Limits for extractable metals	4
7 Lubricant	4
8 Tolerance on nominal capacity	4
9 Graduated scale	5
9.1 Scale	5
9.2 Position of scale	5
10 Barrel	5
10.1 Dimensions	5
10.2 Barrel flanges	5
11 Plunger stopper/plunger assembly	5
11.1 Design	5
11.2 Fit of the plunger stopper/plunger in the barrel	6
11.3 Fiducial line	6
12 Needle	6
12.1 General	6
12.2 Integrated needle	6
12.3 Non-integrated needle	6
12.4 Sharps protection features	6
13 Performance	7
13.1 General	7
13.2 Dead space	7
13.3 Freedom from air and liquid leakage	7
13.4 Auto-disable syringe feature	7
13.5 Performance after shipping	7
14 Packaging	8
14.1 Unit packaging providing sterile barrier	8
14.2 Multiple unit pack	8
14.3 User packaging	8
15 Information supplied by the manufacturer	8
15.1 General	8
15.2 Syringes	8
15.3 Unit packaging providing sterile barrier	8
15.4 User packaging	9
15.5 Storage containers	9
15.6 Transport wrapping	10
Annex A (normative) Method for preparation of extracts	11
Annex B (informative) Test method for forces required to operate piston	12

Annex C (normative) Test method for testing auto-disable syringe feature	14
Bibliography	15

iTeh Standards
(<https://standards.itih.ai>)
Document Preview

[ISO 7886-3:2020](https://standards.itih.ai/catalog/standards/iso/1d7c598a-4ffe-409b-99a2-a45cbe7a87b1/iso-7886-3-2020)

<https://standards.itih.ai/catalog/standards/iso/1d7c598a-4ffe-409b-99a2-a45cbe7a87b1/iso-7886-3-2020>