
Dentistry — Cartridge syringes

Médecine bucco-dentaire — Seringues pour cartouches

iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

[ISO 9997:2020](https://standards.iteh.ai/catalog/standards/iso/7c9ea256-3463-4246-adb9-322b9ad881cf/iso-9997-2020)

<https://standards.iteh.ai/catalog/standards/iso/7c9ea256-3463-4246-adb9-322b9ad881cf/iso-9997-2020>



iTeh Standards
(<https://standards.iteh.ai>)
Document Preview

ISO 9997:2020

<https://standards.iteh.ai/catalog/standards/iso/7c9ea256-3463-4246-adb9-322b9ad881cf/iso-9997-2020>



COPYRIGHT PROTECTED DOCUMENT

© ISO 2020

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
CP 401 • Ch. de Blandonnet 8
CH-1214 Vernier, Geneva
Phone: +41 22 749 01 11
Fax: +41 22 749 09 47
Email: copyright@iso.org
Website: www.iso.org

Published in Switzerland

Contents

Page

Foreword	iv
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Classification	2
5 Requirements	2
5.1 General.....	2
5.2 Materials.....	2
5.3 Dimensions.....	2
5.4 Special requirements.....	3
5.4.1 Loading and cartridge size.....	3
5.4.2 Viewing of contents.....	4
5.4.3 Plunger rod.....	4
5.4.4 Aspirating syringes.....	4
5.5 Resistance to reprocessing.....	5
6 Measurement and test methods	5
6.1 Visual inspection.....	5
6.2 Dimensions.....	5
6.3 Secure loading.....	5
6.4 Aspirating test for cartridge syringes of Type 2.....	5
6.4.1 Reagent.....	5
6.4.2 Procedure.....	5
6.4.3 Observation.....	5
6.5 Aspirating test for cartridge syringes of Type 3.....	5
6.5.1 Reagent.....	5
6.5.2 Procedure.....	5
6.5.3 Observation.....	6
6.6 Plunger rod.....	6
6.6.1 Plunger rod movement.....	6
6.6.2 Plunger rod displacement.....	6
6.7 Resistance to reprocessing.....	6
7 Instructions for use	6
8 Marking	7
8.1 Labelling on the unit pack.....	7
8.2 Marking of syringe.....	7
Annex A (informative) Imperial thread sizes	8
Bibliography	9